
→ SCIENTIFIC ABSTRACTS

18 years of clinical research

## Dear Doctor,

The present collection of scientific abstracts gives you an overview of the existing scientific evidence for mectron PIEZOSURGERY®. Sixteen years ago, Prof. Dr. Tomaso Vercellotti came up with the idea of piezoelectric bone surgery. From the beginning of this revolutionary surgical technique, many universities and specialized surgeons have been involved in the clinical and scientific validation of the piezoelectric bone surgery. Today, a continuously growing number of publications confirm the clinical benefits of mectron PIEZOSURGERY® :

- micrometric cuts
- selective cuts
- maximum intra-operative visibility

With the scientific proof statement of these benefits, mectron PIEZOSURGERY® has defined a new dimension in bone surgery and today is considered worldwide the only evidence based alternative to traditional surgical tools.

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## ⇒ TABLE OF CONTENT

<b>TECHNIQUE ADVANTAGES</b>	<b>Pereira CC, Gealh WC, Nogueira LM, Garcia Junior IR, Okamoto R.</b> Piezosurgery applied to implant dentistry: clinical and biological aspects. J Oral Implantol. 2014 Jul;40 Spec No:401-8.	18
	<b>Zhao J, Huang C.</b> The advanced techniques of dentoalveolar surgery. CHua Xi Kou Qiang Yi Xue Za Zhi. 2014 Jun;32(3):213-6. (chinese)	18
	<b>Bauer SE, Romanos GE.</b> Morphological characteristics of osteotomies using different piezosurgical devices. A scanning electron microscopic evaluation. Implant Dent. 2014 Jun;23(3):334-42.	19
	<b>Patel A.</b> The role of piezosurgery in implant dentistry: Adam Patel details the huge potential of piezosurgery for improving the predictability and ease of dental implant treatment. IDT May 2014	20
	<b>Franco S, Miccoli S, Limongelli L, Tempesta A, Favia G, Maiorano E, Favia G.</b> New Dimensional Staging of Bisphosphonate-Related Osteonecrosis of the Jaw Allowing a Guided Surgical Treatment Protocol: Long-Term Follow-Up of 266 Lesions in Neoplastic and Osteoporotic Patients from the University of Bari. Int J Dent. 2014 Jun;2014:935657.	20
	<b>Gülnahar Y, Hüseyin Kösgen H, Tutar Y.</b> A comparison of piezosurgery and conventional surgery by heat shock protein 70 expression. Int J Oral Maxillofac Surg. 2013 Apr;42(4):508-10.	21
	<b>Claire S, Lea SC, Walmsley AD.</b> Characterisation of bone following ultrasonic cutting. Clin Oral Investig. 2013 Apr;17(3):905-12.	22
	<b>Simonetti M, Facco G, Barberis F, Signorini G, Capurro M, Rebaudi A, Sammartino G.</b> Bone characteristics following osteotomy surgery: an in vitro SEM study comparing traditional Lindemann drill with sonic and ultrasonic instruments. POSEIDO 2013;1(3):187.	22
	<b>Schütz S, Egger J, Kühl S, Filippi A, Lambrecht JT.</b> Intraosseous temperature changes during the use of piezosurgical inserts in vitro. Int J Oral Maxillofac Surg. 2012 Nov;41(11):1338-43.	23
	<b>Itiro A, Lupo G, Carotenuto A, Filipi M, Cocozza E, Marra A.</b> I vantaggi della chirurgia piezoelettrica nella chirurgia orale e nella chirurgia maxillo-facciale. Revisione della letteratura. Minerva Stomatol. 2012 May;61(5):213-24.	24
	<b>Parmar D, Mann M, Walmsley AD, Lea SC.</b> Cutting characteristics of ultrasonic surgical instruments. Clin Oral Implants Res. 2011 Dec;22(12):1385-90.	24

→ TABLE OF CONTENT

<b>Pavlíková G, Foltán R, Burian M, Horká E, Adámek S, Hejčl A, Hanzelka T, Sedý J.</b> Piezosurgery prevents brain tissue damage: an experimental study on a new rat model. Int J Oral Maxillofac Surg. 2011 Aug;40(8):840-4.	25
<b>Pavlíková G, Foltán R, Horká M, Hanzelka T, Borunská H, Sedý J.</b> Piezosurgery in oral and maxillofacial surgery. Int J Oral Maxillofac Surg. 2011 May;40(5):451-7.	26
<b>Maurer P, Kriwalsky MS, Block Veras R, Vogel J, Syrowatka F, Heiss C.</b> Micromorphometrical analysis of conventional osteotomy techniques and ultrasonic osteotomy at the rabbit skull. Clin Oral Implants Res. 2008 Jun;19(6):570-5.	26
<b>Labanca M, Azzola F, Vinci R, Rodella LF</b> Piezoelectric surgery: Twenty years of use. Br J Oral Maxillofac Surg. 2008; 46(4):265-269.	27
<b>Stübinger S, Landes C, Seitz O, Zeilhofer HF, Sader R.</b> Ultraschallbasiertes Knochenschneiden in der Oralchirurgie: eine Übersicht anhand von 60 Patientenfällen. Ultraschall Med. 2008; 29(1):66-71.	27
<b>Peivandi A, Bugnet R, Debize E.</b> Méthode piézoélectrique d'ostéotomie: une nouvelle technique d'aménagement osseux en chirurgie implantaire. Revue Implantologie; 2007 Nov;15-23.	29
<b>Maurer P, Kriwalsky MS, Block Veras R, Brandt J, Heiss C.</b> Auflichtmikroskopische Untersuchungen an der Kaninchenkalotte nach ultraschallgestützter und konventioneller Osteotomie. Biomed Tech (Berl). 2007;52(5):351-5.	29
<b>Blakenburg JJ, Both CJ, Borstlap WA, van Damme PA.</b> Geluidsniveau van de PIEZOSURGERY®: Risico van permanente gehoorschade. Ned Tijdschr Tandheelkd. 2007 Nov;114(11):451-4.	31
<b>Peivandi A, Bugnet R, Debize E, Gleizal A, Dohan DM.</b> Piezoelectric osteotomy: applications in periodontal and implant surgery. Rev Stomatol Chir Maxillofac. 2007; 108(5):431-440.	31
<b>Su YC.</b> Development and clinical application of ultrasonic osteotomy in dentistry. Shanghai kou qiang yi xue = Shanghai journal of stomatology. 2007; 16(1):1-7 [chinese]	32
<b>Sembronio S, Albiero AM, Polini F, Robiony M, Politi M.</b> Intraoral endoscopically assisted treatment of temporomandibular joint ankylosis: Preliminary report. Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2007; 104(1):e7-e10.	33
<b>Ramaglia L, Saviano R, Espedito di Lauro A, Capece G.</b> La guarigione dei tessuti peri-implantari in impianti posizionati in alveoli post-estrattivi di premolari mascellari. Minerva Stomatol. 2006 Apr;55(4):199-207.	33

⇒ TABLE OF CONTENT

<b>Hyvernat P.</b> Le Bistouri Piezoélectrique. Un Outil d'une finesse vraiment adaptée à l'Odontostomatologie. Revue Implantologie, 2006 Feb;5-18.	34
<b>Schlee M, Steigmann M, Bratu E, Garg AK.</b> Piezosurgery: basics and possibilities. Implant Dent. 2006; 15(4):334-340.	35
<b>Stübinger S, Filippi A, Sader R, Zeilhofer HF.</b> Intraoral Piezosurgery: Preliminary Results of a New Technique. J Oral Maxillofac Surg. 2005; 63(9):1283-1287.	35
<b>Schlee M.</b> Ultraschallgestützte Chirurgie - Grundlagen und Möglichkeiten. Z Zahnärztl Impl. 2005; 21(1):48-59.	36
<b>Vercellotti T.</b> Características tecnológicas e indicaciones clínicas de la cirugía ósea piezoelectrica. Revista Mundo Dental. 2005; 26-28.	37
<b>Troiani C, Russo C, Ballarani G, Vercellotti T.</b> Piezoelectric Surgery: A new reality to cut and manage bone. Maxillo Odontostomatologia- International Journal of Maxillo Odontostomatology - S.I.M.O. 2005; 4(1):23-28.	38
<b>Leclercq P, Dohan D</b> De l'intérêt du bistouri ultrasonore en implantologie: technologies, applications cliniques: 2ème partie: applications cliniques. Implantodontie, 2004 Jul-Sep;13(3): 159-165.	38
<b>Leclercq P, Dohan D.</b> De l'intérêt du bistouri ultrasonore en implantologie: technologies, applications cliniques: 1ère partie: technologies. Implantodontie, 2004 Jul-Sep;13(3):151-157.	39
<b>Siervo S, Ruggli-Milic S, Radici M, Siervo P, Jäger K.</b> Piezoelektrische Chirurgie. Eine alternative Methode für die schonende Chirurgie. Schweiz Monatsschr Zahnmed. 2004; 114(4):365-377.	40
<b>Lambrecht JT.</b> Intraorale Piezo-Chirurgie. Schweiz Monatsschr Zahnmed. 2004; 114(1):29-36.	40
<b>Vercellotti T.</b> Technological characteristics and clinical indications of piezoelectric bone surgery. Minerva Stomatol. 2004; 53(5):207-214.	41
<b>Boioli LT, Vercellotti T, Tecucianu JF.</b> La chirurgie piézoélectrique: Une alternative aux techniques classiques de chirurgie osseuse. Inf Dent. 2004; 86(41):2887-2893.	42
<b>Vercellotti T.</b> La Chirurgia Ossea Piezoelettrica. Il Dentista Moderno. 2003; 5:21-55.	44

→ TABLE OF CONTENT

<b>HISTOLOGICAL ADVANTAGES</b>	Sohn DS, Lee JK, An KM, Shin HI. Histomorphometric evaluation of mineralized cancellous allograft in the maxillary sinus augmentation: a 4 case report. Implant Dent. 2009 Apr;18(2):172-81.	44
	Stübinger S, Goethe JW. Bone Healing After PIEZOSURGERY® and its Influence on Clinical Applications. Journal of Oral and Maxillofacial Surgery 2007 Sep;65(9):39.e7-39.e8.	45
	Gleizal A, Li S, Pialat JB, Béziat JL. Transcriptional expression of calvarial bone after treatment with low – intensity ultrasound: An in vitro study. Ultrasound Med Biol. 2006; 32(10):1569-1574.	46
	Vercellotti T, Nevins ML, Kim DM, Nevins M, Wada K, Schenk RK, Fiorellini JP. Osseous Response following Resective Therapy with a Piezosurgery®. Int J Periodontics Restorative Dent. 2005; 25(6):543-549.	47
<b>PERIOSTEUM PREPARATION</b>	Stoetzer M, Felgenträger D, Kampmann A, Schumann P, Rücker M, Gellrich NC, von See C Effects of a new piezoelectric device on periosteal microcirculation after subperiosteal preparation. Microvasc Res. 2014 Jul;94:114-8.	48
	von See C, Gellrich NC, Rücker M, Kokemüller H, Kober H, Stöver E. Investigation of perfusion in osseous vessels in close vicinity to piezoelectric bone cutting. Br J Oral Maxillofac Surg. 2012 Apr;50(3):251-5.	48
<b>IMPLANT SITE PREPARATION</b>	Vercellotti T, Corrias D, Russo C, Carossa S, Schierano G. Studio preliminare della perdita ossea marginale in 156 impianti in titanio posizionati con preparazione a ultrasuoni del sito implantare mediante Piezosurgery® (UISP): follow-up di quattro anni. Minerva Stomatologica – 2014 Apr. vol. 63 suppl. 1 al n. 4 - 232.	49
	da Silva Neto UT, Joly JC, Gehrke SA. Clinical analysis of the stability of dental implants after preparation of the site by conventional drilling or piezosurgery. Br J Oral Maxillofac Surg. 2014 Feb;52(2):149-53.	50
	Gandhi SA, Baker JA, Bairam L, Kim HI, Davis EL, Andreana S. Primary stability comparison using piezoelectric or conventional implant site preparation systems in cancellous bone: a pilot study. Implant Dent. 2014 Feb;23(1):79-84.	51
	Vercellotti T, Stacchi C, Russo C, Rebaudi A, Vincenzi G, Pratella U, Baldi D, Mozzati M, Monaghedu C, Senterini R, Cuneo T, Di Alberti L, Carossa S, Schierano G. Ultrasonic implant site preparation using piezosurgery: a multicenter case series study analyzing 3,579 implants with a 1- to 3-year follow-up. Int J Periodontics Restorative Dent. 2014 Jan-Feb;34(1):11-8.	52
	Bengazi F, Lang NP, Canciani E, Viganò P, Velez JU, Botticelli D. Osseointegration of implants with dendrimers surface characteristics installed conventionally or with Piezosurgery®. A comparative study in the dog. Clin Oral Implants Res. 2014 Jan;25(1):10-5.	52

⇒ TABLE OF CONTENT

<b>Lamazza L, Laurito D, Lollobrigida M, Brugnoletti O, Garreffa G, De Biase A.</b> <b>Identification of possible factors influencing temperatures elevation during implant site preparation with piezoelectric technique.</b> <i>Annali di Stomatologia</i> 2014; V (4):115-122.	53
<b>Stacchi C, Chen ST, Raghoebar GM, Rosen D, Poggio CE, Ronda M, Bacchini M, Di Lenarda R.</b> <b>Malpositioned osseointegrated implants relocated with segmental osteotomies: a retrospective analysis of a multicenter case series with a 1- to 15-year follow-up.</b> <i>Clin Implant Dent Relat Res.</i> 2013 Dec;15(6):836-46.	54
<b>Stacchi C, Vercellotti T, Torelli L, Furlan F, Di Lenarda R.</b> <b>Changes in implant stability using different site preparation techniques: twist drills versus piezosurgery. A single-blinded, randomized, controlled clinical trial.</b> <i>Clin Implant Dent Relat Res.</i> 2013 Apr;15(2):188-97.	55
<b>Baker JA, Vora S, Bairam L, Kim HI, Davis EL, Andreana S.</b> <b>Piezoelectric vs. conventional implant site preparation: ex vivo implant primary stability.</b> <i>Clin Oral Implants Res.</i> 2012 Apr;23(4):433-7.	56
<b>Sentinelli R, Dagnino G.</b> <b>Riabilitazione totale di un mascellare inferiore con protocollo "All on four" mediante preparazione piezoelettrica del sito implantare.</b> <i>Implantologia</i> 2011, 3:21-28.	56
<b>Di Alberti L, Donnini F, Di Alberti C, Camerino M.</b> <b>A comparative study of bone densitometry during osseointegration: piezoelectric surgery versus rotary protocols.</b> <i>Quintessence Int.</i> 2010 Sep;41(8):639-44.	57
<b>Stacchi C, Costantinides F, Biasotto M, Di Lenarda R.</b> <b>Relocation of a malpositioned maxillary implant with piezoelectric osteotomies: a case report.</b> <i>Int J Periodontics Restorative Dent.</i> 2008 Oct;28(5):489-95.	58
<b>Prete G, Martinasso G, Peirone B, Navone R, Manzella C, Muzio G, Russo C, Canuto RA, Schierano G.</b> <b>Cytokines and Growth Factors Involved in the Osseointegration of Oral Titanium Implants Positioned using Piezoelectric Bone Surgery Versus a Drill Technique: A Pilot Study in Minipigs.</b> <i>J Periodontol.</i> 2007; 78(4):716-722.	58
<b>SINUS LIFT TECHNIQUE</b>	
<b>Tatullo M, Marrelli M, Cassetta M, Pacifici A, Stefanelli LV, Scacco S, Dipalma G, Pacifici I, Inchingolo F.</b> <b>Platelet Rich Fibrin (P.R.F.) in reconstructive surgery of atrophied maxillary bones: clinical and histological evaluations.</b> <i>Int J Med Sci.</i> 2012;9(10):872-80.	59
<b>Rodella LF, Labanca M, Boninsegna R, Favero G, Tschaibitscher M, Rezzani R.</b> <b>Anastomosi intraossea nel seno mascellare.</b> <i>Minerva Stomatol.</i> 2010 Jun;59(6):349-54.	60
<b>Tordjman S, Boioli LT, Fayd N.</b> <b>Apport de la Piézochirurgie dans la surélévation du plancher sinusien.</b> <i>Département de Parodontologie de l'UFR de Stomatologie et Chirurgie Maxillo-Faciale. Université de Paris VI - Paris. Revue Implantologie.</i> 2006 Nov: 17-25.	61

→ TABLE OF CONTENT

<b>SINUS LIFT TECHNIQUE BY CRESTAL APPROACH</b>	<b>Vercellotti T, Nevins M, Jensen Ole T.</b> <b>Piezoelectric Bone Surgery for Sinus Bone Grafting.</b> The Sinus Bone Graft, Second Edition. Edited by Ole T. Jensen, Quintessence Books. 2006; 23:273-279.	61
	<b>Vercellotti T.</b> <b>La Chirurgia Piezoelettrica. Tecniche di rialzo del seno mascellare.</b> La chirurgia del seno mascellare e le alternative terapeutiche, Testori T, Weinstein R, Wallace S, Edizioni ACME. 2005; 14:245-255.	62
	<b>Kühl S, Kirmeier R, Platzer S, Bianco N, Jakse N, Payer M.</b> <b>Transcrestal maxillary sinus augmentation: Summers' versus a piezoelectric technique – an experimental cadaver study.</b> Clin. Oral Impl. Res. 2015; 1–4. DOI: 10.1111/cir.12546.	63
	<b>Zhen F, Fang W, Jing S, Zuolin W.</b> <b>The use of a piezoelectric ultrasonic osteotome for internal sinus elevation: a retrospective analysis of clinical results.</b> Int J Oral Maxillofac Implants. 2012 Jul-Aug;27(4):920-6.	63
	<b>Petrucci M, Ceccarelli R, Testori T, Grassi FR.</b> <b>Sinus floor augmentation with a hydropneumatic technique: a retrospective study in 40 patients.</b> Int J Periodontics Restorative Dent. 2012 Apr;32(2):205-10.	64
	<b>Baldi D, Menini M, Pera F, Ravera G, Pera P.</b> <b>Sinus floor elevation using osteotomes or piezoelectric surgery.</b> Int J Oral Maxillofac Surg. 2011 May;40(5):497-503.	65
	<b>Sentineri R, Dagnino G.</b> <b>Sinus Physiolift: a new technique for a less invasive great sinus augmentation with crestal approach.</b> J Osteol Biomat 2011;1:69-75.	66
<b>SINUS LIFT TECHNIQUE BY LATERAL APPROACH</b>	<b>Stacchi C, Vercellotti T, Toschetti A, Speroni S, Salgarello S, Di Lenarda R.</b> <b>Intraoperative Complications during Sinus Floor Elevation Using Two Different Ultrasonic Approaches: A Two-Center, Randomized, Controlled Clinical Trial.</b> Clin Implant Dent Relat Res. 2013 Aug 22.	66
	<b>Cassetta M, Ricci L, Iezzi G, Calasso S, Piattelli A, Perrotti V.</b> <b>Use of piezosurgery during maxillary sinus elevation: clinical results of 40 consecutive cases.</b> Int J Periodontics Restorative Dent. 2012 Dec;32(6):e182-8.	67
	<b>Cortes AR, Cortes DN, Arita ES.</b> <b>Effectiveness of piezoelectric surgery in preparing the lateral window for maxillary sinus augmentation in patients with sinus anatomical variations: a case series.</b> Int J Oral Maxillofac Implants. 2012 Sep-Oct;27(5):1211-5.	68
	<b>Wallace SS, Tarnow DP, Froum SJ, Cho SC, Zadeh HH, Stoupel J, Del Fabbro M, Testori T.</b> <b>Maxillary sinus elevation by lateral window approach: evolution of technology and technique.</b> J Evid Based Dent Pract. 2012 Sep;12(3 Suppl):161-71.	68

⇒ TABLE OF CONTENT

<b>Toscano NJ, Holtzclaw D, Rosen PS.</b> The effect of piezoelectric use on open sinus lift perforation: a retrospective evaluation of 56 consecutively treated cases from private practices. J Periodontol. 2010 Jan;81(1):167-71.	69
<b>Stübinger S, Saldamli B, Seitz O, Sader R, Landes CA.</b> Palatal versus vestibular piezoelectric window osteotomy for maxillary sinus elevation: a comparative clinical study of two surgical techniques. Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2009 May;107(5):648-55.	70
<b>Muñoz-Guerra MF, Naval-Gías L, Capote-Moreno A.</b> Le Fort I osteotomy, bilateral sinus lift, and inlay bone-grafting for reconstruction in the severely atrophic maxilla: a new vision of the sandwich technique, using bone scrapers and PIEZOSURGERY®. J Oral Maxillofac Surg. 2009 Mar;67(3):613-8.	71
<b>Wallace SS, Mazor Z, Froum SJ, Cho SC, Tarnow DP.</b> Schneiderian membrane perforation rate during sinus elevation using piezosurgery: clinical results of 100 consecutive cases. Int J Periodontics Restorative Dent. 2007; 27(5):413-419.	72
<b>Vercellotti T, De Paoli S, Nevins M.</b> The Piezoelectric Bony Window Osteotomy and Sinus Membrane Elevation: Introduction of a New Technique for Simplification of the Sinus Augmentation Procedure. Int J Periodontics Restorative Dent. 2001; 21(6):561-567.	72
<b>RIDGE EXPANSION</b>	
<b>Brugnami F, Caiazzo A, Mehra P.</b> Piezosurgery assisted, flapless split crest surgery for implant site preparation. J Maxillofac Oral Surg. 2014 Mar;13(1):67-72.	73
<b>Kelly A, Flanagan D.</b> Ridge expansion and immediate placement with piezosurgery and screw expanders in atrophic maxillary sites: two case reports. J Oral Implantol. 2013 Feb;39(1):85-90.	74
<b>Anitua E, Begoña I, Orive G.</b> Controlled ridge expansion using a two-stage split-crest technique with ultrasonic bone surgery. Implant Dent. 2012 Jun;21(3):163-70.	75
<b>Horrocks GB.</b> The controlled assisted ridge expansion technique for implant placement in the anterior maxilla: a technical note. Int J Periodontics Restorative Dent. 2010 Oct;30(5):495-501.	75
<b>Bertossi D, Albanese M, Bissolotti G, Bondi V, Nocini P.</b> L'utilizzo di impianti conici in contemporanea all'espansione alveolare piezoelettrica come alternativa alla chirurgia preprotesica con innesti ossei per la riabilitazione dei mascellari atrofici. Quintessenza Internazionale dental Tech – anno 26 – numero 9 bis – Speciale Implantologia 2010.	76

→ TABLE OF CONTENT

<b>Belleggia F, Pozzi A, Rocci M, Barlattani A, Gargari M.</b> Piezoelectric surgery in mandibular split crest technique with immediate implant placement: a case report. Oral Implantol (Rome). 2008 Oct;1(3-4):116-23.	77
<b>Enislidis G, Wittwer G, Ewers R.</b> Preliminary Report on a Staged Ridge Splitting Technique for Implant Placement in the Mandible: A Technical Note. Int J Oral Maxillofac Implants. 2006; 21(3):445-449.	78
<b>Vercellotti T, Russo C, Gianotti S.</b> A New Piezoelectric Ridge Expansion Technique in the Lower Arch - A Case Report. World Dentistry. 2000 online article.	79
<b>Vercellotti T.</b> Piezoelectric Surgery in Implantology: A Case Report - A New Piezoelectric Ridge Expansion Technique. Int J Periodontics Restorative Dent. 2000; 20(4):359-365.	79
<b>ORTHODONTIC MICROSURGERY</b> <b>Yu H, Jiao F, Wang B, Shen SG.</b> Piezoelectric decortication applied in periodontally accelerated osteogenic orthodontics. Craniofac Surg. 2013 Sep;24(5):1750-2.	80
<b>Grenga V, Bovi M.</b> Corticotomy-Enhanced Intrusion of an Overerupted Molar Using Skeletal Anchorage and Ultrasonic Surgery. J Clin Orthod. 2013 Jan;47(1):50-5.	80
<b>Hernández-Alfaro F, Guijarro-Martínez R.</b> Endoscopically assisted tunnel approach for minimally invasive corticotomies: a preliminary report. J Periodontol. 2012 May;83(5):574-80.	81
<b>Vercellotti T, Podestà A.</b> Orthodontic Microsurgery: A New Surgically Guided Technique for Dental Movement. Int J Periodontics Restorative Dent. 2007; 27:325-331.	82
<b>BONE GRAFTING</b> <b>Mouraret S, Houshyar KS, Hunter DJ, Smith AA, Jew OS, Girod S, Helms JA.</b> Cell viability after osteotomy and bone harvesting: comparison of piezoelectric surgery and conventional bur. Int J Oral Maxillofac Surg. 2014 Aug;43(8):966-71.	83
<b>Majewski P.</b> Piezoelectric surgery in autogenous bone block grafts. Int J Periodontics Restorative Dent. 2014 May-Jun;34(3):355-63.	84
<b>Majewski P.</b> Autogenous bone grafts in the esthetic zone: optimizing the procedure using piezosurgery. Int J Periodontics Restorative Dent. 2012 Dec;32(6):e210-7.	84

→ TABLE OF CONTENT

Lakshmiganthan M, Gokulanathan S, Shanmugasundaram N, Daniel R, Ramesh SB. Piezosurgical osteotomy for harvesting intraoral block bone graft. J Pharm Bioallied Sci. 2012 Aug;4(Suppl 2):S165-8.	84
Landes CA, Stübinger S, Laudemann K, Rieger J, Sader R. Bone harvesting at the anterior iliac crest using piezosteotomies versus conventional open harvesting: a pilot study. Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2008; 105(3):e19-e28.	85
Happe A. Use of a piezoelectric surgical device to harvest bone grafts from the mandibular ramus: report of 40 cases. Int J Periodontics Restorative Dent. 2007; 27(3):241-249.	86
Sohn DS, Ahn MR, Lee WH, Yeo DS, Lim SY. Piezoelectric osteotomy for intraoral harvesting of bone blocks. Int J Periodontics Restorative Dent. 2007; 27(2):127-131.	86
Gellrich NC, Held U, Schoen R, Pailing T, Schramm A, Bormann KH. Alveolar zygomatic buttress: A new donor site for limited preimplant augmentation procedures. J Oral Maxillofac Surg. 2007 Feb;65(2):275-80.	87
Patel A, Schofield J. Using Piezosurgery to harvest a block bone from the symphyseal region: a clinical case presentation. Implant Dentistry Today. 2007; 1(4):20-24.	87
Sivolella S, Berengo M, Scarin M, Mella F, Martinelli F. Autogenous particulate bone collected with a piezo-electric surgical device and bone trap: a microbiological and histomorphometric study. Arch Oral Biol. 2006; 51(10):883-891.	89
Stübinger S, Robertson A, Zimmerer SK, Leiggner C, Sader R, Kunz C. Piezoelectric Harvesting of an Autogenous Bone Graft from the Zygomaticomaxillary Region: Case Report. Int J Periodontics Restorative Dent. 2006; 26(5):453-457.	90
Held U, Bormann KH, Schmelzeisen R, Gellrich NC. Augmentation von Alveolarkammdefekten: Autologes Knochentransplantat aus der Crista zygomaticoalveolaris – eine neue Technik. Schweiz Monatsschr Zahnmed. 2005;115(8):692-703.	90
Bader G. Piezo chirurgie et greffes osseuses d'apposition L'information Dentaire, 2005;87(23):1377.	91
Chiriac G, Herten M, Schwarz F, Rothamel D, Becker J. Autogenous bone chips: influence of a new piezoelectric device (Piezosurgery®) on chips morphology, cell viability and differentiation. J Clin Periodontol. 2005; 32(9):994-999.	91

→ TABLE OF CONTENT

	<b>Boioli LT, Etrillard P, Vercellotti T, Tecucianu JF.</b> Piézochirurgie et aménagement osseux préimplantaire. Greffes par apposition de blocs d'os autogène avec prélèvement ramique. Implant. 2005; 11(4):261-274.	92
<b>OSTEOTOMY CLOSE TO NERVES</b>	<b>Bovi M, Manni A, Mavriqi L, Bianco G, Celletti R.</b> The use of piezosurgery to mobilize the mandibular alveolar nerve followed immediately by implant insertion: a case series evaluating neurosensory disturbance. Int J Periodontics Restorative Dent. 2010 Feb;30(1):73-81.	93
	<b>Degerliyurt K, Akar V, Denizci S, Yucel E.</b> Bone lid technique with PIEZOSURGERY® to preserve inferior alveolar nerve. Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2009 Dec;108(6):e1-5.	93
	<b>Schaeren S, Jaquiéry C, Heberer M, Tolnay M, Vercellotti T, Martin I.</b> Assessment of nerve damage using a novel ultrasonic device for bone cutting. Département de Parodontologie de l'UFR de Stomatologie et Chirurgie Maxillo-Faciale. J Oral Maxillofac Surg. 2008; 66(3):593-596.	94
	<b>Sakkas N, Otten JE, Gutwald R, Schmelzeisen R.</b> Transposition of the mental nerve by piezosurgery followed by postoperative neurosensory control: a case report. Br J Oral Maxillofac Surg. 2008; 46(4):270-271.	95
	<b>Tordjman S, Botoli LT.</b> Implants Juxta- canalaire. L'information Dentaire, 2007 May;89(26):1499.	95
	<b>Geha H, Gleizal A, Nimeskern N, Beziat JL.</b> Sensitivity of the Inferior Lip and Chin following Mandibular Bilateral Sagittal Split Osteotomy Using Piezosurgery. Plast Reconstr Surg. 2006; 118(7):1598-1607.	96
	<b>Bovi M.</b> Mobilization of the Inferior Alveolar Nerve with simultaneous implant insertion: A New Technique. A Case Report. Int J Periodontics Restorative Dent. 2005; 25(4):375-383.	97
<b>EXTRACTIONS</b>	<b>Spinato S., Rebaudi A., Bernardello F., Bertoldi C., Zaffe D.</b> Piezosurgical treatment of crestal bone: quantitative comparison of post-extractive socket outcomes with those of traditional treatment. Clin Oral Implants Res. 2015 Jan 30.	98
	<b>Nagori SA, Jose A, Bhutia O, Roychoudhury A.</b> Evaluating success of autotransplantation of embedded/impacted third molars harvested using piezosurgery: a pilot study. Acta Odontol Scand. 2014 Nov;72(8):846-51.	99
	<b>Piersanti L, Dilorenzo M, Monaco G, Marchetti C</b> Piezosurgery or Conventional Rotatory Instruments for Inferior Third Molar Extractions? J Oral Maxillofac Surg. 2014 Sep;72(9):1647-52.	99

⇒ TABLE OF CONTENT

<b>Mozzati M, Gallesio G, Russo A, Staiti G, Mortellaro C.</b> Third-molar extraction with ultrasound bone surgery: a case-control study. Craniofac Surg. 2014 May;25(3):856-9.	100
<b>Koszowski R, Morawiec T, Bubilek-Bogacz A.</b> Use of the piezosurgery technique for cutting bones in the autotransplantation of unerupted third molars. Int J Periodontics Restorative Dent. 2013 Jul-Aug;33(4):477-81.	101
<b>Rullo R, Addabbo F, Papaccio G, D'Aquino R, Festa VM.</b> Piezoelectric device vs. conventional rotative instruments in impacted third molar surgery: relationships between surgical difficulty and postoperative pain with histological evaluations. J Craniomaxillofac Surg. 2013 Mar;41(2):e33-8.	101
<b>Itro A, Lupo G, Marra A, Carotenuto A, Cocozza E, Filipi M, D'Amato S.</b> Confronto tra tecnica di osteotomia piezoelettrica e tecnica con strumenti rotanti nella chirurgia dei terzi molari inclusi. Uno studio clinico. Minerva Stomatol. 2012 Jun;61(6):247-53	102
<b>Guo ZZ, Zhang H, Li Y, Li X, Liu Y, Wang Y, Yuan CX, Liu X.</b> Comparative study of complications among routine method, high speed turbine handpiece and piezosurgery device after extraction of impacted wisdom teeth. Shanghai Kou Qiang Yi Xue. 2012 Apr;21(2):208-10.	103
<b>Gao Y, Jiang A, Li B, Yang L.</b> Comparison of piezosurgery and chisel osteotomy in the extraction of mandibular impacted third molars. Hua Xi Kou Qiang Yi Xue Za Zhi. 2011 Aug;29(4):372-4.	104
<b>Sortino F, Pedullà E, Masoli V.</b> The piezoelectric and rotatory osteotomy technique in impacted third molar surgery: comparison of postoperative recovery. J Oral Maxillofac Surg. 2008 Dec;66(12):2444-8.	105
<b>Sivolella S, Berengo M, Fiorot M, Mazzuchin M.</b> Retrieval of blade implants with piezosurgery: two clinical cases. Minerva Stomatol. 2007; 56(1-2):53-61.	105
<b>Grenga V, M. Bovi.</b> Piezoelectric Surgery for Exposure of Palatally Impacted Canines. J Clin Orthod. 2004; 38(8):446-448.	106
<b>EXPLANTATION</b>	
<b>Marini E, Cisterna V, Messina AM.</b> The removal of a malpositioned implant in the anterior mandible using piezosurgery. Oral Surg Oral Med Oral Pathol Oral Radiol. 2013 May;115(5):e1-5.	107
<b>Moorhouse J, Campbell I.</b> Recovery of a poorly-placed implant using piezosurgery. Implant Dentistry Today – 2012 September.	107
<b>Sammartino G, Trosino O, di Lauro AE, Amato M, Cioffi A.</b> Use of piezosurgery device in management of surgical dental implant complication: a case report. Implant Dent. 2011 Apr;20(2):e1-6.	108

→ TABLE OF CONTENT

<b>ENDO SURGERY</b>	<b>Abella F, de Ribot J, Doria G, Duran-Sindreu F, Roig M.</b> Applications of piezoelectric surgery in endodontic surgery: a literature review. J Endod. 2014 Mar;40(3):325-32.	108
	<b>Del Fabbro M, Tsesis I, Rosano G, Bortolin M, Taschieri S.</b> Scanning electron microscopic analysis of the integrity of the root-end surface after root-end management using a piezoelectric device: a cadaveric study. J Endod. 2010 Oct;36(10):1693-7.	109
<b>PERIODONTAL SURGERY</b>	<b>Seshan H, Konuganti K, Zope S</b> PIEZOSURGERY® in periodontology and oral implantology. J Indian Soc Periodontol. 2009 Sep;13(3):155-6.	110
	<b>Vercellotti T, Pollack AS.</b> A New Bone Surgery Device: Sinus Grafting and Periodontal Surgery. Compend Contin Educ Dent. 2006; 27(5):319-325.	110
<b>DISTRACTION OSTEOGENESIS</b>	<b>Menini I, Zornitta C, Menini G.</b> Distraction Osteogenesis for Implant Site Development Using a Novel Orthodontic Device: A Case Report. Int J Periodontics Restorative Dent. 2008; 28(2):189-196.	111
	<b>Lee HJ, Ahn MR, Sohn DS.</b> Piezoelectric distraction osteogenesis in the atrophic maxillary anterior area: a case report. Implant Dent. 2007 Sep;16(3):227-34.	111
	<b>González-García A, Diniz-Freitas M, Somoza-Martín M, García-García A.</b> Piezoelectric Bone Surgery Applied in Alveolar Distraction Osteogenesis: A Technical Note. Int J Oral Maxillofac Implants. 2007; 22(6):1012-1016.	112
<b>MAXILLO-FACIAL SURGERY</b>	<b>Pappalardo S, Guarnieri R.</b> Randomized clinical study comparing piezosurgery and conventional rotatory surgery in mandibular cyst enucleation. J Craniomaxillofac Surg. 2014 Jul;42(5):e80-5.	112
	<b>D'Amato S, Sgaramella N, Vanore L, Piombino P, Orabona GD, Santagata M.</b> Piezoelectric bone surgery in the treatment of an osteoma associated with an impacted inferior third molar: a case report. Clin Cases Miner Bone Metab. 2014 Jan;11(1):73-6.	113
	<b>Rodriguez JG, Eldibany RM.</b> Vertical splitting of the mandibular body as an alternative to inferior alveolar nerve lateralization. Int J Oral Maxillofac Surg. 2013 Sep;42(9):1060-6.	113
	<b>Bertossi D, Lucchese A, Albanese M, Turra M, Faccioni F, Nocini P, Rodriguez Y Baena R.</b> Piezosurgery versus conventional osteotomy in orthognathic surgery: a paradigm shift in treatment. J Craniofac Surg. 2013 Sep;24(5):1763-6.	114

⇒ TABLE OF CONTENT

- Rana M, Gellrich NC, Rana M, Piffkó J, Kater W.**  
Evaluation of surgically assisted rapid maxillary expansion with piezosurgery versus oscillating saw and chisel osteotomy - a randomized prospective trial.  
*Trials.* 2013 Feb 17;14:49. 115
- Sammartino G, Ricciutello F, Trosino O, Marenzi G, Cioffi A, Mortellaro C.**  
Uso del dispositivo piezochirurgico nella gestione delle complicazioni di chirurgia orale: relazione su caso clinico ed esperienza clinica.  
*Minerva Stomatol.* 2012 May;61(5):225-31. 116
- Wagner ME, Rana M, Traenkenschuh W, Kokemueller H, Eckardt AM, Gellrich NC.**  
Piezoelectric-assisted removal of a benign fibrous histiocytoma of the mandible: an innovative technique for prevention of dentoalveolar nerve injury.  
*Head Face Med.* 2011 Oct 31;7:20. 116
- Robiony M, Polini F.**  
Piezosurgery: a safe method to perform osteotomies in young children affected by hemifacial microsomia.  
*J Craniofac Surg.* 2010 Nov;21(6):1813-5. 117
- Beziat JL, Faghahati S, Ferreira S, Babic B, Gleizal A.**  
Blocage maxillomandibulaire: technique et intérêt dans le clivage sagittal piézoélectrique.  
*Rev Stomatol Chir Maxillofac.* 2009 Nov;110(5):273-7. Epub 2009 Oct 20. 117
- Bader G, Morais D.**  
Apport de la piézochirurgie pour l'avancée des génotubercules dans le syndrome d'apnées obstructives du sommeil.  
*Rev Stomatol Chir Maxillofac.* 2008 Dec;109(6):375-8. 118
- Landes CA, Stübinger S, Ballon A, Sader R**  
Piezoosteotomy in orthognathic surgery versus conventional saw and chisel osteotomy.  
*Oral Maxillofac Surg.* 2008 Sep;12(3):139-47. 119
- Landes CA, Stübinger S, Rieger J, Williger B, Ha TK, Sader R.**  
Critical evaluation of piezoelectric osteotomy in orthognathic surgery: operative technique, blood loss, time requirement, nerve and vessel integrity.  
*J Oral Maxillofac Surg.* 2008 Apr;66(4):657-74. 120
- J. González-Lagunas, J. Mareque**  
Piezosurgery®: its role in TMJ surgery.  
*Ital J Maxillofac Surg.* 2008 December;19(3):119-21. 121
- Beziat JL, Béra, JC, Lavandier B, Gleizal A.**  
Ultrasonic osteotomy as a new technique in craniomaxillofacial surgery.  
*International Journal of Maxillo-facial Surgery,* 2007;36(6):493-500. 121
- Robiony M, Polini F, Costa F, Zerman N, Polit M.**  
Ultrasound bone cutting for surgically assisted rapid maxillary expansion (SARME) under local anaesthesia.  
*Int J Oral Maxillofac Surg.* 2007;36(3):267-9. 122

→ TABLE OF CONTENT

- Gleizal A, Béra JC, Lavandier B., Béziat JL.**  
Piezoelectric osteotomy: a new technique for bone surgery – advantages in craniofacial surgery.  
Childs Nerv Syst . 2007;23(5):509-513. 122
- Cipriano L, Cimmino R, De Paolis G, Guerra F, Pillon A, Caputo M, Izzo P, Trombetta S, Basso L, Izzo L.**  
Asportazione di enostosi mandibolare mediante tecnica piezoelettrica: case report.  
G Chir 2007 May; 28(5):222-6. 123
- Beziat JL, Vercellotti T, Gleizal A.**  
Qu'est-ce que la PIEZOSURGERY®? Intérêt en Chirurgie cranio-maxillofaciale.  
A propos de deux ans d'expérience.  
Revue de Stomatologie et Chir Maxillofaciale, 2007 Apr;108(2):101-107. 124
- Guo ZZ, Liu X, Li Y, Deng YF, Wang Y**  
The use of PIEZOSURGERY® osteotomy in treatment of longstanding maxillary fractures: report of 12 consecutive patients.  
Shangai Kou Qiang Yi Xue. 2007 Feb;16(1):97-9. 125
- Robiony M, Polini F, Costa F, Vercellotti T, Politi M.**  
Piezoelectric Bone Cutting in multipiece maxillary osteotomies. Technical Note.  
J Oral Maxillofac Surg. 2004; 62:759-761. 126
- Eggers G, Klein J, Blank J, Hassfeld S.**  
PIEZOSURGERY®: an ultrasound device for cutting bone and its use and limitations in maxillofacial surgery.  
Br J Oral Maxillofac Surg. 2004 Apr;42(5):451-3. 127

Pereira CC, Gealh WC, Nogueira LM, Garcia Junior IR, Okamoto R.

## Piezosurgery applied to implant dentistry: clinical and biological aspects.

EN J Oral Implantol. 2014 Jul;40 Spec No:401-8.

### → Piezosurgery applied to implant dentistry: clinical and biological aspects.

The piezosurgery is a new and modern technique of bone surgery in implantology. Selective cutting is possible for different ultrasonic frequencies acting only in hard tissues (mineralized), saving noble anatomical structures. With the piezoelectric osteotomy technique, the preparation of receptor sites for implants, obtaining autogenous bone graft (particles and blocks), osteotomy for the alveolar bone crest expansion, maxillary sinus lifting and removal of dental implants can be performed accurately and safely, providing excellent clinical and biological results, especially for the osteocytes viability. The aim of this study was, through literature, present clinical applications of piezosurgery in implant dentistry, as well as their advantages and disadvantages over conventional surgical systems. Moreover, to address the biological aspects related to piezosurgery that differentiates of the bone tissue approaches. Concluding that piezosurgery transform critical operations in simple and fully executable procedures. Effectively, surgeries performed in difficult access areas come to be of less risk to soft and neurovascular tissues.

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Zhao J, Huang C.

## The advanced techniques of dentoalveolar surgery.

CN CHua Xi Kou Qiang Yi Xue Za Zhi. 2014 Jun;32(3):213-6.

### → The advanced techniques of dentoalveolar surgery.

During the recent decade, the advanced medical technology has brought the rapid development in the dento-alveolar surgery, and the medical concepts of humanity, painless, minimally invasive, safe and comfortable are gradually accepted by the patients and doctors. Many advanced techniques and equipments have been used in the dentoalveolar surgery. This paper would like to make a review on the clinical application and experience of the advanced medical technologies in the tooth extraction, such as the general anaesthetic technique, the nitrous oxide sedative technique, the computer-controlled local anesthesia delivery devices, the electrocardiogram monitor, the location of embedded tooth by cone beam CT and minimally invasive surgical technique. Additionally, in this paper we also described the indication, contraindication, advantages and operating essentials of the orthodontic traction technique, ultrasonic bone surgery (piezosurgery) device and the surgery technology for alveolar bone shape and functional preservation in the extraction of tooth. Finally, we introduced the immediate implantation technique used in the dentoalveolar surgery.

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Bauer SE, Romanos GE.

## Morphological characteristics of osteotomies using different piezosurgical devices. A scanning electron microscopic evaluation.

EN Implant Dent. 2014 Jun;23(3):334-42.

⇒ **Morphological characteristics of osteotomies using different piezosurgical devices. A scanning electron microscopic evaluation.**

**Purpose:** The objective of this study was to compare morphological characteristics of osteotomies performed by 6 Piezosurgical devices.

**Materials and Methods:** The 6 Piezosurgical units were: (a) Piezotom, (b) SurgySonic, (c) Piezon Master Surgery, (d) VarioSurg, (e) Surgybione, and (f) Piezosurgery 3. Osteotomies on 9 freshly slaughtered cattle ribs (2 cuts by each unit, per rib) from the cortical (first cut at 5 mm) to the cancellous (second cut at 3 mm) bone layer were performed. The osteotomy margins were compared using scanning electron microscopy analysis. The cutting areas, osteotomy bottoms, and osteotomy margins were analyzed morphologically. Statistical evaluation of the 2 cuts regarding the design of the tips (a-d: tapered tool shanks, f: parallel tool shank) was performed by an unpaired t test.

**Results:** Morphological characteristics were different for each Piezosurgical unit and each examined area. A significant difference ( $P = 0.0209$ ) of the upper width of the first cut between tapered and parallel tips was shown.

**Conclusions:** The morphological characteristics of the produced Piezosurgical osteotomies vary and depend on the Piezosurgical unit and tip.

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Patel A.

## The role of piezosurgery in implant dentistry.

EN IDT May 2014

→ **The role of piezosurgery in implant dentistry:** Adam Patel details the huge potential of piezosurgery for improving the predictability and ease of dental implant treatment.

Piezosurgery was first introduced in 1988. Its development was encouraged by the need for high levels of precision and safety in bone surgery compared with that achieved by standard bur and saw instruments (Landes et al, 2008). Piezosurgery as a technique spread because of its ease of use and safety. The piezoelectric effect occurs when an electric current is passed around a stack of crystals and they start to vibrate at a precise frequency. The piezoelectric instrument produces a modulated ultrasonic frequency of 24 to 29 kHz, and a microvibration amplitude between 60 and 200 mm/sec (Sortino et al, 2008). The amplitude of these microvibrations allows a clean, precise and controlled cut of bony structures without causing destruction of soft tissue (including nerves, blood vessels and oral mucosa) (Eggers et al, 2004). Since its introduction, piezosurgery has established an important role in various aspects of dentistry and dental implantology (see the box: Uses of piezosurgery in dental implantology). This article will discuss the role of piezosurgery in these areas, including its advantages and disadvantages.

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Franco S, Miccoli S, Limongelli L, Tempesta A, Favia G, Maiorano E, Favia G.

## New Dimensional Staging of Bisphosphonate-Related Osteonecrosis of the Jaw Allowing a Guided Surgical Treatment Protocol: Long-Term Follow-Up of 266 Lesions in Neoplastic and Osteoporotic Patients from the University of Bari.

EN Int J Dent. 2014 Jun;2014:935657.

→ **Clinical Study**

**New Dimensional Staging of Bisphosphonate-Related Osteonecrosis of the Jaw Allowing a Guided Surgical Treatment Protocol: Long-Term Follow-Up of 266 Lesions in Neoplastic and Osteoporotic Patients from the University of Bari.**

Bisphosphonate-related osteonecrosis of the jaw (BRONJ) is the most serious side effect in patients receiving bisphosphonates (BPs) for neoplastic disease and osteoporosis. The aim of this study is to propose a new dimensional stage classification, guiding the surgical treatment of BRONJ patients, and to evaluate the success rate of this new management. From 2004 to 2013, 203 neoplastic and osteoporotic patients with 266 BRONJ lesions were referred to the Odontostomatology Unit of the University of Bari. All patients underwent surgery after suspension of BPs therapy and antibiotic treatment. The surgical procedure was complemented by piezosurgery and followed by the application of hyaluronate and amino acids. The newdimensional staging suggests the choice of the surgical approach, and allows the prediction of postoperative complications and soft and hard tissues healing time, guiding the surgical treatment protocol. This protocol could be a successful management strategy for BRONJ, considering the low recurrences rate and the good stabilisation of the surgical sites observed after a long-term follow-up.

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Gülnahar Y, Hüseyin Köşger H, Tutar Y.

## A comparison of piezosurgery and conventional surgery by heat shock protein 70 expression.

EN Int J Oral Maxillofac Surg. 2013 Apr;42(4):508-10.

⇒ A comparison of piezosurgery and conventional surgery by heat shock protein 70 expression.

The effects of mechanical instruments on the viability of cells are essential in terms of regeneration. There is considerable interest in cell repair following damage mediated by dental surgical procedures. Cells can tolerate stress by expressing heat shock protein 70 (Hsp70). During and after surgical tooth removal, oxidative stress can activate Hsp70 expression proportional to the intensity of the stress signal stimulus to cope with stress. This study examined the expression of Hsp70 as a potential biomarker of immediate postoperative stress in patients undergoing two different surgical procedures of different severity. Expression of Hsp70 both at mRNA and at protein level in the conventional group was two-fold higher than that of the piezo group. This suggests that tooth movement by the piezo method causes relatively lower stress in the alveolar bone. Piezosurgery provides relatively low stress to the patients and this may help cell repair after the surgical procedure. Patients undergoing more aggressive surgery using conventional methods showed a significant increase in Hsp70 in the immediate postoperative period. Therefore, Hsp70 induction can be a potential tool as a prognostic surgical marker.

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## → TECHNIQUE ADVANTAGES

Claire S, Lea SC, Walmsley AD.

### Characterisation of bone following ultrasonic cutting.

EN Clin Oral Investig. 2013 Apr;17(3):905-12.

#### ⇒ Characterisation of bone following ultrasonic cutting.

**Objectives:** Ultrasonic surgery is an increasingly popular technique for cutting bone, but little research has investigated how the ultrasonic tip oscillations may affect the cuts they produce in bone. The aim of this investigation was to evaluate the oscillation and cutting characteristics of an ultrasonic surgical device.

**Materials and Methods:** A Piezosurgery 3 (Mectron, Carasco, Italy) ultrasonic cutting system was utilised with an OP3 style tip. The system was operated with the tip in contact with porcine bone samples (loads of 50 to 200 g) mounted at 45° to the vertical insert tip and with a water flow of 57 ml/min. Tip oscillation amplitude was determined using scanning laser vibrometry. Bone surfaces defects were characterised using laser profilometry and scanning electron microscopy.

**Results:** A positive relationship was observed between the magnitude of tip oscillations and the dimensions of defects cut into the bone surface. Overloading the tip led to a reduction in oscillation and hence in the defect produced. A contact load of 150 g provided the greatest depth of cut. Defects produced in the bone came from two clear phases of cutting.

**Conclusions:** The structure of the bone was found to be an important factor in the cut characteristics following piezosurgery.

**Clinical Relevance:** Cutting of bone with ultrasonics is influenced by the load applied and the setting used. Care must be used to prevent the tip from sliding over the bone at low loadings.

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### Bone characteristics following osteotomy surgery: an in vitro SEM study comparing traditional Lindemann drill with sonic and ultrasonic instruments.

EN POSEIDO 2013;1(3):187.

#### ⇒ Bone characteristics following osteotomy surgery: an in vitro SEM study comparing traditional Lindemann drill with sonic and ultrasonic instruments.

**Background and objectives:** Osteotomy surgery is widely used in dental surgery for implant site preparation, bone grafting and GBR. In this study, the characteristics of bone surfaces were examined after bone osteotomy surgery performed with the Lindemann bur, sonic (Komet Sonosurgery) and ultrasonic (Mectron Piezosurgery) instruments.

**Materials and Methods:** Anatomic integrity and osteotomic precision were analyzed using Scanning Electron Microscopy (SEM) to observe vascular canals, microfractures, exfoliations and bone debris on cortical and cancellous surfaces cut with

the 3 types of instruments.

**Results:** The use of ultrasonic instruments resulted in extremely precise cuts and reduced bone damage. The sonic instrument was precise in cortical bone but showed minor signs of bone damage in cancellous bone. Lindemann bur showed less precision and higher bone damage both in cortical and in cancellous bone. In cortical bone, ultrasonic and sonic cuts showed nicely opened bone vascular canals, while Lindemann bur showed many canals closed by abrasions, exfoliation and cracks by dragging attrition. In cancellous bone, ultrasonic cut showed intact trabeculae and trabecular spaces free of debris, while sonic cut showed more debris accumulation in trabecular spaces. Lindemann bur showed huge quantity of bone debris that filled trabecular spaces.

**Discussion and Conclusion:** For all parameters, the ultrasonic cut offered the most precise and atraumatic bone cut. Ultrasonic and sonic instruments both showed more precise and less traumatic results than the Lindemann bur.

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## Intraosseous temperature changes during the use of piezosurgical inserts in vitro.

EN Int J Oral Maxillofac Surg. 2012 Nov;41(11):1338-43.

⇒ **Intraosseous temperature changes during the use of piezosurgical inserts in vitro.**

This study concerns intraosseous temperature changes during the use of piezosurgical inserts. On six fresh pig jaws heated to body temperature (36°C), osteotomies and osteoplasties were performed in vitro with the Piezosurgery 3 device (Mectron, Carasco, Italy) and various inserts. The intraosseous temperature increases were measured at a depth of 3 mm and at a distance of 1 mm from the working site using nickelchromium/nickel temperature sensors. 20°C Ringer's solution was used for cooling in an initial test series and 10°C Ringer's in a second series. The processed bone was examined using digital volume tomography images to determine the ratio of cortical to cancellous bone thickness. Mean temperature increases of 4.4 - 10.9°C were found; maximum temperature peaks were over 47°C for an average of only 8.5 s. The type of piezosurgical insert had a marked influence on intraosseous temperature generation ( $p=0.026$ ); the thickness of the cortical bone and the temperature of the coolant did not. Coolant temperature had an influence on the bone cooling time ( $p=0.013$ ). The results show that correct use of the piezosurgery device does not give rise to prolonged temperature increases over 47°C and hence does not cause any irreversible thermal damage in the bone.

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**IT** Minerva Stomatol. 2012 May;61(5):213-24.

⇒ **Benefits of piezoelectric surgery in oral and maxillofacial surgery. Review of literature.**

Piezoelectric surgery is based on the use of ultrasound for the cutting of bones. It represents an innovative technique, as it offers the maxillofacial surgeon the opportunity of making precise bone cuts without damaging any soft tissue, minimizing the invasiveness of surgical procedure, and the opportunity of working in a field which is almost totally blood-free. It reduces the impact on soft tissues (vessels and nerves) which lie adjacent to the areas of treatment. Compared to traditional methods, it enables optimal healing because it reduces the postsurgery swelling and discomfort. In this article the authors realized a review of the literature.

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**Parmar D, Mann M, Walmsley AD, Lea SC.**

## Cutting characteristics of ultrasonic surgical instruments.

**EN** Clin Oral Implants Res. 2011 Dec;22(12):1385-90.

⇒ **Cutting characteristics of ultrasonic surgical instruments.**

**Objective:** Ultrasonic surgical devices are becoming increasingly popular and work is required to understand the performance of the cutting tips. This experimental study looks to investigate the way in which ultrasonic bone cutting tools oscillate and how this oscillation is modified when contacted against bone surfaces with varying loads. The defects produced in instrumented bone surfaces were measured and related to the tip motion.

**Methods:** An ultrasonic cutting probe was scanned, unloaded, using a scanning laser vibrometer to determine its free oscillation pattern and amplitude. This probe tip was then contacted against bone under various loads to assess the modification in oscillation characteristics. Cuts were performed over a period of 10 s. The cut bone surfaces were assessed using laser profilometry to determine defect depths.

**Results:** The average vibration displacement amplitude at the probe tip, under load, was <12 µm in the longitudinal direction and was greatest for the cortical cutting mode. Elliptical probe motion was successfully mapped out under the range of loads tested. Defect depths of up to 0.36 mm were detected and were greatest when the tip was in contact with the bone with a load of 100 g.

**Conclusions:** This work showed that the nature of the surface being cut may significantly alter the mode shape and magnitude of the probe oscillation. The maximum depth of cut with minimum restraining of tip motion was achieved at 100 g contact load.

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## Piezosurgery prevents brain tissue damage: an experimental study on a new rat model.

EN Int J Oral Maxillofac Surg. 2011 Aug;40(8):840-4.

⇒ **Piezosurgery prevents brain tissue damage: an experimental study on a new rat model.**

Piezosurgery is a promising meticulous system for bone cutting, based on ultrasound microvibrations. It is thought that the impact of piezosurgery on the integrity of soft tissue is generally low, but it has not been examined critically. The authors undertook an experimental study to evaluate the brain tissue response to skull bone removal using piezosurgery compared with a conventional drilling method. In Wistar male rats, a circular bone window was drilled to the parietal bone using piezosurgery on one side and a conventional bone drill on the other side. The behavioural performance of animals was evaluated using the motor BBB test and sensory plantar test. The brains of animals were evaluated by magnetic resonance imaging (MRI) and histology. The results of MRI showed significantly increased depth and width of the brain lesion in the region of conventional drilling compared with the region where piezosurgery was used. Cresylviolet and NF 160 staining confirmed these findings. There was no significant difference in any of the behavioural tests between the two groups. In conclusion, piezosurgery is a safe method for the performance of osteotomy in close relation to soft tissue, including an extremely injury-sensitive tissue such as brain.

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## Piezosurgery in oral and maxillofacial surgery.

EN Int J Oral Maxillofac Surg. 2011 May;40(5):451-7.

### → Piezosurgery in oral and maxillofacial surgery.

This review summarizes current knowledge and experience with piezosurgery, a promising, meticulous and soft tissue-sparing system for bone cutting, based on ultrasonic microvibrations. The main advantages of piezosurgery include soft tissue protection, optimal visibility in the surgical field, decreased blood loss, less vibration and noise, increased comfort for the patient and protection of tooth structure. To date it has been indicated for use in oral and maxillofacial surgery, otorhinolaryngology, neurosurgery, ophthalmology, traumatology and orthopaedics. The main indications in oral surgery are sinus lift, bone graft harvesting, osteogenic distraction, ridge expansion, endodontic surgery, periodontal surgery, inferior alveolar nerve decompression, cyst removal, dental extraction and impacted tooth removal. In conclusion, piezosurgery is a promising technical modality for different aspects of bone surgery with a rapidly increasing number of indications throughout the whole field of surgery.

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**Maurer P, Kriwalsky MS, Block Veras R, Vogel J, Syrowatka F, Heiss C.**

## Micromorphometrical analysis of conventional osteotomy techniques and ultrasonic osteotomy at the rabbit skull.

EN Clin Oral Implants Res. 2008 Jun;19(6):570-5.

### → Micromorphometrical analysis of conventional osteotomy techniques and ultrasonic osteotomy at the rabbit skull.

**Objectives:** The ultrasonic osteotome, which was recently introduced, is an alternative to conventional methods of osteotomy. The aim of the present study was to establish the differences between three osteotomy techniques and to perform a quantitative roughness analysis of the osteotomized bone surfaces.

**Materials and Methods:** Fresh bony samples of standardized size were taken from the rabbit skull. The techniques used were as follows: reciprocate micro-saw, Lindemann bur, ultrasonic osteotome with the two insert tips OT6 (rough) and OT7 (fine). The prepared surfaces were examined by light microscopy, environmental surface electron microscopy (ESEM) and by confocal laser scanning microscopy (CLSM).

**Results:** It was difficult to distinguish between cortical and cancellous bone after using the conventional osteotomy technique. The ultrasonic technique preserved the original structure of the bone. The values observed for superficial roughness were as follows: 3.97 microm (micro-saw), 5.7 microm (Lindemann bur), 2.48 microm (OT7) and 3 microm (OT6). There were statistical differences between the values of the bur and insert tip OT6 ( $P=0.015$ ) as well as between the bur and insert tip OT7 ( $P=0.003$ ).

## → TECHNIQUE ADVANTAGES

**Conclusions:** In the present study micromorphological differences after using various osteotomy techniques could be clearly identified.

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**Labanca M, Azzola F, Vinci R, Rodella LF.**

### **Piezoelectric surgery: Twenty years of use.**

**EN** Br J Oral Maxillofac Surg. 2008; 46(4):265-269

⇒ **Piezoelectric surgery: Twenty years of use.**

The use of ultrasonic vibrations for the cutting of bone was first introduced two decades ago. Piezoelectric surgery is a minimally invasive technique that lessens the risk of damage to surrounding soft tissues and important structures such as nerves, vessels, and mucosa. It also reduces damage to osteocytes and permits good survival of bony cells during harvesting of bone.

Piezoelectric surgery was first used by oral and maxillofacial surgeons for osteotomies, but recently some specific applications in neurosurgery and orthopaedics have been proposed. We review the different applications of piezoelectric surgery.

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### ⇒ Ultrasonic Bone Cutting in Oral Surgery: a Review of 60 Cases.

**Purpose:** Surgical reconstruction of bony defects in the oral cavity can often be challenging since thin and fragile bony structures are especially prone to fracture caused by bulky cutting tips or the application of significant pressure by conventional mechanical instruments. The risk of accidental damage to adjacent soft tissue structures, such as nerves, by a dental drill or saw is also extremely high. The use of modulated ultrasound (piezosurgery) makes it possible to overcome such complications as a result of the precise and minimally invasive surgery technique which is limited to mineralized hard tissue.

**Materials and Methods:** In 60 patients (38 male, 22 female) a piezosurgery device was used for different bone augmentation procedures before dental implant placement. The instrument uses modulated ultrasound (25 - 30 kHz) and the amplitude of the working tip ranges from 60 µm to 200 µm. The device was employed for sinus floor elevation (25), alveolar ridge augmentation using an autogenous block graft (25), alveolar ridge splitting (5) or lateralization of the alveolar nerve (5). Physiological sodium chloride was used as a cooling solution. For all osteotomies mode boosted burst c and pump 5 were used.

**Results:** Piezoelectric osteotomy permitted micrometric selective cutting and a clear surgical site due to the cavitation effect created by the cooling solution and the oscillating tip. No excessive bleeding was encountered. The risk of accidental soft tissue harm, such as perforating the sinus membrane or damaging adjacent nerves, was definitely lower than in the case of a conventional bur. No serious complications were encountered in the postoperative wound healing process after 2, 14, 30 and 90 days. However, the surgical procedures were time-consuming.

**Conclusion:** Piezosurgery is an advantageous osteotomy technique for delicate structures in the oral and maxillofacial region. With respect to osteotomies of thin and fragile bones, the application of ultrasound is superior to other mechanical instruments because of the extremely precise and virtually arbitrary cut geometries, easy handling, efficient bone ablation and minimal accidental damage to adjacent soft tissue structures.

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⇒ **Piezoelectric osteotomy: a new technique for cutting bone in implant surgery.**

Within a few years of ultrasonic technology, a very powerful ultrasonic device has been developed. This device, together with its modified tips, is a tool that allows impressive cutting, providing numerous very practical solutions in implant surgery. This review article describes the principles of piezoelectric osteotomy as well as a survey of the numerous known clinical applications of this instrument in implant surgery.

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→ Light microscopic examination of rabbit skulls following conventional and ultrasonic osteotomy.

**Introduction:** The novel ultrasonic osteotomy technique (PIEZOSURGERY®) is an alternative to conventional osteotomy devices. The aim of the present study was to carry out morphological comparison of the bone surface using conventional osteotomy techniques in comparison to the rather new ultrasonic osteotomy technique by means of a reflected-light microscopic examination.

**Materials and Methods:** Following the sacrifice of 12 rabbits, 24 standardized bone samples were removed from the skull. The osteotomy devices used were a rotating instrument (Lindemann bur), an oscillating micro-saw, and an ultrasonic osteotomy device (PIEZOSURGERY®) with insert tips OT6 and OT7. The times needed for osteotomy were measured. The bone surfaces were examined using a reflected-light microscope with a magnification of 40x and 100x.

**Results:** Osteotomy with PIEZOSURGERY® is significantly more time consuming than osteotomy with conventional methods ( $p<0,05$ ). Following osteotomy with the ultrasonic device, the reflected-light microscopic examinations of the unmodified bone samples revealed typical bone structure of the calvaria, including compacta externa, diploe and compacta interna. On the contrary, following osteotomy with the conventional devices, the diploe structure presented distinct modifications. The cancellous spaces were filled with bone debris, and the cancellous structure was demolished. The samples prepared by the micro-saw technique showed a superficially condensed and grooved surface.

**Conclusion:** In the present study, well-defined differences were observed following osteotomy with conventional devices and osteotomy with the ultrasonic device. The integrity of the bony structure observed after the ultrasonic technique could benefit the bone healing process. Further studies dealing with the bone healing process after using different osteotomy techniques are recommended.

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⇒ Sound levels of the PIEZOSURGERY®: Risk of permanent damage to hearing.

In the past, research has regularly been carried out concerning the sound levels of various drilling devices and the impact these have on those who regularly use these devices. The present research is concerned with the possible permanent damage to hearing which can occur during the use of a newly developed instrument, the PIEZOSURGERY®. Measurements have been performed in a standardized set-up in which use is made of the lower jaws of pigs and a weight scale for measuring various degrees of pressure. The boundary values of the permissible exposure to noise were determined. The values of the PIEZOSURGERY® were compared with conventional drilling machines. It was concluded that using the PIEZOSURGERY® for less than 1.5 hours per day implies no risk of permanent damage to hearing. This means that the risk in the case of PIEZOSURGERY® is somewhat lower than that of conventional drilling devices.

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**Peivandi A, Bugnet R, Debize E, Gleizal A, Dohan DM.**

## Piezoelectric osteotomy: applications in periodontal and implant surgery.

**EN** Rev Stomatol Chir Maxillofac. 2007; 108(5):431-440.

⇒ Piezoelectric osteotomy: applications in periodontal and implant surgery.

Reverse piezoelectricity is currently used in Dentistry for ultrasonic scaling devices. Using this technology, increasingly more powerful ultrasonic surgical lancets were developed within a few years, and these new tools have provided many practical solutions in oral and maxillofacial surgery. This article reviews the basic principles of piezosurgery, and its numerous clinical applications: in oral surgery (atraumatic extractions, wisdom teeth exposure, periapical cyst debridement, pre-prosthetic surgery), more specifically in periodontal (root planning and bone remodeling, crown lengthening) and implant surgery (sinus lift, removal

of fractured implants, bone ridge augmentation, bone graft harvesting (crestal, retromolar, chin or tori)). The recent increased power of these instruments allows emphasizing a much wider range of clinical applications, extending to all surgical fields.

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Su YC.

## Development and clinical application of ultrasonic osteotomy in dentistry.

CN Shanghai kou qiang yi xue = Shanghai journal of stomatology. 2007; 16(1):1-7 [chinese]

⇒ **Development and clinical application of ultrasonic osteotomy in dentistry.**

Ultrasonic osteotomy is a new technique for osseous surgery utilizing an innovative ultrasonic surgical apparatus. Oral surgeons have started using this newly developed ultrasonic bone scalpel when performing precise osteotomies. This technique is known as piezosurgery also and has found wide applications in dental and medical specialities. Numerous protocols including conventional diamond round bur, safe scraper, rongeur pliers and gouge shaped bone chisel have been proposed to simplify osteotomy and osteoplasty in bone surgery before this technique was created. But only ultrasonic osteotomy could reach high levels of precision and safety as compared to those instruments. Ultrasonic osteotomy device is characterized by piezoelectric ultrasonic vibrations of a frequency of 25-29 kHz and a range between 60-200 µm. Due to low-frequency ultrasonic vibration, the device precisely cuts bone without cutting soft tissue even in case of accidental contact. It also maintains a blood-free site during cutting bone. This technique was first used in sinus floor elevation for bony window osteotomy and sinus membrane elevation. After several years, this technique is utilized in more and more dental indications such as ridge expansion, maxillary sinus floor elevation, harvesting of bone graft, extraction and socket debridement for the installation of immediate-loaded implants, extraction of third molars, periodontal surgery, cyst exeresis and inferior alveolar nerve transposition. It's also reported in hand surgery, cranial osteoplasty and segmental maxillary Le Fort osteotomy. Ultrasonic osteotomy is different from the conventional techniques utilized and proved to be a safe and useful technique in bone surgery.

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Sembronio S, Albiero AM, Polini F, Robiony M, Politi M.

## Intraoral endoscopically assisted treatment of temporomandibular joint ankylosis: Preliminary report.

**EN** Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2007; 104(1):e7-e10.

### ⇒ Intraoral endoscopically assisted treatment of temporomandibular joint ankylosis: Preliminary report.

Temporomandibular joint (TMJ) ankylosis is characterized by the formation of a bony or fibrous mass that replaces the normal articulation. To avoid a possible re-ankylosis it is mandatory to perform a radical, complete resection of the bony/fibrous mass. We treated a patient affected by right temporomandibular joint ankylosis performing the osteotomy of the ankylosic mass through a preauricular and intraoral approach under endoscopic control. The use of a piezoelectric instrument, whose microvibrations are adjusted to target only mineralized tissue, reduced the risk of damage to soft tissues. Then a temporalis muscle and fascia flap were used as the interpositional material. Through the endoscope it was easy to check the medial aspect of the resection and suture the flap. At 1-year follow-up the patient had significantly increased maximal mouth opening. No evidence of relapse of the joint ankylosis was shown by radiological studies. Intraoral endoscopic assistance may be useful to make the removal of the ankylosic mass safer, and the anchorage of the temporalis muscle and fascia flap more accurate, reducing the risk of re-ankylosis.

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**IT** Minerva Stomatol. 2006 Apr;55(4):199-207.

### ⇒ Peri-implant tissue healing in implant placed in post-extraction sockets of maxillary premolars.

**Aim:** The aim of this study was the clinical and radiographic evaluation of peri-implant tissues healing associated with two-stage implants performed with one-stage surgery and placed into post-extraction sockets of maxillary premolars.

**Methods:** Ten natural tapered submerged titanium implants with double acid-etched surface were installed in post-extraction sockets of maxillary premolars showing no bone dehiscence or fenestrations. Implants were selected according to the greatest dimension compatible to vertical and horizontal diameters of the post-extraction sockets. Peri-implant defects showing more than 2 mm of distance between the marginal palatal bone and the implant margin, were not treated with a regenerative procedure but with a bone swaging technique by means of piezoelectric surgery. At 16 weeks osteointegration of implants was examined with a clinical and standardized radiographic evaluation.

**Results:** All patients showed good clinical healing without any complication. At 16 weeks all implants, either those with reduced marginal defect or those with modified crestal bone profile, were osteointegrated.

**Conclusion:** The use of two-stage implants performed with one-stage surgery in 4 walls post-extraction sockets of maxillary premolars, seems to allow osteointegration both in reduced peri-implant horizontal defects as well as in defects larger than 2 mm. Furthermore, these procedures seem to simplify the management of soft peri-implant tissues.

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... The piezoelectric scalpel: a delicate tool highly suitable for use in dentistry.

Bone surgery (osteotomy and osteoplasty) has gained in refinement through PIEZOSURGERY®. The principle is based on use of an ultrasonic scalpel with inserts connected to a handpiece identical with that of a scaler, allowing a very fine cut, up to 0.2 mm for the finest inserts. The effect on mineralized tissues is obtained by means of double oscillation of the inserts, a vertical amplitude of 20 to 60 µm and a horizontal amplitude of 60 to 200 µm, by using ultrasonic frequencies that can be modulated from 22 to 30 kHz and simultaneously combining these frequencies. In this frequency range, the second advantage is the absence of any action of the inserts on the soft tissues, which is of particular interest with regard to anatomical structures at risk (neurovascular bundles, mucosa). We present PIEZOSURGERY® by means of five clinical cases showing some of the possibilities of use of the system.

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**Schlee M, Steigmann M, Bratu E, Garg AK.**

## Piezosurgery: basics and possibilities.

**EN** Implant Dent. 2006; 15(4):334-340.

### ⇒ Piezosurgery: basics and possibilities.

Useful in a variety of oral surgery procedures, piezosurgery has therapeutic features that include a micrometric cut (precise and secure action to limit tissue damage, especially to osteocytes), a selective cut (affecting mineralized tissues, but not surrounding soft tissues), and a clear surgical site (the result of the cavitation effect created by an irrigation/cooling solution and oscillating tip). Because the instrument's tip vibrates at different ultrasonic frequencies, since hard and soft tissues are cut at different frequencies, a „selective cut“ enables the clinician to cut hard tissues while sparing fine anatomical structures (e.g., schneiderian membrane, nerve tissue). An oscillating tip drives the cooling-irrigation fluid, making it possible to obtain effective cooling as well as higher visibility (via cavitation effect) compared to conventional surgical instruments (rotating burs and oscillating saws), even in deep spaces. As a result, implantology surgical techniques such as bone harvesting (chips and blocks), crestal bone splitting, and sinus floor elevation can be performed with greater ease and safety.

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**Stübinger S, Filippi A, Sader R, Zeilhofer HF.**

## Intraoral Piezosurgery: Preliminary Results of a New Technique.

**EN** J Oral Maxillofac Surg. 2005; 63(9):1283-1287.

### ⇒ Intraoral Piezosurgery: Preliminary Results of a New Technique.

The piezosurgery instrument, developed in 1988, uses a modulated ultrasonic frequency that permits highly precise and safe cutting of hard tissue. Nerves, vessels, and soft tissue are not injured by the microvibrations (60 to 200 mm/sec), which are optimally adjusted to target only mineralized tissue. The selective and thermally harmless nature of the piezosurgery instrument results in a low bleeding tendency. In addition, the instrument can be used in operations requiring either local or general anaesthesia. The precise nature of the instrument allows exact, clean, and smooth cut geometries during surgery. The difference in time requirement for surgical procedures using the piezosurgery instrument in comparison with the conventional drill is negligible. Postoperatively, excellent wound healing, with no nerve and soft tissue injuries, is observed. It is apparent that the range of application of piezosurgery is not limited to minor operations. Because of its highly selective and accurate nature, with its cutting effect exclusively targeting hard tissue, its use may be extended to more complex oral surgery cases, as well as to other interdisciplinary problems.

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### → Piezosurgery - rudiments and possibilities.

This article describes the therapeutical potential of the Piezoelectric Bone Surgery in implantology. PIEZOSURGERY enables more precise and non traumatic cutting of bone in comparison to conventional methods (micrometric cut). The instrument vibrates with a modulated ultrasonic frequency. Because the vibration frequency of Piezosurgery is optimal for mineralized tissue it does not cut soft tissue. In fact, in order to cut soft tissue other ultrasonic frequencies are necessary. This allows an osteotomy which prevents delicate anatomic soft tissue structures like nerves or the Schneiderian Membrane of being damaged (selective cut). The cooling liquid flows laminar along the tip. Thereby, a more effective cooling also in deeper regions is achieved. Methods like crest expansion, sinus lift elevation, harvesting of bone chips and bone blocks are becoming much more secure and easier.

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### ⇒ Technological characteristics and clinical indications of the piezoelectric bone surgery.

Piezoelectric bone surgery, also known as Piezosurgery, is a new osteotomy and osteoplasty technique using an innovative piezoelectric surgery device. The technique in question was developed by the author in order to meet the need to achieve a higher level of precision and safety compared to manual and motorised instruments normally used in bone surgery. The main feature of the device which has enabled these results to be achieved, known as the Mectron Piezosurgery Device, is its piezoelectric vibration at a frequency of 29kHz and with an amplitude comprised between 60 and 200 µm. Thanks to their features, the microvibrations enable selective cutting of the mineralised structure without causing any damage to the soft tissues even if these are touched by accident. In addition, the micrometric vibrations enable great precision of cutting and at the same time they keep the surgical site clean (bloodless) thanks to the physical phenomenon of cavitation. The absence of macrovibrations makes it easy to grip the device in such a way as to enable greater intra-operative control, thus improving safety and precision of cutting. This means that it can also be used in anatomical areas that are difficult to access. Piezoelectric bone surgery is an innovative technique, and its use differs from traditional techniques. This means that the manual skills required are also different. They can be learnt according to a suitable learning curve.

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Troiani C, Russo C, Ballarani G, Vercellotti T.

## Piezoelectric Surgery: A new reality to cut and manage bone.

**EN** Maxillo Odontostomatologia - International Journal of Maxillo Odontostomatology - S.I.M.O. 2005; 4(1):23-28.

### ⇒ Piezoelectric Surgery: A new reality to cut and manage bone.

Piezosurgery is a new osteotomic and osteoplastic technique that uses ultrasound technology as its functioning principle: this system aims to improve the precision and safety of the traditional instruments normally used in bone surgery. Its special characteristics include: the selective cutting due to the use of low frequency ultrasound of 29.000 Hz (respecting soft tissues, vessels, nerves, mucosa); the precision of the cut guaranteed by the micro vibrations of the insert (40-200 µm); the cavitation phenomenon which allows a blood free and clean site, absence of macro vibrations that guarantees greater control and greater operative safety even in the most complex anatomical zones. The surgeon with minor experience in different surgical techniques will need a reasonable amount of practice to achieve a satisfactory learning curve.

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### ⇒ The ultrasonic cut in implantology; 2nd part: clinical applications

The ultrasonic lancet makes possible to realize with extreme precision and safety a great number of oral and maxillofacial surgeries. However, it is during four quite particular interventions that this tool will present really higher performances than the conventional instrumentation: non-traumatic removal of osteo-integrated implants, symphysar and retromolar bone harvesting and inferior alveolar nerve lateralization. The objective of this article is to present in detail each one of these four applications and to discuss their advantages and disadvantages compared to the former techniques. Indeed, the use of the piezoelectric lancet considerably simplifies these surgical protocols with non-traumatic cut and ultrasonic solid interfaces cleavage, which supports the bone fragments rupture.

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### ⇒ **The ultrasonic cut in implantology; 1st part: technologies**

The ultrasonic lancet is a surgical device able to cut out with precision hard tissue and to facilitate the solid interfaces cleavage. It uses microvibrations of intermediate frequency generated by a piezoelectric transducer and applied to titanium nitride hardened and/or diamond-coated inserts. With its vast range of inserts, it finds many applications in oral and maxillofacial surgery, such as non-traumatic dental avulsions, root planing and bone defect debridement, or cyst removal. It also proposes a simplified protocol for the sinus lift surgery. Lastly, it offers a true revolution in the bone grafting surgery by allowing precise and non-traumatic graft harvesting. This article presents its physical, technological and clinical aspects and discusses its most promising applications. Indeed, although the ultrasonic lancet remains a safe and powerful tool in many circumstances, it is however significant to define a reasoned application field to it.

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### ⇒ Piezoelectric surgery. An alternative method of minimally invasive surgery.

The piezoelectric surgery is an ultrasonic surgery which represents a novel and alternative method to the conventional hard and soft tissue management with rotating instruments in the oral district. The innovation of this technique is mainly due to three major characteristics: a) a micrometric cut with an amplitude between 60 µm and 200 µm, b) a selective cut which works on hard tissues but not on soft tissues, c) a relative blood free surgical field due to the air-water cavitation effect of the ultrasonic device. The final result of these characteristics is a definite clinical advantage with regard to the cut precision, the sparing of vital nervous soft tissues, the better visualisation of the surgical area. The present work shows through a series of case reports the possible use of the piezoelectric surgery in the oral district, presenting observed advantages and disadvantages. The entire clinical work was performed with the Piezosurgery device Mectron (Mectron Medical Technology, 16042 Carasco, Italy).

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### ⇒ Intraoral Piezo-Surgery.

The Piezo “effect” describes a physical interaction, which is used in different fields of clinical medicine. An ultrasound technology results in the power of a surgical cut technique, which is able to distinguish between hard and soft tissue. With the presented instruments it is possible to treat bone mechanically in a piezo-surgical manner without hurting the surrounding soft tissues when touching them.

Four examples:

- Tooth extraction before implantation under minimal grinding of bone,
- Gaining bone for peri-implant transplantation,
- Preparation of Sinus Lift,
- Exposing the inferior alveolar nerve

demonstrate the use of Piezo-Surgery in the intraoral area. This method shows a remarkable addition to the intraoral operation techniques.

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Vercellotti T.

## Technological characteristics and clinical indications of piezoelectric bone surgery.

**EN** Minerva Stomatol. 2004; 53(5):207-214

### ⇒ Technological characteristics and clinical indications of piezoelectric bone surgery.

Piezoelectric bone surgery, also simply known as piezosurgery, is a new technique for osteotomy and osteoplasty utilizing an innovative ultrasonic surgical apparatus. This technique was created and developed in response to the need to reach major levels of precision and safety in bone surgery, as compared to that available by the usual manual and motorized instruments. The instrument which offers these results, known as the Mectron Piezosurgery Device, is characterized by piezoelectric ultrasonic vibrations of a frequency of 29 kHz and a range between 60/200 Hz. Due to their characteristics, the micro-vibrations allow a selective cut of only mineralized structures without damaging soft tissues, which remain undamaged even in case of accidental contact. The micrometric vibration ensures precise cutting action and at the same time maintains a blood-free site because of the physical phenomenon of cavitation. The micrometric vibration makes the instrument manageable and permits major interoperative control with a consequent increase in safety especially in anatomically difficult areas. The absence of macrovibrations makes the instrument more manageable and allows greater interoperative control with a significant increase in the cutting safety in the more difficult anatomical cutting zones. Given its innovative nature, piezosurgery distinguishes itself as being different from the conventional techniques utilised in bone surgery. Consequently, it requires different surgical skill. To master the right surgical skill it is essential to demonstrate an adequate learning curve.

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⇒ **Piezoelectric surgery. An alternative to conventional bone-surgery techniques.**

Piezoelectric bone surgery, known as piezosurgery, is a technique that can be used to carry out osteotomies and osteoplasties thanks to the ultrasonic micro-vibrations of the inserts. This technique has been devised by Dr Vercellotti in order to broaden the indications of bone surgery, to simplify surgical protocols and improve the predictability of the results. At the moment, bone burs and saws are used for this type of surgery. The burs are operated by micromotors and the user must counter the reactive movements induced by the torque of the instrument, thus forcing him to make an additional effort. This lowers the surgeon's sensitivity, above all when structures with different degrees of mineralization are encountered.

Similarly, saws produce macrovibrations that require resistance to be exerted while the cut is being made. Although saws ensure perfect linearity, they do not enable control of the depth of the cut, and this means that it is necessary to use additional instruments such as hammers and chisels to complete the incision.

When using the piezosurgery technique, on the other hand, the effort required to make a cut is very slight. This means that greater precision is achieved, guaranteed by the microvibrations of the insert. The thickness of the cut obtained will be determined solely by the thickness of the insert used, and the surgeon's macromovements will be minimal. The frequencies produce an effective cutting action only on the mineralised tissues, which makes for safer bone surgery, even in the proximity of anatomical structures to be preserved, such as vascular and nerve masses or other soft tissues. It is thus possible to minimise the risk of lesions to these structures, considering that they are not sensitive to the frequencies used. This gives rise to a high-precision technique, with a selective cutting action.

A histological study has enabled the healing mechanisms of the tissues following piezoelectric cutting to be studied. This study, carried out on dogs, showed that it was possible to note the absence of necrosis traces on the surface of the cut, and the presence of nucleated osteocytes, indicative of the only slightly traumatic effect of this technique.

Before using piezoelectric surgery it is essential to acquire the technique and to become familiar with the instruments. The surgeon's movements differ from those typical of other bone surgery techniques, and the cut is generated by the micro-vibrations of the insert. It is necessary to find and acquire a perfect balance between the pressure exerted and the speed of movement of the insert.

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IT Il Dentista Moderno. 2003; 5:21-55

### → Piezoelectric bone surgery.

This article aims at illustrating the origin, development and clinical advantages of piezoelectric surgery, a new, revolutionary bone cutting technique. First, the features of traditional instruments are described highlighting their limitations in terms of intraoperative precision and safety. Piezoelectric surgery, however, goes beyond these limits thanks to the physical characteristics of the mechanical action leading to osteotomy. Piezoelectric cutting can be defined as micrometric and selective in nature for it is not visible to the human eye (the cutting mechanical action is obtained thanks to the modulated, ultrasonic micro-vibration of only 60-200 mm of amplitude); it is selective for piezoelectric ultrasounds used at low frequency allow to cut mineralized tissue with great effectiveness, whilst they are inert on soft tissue. Thanks to these and other physical features, piezoelectric cutting offers a number of clinical advantages, such as cutting precision, safe performance and excellent intraoperative visibility. A number of histological results are also presented from the wealth of multidisciplinary scientific research that has characterized the development and evolution of this young surgical technique from the start. Dental indications are presented, such as, for example, piezoelectric impacted or ankylosed root extraction techniques, and socket debridement for immediate post-extraction implant placement. Lastly, operative protocols for periodontal surgery, ridge expansion and maxillary sinus lift are also detailed.

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## → HISTOLOGICAL ADVANTAGES

Sohn DS, Lee JK, An KM, Shin HI.

## Histomorphometric evaluation of mineralized cancellous allograft in the maxillary sinus augmentation: a 4 case report.

EN Implant Dent. 2009 Apr;18(2):172-81.

### → Histomorphometric evaluation of mineralized cancellous allograft in the maxillary sinus augmentation: a 4 case report.

**Purpose:** The aim of this article is to report the histomorphometric analysis and clinical study of mineralized cancellous bone allografts in maxillary sinus augmentation procedures in 4 cases.

**Materials and Methods:** Sinus bone augmentation, using a piezoelectric device and mineralized cancellous allografts, was performed in 3 patients. After an average of about 8 months of healing, a trephine core sample was gained from the lateral

## ⇒ HISTOLOGICAL ADVANTAGES

windows area for histomorphometric analysis. The mineralized cancellous allografts demonstrated normal healing of hard tissues.

**Results:** Upon reentry for dental implant placement or secondary surgery, the graft materials showed hard bony structures within the maxillary sinus. Histomorphometric analysis of the cores revealed an average vital new bone content of 17.88%.

**Conclusion:** Mineralized cancellous allografts showed favorable bone regeneration histologically and was clinically useful in augmenting bone volume for implant placement in the sinus bone graft.

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Stübinger S, Goethe JW.

## Bone Healing After PIEZOSURGERY® and its Influence on Clinical Applications.

EN Journal of Oral and Maxillofacial Surgery 2007 Sep;65(9):e7-39.e8.

⇒ **Bone Healing After PIEZOSURGERY® and its Influence on Clinical Applications.**

**Statement:** Osteotomies are well-established procedures in oral and maxillofacial surgery for either corrective measures or for gaining access to subjacent hard and soft tissues. Yet surgical reconstruction of bony defects can often be very challenging as thin and fragile bony structures are especially prone to fracture due to bulky cutting tips or massive application of pressure by conventional mechanical instruments. By employing the new technique of PIEZOSURGERY® such complications can be overcome. It was the aim of this animal study to analyse the bone remodelling process after ultrasound osteotomy in the sheep tibia under functional loading after 8 and 12 weeks and to evaluate its impact on clinical applications.

**Materials and Methods:** In 12 sheep the PIEZOSURGERY® device (Mectron, Italy) was used to perform a midshaft osteotomy of the tibia diaphysis. The instrument operates with modulated ultrasound and thereby generates micromovements between 60 and 200 mm/sec. Physiological sodium chloride solution of approx. 4°C was used as an irrigant. For all osteotomies bone saw OT7 with parameters mode boosted burst c and pump 5 were used according to the manufacturer's recommendation. All osteotomies were performed by one surgeon in a latero-lateral movement as well as in a turning guidance around the bone.

**Method of Data Analysis:** After 2 and 3 months specimen were evaluated by radiographic and histologic analysis (toluidin-blue-staining, fuchsin-vitalstaining). A complete, semi-quantitative and histomorphometrical evaluation was undertaken including fluorochromatic stainings indicative for bone remodeling. Data were compared to results gained by analogous experiments performed with other osteotomy techniques.

**Results:** Piezoelectric osteotomy permitted a micrometric, selective cut and a clear surgical site by the cavitation effect created by irrigation/cooling solution and oscillating tip. No excessive bleeding was encountered. The evaluation of the PIEZOSURGERY® specimen proved 8 weeks as well as 12 weeks after surgery an ingrowth of vital bone-forming tissue into the osteotomy gap. The remodelling in the compact bone was undisturbed and the osteotomy gap was filled with new bone. Additional radiological findings supported these findings. The bone fragments were completely healed, the bone marrow cavity restored as well as the external callus formation was subsided undergoing piezoelectric surgery. In general bone healing was faster than known from conventional methods.

**Conclusion:** PIEZOSURGERY® definitely enhances handling of delicate structures in the oral and maxillofacial region. Concerning osteotomies of thin and fragile bones, application of ultrasound is assessed to be superior to other mechanical

## ⇒ HISTOLOGICAL ADVANTAGES

instruments, because of easy handling, efficient bone ablation and minimal accidental harm to adjacent soft tissue structures. As bone healing is not disturbed by the PIEZOSURGERY®, but even seems to be improved, this method will have a major influence on new minimally invasive bone surgery techniques with special regard to biomechanics.

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Gleizal A, Li S, Pialat JB, Béziat JL.

## Transcriptional expression of calvarial bone after treatment with low-intensity ultrasound: An in vitro study.

EN Ultrasound Med Biol. 2006; 32(10):1569-1574

### ⇒ Transcriptional expression of calvarial bone after treatment with low - intensity ultrasound: An in vitro study.

The capacity to reossify a calvarial bone defect is very small in mature animals and in humans greater than 2 y of age. The clinical treatment of injured tissue sites of bones by low-intensity pulsed ultrasound is widespread, but little is known about the precise effects of ultrasound on the fundamental processes that promote repair and regeneration. In this study, we used real-time polymerase chain reaction (RT-PCR) to investigate the expression of osteogenesis-associated genes after stimulation by low-intensity ultrasound in adult mouse osteoblast from the parietal calvaria. The gene associated with the Runx2 pathway had notably higher levels after 1, 2 and 3 days of stimulation. Therefore, low-intensity ultrasound seems to have an effect on the transcriptional gene expression of the calvarial bone in vitro.

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Vercellotti T, Nevins ML, Kim DM, Nevins M, Wada K, Schenk RK, Fiorellini JP.

### Osseous Response following Resective Therapy with a Piezosurgery®.

EN Int J Periodontics Restorative Dent. 2005; 25(6):543-549

⇒ **Osseous Response following Resective Therapy with a Piezosurgery®.**

A piezoelectric instrument vibrating in the ultrasonic frequency range was investigated for its potential use in periodontal resective therapy. The rate of postoperative wound healing (baseline and 14, 28, and 56 days after surgery) in a dog model following surgical ostectomy and osteoplasty was the marker used to compare the efficacy of this instrument (PS) with a commonly used carbide bur (CB) or a diamond bur (DB). The surgical sites treated by CB or DB lost bone, in comparison to baseline measurements, by the 14th day, while the surgical sites treated by PS revealed a gain in the bone level. By day 28, the surgical sites treated by all three instruments demonstrated an increased bone level and regeneration of cementum and periodontal ligament. However, by day 56, the surgical sites treated by CB or DB evidenced a loss of bone, versus a bone gain in the PS-treated sites. Thus, it appears that PS provided more favourable osseous repair and remodeling than CB or DB when surgical ostectomy and osteoplasty procedures were performed. Therefore, PS could be regarded as being efficacious for use in osseous surgery.

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### Effects of a new piezoelectric device on periosteal microcirculation after subperiosteal preparation.

**EN** Microvasc Res. 2014 Jul;94:114-8.

⇒ Effects of a new piezoelectric device on periosteal microcirculation after subperiosteal preparation.

**Introduction:** Subperiosteal preparation using a periosteal elevator leads to disturbances of local periosteal microcirculation. Soft-tissue damage can usually be considerably reduced using piezoelectric technology. For this reason, we investigated the effects of a novel piezoelectric device on local periosteal microcirculation and compared this approach with the conventional preparation of the periosteum using a periosteal elevator.

**Materials and Methods:** A total of 20 Lewis rats were randomly assigned to one of two groups. Subperiosteal preparation was performed using either a piezoelectric device or a periosteal elevator. Intravital microscopy was performed immediately after the procedure as well as three and eight days postoperatively. Statistical analysis of microcirculatory parameters was performed offline using analysis of variance (ANOVA) on ranks ( $p < 0.05$ ).

**Results:** At all time points investigated, intravital microscopy demonstrated significantly higher levels of periosteal perfusion in the group of rats that underwent piezosurgery than in the group of rats that underwent treatment with a periosteal elevator.

**Conclusion:** The use of a piezoelectric device for subperiosteal preparation is associated with better periosteal microcirculation than the use of a conventional periosteal elevator. As a result, piezoelectric devices can be expected to have a positive effect on bone metabolism.

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### Investigation of perfusion in osseous vessels in close vicinity to piezoelectric bone cutting.

**EN** Br J Oral Maxillofac Surg. 2012 Apr;50(3):251-5.

⇒ Investigation of perfusion in osseous vessels in close vicinity to piezoelectric bone cutting.

Cutting bones by piezosurgery leads to failure of perfusion at the site of the osteotomy, the cause of which cannot be identified immediately. Among other things the formation of vascular thrombi by the transmission of oscillations from the piezoelectric unit to the bone may be responsible. We used three output levels of oscillation that were predefined by the system. The outer cortical bone of the calvaria of rats ( $n=24$ ) was removed horizontally and the intraosseous vessels exposed at the surface of the osteotomy. The blood flow was then examined repeatedly using intravital fluorescence microscopy. To calculate the transmission of oscillations to the bone, the spatial oscillation frequency of each calvarium and the contact pressure during removal of bone in vitro ( $n=18$ ) were also examined. After removal of the bone there was constant blood

flow at all three levels of oscillation output. In no case did an individual vessel seem to be occluded. The excitation oscillation of the bone was established at 2000 Hz in all spatial directions, irrespective of the predefined oscillation output. The application of piezosurgery does not cause the formation of vascular thrombi in the bone. This probably results from the oscillation damping properties of bone.

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IT Minerva Stomatologica – 2014 Apr. vol. 63 suppl. 1 al n. 4 - 232.

→ Preliminary study of marginal bone loss of 156 titanium implants placed with ultrasonic implant site preparation with Piezosurgery® (UISP): four years follow-up.

**Aim:** marginal bone level is considered a criterion for implant success. The use of titanium implants may incur in a peri-implant infection and in a marked marginal bone loss. The aim of the present study was to evaluate implant survival and radiographic marginal bone loss around titanium dental implants placed with Ultrasonic Implant Site Preparation (UISP) after 4 years of function.

**Materials and Methods:** A total of 28 totally and partially edentulous patients were rehabilitated with 156 titanium implants (100 Speedy TiUniteTM, Nobel Biocare AB, Agrate Brianza (MI), Italy and 56 Osseotite®, Biomet 3I, Vicenza, Italy) and fixed prosthetic rehabilitations. All implants were placed with Ultrasonic Implant Site Preparation (UISP) with Piezosurgery® using two-stage procedure in native bone. A total of 17 single, 38 partial and 6 full arch rehabilitations were delivered. Marginal bone loss at mesial and distal aspects was measured on intra-oral paralleling digital radiographies, taken to the impression phase and after 4 years of function, as the difference between baseline (distance from implant/transfert junction [IT] to the first bone contact [To] and bone level after 4 years of function (distance from implant/abutment junction [IA] to the bone contact [T1]). The formula used was: Periimplant bone loss ( $To/T_1$ ) =  $IA/T_1 - IT/To$ . Marginal bone loss was related to different variables: gender, jaws, bone quality, implant surface, connection and diameter and type of rehabilitation. Statistical analysis was carried out using SAS 9.0 a P-value less than 0,05 was considered the threshold for statistical significance.

**Result:** No operative complications were reported. Three implants (1,92%) were removed at the second surgery stage, with a overall Survival Rate (SR) of 98,08% after four years of function. The mean marginal bone loss was  $0,52 \pm 0,33$  mm ( $0,51 \pm 0,35$  at mesial aspect and  $0,53 \pm 0,35$  mm at distal aspect). Maxilla ( $0,52 \pm 0,32$  mm) and mandible ( $0,52 \pm 0,35$  mm) was almost the same. Statistical analysis showed a significantly greater bone loss in male as compared with female both in maxilla and in mandible full arch than maxilla and in maxilla Osseotite® surface than mandible.

**Conclusions:** Ultrasonic Implant Site Preparation with Piezosurgery® (UISP) can be considered a valid technique for implant placement, showing a Survival Rate (SR) and a mean marginal bone loss comparable to drilling surgery technique.

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da Silva Neto UT, Joly JC, Gehrke SA.

## Clinical analysis of the stability of dental implants after preparation of the site by conventional drilling or piezosurgery.

EN Br J Oral Maxillofac Surg. 2014 Feb;52(2):149-53.

⇒ Clinical analysis of the stability of dental implants after preparation of the site by conventional drilling or piezosurgery.

We used resonance frequency analysis to evaluate the implant stability quotient (ISQ) of dental implants that were installed in sites prepared by either conventional drilling or piezoelectric tips. We studied 30 patients with bilateral edentulous areas in the maxillary premolar region who were randomised to have the implant inserted with conventional drilling, or with piezoelectric surgery. The stability of each implant was measured by resonance frequency analysis immediately after placement to assess the immediate stability (time 1) and again at 90 days (time 2) and 150 days (time 3). In the conventional group the mean (SD) ISQ for time 1 was 69.1 (6.1) (95% CI 52.4-77.3); for time 2, 70.7 (5.7) (95% CI 60.4-82.8); and for time 3, 71.7 (4.5) (95% CI 64.2-79.2). In the piezosurgery group the corresponding values were: 77.5 (4.6) (95% CI 71.1-84.3) for time 1, 77.0 (4.2) (95% CI 69.7-85.2) for time 2, and 79.1 (3.1) (95% CI 74.5-87.3) for time 3. The results showed significant increases in the ISQ values for the piezosurgery group at each time point ( $p=0.04$ ). The stability of implants placed using the piezoelectric method was greater than that of implants placed using the conventional technique.

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**Gandhi SA, Baker JA, Bairam L, Kim HI, Davis EL, Andreana S.**

## Primary stability comparison using piezoelectric or conventional implant site preparation systems in cancellous bone: a pilot study.

**EN** Implant Dent. 2014 Feb;23(1):79-84.

→ Primary stability comparison using piezoelectric or conventional implant site preparation systems in cancellous bone: a pilot study.

**Purpose:** This study compares implant primary stability achieved in cancellous bone after placement in piezoelectric prepared sites versus conventionally drilled sites.

**Materials and Methods:** Four bovine ribs were randomly assigned and placed in a water bath at 36.5 °C. Five sites per rib (total n = 20 sites) were prepared using piezoelectric system (test) or conventional drills (control) with twenty 10 x 3.6-mm Implantium implants placed. Using Osstell Mentor quantitative analysis, 5 resonance frequency analysis [implant stability quotient {ISQ}] values per implant were recorded at 5 locations for a total of 100 measurements.

**Results:** Independent t test analysis indicated significant difference in primary stabilities between groups:  $t (17) = 2.637$ ,  $P = 0.17$ , with equal variance assumption satisfied ( $P = 0.196$ ). Examination of means indicated a higher mean ISQ for piezoelectric than for conventional: 58.9 (+8.55) versus 49.2 (+7.33), respectively. Analysis of variance indicated a significant difference in mean ISQ value by rib. Tukey test indicated significantly higher ISQ values for rib A (test) than ribs B, C (control), and D (test).

**Conclusion:** Implant site preparation using the piezoelectric system gives higher implant primary stability in cancellous bone. However, variations in quality across bones may have affected the results.

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### Ultrasonic implant site preparation using piezosurgery: a multicenter case series study analyzing 3,579 implants with a 1- to 3-year follow-up.

EN Int J Periodontics Restorative Dent. 2014 Jan-Feb;34(1):11-8.

#### ⇒ Ultrasonic implant site preparation using piezosurgery: a multicenter case series study analyzing 3,579 implants with a 1- to 3-year follow-up.

This multicenter case series introduces an innovative ultrasonic implant site preparation (UISP) technique as an alternative to the use of traditional rotary instruments. A total of 3,579 implants were inserted in 1,885 subjects, and the sites were prepared using a specific ultrasonic device with a 1- to 3-year follow-up. No surgical complications related to the UISP protocol were reported for any of the implant sites. Seventy-eight implants (59 maxillary, 19 mandibular) failed within 5 months of insertion, for an overall osseointegration percentage of 97.82% (97.14% maxilla, 98.75% mandible). Three maxillary implants failed after 3 years of loading, with an overall implant survival rate of 97.74% (96.99% maxilla, 98.75% mandible).

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Bengazi F, Lang NP, Canciani E, Viganò P, Velez JU, Botticelli D.

### Osseointegration of implants with dendrimers surface characteristics installed conventionally or with Piezosurgery®. A comparative study in the dog.

EN Clin Oral Implants Res. 2014 Jan;25(1):10-5.

#### ⇒ Osseointegration of implants with dendrimers surface characteristics installed conventionally or with Piezosurgery®. A comparative study in the dog.

**Aim:** The first aim of the present experiment was to compare bone healing at implants installed in recipient sites prepared with conventional drills or a piezoelectric device. The second aim was to compare implant osseointegration onto surfaces with and without dendrimers coatings.

**Materials and Methods:** Six Beagles dogs were used in this study. Five implants with two different surfaces, three with a ZirTi® surface (zirconia sand blasted, acid etched), and two with a ZirTi®-modified surface with dendrimers of phosphoserine and polylysine were installed in the right side of the mandible. In the most anterior region (P2, P3), two recipient sites were prepared with drills, and one implant ZirTi® surface and one coated with dendrimers implants were installed at random. In the posterior region (P4 and M1), three recipient sites were randomly prepared: two sites with a Piezosurgery® instrument and one site with drill and two ZirTi® surface and one coated with dendrimers implants installed. Three months after the surgery, the animals were sacrificed for histological analysis.

**Results:** No complications occurred during the healing period. Three implants were found not integrated and were excluded from analysis. However, n = 6 was obtained. The distance IS-B at the buccal aspect was  $2.2 \pm 0.8$  and  $1.8 \pm 0.5$  mm, while IS-C was  $1.5 \pm 0.9$  and  $1.4 \pm 0.6$  mm at the Piezosurgery® and drill groups, respectively. Similar values were obtained between the dendrimers-coated and ZirTi® surface implants. The BIC% values were higher at the drill (72%) compared to the Piezosurgery® (67%) sites. The BIC% were also found to be higher at the ZirTi® (74%) compared to the dendrimers-coated

(65%) implants, the difference being statistically significant.

**Conclusion:** This study has revealed that oral implants may osseointegrate equally well irrespective of whether their bed was prepared utilizing conventional drills with abundant cooling or Piezosurgery®. Moreover, the surface coating of implants with dendrimers phosphoserine and polylysine did not improve osseointegration.

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Lamazza L, Laurito D, Lollobrigida M, Brugnoletti O, Garreffa G, De Biase A.

## Identification of possible factors influencing temperatures elevation during implant site preparation with piezoelectric technique.

EN Annali di Stomatologia 2014; V (4): 115-122.

⇒ **Identification of possible factors influencing temperatures elevation during implant site preparation with piezoelectric technique.**

**Background:** Overheating during implant site preparation negatively affects the osseointegration process as well the final outcome of implant rehabilitations. Piezoelectric techniques seem to provide to a gentle implant preparation although few scientific reports have investigated the heat generation and its underlying factors.

**Purpose:** To investigate, through a proper methodological approach, the main factors influencing temperature rise during piezoelectric implant site preparation.

**Materials and Methods:** Different piezoelectric tips (IM1s, IM2, P2-3, IM3, Mectron Medical Technology, Carasco, Italy) have been tested. The experimental set-up consisted in a mechanical positioning device equipped with a load cell and a fluoroptic thermometer.

**Results:** The first tip of the sequence (IM1s) generated the highest temperature increasing ( $\Delta T$ ). The diamond tips (IM1s and P2-3) determined higher  $\Delta T$  values than the smooth tips (IM2 and IM3). Further tests with IM1s suggested that the temperature elevation during the first thirty seconds may be predictive of the maximal temperature as well as of the overall thermal impact.

**Conclusions:** Working load, working movements management and irrigation resulted to be the main factors influencing temperature rise during piezoelectric implant site preparation. Cortical thickness, irrigant temperature and clogging effect may also synergically contribute to the heat generation.

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## Malpositioned osseointegrated implants relocated with segmental osteotomies: a retrospective analysis of a multicenter case series with a 1- to 15-year follow-up.

EN Clin Implant Dent Relat Res. 2013 Dec;15(6):836-46.

⇒ Malpositioned osseointegrated implants relocated with segmental osteotomies: a retrospective analysis of a multicenter case series with a 1- to 15-year follow-up.

**Purpose:** This multicenter case series evaluates retrospectively the clinical outcomes of malpositioned implants surgically relocated in a more convenient position by segmental osteotomies.

**Materials and Methods:** Authors who published, on indexed journals or books, works about malpositioned implant correction by segmental osteotomies were contacted. Five centers, out of 11 selected, accepted to participate in this study. The dental records of patients who underwent implant relocation procedures were reviewed. Implant survival rates were analyzed and a blinded assessor examined clinical photos and periapical radiographs to evaluate esthetic outcome (pink esthetic score [PES]) and changes of marginal bone level over time. Patients were requested to fill a verbal rating scale form about discomfort, compliance, and satisfaction related to the procedure.

**Results:** Fifteen malpositioned implants relocated by segmental osteotomies were followed for a period ranging from 1 to 15 years (mean  $6.0 \pm 3.9$  years). The overall implant survival rate from baseline to the last follow-up visit was 100%. The mean marginal bone loss was 0.36 mm at the 12-month follow-up visit and no relevant further changes were observed at the following examinations. Significant esthetic improvement was recorded at 1-year examination with PES evaluation ( $p < .0001$ ). Patients' feedback described this procedure as not excessively invasive and uncomfortable, reporting a high final satisfaction rate.

**Conclusions:** The present study suggests that implant relocation with segmental osteotomies could be an effective alternative method to correct the position of unrestorable malpositioned implants in a single-stage surgery.

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## Changes in implant stability using different site preparation techniques: twist drills versus piezosurgery. A single-blinded, randomized, controlled clinical trial.

**EN** Clin Implant Dent Relat Res. 2013 Apr;15(2):188-97.

⇒ Changes in implant stability using different site preparation techniques: twist drills versus piezosurgery.

A single-blinded, randomized, controlled clinical trial.

**Purpose:** The objective of the present investigation was to longitudinally monitor stability changes of implants inserted using traditional rotary instruments or piezoelectric inserts, and to follow their variations during the first 90 days of healing.

**Materials and Methods:** A randomized, controlled trial was conducted on 20 patients. Each patient received two identical, adjacent implants in the upper premolar area: the test site was prepared with piezosurgery, and the control site was prepared using twist drills. Resonance frequency analysis measurements were taken by a blinded operator on the day of surgery and after 7, 14, 21, 28, 42, 56, and 90 days.

**Results:** At 90 days, 39 out of 40 implants were osseointegrated (one failure in the control group). Both groups showed an initial decrease in mean implant stability quotient (ISQ) values: a shift in implant stability to increasing ISQ values occurred after 14 days in the test group and after 21 days in the control group. The lowest mean ISQ value was recorded at 14 days for test implants (97.3% of the primary stability) and at 21 days for the control implants (90.8% of the primary stability). ISQ variations with respect to primary stability differed significantly between the two groups during the entire period of observation: from day 14 to day 42, in particular, the differences were extremely significant ( $p < .0001$ ). All 39 implants were in function successfully at the visit scheduled 1 year after insertion.

**Conclusions:** The findings from this study suggest that ultrasonic implant site preparation results in a limited decrease of ISQ values and in an earlier shifting from a decreasing to an increasing stability pattern, when compared with the traditional drilling technique. From a clinical point of view, implants inserted with the piezoelectric technique demonstrated a short-term clinical success similar to those inserted using twist drills.

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Baker JA, Vora S, Bairam L, Kim HI, Davis EL, Andreana S.

## Piezoelectric vs. conventional implant site preparation: ex vivo implant primary stability.

EN Clin Oral Implants Res. 2012 Apr;23(4):433-7.

### ⇒ Piezoelectric vs. conventional implant site preparation: ex vivo implant primary stability.

**Objectives:** This study aims to determine differences in primary stability between implants placed in cortical bone following Piezoelectric or conventional site preparation, as assessed by resonance frequency analysis (RFA) and reverse torque testing (RTT).

**Materials and Methods:** Four fresh bovine ribs were acquired and surgical guides fabricated with five sites per rib ( $n = 20$ ), for proper site preparation. Implant sites were prepared via conventional drilling technique as per manufacturer's instruction (Implantium) or via Piezoelectric (Mectron) implant site preparation using the Implant Prep kit. Twenty 10 mm long, 3.6 mm diameter Implantium implants were placed with 35 Ncm torque; 10 implants per preparation method. RFA was assessed via the Osstell Mentor. Five values were taken per implant. All implants where subjected to a reverse torque in increasing increments of 5 until 50 Ncm force was reached.

**Results:** The five RFA values per site were averaged and plotted by placement technique. A paired t-test statistical analysis was run. The average RFA values showed no statistical significance between the 10 test (RFA =  $69.04 \pm 5.11$ ) and 10 control (RFA =  $70.94 \pm 6.41$ ) sites ( $P > 0.05$ ). All implants in both groups withstood RTT up to 50 Ncm force without movement and thus showed no statistical differences.

**Conclusion:** Results of this ex vivo study imply that the Piezoelectric implant site preparation affords similar primary implant stability in comparison to conventional rotary instrumentation in cortical bone.

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IT Implantologia 2011; 3:21-28.

⇒ Total rehabilitation of the jaw with the “all on four” protocol and piezoelectric preparation of the implant site.

**Purpose:** The purpose of this paper is to present a new surgical technique that, through the use of piezoelectric technology, allows an implant-prosthetic rehabilitation with immediate function on four implants “All on four” reducing the risk of injury to the inferior alveolar nerve.

**Materials and Methods:** The authors present a case of bilateral partial edentulism with severe bone resorption. The patient’s requirement to be rehabilitated with a fixed solution did opt for an extraction of all the residual teeth and a surgical protocol for immediate function with four implants, of which the distal inclined distal-mesial. The preparation of the implant sites was performed with special piezoelectric long inserts.

**Results:** The clinical and radiographic controls four months after prosthetic loading showed no signs of failed osseointegration.

**Conclusion:** Although this work considers a single case, the simplicity of operation and security related to the use of piezoelectric inserts for implant site preparation, can achieve immediate loading protocols with a reduced risk of nerve damage caused by accidental contact respect using traditional burs.

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## A comparative study of bone densitometry during osseointegration: piezoelectric surgery versus rotary protocols.

EN Quintessence Int. 2010 Sep;41(8):639-44.

⇒ A comparative study of bone densitometry during osseointegration: piezoelectric surgery versus rotary protocols.

**Objectives:** To date, there have been no studies on the outcome of osseointegration of alveolar bone around dental implants inserted with piezoelectric osteotomy versus conventional osteotomy. The aim of this study was to compare the radiographic differences, through evaluation of peri-implant bone density, between implant insertion using traditional surgical technique and piezoelectric technique.

**Method and Materials:** Forty patients were selected whose treatment consisted of a minimum of two implants placed in nonpathologic native bone. A single type of implant surface (SLA) was chosen. The implants were placed following the manufacturer protocol for traditional surgical technique and piezoelectric technique. Radiographs were taken following surgery and 30, 60, and 90 days after surgery. The bone density was studied with the densitometry application.

**Results and conclusion:** All patients completed the study period with success. Despite a limited number of treated

patients, the results of this pilot study demonstrated that (1) piezoelectric implant site preparation promotes better bone density and osteogenesis, and (2) the piezoelectric technique is predictable, with a 100% success rate in this study.

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## Relocation of a malpositioned maxillary implant with piezoelectric osteotomies: a case report.

EN Int J Periodontics Restorative Dent. 2008 Oct;28(5):489-95.

⇒ Relocation of a malpositioned maxillary implant with piezoelectric osteotomies: a case report.

Implant relocation is a new surgical technique for correcting the alignment of malpositioned implants by mobilizing them with the surrounding bone until the desired position is achieved. In this case report, a 25-year-old woman was treated for the malposition of an implant in the maxillary left canine site. The use of a piezoelectric scalpel permits narrow, precise, and safe osteotomies, thus preventing involvement of the soft tissue and producing better healing potential compared to burs or saws. The results suggest that inadequately axially inclined implants can be successfully reconfigured using segmental piezoelectric osteotomies.

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## Cytokines and Growth Factors Involved in the Osseointegration of Oral Titanium Implants Positioned using Piezoelectric Bone Surgery Versus a Drill Technique: A Pilot Study in Minipigs.

EN J Periodontol. 2007; 78(4):716-722.

⇒ Cytokines and Growth Factors Involved in the Osseointegration of Oral Titanium Implants Positioned using Piezoelectric Bone Surgery Versus a Drill Technique: A Pilot Study in Minipigs.

**Background:** Most dental implants are positioned using a drilling surgery technique. However, dentistry recently experienced the implementation of piezoelectric surgery. This technique was introduced to overcome some of the limitations involving rotating instruments in bone surgery. This study used biomolecular and histologic analyses to compare the osseointegration of porous implants positioned using traditional drills versus the piezoelectric bone surgery technique.

**Methods:** Porous titanium implants were inserted into minipig tibias. Histomorphology and levels of bone morphogenetic protein (BMP)-4, transforming growth factor (TGF)-beta2, tumor necrosis factor alpha, and interleukin-1beta and -10 were evaluated in the peri-implant osseous samples.

**Results:** Histomorphological analyses demonstrated that more inflammatory cells were present in samples from drilled sites. Also, neo-osteogenesis was consistently more active in bone samples from the implant sites that were prepared using piezoelectric bone surgery. Moreover, bone around the implants treated with the piezoelectric bone surgery technique showed an earlier increase in BMP-4 and TGF-beta2 proteins as well as a reduction in proinflammatory cytokines.

**Conclusion:** Piezoelectric bone surgery appears to be more efficient in the first phases of bone healing; it induced an earlier increase in BMPs, controlled the inflammatory process better, and stimulated bone remodeling as early as 56 days post-treatment.

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**Platelet Rich Fibrin (P.R.F.) in reconstructive surgery of atrophied maxillary bones: clinical and histological evaluations.**

EN Int J Med Sci. 2012;9(10):872-80.

⇨ Platelet Rich Fibrin (P.R.F.) in reconstructive surgery of atrophied maxillary bones: clinical and histological evaluations.

**Introduction:** Maxillary bone losses often require additional regenerative procedures: as a supplement to the procedures of tissue regeneration, a platelet concentrate called PRF (Platelet Rich Fibrin) was tested for the first time in France by Dr. Choukroun. Aim of the present study is to investigate, clinically and histologically, the potential use of PRF, associated with deproteinized bovine bone (Bio-Oss), as grafting materials in pre-implantology sinus grafting of severe maxillary atrophy, in comparison with a control group, in which only deproteinized bovine bone (Bio-Oss) was used as reconstructive material.

**Materials and Methods:** 60 patients were recruited using the cluster-sampling method; inclusion criteria were maxillary atrophy with residual ridge < 5mm. The major atrophies in selected patients involved sinus-lift, with a second-look reopening for the implant insertion phase. The used grafting materials were: a) Bio-Oss and b) amorphous and membranous PRF

## ⇒ SINUS LIFT TECHNIQUE

together with Bio-Oss. We performed all operations by means of piezosurgery in order to reduce trauma and to optimize the design of the operculum on the cortical bone. The reopening of the surgical area was scheduled at 3 different times.

**Results:** 72 sinus lifts were performed with subsequent implants insertions. We want to underline how the histological results proved that the samples collected after 106 days (Early protocol) with the adding of PRF were constituted by lamellar bone tissue with an interposed stroma that appeared relaxed and richly vascularized.

**Conclusion:** The use of PRF and piezosurgery reduced the healing time, compared to the 150 days described in literature, favoring optimal bone regeneration. At 106 days, it is already possible to achieve good primary stability of endosseous implants, though lacking of functional loading.

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**IT** Minerva Stomatol. 2010 Jun;59(6):349-54.

### ⇒ Intraosseous anastomosis in the maxillary sinus.

**Aim:** During sinus-lift surgery, certain intra-osseous vessels may be accidentally cut and this cause bleeding complications in approximately 20% of osteotomies. Therefore, understanding vascular details of the maxilla is very important for the surgeon. Here, we have given an anatomical overview of maxillary sinus vascularization through anatomical dissection. We have analyzed the distribution, localization and distance from the alveolar ridge of intraosseous branches of the maxillary artery found during sinus lift surgery.

**Methods:** Fifty-six maxillary bone doors were made bilaterally in twenty-eight unfixed cadavers; the doors were made between the first molar and the second molar (24 doors) or between the first and the second premolar (32 doors).

**Results:** Intraosseous arteries were found in 37 maxillary bones (66%). The average height of the artery from the alveolar crest was 13+/-3.2 mm in the distal doors and 18 +/- 6.1 mm in the mesial doors. Generally, the intraosseous maxillary branches ran caudo-rostrally; but in five maxillae, we found two parallel arteries, while in three cases the maxillary artery ran vertically. No differences were found between the left and right side.

**Conclusion:** The risk of vascular damage in sinus floor elevation surgery is a real problem for the oral surgeon. Detailed anatomical knowledge about sinus vascularization is very important to reduce the risk of vascular damage and bleeding. In addition the visualization of sinus anastomosis by radiology and less invasive surgery, such as piezosurgery, could be helpful.

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#### ⇒ Contribution of Piezosurgery to sinus elevation.

Piezoelectric surgery is a new technique providing an interesting alternative to traditional bone surgery techniques. The insert performs a selective cutting action on mineralised tissues but is ineffective on the soft tissues. It can be used for osteotomies and for osteoplastic surgery exploiting the microvibrations of the insert. Because of the small size of the insert, it is possible to obtain cuts of sub-millimetric width. Furthermore, the cooling of the insert by the physiological liquid contributes towards achieving a bloodless surgical field. The advantages of piezosurgery by comparison with the techniques currently used for sinus elevation are presented below.

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#### ⇒ Piezoelectric Bone Surgery for Sinus Bone Grafting.

Over the last half century, applications of ultrasonic instrumentation in bone surgery largely failed, particularly in light of reports of their substantial ineffectiveness or even negative influence on the healing process. Studies comparing different

types of instruments (ultrasonic instrument, bur, and scalpel) have obtained conflicting results based on a variety of experimental protocols. Several authors explored the possibility of using ultrasonics in maxillary sinus surgery, but the effort was curtailed because of insurmountable difficulties. Significant limiting factors include the failure of traditional ultrasonic instruments to cut bone structures that are more than 1 mm thick and a tendency for excessive overheating leading to necrosis. However, recent scientific-technologic research of bone's response to mechanical stimuli, particularly ultrasonic vibration, has led to a new and highly effective technique for osteotomy and osteoplasty called piezoelectric bone surgery. The revolutionary properties of piezo-electric surgery have simplified many common osseous surgical procedures, including sinus bone grafting.

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⇒ **Piezoelectric Surgery. Maxillary sinus lift techniques.**

Bone surgery, the one used for the maxillary sinus lift in Implantology in particular, has to follow two main purposes. The first purpose consists in the predictability of the results which are wished to be obtained by means of surgery so that a suitable therapeutic valence could be reached, according to the clinical needs of each patient. The second purpose consists in reaching the goals just described before with the less morbidity as possible during the post-operative course, both in the sinus region where bone grafting is performed and in the potential donor site from which the graft is harvested (Chanavaz 1990; Sinus Graft Consensus Conference 1996). The purpose of this chapter is to introduce the new Piezoelectric Surgery techniques dedicated to sinus lift pointing out the surgeon intraoperative safety and the patient lowered morbidity.

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## → SINUS LIFT TECHNIQUE BY CRESTAL APPROACH

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### Transcrestal maxillary sinus augmentation: Summers' versus a piezoelectric technique – an experimental cadaver study.

**EN** Clin. Oral Impl. Res. 00, 2015, 1–4. DOI: 10.1111/cir.12546.

→ Transcrestal maxillary sinus augmentation: Summers' versus a piezoelectric technique – an experimental cadaver study.

**Objectives:** Sinus floor augmentation using transalveolar techniques is a successful and predictable procedure. The aim of the study was to compare the performance of conventional hand instruments using mallets and osteotomes with that of piezoelectric-hydrodynamic devices for maxillary sinus floor elevation.

**Material and Methods:** In 17 undamaged cadaver heads on randomly allocated sites, Schneiderian membrane elevation was carried out transcrestally using piezosurgery and a hydrodynamic device or by conventional hand instrumentation. After simulation of sinus augmentation by the use of a radiopaque impression material, a post-operative CT scan was carried out and volumes were determined. Statistic significant differences between the two methods were evaluated by nonparametric Mann-Whitney U-test with  $P < 0.05$ .

**Results:** A mean graft volume of  $0.29 \pm 0.18 \text{ cm}^3$  ( $0.07$ – $0.60 \text{ cm}^3$ ) was measured for the Summers' technique compared to  $0.39 \pm 0.32 \text{ cm}^3$  ( $0.05$ – $1.04 \text{ cm}^3$ ) for the Sinus Physiolift® technique. There is no statistically significant difference with regard to trauma to the Schneiderian membrane or augmented volume.

**Conclusions:** Both techniques generate expedient augmentation volume in the posterior atrophic maxilla. The piezoelectric technique can be recommended as an alternative tool to graft the floor of human maxillary sinuses.

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Zhen F, Fang W, Jing S, Zuolin W.

### The use of a piezoelectric ultrasonic osteotome for internal sinus elevation: a retrospective analysis of clinical results.

**EN** Int J Oral Maxillofac Implants. 2012 Jul-Aug;27(4):920-6.

→ The use of a piezoelectric ultrasonic osteotome for internal sinus elevation: a retrospective analysis of clinical results.

**Purpose:** To explore the possibility and evaluate the clinical outcome of accomplishing maxillary internal sinus floor augmentation through the use of a piezoelectric osteotome in conjunction with dental implant placement and to discuss this technique.

**Material and Methods:** Patients attending the outpatient clinic of the Department of Dental Implantology, Hospital of Stomatology, Tongji University, between July 2007 and September 2009, who had insufficient bone volume to harbor endosseous implants at least 8 mm long in the lateral/posterior maxilla because of sinus pneumatization were enrolled in the study. Sinus augmentations were accomplished with a piezoelectric osteotome, followed by implant placement, either

## ⇒ SINUS LIFT TECHNIQUE BY CRESTAL APPROACH

immediately or delayed (6 months after augmentation, if the residual ridge height was less than 4 mm).

**Results:** Thirty patients with 36 maxillary molar sites with insufficient alveolar bone height as a result of pneumatization of the sinus were included in this study. The residual vertical bone height ranged from 2 to 8 mm. Twenty-eight implants were placed into 24 patients immediately after sinus augmentation. Another eight implants were placed into 6 patients 6 months after sinus augmentation. Only one sinus membrane perforated (failure rate: 2.78%). Only one implant was lost during the observation period. No other implant mobility or rapid bone loss was seen during a follow-up period of 5 to 27 months.

**Conclusions:** Application of a piezoelectric osteotome for internal sinus elevation simplified manipulation of the membrane and greatly reduced the chance of perforation. The pressure gradient between sinus and implant cavity was helpful in accomplishing this technique.

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Petrucci M, Ceccarelli R, Testori T, Grassi FR.

## Sinus floor augmentation with a hydropneumatic technique: a retrospective study in 40 patients.

EN Int J Periodontics Restorative Dent. 2012 Apr;32(2):205-10.

⇒ **Sinus floor augmentation with a hydropneumatic technique: a retrospective study in 40 patients.**

The use of a hydropneumatic balloon for the elevation of the sinus membrane is a new technique for sinus floor augmentation procedures. Few cases using such a technique are reported in the English medical literature. This report describes 40 patients who were treated with this technique and studied retrospectively. Forty consecutive patients with an alveolar crest-sinus floor distance (bone height)  $\leq 12$  mm were enrolled. Under microscopy (40x) and using piezosurgical instruments, hydropneumatic sinus membrane elevation was performed, and a calcium sulphate solution was injected under the elevated antral membrane using a syringe. In the same surgical session, 4.00- to 6.50-mm-diameter implants were placed. Bone height at 12 months, complications related to the surgical technique, and implant failure were all recorded. Bone height at 12 months was  $14.66 \pm 1.48$  mm, with a sinus membrane elevation of  $9.01 \pm 3.01$  mm. Fifty-six implants were placed, and no failures were observed after 1 year. One macrolaceration and two microlacerations were the only complications related to the technique. Minimal invasiveness and reduced trauma characterize this new approach. In fact, gradual balloon inflation provides a controlled and atraumatic preparation of the sinus floor membrane. Piezoelectric instruments and microscopy make this technique predictable and safe. The relatively short learning curve of this approach for sinus floor elevation allows for its use in private practice.

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**EN** Int J Oral Maxillofac Surg. 2011 May;40(5):497-503.

### ⇒ Sinus floor elevation using osteotomes or piezoelectric surgery.

The aim of this paper is to describe a technique for sinus floor augmentation with a 1-step crestal approach where the residual bone is  $\leq 7.5$  mm. 36 implants were installed in 25 patients in the atrophic posterior maxilla immediately after sinus floor elevation. Sinus floor elevation was performed with a crestal approach using either osteotomes and burs or piezosurgery. Standardized intraoral radiographs were taken prior to surgery and 1 year after surgery. The mean residual bone height was 5.61 mm (range 3-7.5 mm). The mean gain of sinus elevation was 6.78 mm (range 3.5-10 mm) at 1 year after surgery. Two patients dropped out of the study. Of the 23 patients completing the study, one implant failed, whilst the remaining 33 implants were stable 12 months after surgery (cumulative survival rate 97%). A statistically significantly higher bone height was achieved with tapered implants compared with cylindrical implants ( $P < 0.05$ ). No statistically significant differences were found in bone level using osteotomes or piezosurgery. Piezosurgery was considered to provide less discomfort for the patient and greater convenience for the surgeon.

### ⇒ Rialzo del pavimento del seno mascellare mediante osteotomi o chirurgia piezoelettrica.

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Sentineri R, Dagnino G.

## Sinus Physiolift: a new technique for a less invasive great sinus augmentation with crestal approach

EN J Osteol Biomat 2011;1:69-75.

→ Sinus Physiolift: a new technique for a less invasive great sinus augmentation with crestal approach.

**Aim:** The purpose of this article is to present an innovative surgical technique that produces a big maxillary sinus lift by the crestal approach through the use of hydrodynamic pressure for detaching the Schneiderian membrane.

**Materials and Methods:** Specific hollow elevators were designed, which due to their specific shape enable a closed system and exploits the Pascal's principle of the incompressibility of liquids. With a micrometric device the physiological liquid was injected into the sub-Schneiderian space in order to detach the membrane.

**Results:** This Sinus Physiolift technique uses piezoelectric surgery to reduce the percentage of perforations of the sinus membrane compared to traditional drills and osteotomes.

**Conclusions:** This technique significantly reduces the risk of perforation of the membrane when compared to previous methods, but the most important benefit is a much less debilitating postoperative phase for the patient. In addition, the simplicity of the procedure reduces stress for the patient as well as the surgeon's discomfort.

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## → SINUS LIFT TECHNIQUE BY LATERAL APPROACH

Stacchi C, Vercellotti T, Toschetti A, Speroni S, Salgarello S, Di Lenarda R.

## Intraoperative Complications during Sinus Floor Elevation Using Two Different Ultrasonic Approaches: A Two-Center, Randomized, Controlled Clinical Trial.

EN Clin Implant Dent Relat Res. 2013 Aug 22.

→ Intraoperative Complications during Sinus Floor Elevation Using Two Different Ultrasonic Approaches: A Two-Center, Randomized, Controlled Clinical Trial.

**Purpose:** The aim of this study was to assess the prevalence of intraoperative complications during maxillary sinus elevation with lateral approach using a piezoelectric device with two different surgical techniques.

**Material and Methods:** Antrostomies were randomly performed by outlining a window (group A, 36 patients) or by eroding the cortical wall with a grinding insert until the membrane was visible under a thin layer of bone, before outlining the window (group B, 36 patients). Occurrence of membrane perforation, laceration of vascular branches, and surgical time was recorded.

**Results:** Seventy-two patients underwent sinus floor elevation: four perforations (11.1%) were observed in group A (two occurred during elevation with hand instruments) and zero perforations in group B ( $p < .05$ ). No evidence of vascular

## → SINUS LIFT TECHNIQUE BY LATERAL APPROACH

lacerations was registered in both groups. A clinically insignificant but statistically shorter surgical time was recorded in group A ( $9.2 \pm 3.7$  minutes) than in group B ( $13.3 \pm 2.4$  minutes;  $p < .05$ ).

**Conclusions:** Within the limits of the present study, it may be concluded that ultrasonic erosion of the lateral wall of the sinus is a more predictable technique than piezoelectric outlining of a bone window in preventing from accidental perforations of Schneiderian membrane during sinus augmentation procedures.

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### Use of piezosurgery during maxillary sinus elevation: clinical results of 40 consecutive cases.

EN Int J Periodontics Restorative Dent. 2012 Dec;32(6):e182-8.

⇒ **Use of piezosurgery during maxillary sinus elevation: clinical results of 40 consecutive cases.**

The aim of this study was to evaluate the performance of piezoelectric devices during sinus elevation to determine the percentage of sinus membrane perforation and the time required to perform the antrostomy and elevation of the membrane. A total of 35 patients and 40 grafted sinuses were included. The parameters recorded were bony window length and height, bone thickness, osteotomy area, operative time, and number of perforations. Seven (17.5%) membrane perforations were observed, which were repaired with resorbable membranes. The mean length, height, and thickness of the osteotomy were  $13.8 \pm 2.9$  mm,  $6.9 \pm 1.4$  mm, and  $1.4 \pm 0.4$  mm, respectively. The mean osteotomy area was  $96.8 \pm 32.2$  mm<sup>2</sup>, and the mean operative time was  $10.3 \pm 2.1$  minutes. This study demonstrated that a piezoelectric device could be an attractive alternative for successful sinus augmentation.

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## ⇒ SINUS LIFT TECHNIQUE BY LATERAL APPROACH

Cortes AR, Cortes DN, Arita ES.

### Effectiveness of piezoelectric surgery in preparing the lateral window for maxillary sinus augmentation in patients with sinus anatomical variations: a case series.

EN Int J Oral Maxillofac Implants. 2012 Sep-Oct;27(5):1211-5.

⇒ Effectiveness of piezoelectric surgery in preparing the lateral window for maxillary sinus augmentation in patients with sinus anatomical variations: a case series.

**Purpose:** The present article sought to evaluate the effectiveness of a piezoelectric surgical unit for maxillary sinus augmentation surgeries in avoiding perforation of the sinus membrane and other possible procedural complications in patients with anatomical variations of the sinus.

**Material and Methods:** Twenty-five patients presenting sinus anatomical variations, who were indicated for a total of 40 sinus grafting procedures performed by the lateral window approach with a piezoelectric device, were analyzed. After 6 months of healing, implants were placed. Information collected included clinical and computed tomographic information on anatomical variations in the sinus bone walls, in the size of the sinus, and in the thickness of the sinus membrane. Occurrence of sinus membrane perforation and computed tomographic measurements of the amount of bone height gained with the grafting procedures were also recorded.

**Results:** Only two patients presented a small perforation (less than 5 mm in diameter) of the sinus membrane, which occurred only after osteotomies of the lateral windows and did not compromise the surgical outcome. No implants were lost during a mean follow-up period of 19 months.

**Conclusion:** The use of piezoelectric surgery allowed for the accomplishment of all rehabilitation treatments within the follow-up period of this study.

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### Maxillary sinus elevation by lateral window approach: evolution of technology and technique.

EN J Evid Based Dent Pract. 2012 Sep;12(3 Suppl):161-71.

⇒ Maxillary sinus elevation by lateral window approach: evolution of technology and technique.

**Context:** The maxillary sinus elevation procedure has become an important pre-prosthetic surgical procedure for the creation of bone volume in the edentulous posterior maxilla for the placement of dental implants. Research and clinical

## → SINUS LIFT TECHNIQUE BY LATERAL APPROACH

experience over the past 30 years has increased the predictability of this procedure as well as reduced patient morbidity.

**Evidence acquisition:** Data on grafting materials and implant survival rates comes from 10 published evidence-based reviews that include all relevant published data from 1980 to 2012. Supporting clinical material comes from the experience of the authors.

**Evidence synthesis:** The evidence-based reviews report and compare the implant survival rates utilizing various grafting materials, implant surfaces, and the use or non-use of barrier membranes over the lateral window. Clinical studies report on complication rates utilizing piezoelectric surgery and compare them to complication rates with rotary instrumentation.

**Conclusions:** The conclusions of all the evidence-based reviews indicate that the utilization of bone replacement grafts, rough-surfaced implants, and barrier membranes result in the most positive outcomes when considering implant survival. Further, the utilization of piezoelectric surgery, rather than rotary diamond burs, for lateral window preparation and membrane separation leads to a dramatic reduction in the occurrence of the intraoperative complications of bleeding and membrane perforation.

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Toscano NJ, Holtzclaw D, Rosen PS.

### The effect of piezoelectric use on open sinus lift perforation: a retrospective evaluation of 56 consecutively treated cases from private practices.

EN J Periodontol. 2010 Jan;81(1):167-71.

#### → The effect of piezoelectric use on open sinus lift perforation: a retrospective evaluation of 56 consecutively treated cases from private practices.

**Background:** The lateral window approach to maxillary sinus augmentation is a well-accepted treatment option in implant dentistry. The most frequent complication reported with traditional techniques has been the perforation of the Schneiderian membrane, with perforation rates ranging from 11% to 56%. The purpose of this retrospective, consecutive case series from two private practices was to report on the rate of Schneiderian membrane perforations and arterial lacerations when a piezoelectric surgical unit was used in conjunction with hand instrumentation to perform lateral window sinus elevations.

**Methods:** Clinical data (Schneiderian membrane perforation, Underwood septa, and laceration of the lateral arterial blood supply to the maxillary sinus) were obtained retrospectively from two private practices and pooled for analysis. The information was collated after an exhaustive chart review. Fifty-six consecutively treated lateral window sinus lifts were performed on 50 partially or completely edentate patients.

**Results:** Zero perforations of the Schneiderian membrane occurred during the piezoelectric preparation of the lateral antrostomies, whereas two perforations were noted during subsequent membrane elevations using hand instrumentation. In both instances, membrane perforations were associated with sinus septa. The overall sinus perforation rate was 3.6%. Arterial branches of the posterior superior alveolar artery were encountered in 35 cases, and there were zero instances of arterial laceration.

## ⇒ SINUS LIFT TECHNIQUE BY LATERAL APPROACH

**Conclusions:** This retrospective case series from clinical private practices confirmed that a lateral window approach to sinus elevation incorporating piezoelectric technology in conjunction with hand instrumentation was an effective means to achieve sinus elevation while minimizing the potential for intraoperative complications. Further prospective and randomized controlled studies are warranted to qualify these observations.

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## Palatal versus vestibular piezoelectric window osteotomy for maxillary sinus elevation: a comparative clinical study of two surgical techniques.

EN Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2009 May;107(5):648-55.

⇒ **Palatal versus vestibular piezoelectric window osteotomy for maxillary sinus elevation: a comparative clinical study of two surgical techniques.**

**Objectives:** The goal of this study was to compare the surgical advantages and disadvantages of a new palatal access osteotomy for sinus elevation with a conventional lateral approach.

**Study design:** In 32 patients, either a palatal (n = 16) or a lateral (n = 16) osteotomy to the maxillary sinus was performed under local anesthesia. The palatal access included a circular paramarginal incision and elevation of a palatal mucosal flap based on a median pedicle. The lateral access was performed by vestibular standard incision and development of a mucoperiosteal flap with a vestibular and superior basis. For all osteotomies a piezoelectric device was used. The sinus cavity was augmented with synthetic nanostructured hydroxyapatite graft material.

**Results:** Intraoperative complications during both procedures were minimal and wound healing was uneventful. Membrane perforation occurred in 19% of the palatal group and in 19% of the lateral group. Soft tissue management of the palatal technique was superior to that of the lateral approach, because the vestibular anatomy was not altered and consequently no disharmonious soft tissue scarring and no postoperative swelling occurred.

**Conclusion:** The palatal approach permitted higher postoperative comfort, especially for edentulous patients, because full dentures could be incorporated directly after surgery with almost perfect fit.

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**Muñoz-Guerra MF, Naval-Gías L, Capote-Moreno A.**

### Le Fort I osteotomy, bilateral sinus lift, and inlay bone-grafting for reconstruction in the severely atrophic maxilla: a new vision of the sandwich technique, using bone scrapers and PIEZOSURGERY®.

**EN** J Oral Maxillofac Surg. 2009 Mar;67(3):613-8.

⇒ **Le Fort I osteotomy, bilateral sinus lift, and inlay bone-grafting for reconstruction in the severely atrophic maxilla: a new vision of the sandwich technique, using bone scrapers and PIEZOSURGERY®.**

Severe atrophy of the edentulous maxilla and progressive pneumatisation of the maxillary sinus can compromise the insertion of dental implants. In this context, ideal implant positioning is limited by inadequate height, width, and quality of the bone. Le Fort I osteotomy and interpositional bone graft is an excellent treatment concept for the dental rehabilitation of patients with atrophied maxilla and reversed intermaxillary relationship. In this report, we indicate the transcendent aspect of elevation and preservation of maxillary sinus and nasal mucosa, modifying the sandwich technique by the use of bone scrapers and PIEZOSURGERY®. The procedure is described including a 1-stage approach using cortico-cancellous bone blocks through which implants are placed. In the extremely atrophied alveolar process of the maxilla, this technique provides the desired gain of bone, allows for the ideal placement of dental implants, and improves any discrepancy between the upper and lower arches.

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## ⇒ SINUS LIFT TECHNIQUE BY LATERAL APPROACH

Wallace SS, Mazor Z, Froum SJ, Cho SC, Tarnow DP.

### Schneiderian membrane perforation rate during sinus elevation using piezosurgery: clinical results of 100 consecutive cases.

EN Int J Periodontics Restorative Dent. 2007; 27(5):413-419

⇒ Schneiderian membrane perforation rate during sinus elevation using piezosurgery: clinical results of 100 consecutive cases.  
The lateral window sinus elevation procedure has become a routine and highly successful preprosthetic procedure that is used to increase bone volume in the posterior maxilla for the placement of dental implants. Many surgical techniques have been proposed that provide access to the maxillary sinus through the lateral wall to allow for elevation of the sinus membrane. Among these are the multiple variations of the hinge and complete osteotomy techniques, which make use of rotary cutting instruments for the antrostomy. The most common intraoperative complication with these surgical approaches is perforation of the schneiderian membrane, with perforation rates of 14% to 56% reported in the literature. In most instances, perforation occurs either while using rotary instruments to make the window or when using hand instruments to gain initial access to begin the elevation of the membrane from the sinus walls. This article presents an alternative approach that uses a piezoelectric instrument for the sinus elevation procedure. Although new to the United States, this approach has been used successfully in Europe for many years. The membrane perforation rate in this series of 100 consecutive cases using the piezoelectric technique has been reduced from the average reported rate of 30% with rotary instrumentation to 7%. Furthermore, all perforations with the piezoelectric technique occurred during the hand instrumentation phase and not with the piezoelectric inserts.

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Vercellotti T, De Paoli S, Nevins M.

### The Piezoelectric Bony Window Osteotomy and Sinus Membrane Elevation: Introduction of a New Technique for Simplification of the Sinus Augmentation Procedure.

EN Int J Periodontics Restorative Dent. 2001; 21(6):561-567.

⇒ The Piezoelectric Bony Window Osteotomy and Sinus Membrane Elevation: Introduction of a New Technique for Simplification of the Sinus Augmentation Procedure.

All of the surgical techniques to elevate the maxillary sinus present the possibility of perforating the schneiderian membrane. This complication can occur during the osteotomy, which is performed with burs, or during the elevation of the membrane using manual elevators. The purpose of this article is to present a new surgical technique that radically simplifies maxillary sinus surgery, thus avoiding perforating the membrane. The piezoelectric bony window osteotomy easily cuts

mineralized tissue without damaging the soft tissue, and the piezoelectric sinus membrane elevation separates the schneiderian membrane without causing perforations. The elevation of the membrane from the sinus floor is performed using both piezoelectric elevators and the force of a physiologic solution subjected to piezoelectric cavitation. Twenty-one piezoelectric bony window osteotomy and piezoelectric sinus membrane elevations were performed on 15 patients using the appropriate surgical device (Mectron Piezosurgery System). Only one perforation occurred during the osteotomy at the site of an underwood septa, resulting in a 95% success rate. The average length of the window was 14 mm; its height was 6 mm, and its thickness was 1.4 mm. The average time necessary for the piezoelectric bony window osteotomy was approximately 3 minutes, while the piezoelectric sinus membrane elevation required approximately 5 minutes.

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## → RIDGE EXPANSION

Brugnami F, Caiazzo A, Mehra P.

### Piezosurgery assisted, flapless split crest surgery for implant site preparation.

EN J Maxillofac Oral Surg. 2014 Mar;13(1):67-72.

→ Piezosurgery-assisted, flapless split crest surgery for implant site preparation.

**Introduction:** Bucco-lingual resorption of the alveolar ridge can, at times, be predictably corrected at the time of implant placement. Among the different options available to achieve this are a group of surgical techniques described as split crest or split ridge procedures. Most of these procedures require the use of a mallet and some type of chisels and/or osteotomes; they are very technique-sensitive and can be uncomfortable for patients. Recently, alternative tools to split the crest have been presented, and these include the newer bone expanders and the piezoelectric scalpel. A flapless approach to implant dentistry has become popular with the aim to alleviate post treatment side effects, accelerate healing and avoid bone resorption caused by flap elevation.

**Methods:** We present a technique combining the use of a piezoelectric scalpel and a tapered bone expander in a flapless fashion as a novel way to perform split crest procedures with an aim to optimize outcomes and acceptability by patients.

**Results:** All implants were successfully placed and the resorbed ridge expanded in the same setting. Findings were confirmed by postoperative cone beam cat scan (CBCT) evaluation.

**Conclusions:** This new technique is a predictable approach for split crest procedures and has high acceptability by patients and is technically simple for surgeons.

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**Kelly A, Flanagan D.**

## Ridge expansion and immediate placement with piezosurgery and screw expanders in atrophic maxillary sites: two case reports.

**EN** J Oral Implantol. 2013 Feb;39(1):85-90.

### ⇒ Ridge expansion and immediate placement with piezosurgery and screw expanders in atrophic maxillary sites: two case reports.

Endosseous dental implants may require bone augmentation before implant placement. Herein is described an approach to edentulous ridge expansion with the use of piezosurgery and immediate placement of implants. This may allow for a shortened treatment time and the elimination of donor-site morbidity. Two cases are reported. This technique uses a piezoelectric device to cut the crestal and proximal facial cortices. Space is then created with motorized osteotomes to widen the split ridge. This technique allows for expansion of narrow, anatomically limiting, atrophic ridges, creating space for immediate implant placement. The facial and lingual cortices provide support with vital osteocytes for osteogenesis. The 2 patients presented had adequate bone height for implant placement but narrow edentulous ridges. In patient 1 at site #11, the ridge crest was 3.12 mm thick and was expanded to accept a 4.3 mm x 13 mm implant. The resulting ridge width was 8.88 mm, which was verified using cone beam computerized tomography (CBCT). In patient 2 at site #8 and site #9, the narrow ridge was expanded using the same technique to accept 2 adjacent 3.5 mm x 14 mm implants. The implants were restored to a functional and esthetic outcome..

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Anitua E, Begoña L, Orive G.

## Controlled ridge expansion using a two-stage split-crest technique with ultrasonic bone surgery.

EN Implant Dent. 2012 Jun;21(3):163-70.

→ Controlled ridge expansion using a two-stage split-crest technique with ultrasonic bone surgery.

**Purpose:** The aim of this study is to evaluate the two-stage split-crest technique with ultrasonic bone surgery for implant placement in patients with very narrow ridges, and to determine the status of soft and hard tissues and implant success rate.

**Material and Methods:** Six patients received nine implants (BTI Biotechnology Institute, Vitoria, Spain) after a two-stage split-crest technique. Plasma rich in growth factors-Endoret was used during first and second stages to promote tissue regeneration. All implants were treated with plasma rich in growth factors-Endoret to promote osseointegration. Patients were recalled for clinical evaluation at least 6 months after implant loading.

**Results:** The status of soft tissues was very good showing adequate plaque index, bleeding index, and probing depth values. Success rate of implants at the end of follow-up (between 15 and 22 months after first surgery) was 100%. Bone ridge was measured and compared at final examination showing a mean ridge expansion of 5.60 mm (SD = 1.9) at apical and 7.33 mm (SD = 1.73) at occlusal.

**Conclusions:** Two-stage split-crest with ultrasonic bone surgery can be considered a safe and predictable bone expansion technique for very narrow ridges.

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Horrocks GB.

## The controlled assisted ridge expansion technique for implant placement in the anterior maxilla: a technical note.

EN Int J Periodontics Restorative Dent. 2010 Oct;30(5):495-501.

→ The controlled assisted ridge expansion technique for implant placement in the anterior maxilla: a technical note.

This technical report describes the simple and valuable technique of controlled assisted ridge expansion (CARE) for the esthetic zone of the anterior maxilla. Atrophic bone and narrow alveolar crests make implant bed preparation difficult. The highly predictable surgical technique of CARE for alveolar expansion uses periosteal containment and subperiosteal intrabony vertical releasing incisions efficiently. This allows for the expansion of narrow, anatomically limiting edentulous ridges. The precise and conservative use of Piezosurgery and assisted manipulation of the labialized bony pedicle by horizontal spreaders and rotary osteocondensers affords a highly predictable surgical procedure, which allows implants to be placed at the time of surgery, reducing patient treatment time, morbidity, functional losses, and overall cost, as well as increasing the native bone volume in the atrophic maxillary ridge.

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→ The use of tapered implants simultaneously with the piezoelectric alveolar expansion as an alternative to preprosthetic surgery with bone grafts for the rehabilitation of atrophic maxilla.

**Objective:** To evaluate the possibility of the combined technique of piezoelectric expansion with the use of tapered implants for the rehabilitation of atrophic maxilla.

**Methods:** We considered 10 patients (aged 45-58) suffering from class 4-5 atrophy of the upper jaw, according to the classification of Cawood and Howell. The patients were subjected to alveolar expansion using the piezoelectric technique and simultaneous implant surgery with the use of tapered implants. The postoperative radiological tests were a standard orthopanoramic and a TC dental scan. Standard antibiotic therapy was administered to each patient. The piezoelectric bone expansion procedure has been associated to the use of alloplastic bone substitutes to fill spaces in the bone and collagen membranes in order to improve the thickness of the alveolar bone and to allow for a correct implant. Each patient was subjected to the insertion of between 3 to 6 tapered implants with a diameter of 3.75 and a length of 10-11.5 mm. The healing screws were inserted 8 months after the primary procedure.

**Results:** the patients were monitored clinically and the changes of the base of the alveolar bone were evaluated by means of intraoral radiography with centring device and millimetre grid measurements, at the time of the primary surgery and subsequent checks up to a year after the procedure. The average survival of implants was 96.5% (n=65). The average resorption of the base of the alveolar bone was 1.19 (1.01) mm one year after the procedure. The use of tapered implants in procedures of alveolar expansion seems to be indicated for the feature of reducing the compressive forces between the fixture and both the cortical and trabecular bone.

**Discussion:** the use of tapered implants proved to be ideal in patients suffering from severe maxillary atrophy. The conical design of the fixture allows for the reduction of the interfacial compressive forces between the implant and the bone, thereby improving primary and secondary stability when used simultaneously with alveolar expansion procedures.

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**Belleggia F, Pozzi A, Rocci M, Barlattani A, Gargari M.**

## Piezoelectric surgery in mandibular split crest technique with immediate implant placement: a case report.

**EN** Oral Implantol (Rome). 2008 Oct;1(3-4):116-23.

### ⇒ Piezoelectric surgery in mandibular split crest technique with immediate implant placement: a case report.

Implant supported rehabilitation of thin edentulous ridges with horizontal atrophy necessitates a regenerative approach. Within the procedures for horizontal bone defects augmentation, ridge expansion techniques permit dislocation of the buccal bone plate in a labial direction and simultaneous implant insertion in single-stage surgery, abbreviating overall treatment time. The piezoelectric ridge expansion technique permits to obtain the expansion of very mineralized bone crests without excessive traumas or the risk of ridge fractures. The case reported shows an implant treatment for partial edentulous lower arch rehabilitation. A full-split thickness flap was raised. The mucoperiosteal reflection permitted to identify alveolar crest contour where osteotomies had to be performed. Split thickness dissection allowed periosteal blood supply to be maintained on the buccal bone plate. After horizontal and vertical osteotomies were performed with OT7 piezoelectric microsaw (Piezosurgery, Mectron), a single-bevel scalpel was used to move the buccal bone plate to the labial. Two Straumann TE 3,3/4,8 mm wide implants were inserted in the lower right premolar area, and 1 Straumann 4,8 mm Wide Neck implant was inserted to replace lower right first molar. The amount of bone expansion was equal to the cervical diameter of the placed implants (4,8 mm) and residual bone gap was packed with Bio-Oss granules (Geistlich). Healing was uneventful and 3 months later final restorations with implant-supported porcelain-fused-to-metal crowns were cemented.

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Enislidis G, Wittwer G, Ewers R.

## Preliminary Report on a Staged Ridge Splitting Technique for Implant Placement in the Mandible: A Technical Note.

EN Int J Oral Maxillofac Implants. 2006; 21(3):445-449

→ Preliminary Report on a Staged Ridge Splitting Technique for Implant Placement in the Mandible: A Technical Note.

**Purpose:** Narrow edentulous alveolar ridges less than 5 mm wide require horizontal augmentation for the placement of screw-type dental implants. A staged approach to ridge splitting in the mandible to decrease the risk of malfracture during osteotomy is presented.

**Materials and Methods:** Five consecutive patients with 6 long-span edentulous areas of the mandibular ridge were included in this study. After corticotomy of a rectangular buccal segment and a 40-day healing period, the mandibular ridge was split, leaving the buccal periosteum attached to the lateralized segment. Seventeen dental implants were placed, and the gap between the implants and the bone filled with a mixture of venous blood and a porous algae-derived hydroxyapatite.

**Results:** All buccal segments fractured as planned at the basal corticotomy during ridge splitting. After 6 months, all implants were stable and surrounded by bone; prosthetic loading with fixed partial dentures was successful in all cases.

**Discussion:** In the mandible, greenstick fracture during widening with osteotomes has not been controllable to date because of cortical thickness of the bone; the risk of malfracture during single-stage ridge splitting was high. With this approach, the location of the greenstick fracture is predetermined, and the perfusion for the buccal segment remains intact, although vascularization shifts from internal perfusion from spongy bone after the first intervention to external perfusion from the periosteum after the second intervention. The buccal cortical segment remains a pedicled graft after ridge splitting.

**Conclusion:** The preliminary results of this report indicate that staged ridge splitting can be a safe technique which overcomes the problems associated with single-stage ridge expansion/ridge splitting procedures without causing significant delay in treatment.

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Vercellotti T, Russo C, Gianotti S.

## A New Piezoelectric Ridge Expansion Technique in the Lower Arch - A Case Report.

**EN** World Dentistry. 2000; online article

### ⇒ A New Piezoelectric Ridge Expansion Technique in the Lower Arch - A Case Report.

This article presents a new piezoelectric ridge expansion technique that permits the use of the bone flap procedure even in areas of highly mineralized bone, such as in the anterior region of the mandible, thus avoiding the risk of fracture. The new piezoelectric ridge expansion technique makes use of a specifically designed piezoelectric device with bone scalpels that vibrate at a variable ultrasound frequency, cutting bone easily, precisely, and non-traumatically. The case reported is that of a patient with an edentulous mandibular ridge defect with thickness reduced to 3 mm. The bone flap expansion technique was performed, followed by a piezoelectric ridge osteotomy to cut the ridge and obtain separation of the buccal and palatal osseous flaps, allowing the positioning of the implants in a single stage surgery. The defect obtained by the distraction of the two cortical walls was filled with bioactive glass mixed with platelet concentrate (PC) and covered with a platelet rich plasma (PRP) membrane. When the site was reopened after six months, the ridge bone appeared to be well stabilized with a thickness of 5 mm, an increase of 2 mm; the implants also appeared to be osteointegrated.

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Vercellotti T.

## Piezoelectric Surgery in Implantology: A Case Report - A New Piezoelectric Ridge Expansion Technique.

**EN** Int J Periodontics Restorative Dent. 2000; 20(4):359-365

### ⇒ Piezoelectric Surgery in Implantology: A Case Report - A New Piezoelectric Ridge Expansion Technique.

The purpose of this preliminary article is to present a new surgical technique that, thanks to the use of modulated-frequency piezoelectric energy scalpels, permits the expansion of the ridge and the placement of implants in single-stage surgery in positions that were not previously possible with any other method. The technique involves the separation of the vestibular osseous flap from the palatal flap and the immediate positioning of the implant between the 2 cortical walls. The case report illustrates the ridge expansion and positioning of implants step by step in bone of quality 1 to 2 with only 2 to 3 mm of thickness that is maintained for its entire height. To obtain rapid healing, the expansion space that was created for the positioning of the implant was filled, following the concepts of tissue engineering, with bioactive glass synthetic bone graft material as an osteoconductive factor and autogenous platelet-rich plasma as an osteoinductive factor. The site was covered with a platelet-rich plasma membrane. A careful evaluation of the site when reopened after 3 months revealed that the ridge was mineralized and stabilized at a thickness of 5 mm and the implants were osseointegrated.

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## → ORTHODONTIC MICROSURGERY

**Yu H, Jiao F, Wang B, Shen SG.**

### Piezoelectric decortication applied in periodontally accelerated osteogenic orthodontics.

**EN** Craniofac Surg. 2013 Sep;24(5):1750-2

#### → Piezoelectric decortication applied in periodontally accelerated osteogenic orthodontics.

Our aim was to evaluate the application of piezoelectric decortication in periodontally accelerated osteogenic orthodontics (PAOO). One hundred fifty-six patients with severe skeletal malocclusions were enrolled in this study. Ultrasonic decortications were performed in 187 labial or lingual PAOO of the maxillary and mandibular anterior teeth. Orthodontic decompensation started from the fifth day after operation. All patients healed uneventfully and no severe periodontic complications were recorded. Rapid teeth movement and relatively short treatment duration were realized. Alveolar fenestration and bony dehiscence was successfully addressed. With physical and mechanical properties of absence of macrovibration, ease of use and control, piezosurgery showed its great values in PAOO.

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**Grenga V, Bovi M.**

### Corticotomy-Enhanced Intrusion of an Overerupted Molar Using Skeletal Anchorage and Ultrasonic Surgery.

**EN** J Clin Orthod. 2013 Jan;47(1):50-5.

#### → Corticotomy-Enhanced Intrusion of an Overerupted Molar Using Skeletal Anchorage and Ultrasonic Surgery.

Overerupted molars due to the loss of antagonists are a common clinical finding in adults. Because intrusion of an overerupted molar to permit prosthetic replacement is difficult to achieve without adverse side effects on the neighboring teeth, skeletal anchorage is now commonly used. Several surgical techniques have been proposed to allow faster

orthodontic movement in these cases. Wilcko and colleagues developed an accelerated osteogenic technique involving full-thickness flap reflection, selective alveolar decortication, ostectomy, and bone grafting. Other authors have recommended selective alveolar corticotomies involving full-thickness flaps and vertical and horizontal corticotomies limited to the cortical bone. Kim and colleagues proposed a minimally invasive technique called corticision to induce accelerated tooth movement by stimulating osteoblasts and bending the surgically separated alveolar bone. In this technique, a reinforced scalpel is used as a thin chisel to separate the interproximal cortices transmucosally without reflecting a flap. Vercellotti and Podesta recently introduced a new approach to bone surgery involving flap elevation and vertical and horizontal osteotomies performed with an ultrasonic piezosurgical device. Piezocision, an evolution of this technique developed by Dibart and colleagues, is a minimally invasive procedure that combines microincisions and localized piezoelectric surgery without flap elevation. In a previous article, we described the use of the Mectron Piezosurgery® unit for exposure of a palatally impacted canine. This device is a powerful piezoelectric handpiece with a functional frequency of 25-29 kHz and a digital modulation capability of 30Hz. It improves surgical control by performing a precise, selective cut, using ultrasonic waves transmitted to the hard tissues through an insert mounted on the handpiece. The inserts vibrate linearly between 60 and 210 microns, providing the handpiece with power in excess of 5W. A high-powered pump emits a biocompatible 9% sodium chloride solution for cooling. This article shows how intrusion of an overerupted molar with skeletal anchorage can be enhanced by selective alveolar corticotomies using ultrasonic bone surgery.

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Hernández-Alfaro F, Guijarro-Martínez R.

## Endoscopically assisted tunnel approach for minimally invasive corticotomies: a preliminary report.

EN J Periodontol. 2012 May;83(5):574-80.

⇒ **Endoscopically assisted tunnel approach for minimally invasive corticotomies: a preliminary report.**

**Background:** The dental community has expressed low acceptance of traditional corticotomy techniques for corticotomy-facilitated orthodontics. These procedures are time consuming, entail substantial postoperative morbidity and periodontal risks, and are often perceived as highly invasive.

**Methods:** A total of 114 interdental sites were treated in nine consecutive patients. Under local anesthesia, a tunnel approach requiring one to three vertical incisions per arch (depending on the targeted teeth) was used. Piezosurgical corticotomies and elective bone augmentation procedures were performed under endoscopic assistance. Postoperative

cone-beam computerized tomography evaluation was used to confirm adequate corticotomy depth.

**Results:** Procedures were completed in a mean time of 26 minutes. Follow-up evaluations revealed no loss of tooth vitality, no changes in periodontal probing depth, good preservation of the papillae, and no gingival recession. No evidence of crestal bone height reduction or apical root resorption was detected.

**Conclusions:** The tunnel approach minimizes soft-tissue debridement and permits effective cortical cuts. The combination of piezosurgery technique with endoscopic assistance provides a quick, reliable means to design and perform these corticotomies while maximizing root integrity preservation. Moreover, the sites needing bone augmentation are selected under direct vision. Compared to traditional corticotomies, this procedure has manifest advantages in surgical time, technical complexity, patient morbidity, and periodontium preservation.

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Vercellotti T, Podestà A.

## Orthodontic Microsurgery: A New Surgically Guided Technique for Dental Movement.

EN Int J Periodontics Restorative Dent. 2007; 27:325-331

⇒ **Orthodontic Microsurgery: A New Surgically Guided Technique for Dental Movement.**

Eight patients with malocclusions were treated with a new orthodontic-surgical technique that reduces the duration of treatment compared to conventional techniques. The monocortical tooth dislocation and ligament distraction (MTDLD) technique combines two different dental movements that work separately but simultaneously on opposite root surfaces. On the root surface corresponding to the direction of movement, vertical and horizontal microsurgical corticotomies are performed around each tooth root with a piezosurgical microsaw to eliminate cortical bone resistance. The immediate application of strong biomechanical forces produces rapid dislocation of the root and the cortical bone together. On the root surface opposite the direction of movement, the force of dislocation produces rapid distraction of ligament fibers. During the osteogenic process that follows, application of normal orthodontic biomechanics achieves the final tooth movement. All eight patients underwent periodontal and radiologic examinations for more than 1 year after treatment. No periodontal defects were observed in any of the patients, including one with a severe malocclusion and a thin periodontal tissue biotype. Compared to traditional orthodontic therapy, the average treatment time with the MTDLD technique in the mandible and maxilla was reduced by 60% and 70%, respectively.

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## → BONE GRAFTING

**Mouraret S, Houshyar KS, Hunter DJ, Smith AA, Jew OS, Girod S, Helms JA.**

### **Cell viability after osteotomy and bone harvesting: comparison of piezoelectric surgery and conventional bur.**

**EN** Int J Oral Maxillofac Surg. 2014 Aug;43(8):966-71.

#### **→ Cell viability after osteotomy and bone harvesting: comparison of piezoelectric surgery and conventional bur.**

The aim of this study was to evaluate and compare the influence of a piezoelectric device versus a conventional bur on osteocyte viability and osteoblast and osteoclast activity using an in vivo mouse model. Osteotomies were created and bone grafts were harvested using either a conventional bur or a piezoelectric device; the resulting injuries and bone grafts were evaluated over an extended time-course using molecular and cellular assays for cell death (TUNEL assay), cell viability (4',6-diamidino-2-phenylindole (DAPI) staining), the onset of mineralization (alkaline phosphatase activity), and bone remodelling (tartrate-resistant acid phosphatase activity). Osteotomies created with a piezoelectric device showed greater osteocyte viability and reduced cell death. Bone grafts harvested with a piezoelectric device exhibited greater short-term cell viability than those harvested with a bur, and exhibited slightly more new bone deposition and bone remodelling. The difference in response of osteocytes, osteoblasts, and osteoclasts to bone cutting via a bur and via a piezoelectric device is negligible in vivo. Given the improved visibility and the margin of safety afforded by a piezoelectric device, they are the instrument of choice when cutting or harvesting bone to preserve soft tissue.

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Majewski P.

## Piezoelectric surgery in autogenous bone block grafts.

EN Int J Periodontics Restorative Dent. 2014 May-Jun;34(3):355-63.

### ⇒ Piezoelectric surgery in autogenous bone block grafts.

This article describes alveolar ridge reconstruction in the esthetic zone using autogenous bone blocks harvested from the chin, taking into account the way the bone block is harvested, stabilized, and contoured in the recipient site. The 38 procedures were divided into two groups: group 1, using piezoelectric surgery, and group 2, using rotary instruments. The piezoelectric surgery technique made it possible to introduce surgical modifications. An observation of bone regeneration and follow-up clinical observations 5 to 7 years after the procedure revealed that the piezoelectric surgery technique provides better and more predictable clinical results for bone regeneration.

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Majewski P.

## Autogenous bone grafts in the esthetic zone: optimizing the procedure using piezosurgery.

EN Int J Periodontics Restorative Dent. 2012 Dec;32(6):e210-7.

### ⇒ Autogenous bone grafts in the esthetic zone: optimizing the procedure using piezosurgery.

Soft and hard tissue defects pose a therapeutic challenge in modern implant dentistry. There are a multitude of surgical techniques available, and it is necessary to match the problem with the solution. This report describes the reconstruction of the alveolar ridge in the esthetic zone with the help of autogenous bone blocks harvested from the chin that were shaped to fit and stabilized at the recipient site. The procedures were performed using Piezosurgery, which made it possible to introduce surgical modifications and had a significant impact on the accuracy of the procedure. An observation period of 2 to 7 years showed positive stable results for treatment in terms of function and esthetics.

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Lakshmiganthan M, Gokulanathan S, Shanmugasundaram N, Daniel R, Ramesh SB.

## Piezosurgical osteotomy for harvesting intraoral block bone graft.

EN J Pharm Bioallied Sci. 2012 Aug;4(Suppl 2):S165-8.

### ⇒ Piezosurgical osteotomy for harvesting intraoral block bone graft.

The use of ultrasonic vibrations for the cutting of bone was first introduced two decades ago. Piezoelectric surgery is a minimally invasive technique that lessens the risk of damage to surrounding soft tissues and important structures such as

nerves, vessels, and mucosa. It also reduces damage to osteocytes and permits good survival of bony cells during harvesting of bone. Grafting with intraoral bone blocks is a good way to reconstruct severe horizontal and vertical bone resorption in future implants sites. The piezosurgery system creates an effective osteotomy with minimal or no trauma to soft tissue in contrast to conventional surgical burs or saws and minimizes a patient's psychological stress and fear during osteotomy under local anesthesia. The purpose of this article is to describe the harvesting of intraoral bone blocks using the piezoelectric surgery device.

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## Bone harvesting at the anterior iliac crest using piezoosteotomy versus conventional open harvesting: a pilot study.

**EN** Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2008; 105(3):e19-e28

⇒ **Bone harvesting at the anterior iliac crest using piezoosteotomy versus conventional open harvesting: a pilot study.**

**Objective:** The objective of this study was to evaluate piezoosteotomy for bone harvesting at the iliac crest in a pilot collective.

**Study design:** Thirteen patients underwent anterior iliac crest bone graft harvesting using piezoosteotomy. These patients were compared to an age- and sex-matched retrospective cohort of 13 patients who underwent bone graft harvesting at the anterior iliac crest using conventional chisels and saws in an otherwise identical protocol.

**Results:** Harvested bone volumes and operation times were comparable; hospitalization time was briefer and pain levels at the first 2 postoperative days less in the piezoosteotomy group; pain medication requirement was comparable. Spearman's-rho correlation showed a strong tie between pain medication requirement and harvested bone volume in the piezoosteotomy group and between harvested bone volume and operation time in the conventional surgery collective.

**Conclusions:** The correlation between required pain medication and harvested bone volume together with lower pain levels in the piezoosteotomy group indicate pain more exclusively related to local osseous damage than to correlated soft tissue lacerations. Conventional technique shows a correlation between operation time and harvested bone volume indicating a more time-effective procedure, although total operation time was comparable. Piezoosteotomy appears to have potential use in iliac crest bone harvesting.

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**Happe A.**

## Use of a piezoelectric surgical device to harvest bone grafts from the mandibular ramus: report of 40 cases.

EN Int J Periodontics Restorative Dent. 2007; 27(3):241-249

### ⇒ Use of a piezoelectric surgical device to harvest bone grafts from the mandibular ramus: report of 40 cases.

Forty patients (12 men, 28 women; mean age 44.9 +/- 13.95 years) with inadequate bone volume for implant placement in either the maxilla or the mandible were scheduled for bone-grafting procedures. Bone grafts were harvested with a piezoelectric surgical device (Piezosurgery, Mectron) from 45 donor sites and grafted into recipient sites prior to implant placement. After healing of the grafted sites, 109 implants were placed. Clinical evaluations were performed to assess the size and quality of the obtained grafts, complications at the donor and recipient sites, integration of the grafts, bone quality, resorption of the grafts, and ability to properly place implants. Piezosurgery allowed precise, clean, and smooth cutting with excellent visibility. Mean graft size was 1.15 cm<sup>3</sup> (SD 0.5), with a maximum of 2.4 cm<sup>3</sup>. The quality of the obtained grafts was mainly cortical (71%), 42 of the 45 donor sites healed uneventfully (93.33%), and 50 of the 52 grafted sites healed uneventfully (96.15%).

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**Sohn DS, Ahn MR, Lee WH, Yeo DS, Lim SY.**

## Piezoelectric osteotomy for intraoral harvesting of bone blocks.

EN Int J Periodontics Restorative Dent. 2007; 27(2):127-131

### ⇒ Piezoelectric osteotomy for intraoral harvesting of bone blocks.

Grafting with intraoral bone blocks is a good way to reconstruct severe horizontal and vertical bone resorption in future implant sites. The Piezosurgery System (Mectron) creates an effective osteotomy with minimal or no trauma to soft tissue, in contrast to conventional surgical burs or saws. In addition, piezoelectric surgery produces less vibration and noise because it uses microvibration, in contrast to the macrovibration and extreme noise that occur with a surgical saw or bur. Micro-vibration and reduced noise minimize a patient's psychologic stress and fear during osteotomy under local anesthesia. The purpose of this article is to describe the harvesting of intraoral bone blocks using the piezoelectric surgery device.

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**Gellrich NC, Held U, Schoen R, Pailing T, Schramm A, Bormann KH.**

## Alveolar zygomatic buttress: A new donor site for limited preimplant augmentation procedures.

**EN** J Oral Maxillofac Surg. 2007 Feb;65(2):275-80.

→ **Alveolar zygomatic buttress: A new donor site for limited preimplant augmentation procedures.**

**Purpose:** Augmentation of anterior maxillary alveolar bone defects before placement of endosseous implants still poses a clinical challenge. Previous traumatic loss of upper anterior teeth is often concomitant with significant bone loss. The aim of augmentation is to reconstitute as far as possible the original hard tissue contour. This is a prerequisite for later optimal positioning of the dental implants in accordance with prosthodontic and functional principles.

**Materials and Methods:** This article describes a new method for the reconstruction of small anterior maxillary alveolar bone defects using donor bone from the zygomatic buttress region.

**Results:** This region provides harvested bone with a natural convex shape, which is ideally suited for the anterior alveolar process region. Later fibrous tissue transplantation to provide more anterior projection to the anterior dental area is also avoided.

**Conclusion:** Bone harvesting was done with the piezoelectric device, which avoids trauma to the mucosal membrane of the maxillary sinus.

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**Patel A, Schofield J.**

## Using Piezosurgery to harvest a block bone from the symphyseal region: a clinical case presentation.

**EN** Implant Dentistry Today. 2007; 1(4):20-24

→ **Using Piezosurgery to harvest a block bone from the symphyseal region: a clinical case presentation.**

The placement of titanium implants in sites that are deficient of alveolar ridge width has always been unpredictable due to the lack of bone around the implant. Implants placed into the available bone can often lead to a compromised aesthetic result due to poor positioning and direction of the implants, especially in the aesthetic zone. A number of surgical techniques have been proposed over the years to create adequate bone volume. These include the use of xenograft, alloplast

and allograft materials with the use of resorbable or non-resorbable membranes (Garg, 1999). It has been shown that autogenous bone grafts are still the most predictable method to achieve good bone quality and quantity for the implant site (Garg, 1999).

Donor sites for autogenous bone include extraoral and intraoral sites. Extraoral sites allow for harvesting of large amounts of bone. These sites involve complex surgical procedures that usually require general anaesthesia with an increase in morbidity at the donor sites (Chiapasco et al, 1999, Sindet-Pedersen and Enemark, 1990). It has been shown that in the facial skeleton, membranous bone (such as that grafted from the mandible) undergoes less resorption than endochondral bone, such as the iliac crest (Zins and Whitacker, 1979).

The main advantage of an intraoral donor site is that harvesting of bone can be performed under local anaesthesia. There are two main intraoral donor sites for autologous bone grafting available to the implantologist (Triplett and Silhow, 1998). The choice of site is based on the size of the defect to be augmented and the type and quantity of bone needed and the access to the donor site (Raghoebar et al, 2001). Only a certain amount of bone can be harvested from the retromolar area. It has been observed that the amount of bone volume obtained from the retromolar area is half of what can be harvested from the symphyseal region (Misch 1997). Retromolar grafts are usually used for single tooth augmentations.

Following harvesting of symphysis bone grafts, some long-term neurosensory disturbances have been described (Nkene et al, 2001, Raghoebar et al, 2001). In a prospective study evaluating intraoral donor sites for bone grafting, it was reported that 29% of symphysis graft patients had altered sensitivity of their mandibular incisor teeth (Misch 1997). All findings were resolved within six months. The majority of studies illustrate a continuous improvement of tooth sensitivity over time. The authors have audited 94 consecutive symphysis grafts performed over the last 10 years. All patients have reported some altered sensation in their lower incisors. The altered sensation lasted between three and six months in all but two cases. Of these two cases, one patient had altered sensation for 11 months. The other patient, on whom a very large chin graft was harvested, had some discolouration of her lower incisors. Despite this, sensitivity testing revealed normal results. The discolouration disappeared after four months.

In contrast, surgical access to the retromolar area might be compromised because of anatomy, i.e. the location of the mandibular canal and width of the ramus. This donor site appears to be associated with less permanent neurosensory disturbances (Nkenke et al, 2002).

When using conventional bone cutting instruments such as reciprocating microsaws or trephines, there is a risk of postoperative necrosis. Comparative histological studies have shown necrosis on cut bone caused by surgical burs, trephines and microsaws (Aro et al, 1981). Over the past five years, piezosurgery has been introduced as a new technique for osteotomy and osteoplasty using ultrasonic vibration. The piezosurgery blades allow for maximum intraoperative precision and minimal tissue damage, and the selective frequency of the scalpel minimises the risk to adjacent soft tissues. Histological studies have also shown minimal cellular damage to the resected bone margins, bony matrix and underlying marrow spaces. A recent study concluded that piezosurgery provided more favourable osseous repair and remodelling than with conventional surgical burs with surgical ostectomy and osteoplasty procedures (Vercellotti et al, 2005).

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**Sivolella S, Berengo M, Scarin M, Mella F, Martinelli F.**

## **Autogenous particulate bone collected with a piezo-electric surgical device and bone trap: a microbiological and histomorphometric study.**

**EN** Arch Oral Biol. 2006; 51(10):883-891

### **→ Autogenous particulate bone collected with a piezo-electric surgical device and bone trap: a microbiological and histomorphometric study.**

The aims of this study were to determine the microbiological and particle size characteristics of particulate bone collected with a piezosurgery device and bone trap, and to reduce bacterial contamination after treatment of debris with rifamycin SV. Samples were taken from 10 patients who underwent surgical extraction of their third lower molars. The ostectomy was performed with a piezosurgery device, and the debris was collected with a surgical aspiration set equipped with a bone trap. Two aliquots were taken from each sample, one of which was treated with rifamycin SV. The second aliquot, used as a control, was treated with a physiological solution. In the samples immersed in antibiotic solution, there was a statistically significant ( $P < 0.005$ ) reduction in bacterial contamination. The stringent protocol followed in this study has proved valid for collection of material, and treatment with rifamycin SV was found to reduce bacterial contamination in collected material.

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## Piezoelectric Harvesting of an Autogenous Bone Graft from the Zygomaticomaxillary Region: Case Report.

EN Int J Periodontics Restorative Dent. 2006; 26(5):453-457

### ⇒ Piezoelectric Harvesting of an Autogenous Bone Graft from the Zygomaticomaxillary Region: Case Report.

The zygomaticomaxillary region offers a large amount of cortical bone that can be obtained simply and safely using the precise and selective cutting properties of a piezosurgical device. A block from this area fits nicely into anterior or premolar maxillary recipient sites and is thus the ideal choice, as no secondary surgical field is needed. As in conventional sinus bone graft procedures, the complication rate is minimal and after a 5-month healing period, the augmented region can be used for stable and esthetic oral implant placement.

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### ⇒ Augmentation of alveolar ridge defects: autologous bone transplant from the zygomatic alveolar crest – a new technique.

The insertion of dental implants with simultaneous loss of alveolar bone still represents a challenge. A traumatic tooth loss within the range of the front of the maxillary is frequently connected with a clear substance defect. The goal of a pre-implantologically treatment must be therefore the hard tissue reconstruction of the defect. Thus the implant camp is optimized and the prosthetically necessary implant position is guaranteed. In this case report the reconstruction of a vertical and horizontal bone loss with a transplant from the Crista zygomatico alveolaris is represented. The withdrawal of the bone transplant took place under indulgence of the vestibular membrane of the maxillary sinus by means of Piezosurgery®. After an appropriate healing time the single tooth gap was supplied with an ITI implant of the company Straumann (Freiburg, Germany).

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⇒ **Piezosurgery and bone grafting.**

Piezosurgery significantly simplifies graft harvesting from the ramus, allowing a bone graft to be taken in a site where the postoperative course is always straightforward without a risk of neurological sequela and with mucosal healing of a quality not obtained with harvesting from the symphysis.

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**Chiriac G, Herten M, Schwarz F, Rothamel D, Becker J.**

### Autogenous bone chips: influence of a new piezoelectric device (Piezosurgery®) on chips morphology, cell viability and differentiation.

**EN** J Clin Periodontol. 2005; 32(9):994-999

⇒ **Autogenous bone chips: influence of a new piezoelectric device (Piezosurgery®) on chips morphology, cell viability and differentiation.**

**Aim:** The aim of the present study was to investigate the influence of a new piezoelectric device, designed for harvesting autogenous bone chips from intra-oral sites, on chip morphology, cell viability and differentiation.

**Methods:** A total of 69 samples of cortical bone chips were randomly gained by either (1) a piezoelectric device (PS), or (2) conventional rotating drills (RD). Shape and size of the bone chips were compared by means of morphometrical analysis. Outgrowing osteoblasts were identified by means of alkaline phosphatase activity (AP), immunhistochemical staining for osteocalcin (OC) synthesis and reverse transcriptase-polymerase chain reaction phenotyping.

**Results:** In 88.9% of the RD and 87.9% of the PS specimens, an outgrowth of adherent cells nearby the bone chips was observed after 6-19 days. Confluence of cells was reached after 4 weeks. Positive staining for AP and OC identified the cells as osteoblasts. The morphometrical analysis revealed a statistically significant more voluminous size of the particles collected with PS than RD.

**Conclusion:** Within the limits of the present study, it may be concluded that both the harvesting methods are not different

from each other concerning their detrimental effect on viability and differentiation of cells growing out of autogenous bone chips derived from intra-oral cortical sites.

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⇒ **Piezosurgery and pre-implant preparation of bone - Grafting by adding blocks of autogenous bone with portions of harvested bone.**

Piezoelectric surgery is a new technique that is an interesting alternative to traditional bone surgery techniques. The insert performs a selective cutting action on mineralised tissues, while it is ineffective on the soft tissues. It can be used to carry out osteotomies and osteoplasties thanks to the micro-vibrations of the active part of the inserts. Because of the small size of the insert, it is possible to make sub-millimetric cuts. Furthermore, the fact that the insert is cooled by physiological liquid contributes towards keeping the surgical field bloodless. Two clinical cases of rehabilitation of vestibulo-palatal bone defects are presented below. Two autogenous bone grafts, taken from the mandibular retromolar triangle, were carried out using the piezosurgery technique to reconstruct the bone and create a volume sufficient to receive an implant.

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## → BONE GRAFTING / → OSTEOTOMY CLOSE TO NERVES

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## → OSTEOTOMY CLOSE TO NERVES

**Bovi M, Manni A, Mavriqi L, Bianco G, Celletti R.**

### The use of piezosurgery to mobilize the mandibular alveolar nerve followed immediately by implant insertion: a case series evaluating neurosensory disturbance.

**EN** Int J Periodontics Restorative Dent. 2010 Feb;30(1):73-81.

→ The use of piezosurgery to mobilize the mandibular alveolar nerve followed immediately by implant insertion: a case series evaluating neurosensory disturbance.

One of the therapeutic options proposed for reconstruction of the atrophic posterior mandible is inferior alveolar nerve (IAN) mobilization with simultaneous implant placement. However, studies on the functionality of this neurovascular bundle after its mobilization have shown mixed results. This variability can be attributed both to the test methodology, which typically requires subjective answers from patients, and to the surgical procedure itself, which is highly dependent on operator technique. This article reports on a series of 10 cases of IAN mobilization using a device specifically engineered to simplify bone surgery. This device enables the oral surgeon to avoid overstretching the nerve by creating a smaller bone window and using an apicocoronal inclination of instruments to capture the neurovascular bundle. Evaluation by means of neurosurgery function tests over a 36-month period found that all patients had a return to normal sensation after a brief period of neurosensory disturbance. Subjective responses to a patient questionnaire confirmed these findings. The implant success rate was 100%.

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**Degerliyurt K, Akar V, Denizci S, Yucel E.**

### Bone lid technique with PIEZOSURGERY® to preserve inferior alveolar nerve.

**EN** Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2009 Dec;108(6):e1-5.

→ Bone lid technique with PIEZOSURGERY® to preserve inferior alveolar nerve.

The inferior alveolar nerve is at risk during surgical removal of impacted mandibular teeth and excision of benign tumors of

## ⇒ OSTEOTOMY CLOSE TO NERVES

mandible. Manual and/or mechanical instruments used in the close proximity of delicate structures do not allow the control of the cutting depth and can damage these structures by accidental contact. Piezoelectric surgery is a new and innovative bone surgery technique using the microvibrations of special scalpels at ultrasonic frequency so therefore soft tissue will not be damaged even upon accidental contact with the cutting tip. This article presents an alternative technique that uses piezoelectricity to minimize trauma to the inferior alveolar nerve, vascular tissues, or surrounding dental tissues.

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**Schaeren S, Jaquiéry C, Heberer M, Tolnay M, Vercellotti T, Martin I.**

### **Assessment of nerve damage using a novel ultrasonic device for bone cutting.**

**EN** J Oral Maxillofac Surg. 2008; 66(3):593-596

### ⇒ **Assessment of nerve damage using a novel ultrasonic device for bone cutting.**

In oral, maxillofacial, and spinal surgery, osteotomy often must be performed in close vicinity to nerve tissue, with the potential risk of transient or permanent neurologic injuries (eg, to trigeminal nerve branches). Traditional tools, such as rotating burs and oscillating saws, are highly effective in cutting bone tissue but are not selective for bone, and thus can produce significant harm to surrounding soft tissues, especially nerves.

Piezosurgery (Mebiotec, Sestri Levante, Italy), a device for bone-cutting based on low-frequency (25 to 29 kHz) ultrasonic vibrations, has recently been introduced in oral and maxillofacial surgery. This device can improve the control and precision of osteotomy and improve bone healing due to reduced local trauma. A recent in vitro study of the use of Piezosurgery for transposition of the inferior alveolar nerve in the mandibles of cadaver sheep showed that this technique caused roughening of the epineurium without affecting deeper structures and induced less injury than a conventional rotary bur. A pilot clinical study found that Piezosurgery reduced neurosensory disturbances in orthognathic surgery of the mandible. To the best of our knowledge, no experimental in vivo study has yet demonstrated that bone cutting using Piezosurgery prevents damage to soft tissue, especially nerves. The aim of the present study was to assess the potential damage of Piezosurgery to a peripheral nerve on direct contact, in 2 possible scenarios.

The first scenario corresponds to the situation in which the surgeon does not immediately realize the contact with the nerve and continues to apply the same force used for cutting bone (working force) for an estimated reaction time of 5 seconds. The second scenario corresponds to an accidental slip of the device on the bone tissue, causing a direct hit on the nerve with a higher force (peak force), estimated to be twice the working force, but for a shorter time (1 second). The rat sciatic nerve was selected for the study because it is an established model for studying sensory and motor nerve damage and recovery and also is of equivalent size of branches of human trigeminal nerve.

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Sakkas N, Otten JE, Gutwald R, Schmelzeisen R.

### Transposition of the mental nerve by piezosurgery followed by postoperative neurosensory control: a case report.

**EN** Br J Oral Maxillofac Surg. 2008; 46(4):270-271

#### ⇒ Transposition of the mental nerve by piezosurgery followed by postoperative neurosensory control: a case report.

Transposition of the mental nerve is a preprosthetic procedure that is effective for patients with hyperesthesia caused by the effect of a dental prosthesis on the alveolar ridge. We present the case of a 74-year-old woman with pain and hyperesthesia of the right inferior alveolar nerve caused by a dental prosthesis. Caudal transposition of the right mental nerve by piezosurgery resulted in postoperative neurosensory controls of the lower lip showing normal nerve function 2 months later.

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#### ⇒ Juxta-canal implants.

Implantology forms part of practitioners' therapeutic arsenal and can be integrated in treatment plans for oral rehabilitation. It obviates the need for removable or fixed prostheses, which sometimes necessitate preparation of healthy teeth. Nevertheless, implant placement depends on the presence of sufficient bone volume at the implant site. If necessary, prior performance of a bone graft may have to be planned (1). Another solution is maximum use of the available bone. The presented clinical case describes an operative technique that allows use of the entire bone height available above the

inferior alveolar nerve in the posterior area. This technique employs piezoelectric bone surgery, developed by Vercellotti (11), an interesting alternative to the classical techniques of bone surgery (2). Piezo surgery has been conceived to simplify surgical protocols and improve the predictability of the results and the safety of the surgical procedure.

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Geha H, Gleizal A, Nimeskern N, Beziat JL.

## Sensitivity of the Inferior Lip and Chin following Mandibular Bilateral Sagittal Split Osteotomy Using Piezosurgery.

EN Plast Reconstr Surg. 2006;118(7):1598-1607

⇒ **Sensitivity of the Inferior Lip and Chin following Mandibular Bilateral Sagittal Split Osteotomy Using Piezosurgery.**

**Background:** Bimaxillary osteotomy, including bilateral sagittal split osteotomy, is the most commonly performed orthognathic surgical intervention in Europe and the United States. Neurosensory perturbation in the territory of the inferior alveolar nerve is a reported adverse effect of bilateral sagittal split osteotomy. Piezosurgery is a relatively new technique that allows bone to be cut while preserving soft tissues, including nerves. The purpose of this study was to assess inferior alveolar nerve function through clinical neurosensory testing after bilateral sagittal split osteotomy using Piezosurgery.

**Methods:** Between February and September of 2004, 20 patients (40 sides) presenting with dento-skeletal deformities underwent bimaxillary osteotomy, including bilateral sagittal split osteotomy. The Mectron Piezosurgery device was used to perform all sagittal splits, with distraction being performed between the two bone valves. The inferior alveolar nerve was evaluated both objectively with clinical neurosensory testing, including pin-prick sensation, light touch sensation, and two-point discrimination tests, and subjectively at the following time points: preoperatively; at days 5, 7, and 10, postoperatively; and at the second month post-operatively.

**Results:** The anatomical integrity of the inferior alveolar nerve was respected in all cases. Observed normal results for the different tests at 10 days were 90, 82, and 70 percent, respectively, for pin-prick sensation, light touch sensation, and two-point discrimination. By computing scores for different clinical neurosensory tests, the authors observed between 75 and 80 percent complete neurosensory recuperation as early as the second postoperative month.

**Conclusions:** Piezosurgery used for bilateral sagittal split osteotomy allows prompt recovery of inferior alveolar nerve neurosensory function within 2 months. No comparison is possible with the results using the standard technique for bilateral sagittal split osteotomy.

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**Bovi M.**

## Mobilization of the Inferior Alveolar Nerve with simultaneous implant insertion: A New Technique. A Case Report.

**EN** Int J Periodontics Restorative Dent. 2005; 25(4):375-383

### ⇒ Mobilization of the Inferior Alveolar Nerve with simultaneous implant insertion: A New Technique. A Case Report.

For reconstruction of an atrophied posterior mandible, different therapeutic options have been proposed, such as autologous bone grafting, guided bone regeneration for vertical ridge augmentation, and inferior alveolar nerve (IAN) mobilization with simultaneous implant placement. The possible dehiscence of soft tissues covering the surgical zone makes the first and second techniques unpredictable. Moreover, two surgical sites are necessary and a long treatment time is required (about 12 months). With IAN mobilization, only one surgical intervention is required and the total treatment time is shorter (about 6 months). However, this technique risks irreversible damage to the IAN, with consequent functional alterations. Current studies have shown extreme variability in the examination of functionality of the neurovascular bundle after its mobilization. This variability can be attributed both to the methodology used for the tests, which evoke subjective answers from the patient, and to the surgical procedure, which is highly dependent on operator technique. Nerve damage can be the result of an overstretched mucoperiosteal flap in the premolar area to achieve optimal visibility of the surgical zone. This article reports a case in which a new surgical approach to IAN mobilization is performed using a specifically engineered device for simplified bone surgery (Mectron Piezosurgery). This device enables the surgeon to cut hard tissue without injuring the soft tissues. Therefore, there is a lower risk of damaging the IAN, and it is possible to reduce over-stretching of the mental nerve by creating a smaller bone window and using an apicocoronal inclination of instruments to capture the neurovascular bundle.

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Spinato S., Rebaudi A., Bernardello F., Bertoldi C., Zaffe D.

## Piezosurgical treatment of crestal bone: quantitative comparison of post-extractive socket outcomes with those of traditional treatment.

EN Clin Oral Implants Res. 2015 Jan 30.

### ⇒ Piezosurgical treatment of crestal bone: quantitative comparison of post-extractive socket outcomes with those of traditional treatment.

**Objective:** The study aimed to quantitatively compare, for the first time, the clinical outcomes of crestal bone volume resorption in sockets undergoing Traditional Extraction Technique (TET) or Piezosurgical Extraction Technique (PET), also considering the influence of buccal plate thickness.

**Material and Methods:** In this prospective study, 19 sockets were randomly treated with TET, and 18 sockets were randomly treated with PET. Furthermore, patients were split into subgroup A, with buccal bone plate thickness (BPT) $\leq$ 1 mm, and subgroup B, with BPT $>$ 1 mm. Buccal (BCH) and palatal (PCH) cortex height, bucco-palatal ridge (BPR) width were monitored at tooth extraction and after the 4-month post-extractive period of natural healing.

**Results:** After 4 months, BCH, PCH and BPR width decreased more in the TET than in the PET group, but only the BPR decrease was statistically significant ( $P = 0.034$ ) after ANOVA test. In both TET and PET groups, all B subgroup patients showed a lower decrease than A subgroup patients for both BCH, PCH and BPR, statistically significant for PCH ( $P=0.019$ ) and BPR ( $P<0.001$ ) of TET group, and BPR ( $P=0.002$ ) of PET group, after ANOVA. Both A and B subgroups of PET showed a statistically significant lower decrease than the corresponding subgroups of TET, comparing A ( $P=0.005$ ) and B ( $P=0.037$ ) subgroups for BPR, after ANOVA.

**Conclusions:** With both thin and thick buccal plates, the piezosurgical extraction technique of teeth significantly decreases the horizontal resorption of the hard tissue ridge, but not the vertical resorption. Moreover, buccal plate thickness seems to be a key factor in post-extractive bone resorption: the thinner the buccal plate the greater the horizontal crestal bone loss.

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## → EXTRactions

Nagori SA, Jose A, Bhutia O, Roychoudhury A.

### Evaluating success of autotransplantation of embedded/impacted third molars harvested using piezosurgery: a pilot study.

**EN** Acta Odontol Scand. 2014 Nov;72(8):846-51.

⇒ Evaluating success of autotransplantation of embedded/impacted third molars harvested using piezosurgery: a pilot study.

**Objective:** To evaluate the success of autogenous transplantation of embedded/impacted third molars harvested using piezosurgery.

**Material and Methods:** This prospective pilot study enrolled 20 healthy patients with non-restorable first/second molars and a caries-free retrievable embedded/impacted third molar. Piezosurgery was used for removing inter-radicular bone at the recipient socket as well as for bone removal around the donor teeth.

**Results:** After an average follow-up of 16.4 months ( $SD = 1.9$ ), 18 cases were successful with formation of periodontal ligament around the teeth. One tooth was lost due to infection at 1 month. One patient was lost to follow-up. There was no root resorption or ankylosis in any of the cases. In six donor teeth with complete root formation, root canal treatment was carried out. All the remaining teeth responded positively with vitality testing.

**Conclusion:** Piezosurgery is an effective device if embedded/impacted third molars are to be harvested for successful autogenous transplantation.

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Piersanti L, Dilorenzo M, Monaco G, Marchetti C

### Piezosurgery or Conventional Rotatory Instruments for Inferior Third Molar Extractions?

**EN** J Oral Maxillofac Surg. 2014 Sep;72(9):1647-52.

⇒ Piezosurgery or Conventional Rotatory Instruments for Inferior Third Molar Extractions?

**Purpose:** The purpose of this study was to compare the discomfort and surgical outcomes of a piezosurgery device with those of rotatory instruments in lower third molar extraction.

**Material and Methods:** A split-mouth, randomized, unblinded clinical study was designed; the 2 molars had to have the same extraction difficulty score. The test side was extracted using a piezosurgery technique and the control side was extracted using a conventional handpiece. The primary endpoint was patient discomfort evaluated with the Postoperative Symptom Severity (PoSSe) scale, which was administered to each patient; secondary endpoints were pain, trismus, swelling, and surgical time evaluation. Pairedsamples t test and repeated-measures analysis of variance were used to compare outcomes within patients.

**Results:** Ten consecutive patients (6 female, 4 male; mean age, 22.4 - 2.3 yr) were recruited. The total score on the PoSSe scale was significantly lower for piezosurgery compared with the conventional rotating handpiece (24.7 - 10.3 vs 36.0 - 7.6;  $t = -4.27$ ;  $P = .002$ ). Moreover, postoperative swelling 1 week after surgery was significantly lower for piezosurgery than for the conventional rotating handpiece (2.75 - 0.23 vs 3.1 - 0.39 cm;  $t = -2.63$ ;  $P = .027$ ).

**Conclusions:** Piezosurgery was associated with less postoperative discomfort and yielded better results for swelling.

## ⇒ EXTRactions

Piezosurgery seems to be a good technique in daily surgical practice, especially if applied in the critical steps in which safety and respect for soft tissue, bone, and nerves are necessary.

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### Third-molar extraction with ultrasound bone surgery: a case-control study.

**EN** Craniofac Surg. 2014 May;25(3):856-9.

#### ⇒ Third-molar extraction with ultrasound bone surgery: a case-control study.

**Purpose:** The aim of this case-control study was to evaluate the postoperative period and healing between 2 surgical methods (traditional and ultrasound bone surgery) that are used for mandibular third-molar extraction.

**Patients and Methods:** Fifteen patients with impaction of both of the lower third molars and indications for their extractions were used in this study. Bilateral-mandibular third-molar extractions were performed at the same surgical time: traditional surgery with burrs was used on 1 side (control site), and ultrasound surgery was used on the other side (test [T] site). After surgery, the patients were examined at 7 and 14 days and at 1 and 3 months to evaluate tissue healing. The following was assessed at every follow-up: pain, trismus, swelling, and alveolar bone level.

**Results:** The study included 15 patients, and 30 mandibular third-molar extractions were performed. We found only 1 post-operative complication: 1 patient had alveolitis in the control site. Complete recoveries without any complications were reported in all of the patients at the T sites.

**Conclusions:** Complete recoveries without any complication were reported in all patients at the T sites. The only disadvantage of the piezoelectric technique was the length of operation time, which was increased by approximately 8 minutes; however, this effect was offset by reducing the morbidity.

**Clinical Relevance:** Our preliminary study showed that Piezosurgery is an excellent tool for reducing the risk of complications and improving the postoperative period.

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## Use of the piezosurgery technique for cutting bones in the autotransplantation of unerupted third molars.

**EN** Int J Periodontics Restorative Dent. 2013 Jul-Aug;33(4):477-81.

### ⇒ Use of the piezosurgery technique for cutting bones in the autotransplantation of unerupted third molars.

Autotransplantation is a well-known method used in oral surgery. However, risk of failure, most commonly resulting from root resorption of the transplanted tooth or ankylosis, is quite high. Piezosurgery with specific device tip vibration frequencies enables selective tissue cutting, and therefore, tooth buds or teeth can easily be removed from bones with little injury to periodontal fibers or bud follicles.

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Rullo R, Addabbo F, Papaccio G, D'Aquino R, Festa VM.

## Piezoelectric device vs. conventional rotative instruments in impacted third molar surgery: relationships between surgical difficulty and postoperative pain with histological evaluations.

**EN** J Craniomaxillofac Surg. 2013 Mar;41(2):e33-8.

### ⇒ Piezoelectric device vs. conventional rotative instruments in impacted third molar surgery: relationships between surgical difficulty and postoperative pain with histological evaluations.

**Purpose:** To investigate and compare the influence of surgical difficulty on postoperative pain after treatment of impacted mandibular third molars by rotatory osteotomy or Piezoelectric surgery.

**Material and Methods:** A prospective, randomized, split-mouth study was performed of 52 patients with bilateral and symmetrically oriented impacted mandibular third molars, who were surgically treated using a burr (Group A) on one random side of the lower jaw and a Piezoelectric device (Group B) on the contralateral side. Surgical difficulty was evaluated using a modified version of the Parant scale to categorize "simple extractions" and "complex extractions". Primary outcome parameters were the comparison of the postoperative pain evaluation rated on the Visual Analogue Scale from day 0 to day 6 postsurgery, and the assessment of differences in surgery time between the groups. Bone biopsies were taken during surgery to assess differences in bone tissue damage levels between the two different techniques.

**Results:** In "complex extractions" lower pain evaluation and significantly shorter surgery times were recorded when rotatory instruments were used. In "simple extractions", similar surgery times were observed for both techniques, but pain was greatest on the day of surgery when the burr was used. Bone heat osteonecrosis was observed only in the rotatory group and a high level of alkaline phosphatase was noted only in the Piezoelectric group.

**Conclusion:** Pain after extraction of a mandibular third molar increases with increased surgical difficulty and especially in longer interventions. The integrity of the bony structure observed after the ultrasonic technique may favour the bone healing process.

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**IT** Minerva Stomatol. 2012 Jun;61(6):247-53

⇒ **The piezoelectric osteotomy technique compared to the one with rotary instruments in the surgery of included third molars. A clinical study.**

**Aim:** The purpose of this study is to compare surgical times, the saving of bonytissues, the postoperative course of the extraction of included third molars in patients treated with piezoelectric surgery compared to those treated with conventional rotary instruments.

**Methods:** The study included 140 patients with upper and lower third molars included. 70 of them were treated with osteotomy using conventional rotary instruments (Group 1) and 70 more with osteotomy through piezosurgery Mectron® (GROUP 2). The treatment protocol was the same for both groups. The time of surgery and the saving of bonytissues in the intra-operative have been estimated; moreover, 24-48-72 hours and 7 days after the surgery, two other parameters have been estimated: facial swelling and trismus. The swelling, trismus and the saving of bony tissues were assessed by a surgery gauge.

**Results:** The average surgical time was 15 minutes in Group 1 and 20 minutes in Group 2. The savings of bonytissues was greater in Group 2 than in group 1 by 2.7 mm. The average facial swelling was 6.23 mm for group 1 and 2.86 mm for group 2 24 hours after the surgery, of 5.22 mm for group 1 and 1.76 mm for group 2 48 hours after the intervention, of 3.75 mm for Group 1 and 0.85 for group 2 72 hours after the surgery, of 0.86 mm for group 1 and 0.12 for group 2. The average trismus was of 14.76 mm in Group 1 and 11.15 mm in Group 2. Statistical analysis showed a significant reduction of swelling, trismus and a bigger saving of bony tissues in Group 2 and a reduction in surgical time in Group 1.

**Conclusion:** The osteotomy technique with a piezoelectric instrument has produced a significant reduction in facial swelling, trismus and a bigger saving of bony tissues than with conventional rotary instruments. A slight lengthening in terms of surgical time has been recorded, compared to the use of conventional rotary instruments.

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Guo ZZ, Zhang H, Li Y, Li X, Liu Y, Wang Y, Yuan CX, Liu X.

## Comparative study of complications among routine method, high speed turbine handpiece and piezosurgery device after extraction of impacted wisdom teeth.

**CN** Shanghai Kou Qiang Yi Xue. 2012 Apr;21(2):208-10.

⇒ Comparative study of complications among routine method, high speed turbine handpiece and piezosurgery device after extraction of impacted wisdom teeth.

**Purpose:** To investigate complications in extraction of complicated impacted wisdom teeth whose root apex near to the inferior alveolar nerve(IAN) by using routine method(chisels),high speed turbine handpiece and piezosurgery device respectively.

**Methods:** Three hundred qualified patients with impacted wisdom teeth were divided into three groups randomly, one hundred patients in group A were extracted by routine method, one hundred patients in group B were extracted by high speed turbine handpiece, and one hundred patients in group C were extracted by piezosurgery device. The operation time, postoperative pain duration,dry socket and IAN injury were compared between each two groups. All statistical analysis was performed using SPSS 13.0 software package. Differences between groups were compared using a paired t test (quantitative data) or Chi-square test (qualitative data).

**Results:** The operation time in group A was  $(14.12 \pm 0.12)$ min,  $(7.22 \pm 0.15)$ min in group B,  $(25.23 \pm 0.32)$ min in group C; Significant difference was found between group A and group B( $P < 0.05$ ), group B and group C( $P < 0.05$ ),group A and group C( $P < 0.05$ ). Postoperative pain duration was  $(62.15 \pm 1.51)$ h in group A,  $(48.23 \pm 1.23)$ h in group B,  $(14.34 \pm 0.80)$ h in group C; Significant difference was found between group A and group B( $P < 0.05$ ), group B and group C( $P < 0.05$ ),group A and group C( $P < 0.05$ ). 9 patients developed dry socket in group A, 2 in group B, and 1 in group C; Significant difference was found between group A and group B( $P < 0.05$ ), group A and group C( $P < 0.05$ ). Six patients had IAN injury in group A, 2 in group B, 0 in group C. Significant difference was found between group A and group C.

**Conclusions:** Compared with routine method, high speed turbine is better in extraction of impacted wisdom teeth, which can shorten operation time, lessen postoperative complications. Although there was longer operation time compared with group A and B, piezosurgery device is more effective in reducing postoperative complications.

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Gao Y, Jiang A, Li B, Yang L.

## Comparison of piezosurgery and chisel osteotomy in the extraction of mandibular impacted third molars.

CN Hua Xi Kou Qiang Yi Xue Za Zhi. 2011 Aug;29(4):372-4.

⇒ Comparison of piezosurgery and chisel osteotomy in the extraction of mandibular impacted third molars.

**Objective:** To study the effect of piezosurgery in the extraction of mandibular impacted third molars.

**Methods:** 228 mandibular impacted third molars with relative difficulty for extraction according to the scheme of degree of difficulty for extraction were included in the study, and were divided into two groups (114 teeth each group): Test group (treated by piezosurgery), control group (treated by traditional chisel osteotomy). The surgery time and pain, restriction of mouth opening and facial swelling after surgery, were evaluated in both groups.

**Results:** The average surgery time was  $(16 \pm 5,2)$  min in test group and  $(30 \pm 8,7)$  min in control group. The surgery time, pain, facial swelling rate, and restriction of mouth opening were lower than control group ( $P < 0,05$ ).

**Conclusion:** Compared with traditional chisel osteotomy, the average surgery time with piezosurgery was significantly shorter and the complications are obviously reduced.

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Sortino F, Pedullà E, Masoli V.

## The piezoelectric and rotatory osteotomy technique in impacted third molar surgery: comparison of postoperative recovery.

EN J Oral Maxillofac Surg. 2008 Dec;66(12):2444-8.

→ The piezoelectric and rotatory osteotomy technique in impacted third molar surgery: comparison of postoperative recovery.

**Purpose:** The aim of this study was the comparison of postoperative outcome in mandibular impacted third molars treated by piezoelectric surgery or by rotatory osteotomy technique.

**Patients and methods:** One hundred patients with impacted mandibular third molars were included in the study. Fifty patients were treated by rotatory osteotomy technique (group A) and 50 patients were treated by piezoelectric osteotomy technique (group B). Therapeutic protocol was the same for both groups. Twenty-four hours after surgery, 2 different parameters, facial swelling and trismus, were evaluated in both groups. A pair of compasses were used for the evaluation of facial swelling and trismus was evaluated.

**Results:** The average surgery time was 17 minutes in group A, 23 minutes in group B; the mean facial swelling was 7.04 mm in group A, 4.22 mm in group B; trismus was 16.76 mm in group A, 12.52 mm in group B. Statistical analysis showed a significant reduction ( $P < .05$ ) of postoperative facial swelling and trismus in group B; however, in this group, a statistically significant increased ( $P < .05$  vs group A) surgery time was required.

**Conclusion:** The piezoelectric osteotomy technique produced a reduced amount of facial swelling and trismus 24 hours after surgery, but a longer surgery time was required when compared with the rotatory osteotomy technique.

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Sivolella S, Berengo M, Fiorot M, Mazzuchin M.

## Retrieval of blade implants with piezosurgery: two clinical cases.

EN Minerva Stomatol. 2007; 56(1-2):53-61

→ Retrieval of blade implants with piezosurgery: two clinical cases.

In this work an ultrasound device was used to perform an osteotomy for the removal of blade implants in order to save as much bone tissue as possible, so that root form implants might later be inserted. Two patients underwent surgery for the removal of two blade implants (one maxillary, the other mandibular) that were no longer functional. The peri-implant osteotomy was carried out with a piezoelectric surgery device. The instrument demonstrated to be effective and precise during osteotomy, providing an extremely thin cutting line. During the course of the operation and at controls after 7 and 30 days, patients did not show any relevant complications and both still had sufficient alveolar bone to be treated with root form implants. The piezosurgery device proved to be an effective instrument in interventions requiring a significant saving of bone tissue, extreme precision in cutting, and respect of soft tissues.

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**Grenga V, M. Bovi.**

## Piezoelectric Surgery for Exposure of Palatally Impacted Canines.

EN J Clin Orthod. 2004; 38(8):446-448

⇒ Piezoelectric Surgery for Exposure of Palatally Impacted Canines.

A precise surgical cut is extremely important when the crown of a palatally impacted maxillary canine is close to the roots or crowns of the central and lateral incisors, as often occurs. Traditional burs do not distinguish between the mineralization and hardness of the bone, the radicular cementum and the enamel. Piezoelectric surgery is precise enough to account for such differences, thus avoiding damage to the adjacent teeth.

The EX1 insert for the Mectron Piezosurgery handpiece provides a good demarcation of the exposed crown. It also removes the follicle faster than manual instruments can, with no loss of cervical bone or consequent gingival recession. The piezoelectric instrument controls bleeding during the surgical procedure, ensuring a dry field for bonding to the impacted tooth and eliminating the need for gauze or electrocoagulation. Total surgical time is greatly reduced with this method. The piezoelectric device is quieter than conventional handpieces, and it avoids the trauma associated with a manual hammer and chisel. Therefore, it improves both the comfort and the cooperation of the patient.

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**Marini E, Cisterna V, Messina AM.**

## The removal of a malpositioned implant in the anterior mandible using piezosurgery.

**EN** Oral Surg Oral Med Oral Pathol Oral Radiol. 2013 May;115(5):e1-5.

### ⇒ The removal of a malpositioned implant in the anterior mandible using piezosurgery.

In oral, cranio, and maxillofacial surgery, a close relationship among the bone, nerves, and blood vessels can be regularly observed. Surgical procedures for the removal of dental implants have the potential to cause vascular injury and bleeding in the floor of the mouth and internal anterior region of the mandible. Furthermore, conventional osteotomy techniques always require extensive protection of adjacent soft tissue because cutting is not limited to bone and could easily affect other tissues when applied improperly. We report the removal by means of piezosurgery of a malpositioned osseointegrated implant that had previously caused a sublingual hematoma during its insertion. The postoperative course was uneventful, no bleeding, infection, or hematoma formation was noted and the patient reported 100% resolution of all symptoms.

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**Moorhouse J, Campbell I.**

## Recovery of a poorly-placed implant using piezosurgery.

**EN** Implant Dentistry Today – 2012 September.

### ⇒ Recovery of a poorly-placed implant using piezosurgery.

High volume implant referral centres are seeing more and more instances of poorly placed implants with often dramatic restorative solutions, where the technician attempts to make the best of a bad situation. Unfortunately, these poor aesthetics go hand-in-hand with periodontal and peri-implant pathology and an unhappy patient. Often the only solution to remove the implant and start again, which presents its own problems because of the concomitant massive loss of bone. The aim of this article is to use a single case to describe the clinical approach to rescuing a failed implant. It will consider the advantages of the Mectron Piezosurgery unit in removing the implant with minimal bone loss, and then in cortical bone grafting to restore the area.

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## Use of piezosurgery device in management of surgical dental implant complication: a case report.

EN Implant Dent. 2011 Apr;20(2):e1-6.

→ Use of piezosurgery device in management of surgical dental implant complication: a case report.

An implant displacement into the maxillary sinus could be a complication of implant surgery in the upper jaw. In such cases, implant removal is needed to avoid the occurrence of sinus pathologies. Piezosurgery techniques could assure safer management of such complications because of clear surgical visibility and a selective ability to cut. The aim of this report is to present a case of an implant-related oral surgery complication that was resolved by means of a piezosurgery technique.

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## → ENDO SURGERY

Abella F, de Ribot J, Doria G, Duran-Sindreu F, Roig M.

## Applications of piezoelectric surgery in endodontic surgery: a literature review.

EN J Endod. 2014 Mar;40(3):325-32.

→ Applications of piezoelectric surgery in endodontic surgery: a literature review.

**Introduction:** Piezosurgery (piezoelectric bone surgery) devices were developed to cut bone atraumatically using ultrasonic vibrations and to provide an alternative to the mechanical and electrical instruments used in conventional oral surgery. Indications for piezosurgery are increasing in oral and maxillofacial surgery, as in other disciplines, such as endodontic surgery. Key features of piezosurgery instruments include their ability to selectively cut bone without damaging adjacent soft tissue, to provide a clear operative field, and to cut without generating heat. Although piezosurgery instruments can be used at most stages of endodontic surgery (osteotomy, root-end resection, and root-end preparation), no published data are available on the effect of piezosurgery on the outcomes of endodontic surgery. To our knowledge, no study has evaluated the effect of piezosurgery on root-end resection, and only 1 has investigated root-end morphology after retrograde cavity preparation using piezosurgery.

**Methods:** We conducted a search of the PubMed and Cochrane databases using appropriate terms and keywords related to the use and applications of piezoelectric surgery in endodontic surgery. A hand search also was conducted of issues published in the preceding 2 years of several journals. Two independent reviewers obtained and analyzed the full texts of the selected articles.

**Results:** A total of 121 articles published between January 2000 and December 2013 were identified. This review summarizes the operating principles of piezoelectric devices and outlines the applications of piezosurgery in endodontic surgery using clinical examples.

**Conclusions:** Piezosurgery is a promising technical modality with applications in several aspects of endodontic surgery, but further studies are necessary to determine the influence of piezosurgery on root-end resection and root-end preparation.

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**Del Fabbro M, Tsesis I, Rosano G, Bortolin M, Taschieri S.**

## Scanning electron microscopic analysis of the integrity of the root-end surface after root-end management using a piezoelectric device: a cadaveric study.

**EN** J Endod. 2010 Oct;36(10):1693-7.

⇒ Scanning electron microscopic analysis of the integrity of the root-end surface after root-end management using a piezoelectric device: a cadaveric study.

**Introduction:** The introduction of piezoelectric instruments in endodontic surgery allowed clinicians to manage the bone tissue with precision reducing the risk of damaging soft tissues. Although such instruments can also be used to prepare root-end cavity, few information is available on the effect of piezosurgery on dentine. This cadaveric study investigated the integrity of root apices after root-end cavity preparation by piezoelectric instruments at different device settings.

**Methods:** Fifty maxillary anterior teeth underwent endodontic treatment and apical resection in situ. Ten teeth were used as control. Retrograde cavities were prepared in 40 teeth (10 per group). In three groups, the piezoelectric device was set at constant vibration mode (CV), and three power levels were used. In another experimental group ( $n = 10$ ), vibration + pulsation mode (VP) was selected, and low power was used. Each root was duplicated and analyzed by scanning electron microscopy for the presence of cracks and marginal chipping.

**Results:** The number of cracks was significantly higher in the VP group, whereas no significant difference was found among CV groups. No difference between groups was found for crack type. The VP group showed a significantly poorer quality of cavity margin respect to the CV groups.

**Conclusion:** Constant vibration mode is recommended for retrograde preparation with piezosurgery.

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## ➡ PERIODONTAL SURGERY

Seshan H, Konuganti K, Zope S.

### PIEZOSURGERY® in periodontology and oral implantology.

EN J Indian Soc Periodontol. 2009 Sep;13(3):155-6.

#### ⇒ Piezosurgery® in periodontology and oral implantology.

Periodontitis is a chronic inflammatory disease of the tooth-supporting structures. The treatment of this condition is largely based on the removal of local factors and restoration of the bony architecture. Moreover, in the era of modern dentistry, successful implant therapy often requires sound osseous support. Traditionally, osseous surgery has been performed by either manual or motor-driven instruments. However, both these methods have their own advantages and disadvantages. Recently, a novel surgical approach using piezoelectric device has been introduced in the field of periodontology and oral implantology. This article discusses about the wide range of application of this novel technique in periodontology.

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Vercellotti T, Pollack AS.

### A New Bone Surgery Device: Sinus Grafting and Periodontal Surgery.

EN Compend Contin Educ Dent. 2006; 27(5):319-325

#### ⇒ A New Bone Surgery Device: Sinus Grafting and Periodontal Surgery.

The Piezosurgery device is a new instrument developed specifically for osseous surgery that has applications in various dental and medical surgical specialties. Using low-frequency piezoelectric ultrasonic vibration, the Piezosurgery device precisely cuts bone without cutting soft tissue. This article presents 3 cases illustrating the use of the device in sinus grafting and periodontal osseous surgery.

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Menini I, Zornitta C, Menini G.

## Distraction Osteogenesis for Implant Site Development Using a Novel Orthodontic Device: A Case Report.

EN Int J Periodontics Restorative Dent. 2008; 28(2):189-196

### ⇒ Distraction Osteogenesis for Implant Site Development Using a Novel Orthodontic Device: A Case Report.

The case report describes the distraction osteogenesis of a small bone fragment supporting a maxillary central incisor that had become ankylosed and infraoccluded after trauma. After initial orthodontic treatment, surgery was performed, entailing a trapezoidal osteotomy with minimal periosteal elevation through a horizontal incision. The osteotomy was performed using a piezoelectric surgical device (Piezosurgery, Mectron) with the thinnest inserts. On this occasion, a device was applied, anchored to the tooth, and activated by a micrometer screw. On completion of distraction, the fragment was retained for 5 months. The ankylosed tooth was then removed, and an implant was inserted and provided with an immediate provisional prosthesis without occlusal loading. Five months later, a definitive prosthesis with occlusal loading was applied and is still in place 1 year later without any problems.

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The reconstruction of a maxillary anterior dento-alveolar defect in patients with trauma has been a challenge for surgeons. Extensive loss of bone and teeth in the anterior maxilla presents a complex problem for reconstruction. This is owing to the difficulty in achieving complete closure using over-lying soft tissue. Tension-free sutures cannot be used after a large bone graft because the overlying soft tissue on severe bone defects of the anterior maxilla is often deficient and is attached to the underlying atrophic bone by scarring. Distraction osteogenesis provides a method to regain both hard tissue and soft tissue without any grafting. We describe a patient who had severe maxillary anterior bony defects that were restored by means of piezoelectric distraction osteogenesis, followed by dental implant placement. Clinical, radiological, and histological results showed that the reconstruction was successful.

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EN Int J Oral Maxillofac Implants. 2007; 22(6):1012-1016

### → Piezoelectric Bone Surgery Applied in Alveolar Distraction Osteogenesis: A Technical Note.

Piezoelectric bone surgery is based on ultrasonic vibration of a device functioning as an osteotome. This device allows precise cuts to be made in bone structures without provoking lesions of adjacent soft tissues and at the same time offering excellent visibility within the surgical field. The use of this technique in alveolar distraction osteogenesis is described. Piezoelectric surgery appears to make the cutting of the transport segment easier and safer.

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## → MAXILLO-FACIAL SURGERY

Pappalardo S, Guarneri R.

## Randomized clinical study comparing piezosurgery and conventional rotatory surgery in mandibular cyst enucleation.

EN J Craniomaxillofac Surg. 2014 Jul;42(5):e80-5.

### → Randomized clinical study comparing piezosurgery and conventional rotatory surgery in mandibular cyst enucleation.

**Purpose:** The aim of the present study is to compare piezoelectric surgery to the conventional rotatory surgery in mandibular cyst enucleation, and to determine the 2 method's suitability and the postoperative outcomes.

**Material and Methods:** Eighty patients were included in the study. 35 male and 45 female, which showed cystic mandibular lesions, compared with the inferior alveolar nerve or the mental nerve. The patients were randomly divided into two groups. In the test group, cystectomy was performed using conventional rotatory instrumentation (rotatorygroup), and in the other one piezoelectric surgery (piezo-group). The swelling was documented 24/48/72 h and 1 week post-surgery and the patients recorded their subjective postoperative pain daily for 7 days using a Visual Analog Scale (VAS).

**Results:** Patients treated with piezoelectric technique have presented a lower VAS, minor swelling and less recovery time compared to the rotatory-group. No lesions of the mandible nerve were detected with piezosurgery whereas surgery with rotary instruments resulted in 8% hypesthesia at least up to one week.

**Conclusion:** The results of this study suggest that piezosurgery may be considered effective in cyst enucleation compared to traditional procedures with burs, since it grants the patients significantly less post-surgical pain and swelling.

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**D'Amato S, Sgaramella N, Vanore L, Piombino P, Orabona GD, Santagata M.**

### Piezoelectric bone surgery in the treatment of an osteoma associated with an impacted inferior third molar: a case report.

**EN** Clin Cases Miner Bone Metab. 2014 Jan;11(1):73-6.

⇒ **Piezoelectric bone surgery in the treatment of an osteoma associated with an impacted inferior third molar: a case report.**

Operative removal of impacted mandibular third molars is a common and not riskless surgical procedure. We present an emblematic case of an osteoma closely associated with an impacted third left mandibular molar treated by Mectron Piezosurgery medical ultrasonic device.

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**Rodriguez JG, Eldibany RM.**

### Vertical splitting of the mandibular body as an alternative to inferior alveolar nerve lateralization.

**EN** Int J Oral Maxillofac Surg. 2013 Sep;42(9):1060-6.

⇒ **Vertical splitting of the mandibular body as an alternative to inferior alveolar nerve lateralization.**

The aim of this study was to present and evaluate a modified technique to inferior alveolar nerve lateralization (IANL) that allows the placement of longer implants in the posterior mandibular region. One hundred and forty-three consecutive patients were enrolled in this study; these patients had between 1.8 and 8mm residual crestal height above the mandibular canal. Vertical splitting of the mandibular body was performed using piezoelectric surgery followed by bone expansion and insertion of special conical implants of 10 and/or 12mm in length. Two hundred and sixty-nine osteotomies were performed and 636 implants were inserted, with a survival rate of 99% at the end of 12 months. Immediately postoperative there was an alteration of sensation in the lip/chin area in 8.5% of cases; 4.1% regained full sensation within 10-14 days, 2.6% after 8 weeks, and 0.7% had persistent paresthesia that did not affect their daily activities. Progressively increasing pain and numbness was present in 1.1%; the implants were removed 6 months postoperatively. This is a relatively simple procedure that has no limitations in clinical situations with minimal bone height. It allows for greater implant stability, and the risk of neurological disturbance is minimal.

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## Piezosurgery versus conventional osteotomy in orthognathic surgery: a paradigm shift in treatment.

**EN** J Craniofac Surg. 2013 Sep;24(5):1763-6.

### ⇒ Piezosurgery versus conventional osteotomy in orthognathic surgery: a paradigm shift in treatment.

The aim of the study was to compare in a randomized controlled clinical trial the use of the piezoelectric osteotomy as an alternative to the conventional approach in terms of surgery time, intraoperative blood loss, cut quality, nerve injury, and costs. One hundred ten patients who had orthognathic surgery procedures with bimaxillary osteotomy were divided into 2 groups: group A was treated with a piezosurgery device, and group B, with a reciprocating saw and bur. The piezosurgical bone osteotomy permitted individualized cut designs. The surgical time in group A was reduced, with a mean for the mandibular osteotomy (1 side) between 3 minutes 31 seconds and 5 minutes 2 seconds, whereas in group B, the surgical time was between 7 minutes 23 seconds and 10 minutes 22 seconds. The surgical time in group A for the Le Fort I osteotomy was between 5 minutes 17 seconds and 7 minutes 55 seconds in group A and between 8 minutes 38 seconds and 15 minutes 11 seconds in group B. All patients in group A had a low blood loss (<300 mL) versus patients of group B who had a medium to high blood loss (medium loss: 400 mL, high loss: >500 mL). Inferior alveolar nerve sensation was retained in 98.2% of group A versus 92.7% in group B at 6 months postoperative testing. Piezoelectric osteotomy reduced surgical time, blood loss, and inferior alveolar nerve injury in bimaxillary osteotomy. Absence of macrovibrations makes the instrument more manageable and easy to use and allows greater intraoperative control with higher safety in cutting in difficult anatomical regions.

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Rana M, Gellrich NC, Rana M, Piffkó J, Kater W.

## Evaluation of surgically assisted rapid maxillary expansion with piezosurgery versus oscillating saw and chisel osteotomy - a randomized prospective trial.

EN Trials. 2013 Feb 17;14:49.

### ⇒ Evaluation of surgically assisted rapid maxillary expansion with piezosurgery versus oscillating saw and chisel osteotomy - a randomized prospective trial.

**Background:** Ultrasonic bone-cutting surgery has been introduced as a feasible alternative to the conventional sharp instruments used in craniomaxillofacial surgery because of its precision and safety. The piezosurgery medical device allows the efficient cutting of mineralized tissues with minimal trauma to soft tissues. Piezoelectric osteotome has found its role in surgically assisted rapid maxillary expansion (SARME), a procedure well established to correct transverse maxillary discrepancies. The advantages include minimal risk to critical anatomic structures. The purpose of this clinical comparative study (CIS 2007-237-M) was to present the advantages of the piezoelectric cut as a minimally invasive device in surgically assisted, rapid maxillary expansion by protecting the maxillary sinus mucosal lining.

**Methods:** Thirty patients (18 females and 12 males) at the age of 18 to 54 underwent a surgically assisted palatal expansion of the maxilla with a combined orthodontic and surgical approach. The patients were randomly divided into two separate treatment groups. While Group 1 received conventional surgery using an oscillating saw, Group 2 was treated with piezosurgery. The following parameters were examined: blood pressure, blood values, required medication, bleeding level in the maxillary sinus, duration of inpatient stay, duration of surgery and height of body temperature.

**Results:** The results displayed no statistically significant differences between the two groups regarding laboratory blood values and inpatient stay. The duration of surgery revealed a significant discrepancy. Deploying piezosurgery took the surgeon an average of 10 minutes longer than working with a conventional-saw technique. However, the observation of the bleeding level in the paranasal sinus presented a major and statistically significant advantage of piezosurgery: on average the bleeding level was one category above the one of the remaining patients.

**Conclusion:** This method of piezoelectric surgery with all its advantages is going to replace many conventional operating procedures in oral and maxillofacial surgery.

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IT Minerva Stomatol. 2012 May;61(5):225-31.

### ⇒ Use of piezosurgery device in management of oral surgery complications: clinical case and clinical experience report.

The root displacement into the maxillary sinus could be a complication of oral surgery in the upper jaw. In these cases, the root removal is needed in order to avoid the occurrence of sinus pathologies. Piezosurgery techniques could assure a safer management of such complications, because of the clear surgical visibility and the selective ability of cut. The aim of this article is to present a case of oral surgery complication (root displacement in the right maxillary sinus), in which piezosurgery technique helped for a correct and safe clinical management, allowing to reduce the soft tissue damage.

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Wagner ME, Rana M, Traenkenschuh W, Kokemueller H, Eckardt AM, Gellrich NC.

## Piezoelectric-assisted removal of a benign fibrous histiocytoma of the mandible: an innovative technique for prevention of dentoalveolar nerve injury.

EN Head Face Med. 2011 Oct 31;7:20.

### ⇒ Piezoelectric-assisted removal of a benign fibrous histiocytoma of the mandible: an innovative technique for prevention of dentoalveolar nerve injury.

In this article, we present our experience with a piezoelectric-assisted surgical device by resection of a benign fibrous histiocytoma of the mandible. A 41 year-old male was admitted to our hospital because of slowly progressive right buccal swelling. After further radiographic diagnosis surgical removal of the yellowish-white mass was performed. Histologic analysis showed proliferating histiocytic cells with foamy, granular cytoplasm and no signs of malignancy. The tumor was positive for CD68 and vimentin in immunohistochemical staining. Therefore the tumor was diagnosed as primary benign fibrous histiocytoma. This work provides a new treatment device for benign mandibular tumour disease. By using a novel piezoelectric-assisted cutting device, protection of the dentoalveolar nerve could be achieved.

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Robiony M, Polini F.

## Piezosurgery: a safe method to perform osteotomies in young children affected by hemifacial microsomia.

EN J Craniofac Surg. 2010 Nov;21(6):1813-5.

⇒ **Piezosurgery: a safe method to perform osteotomies in young children affected by hemifacial microsomia.**

Piezosurgery is a new innovating technique used to perform safe and effective osteotomies using piezoelectric ultrasonic vibrations. It was conceived by Tomaso Vercellotti, MD, DDS, and it was first reported for preprosthetic surgery, as alveolar crest expansion, sinus grafting, and, more recently, alveolar distraction. The purpose of this report was to introduce and justify the use of Piezosurgery to perform osteotomies during mandibular distraction in cases of hemifacial microsomia.

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⇒ **Intermaxillary fixation: technique and benefit for piezosurgical sagittal split osteotomy**

**Introduction:** The aim of this study was to assess piezosurgical sagittal split osteotomy with peroperative inter maxillary fixation.

**Material and method:** We studied 25 bimaxillary osteotomies, 50 sagittal split osteotomies performed with this technique. It included both maxillomandibular fixation during all the split osteotomy and performing split osteotomy in five steps. For each case, we noted the type of dysmorphia, the size of split osteotomy and the time required for surgery, along with common data such as sex, age, etc. The data was compared to results of a previous series of patients also operated with PIEZOSURGERY® but without peroperative maxillomandibular fixation.

**Results:** Using peroperative maxillomandibular fixation during piezosurgical bilateral sagittal osteotomy decreases the

length of surgery by 33%, allows 9 times out of 10 for complete splitting, including the basilar edge, has no adverse effect especially on orthodontic material.

**Discussion:** PIEZOSURGERY® is a great progress for orthognathic surgery because of its precision and ability to preserve soft tissues. But it requires modification of the usual technique for mechanical section. Using peroperative inter maxillary fixation during ultrasonic splitting is a remarkably effective and easy technical modification.

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⇒ **PIEZOSURGERY® for genioglossal advancement in treatment of obstructive sleep apnea syndrome.**

The tongue is a frequent cause of upper airway obstruction in obstructive sleep apnea syndrome. The genioglossi muscles are the front-most muscles of the tongue, attached to the geniotubercle on the inner cortical of the mandibular symphysis. The collapsibility of the tongue can be managed with genioglossus advancement and the obstructive sleep apnea syndrome decreased. Piezo surgery can decrease the risk of apical dental lesion. The perfect depth control during osteotomy and absence of soft tissue lesion considerably limit the risks of oral floor lesion.

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Landes CA, Stübinger S, Ballon A, Sader R.

## Piezosteotomia in orthognathic surgery versus conventional saw and chisel osteotomy.

**EN** Oral Maxillofac Surg. 2008 Sep;12(3):139-47.

### ⇒ Piezosteotomia in orthognathic surgery versus conventional saw and chisel osteotomy.

**Introduction:** Piezosteotomia was assessed as alternative osteotomy method in orthognathic surgery regarding handling, time requirement, nerve and vessel impairment.

**Materials and Methods:** In this comparative clinical experience, 90 patient's orthognathic surgery procedures were performed in typical distribution prospectively by piezosteotomia: 34 (38%) mono-segment, 47 (52%) segmented LeFort I osteotomies, 94 (52%) sagittal split osteotomies, 11 (12%) symphyseal, and 4 (2%) mandibular body osteotomies. As controls served 90 retrospective patients with conventional saw and chisel osteotomy: 58 (64%) monosegment, 27 (30%) segmented LeFort I osteotomies, 130 (72%) sagittal split, and 4 (4%) symphyseal osteotomies.

**Results and discussion:** Piezosteotomies were individually designed to interdigitate the jaw segments after repositioning. The pterygomaxillary suture weakened angulated tools; auxiliary chisels were required in 100% of cases for the nasal septum and lateral nasal walls, in 33% for pterygoid processes. The dorsal maxilla as the pterygoid process were easily reduced; 15% mandibular osteotomies required sawing, while the lingual dorsal osteotomy was performed by manual feedback due to limited visibility. Bloodloss decreased from average  $537 \pm 208$  ml vs.  $772 \pm 338$  ml ( $p = 0.0001$ ). Operation time remained unchanged:  $223 \pm 70$  min vs.  $238 \pm 60$  min ( $p = 0.2$ ) for a conventional bimaxillary procedure. Clinical courses and reossification were unobtrusive. Alveolar inferior nerve sensitivity was retained in 98% of the piezosteotomia collective versus 84% of controls ( $p = 0.0001$ ) at 3 months postoperative testing.

**Conclusion:** Piezoelectric osteotomy did not prolong the operation and reduced blood loss as alveolar nerve impairment. A few patients required additional sawing or chisel. Piezoelectric screw insertion as complex osteotomies may be initiated to simplify the procedure and increase segment interdigitation after repositioning as to minimize the osteofixation time and dimensions.

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## Critical evaluation of piezoelectric osteotomy in orthognathic surgery: operative technique, blood loss, time requirement, nerve and vessel integrity.

EN J Oral Maxillofac Surg. 2008 Apr;66(4):657-74.

⇒ **Critical evaluation of piezoelectric osteotomy in orthognathic surgery: operative technique, blood loss, time requirement, nerve and vessel integrity.**

**Purpose:** Piezo-osteotomy feasibility as a substitute for the conventional saw in orthognathic surgery was evaluated regarding operative technique, blood loss, time requirement, and nerve and vessel integrity.

**Patients and methods:** Fifty patients had orthognathic surgery procedures in typical distribution using piezosurgical osteotomy: 22 (44%) monosegment, 26 (52%) segmented Le Fort I osteotomies; 48 (48%) sagittal split osteotomies, 6 (12%) symphyseal, and 4 (4%) mandibular body osteotomies. Controls were 86 patients with conventional saw and chisel osteotomies: 57 (66%) monosegment, 25 (29%) segmented Le Fort I osteotomies, 126 (73%) sagittal split, and 4 (5%) symphyseal osteotomies.

**Results:** Piezosurgical bone osteotomy permitted individualized cut designs, enabling segment inter-digititation after repositioning. Angulated tools weakened the pterygomaxillary suture; auxiliary chisels were required in 100% of cases for the nasal septum, and lateral nasal walls as 46% pterygoid processes. After downfracture, the dorsal maxillary sinus wall and pterygoid processes were easily reduced. Hemorrhage was successfully avoided with average blood loss of 541 +/- 150 mL versus 773 +/- 344 mL ( $P = .001$ ) for a conventional bimaxillary procedure. Sagittal mandibular osteotomy required considerable time (auxiliary saw in 13%); the lingual dorsal osteotomy was mostly performed tactile. Time investment remained unchanged: 227 +/- 73 minutes per bimaxillary standard osteotomy versus 238 +/- 61 minutes ( $P = .5$ ); clinical courses and reossification were unobtrusive. Alveolar inferior nerve sensitivity was retained in 95% of the study collective versus 85% in the controls ( $P = .0003$ ) at 3 months postoperative testing.

**Conclusions:** Piezoelectric osteotomy reduced blood loss and inferior alveolar nerve injury at no extra time investment. Single cases require auxiliary chiseling or sawing. Piezoelectric drilling for screw insertion and complex osteotomy designs may be developed to maintain bone contact or interdigititation after repositioning and minimize need for osteofixation.

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J. González-Lagunas, J. Mareque

## Piezosurgery®: its role in TMJ surgery.

EN Ital J Maxillofac Surg. 2008 December;19(3):119-21.

### ⇒ Piezosurgery®: its role in TMJ surgery.

This article presents the authors' experience with the use of a piezoelectric device to perform different osteotomies during open temporomandibular (TMJ) surgery. In the last six months six patients underwent TMJ open surgery with a piezoelectric device in order to perform osteotomies both in the condyle and in the articular eminence. This surgical instrument allows to trace an accurate osteotomy line, without injuring the adjacent soft tissues.

Piezosurgery® has gained popularity during the last few years due to the accurate and safe cutting of bone achieved without the risk of injuring the neighboring soft tissues. Most of the applications in oral and maxillofacial surgery have included its use incases when bone cutting was performed near nerves or other weak structures as the sinus membrane. The aim of this paper is to present authors' experience in six cases of TMJ surgery, where piezosurgery allows a safe bone cutting in the vicinity of important vascular structures.

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Beziat JL, Béra JC, Lavandier B, Gleizal A.

## Ultrasonic osteotomy as a new technique in craniomaxillofacial surgery.

EN International Journal of Maxillo-facial Surgery, 2007;36(6):493-500.

### ⇒ Ultrasonic osteotomy as a new technique in craniomaxillofacial surgery.

Ultrasound osteotomy is a new surgical technique used in dentistry to section hard tissues without damaging adjacent soft tissues. It was hypothesized that this could also be useful in craniofacial and orthognathic surgery. An ultrasonic device was employed in the following craniofacial surgical procedures: 144 Le Fort I osteotomies, 140 palatal expansions after Le Fort I osteotomies and 140 bilateral sagittal osteotomies; 2 Le Fort III osteotomies for treatment of Crouzon syndrome in two patients; 12 cases of unicortical calvarial bone grafting; removal of superior orbital roof in 25 cases of craniofaciostenosis; removal of external wall of the orbit in 10 cases of orbital cavity tumour; removal of anterior and posterior walls of the frontal sinuses in four cases of orbital cavity tumour. Integrity of soft tissues and surgical time were evaluated. Functional results were good without any soft-tissue damage being observed, but the overall operative time was increased. Ultrasound osteotomy is a new technical procedure that is advantageous for bone cutting in multiple situations, with minimal to no damage in adjacent soft tissues such as brain, palatal mucosa and the inferior alveolar nerve.

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**Robiony M, Polini F, Costa F, Zerman N, Politi M.**

## **Ultrasound bone cutting for surgically assisted rapid maxillary expansion (SARME) under local anaesthesia.**

**EN** Int J Oral Maxillofac Surg. 2007;36(3):267-9.

⇒ **Ultrasound bone cutting for surgically assisted rapid maxillary expansion (SARME) under local anaesthesia.**

Surgically assisted rapid maxillary expansion (SARME) is a well-established therapy for correction of maxillary transverse deficiency in adults, when consolidation of sutures has just been completed. It can be performed either under general or under local anesthesia and it can be accomplished with many surgical techniques. One of the most critical steps of SARME is the detachment of the pterygo-maxillary junction, due to the risks connected to such procedure. When required to obtain specific expansion patterns, the pterygo-maxillary separation has been suggested until now only for interventions under general anesthesia, due to the dangerousness and the rawness of this surgical step in awake patients. The authors introduce the use of an ultrasonic bone-cutting device to perform all osteotomic steps of SARME under local anesthesia on an outpatient basis, including pterygo-maxillary detachment. This ultrasonic device is unique in that the osteotomic action occurs only when the tool is employed on mineralized tissues, while it stops on soft tissues. It works in a linear pattern of vibration and it allows precise osteotomies without producing any heat damage to osteotomic surfaces and without any dangerous hammer-related stroke. Due to its precision and safety, this device named Piezosurgery®, allows patients to undergo all the steps of SARME under local anesthesia, also without hospitalization.

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**Gleizal A, Béra JC, Lavandier B, Béziat JL.**

## **Piezoelectric osteotomy: a new technique for bone surgery – advantages in craniofacial surgery.**

**EN** Childs Nerv Syst . 2007;23(5):509-513.

⇒ **Piezoelectric osteotomy: a new technique for bone surgery – advantages in craniofacial surgery.**

**Introduction:** Ultrasonic bone cutting is a new surgical technique used in dentistry to section hard tissues without damaging adjacent soft tissues. We hypothesized that such a device could also be useful in craniofacial surgery, particularly during the removal of the superior orbital roof during cranio-faciostenosis surgery.

**Materials and Methods:** An ultrasonic device was employed in different craniofacial surgical procedures: 1. to remove the superior orbital roof in 30 cases of craniofaciostenosis, 2. to perform a Le Fort III osteotomy for the treatment of Crouzon syndrome in two patients, 3. to cut the parietal and frontal bone in 30 cases of craniofaciostenosis. The integrity of soft tissues and surgical time was evaluated.

**Results:** Functional results were good without any soft tissue damage appreciated. The overall operative time, however, was increased.

**Conclusions:** PIEZOSURGERY® is a new technical procedure, which can be advantageous for bone cutting in multiple situations with minimal to no damage in adjacent soft tissues.

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■ G Chir 2007 May; 28(5):222-6.

⇒ **Mandibular enostosis removal with PIEZOSURGERY®: case report.**

The bone surgery has always used manual and rotary instruments. The biomedical engineering with ultrasound working principle has given a new surgery instruments, which allow a selective cutting action of bone tissue and the protection of soft tissue. Our case shows an uncommon clinical lesion surgically dangerous for the narrow adjoining of important anatomical structures as the lower alveolar artery and the inferior alveolar nerve. The clinical result and recovery time go toward a smaller traumatic situation of this methodology of the cutting of bone tissue.

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⇒ What is PIEZOSURGERY®, Two-years experience in craniomaxillofacial surgery

**Introduction:** PIEZOSURGERY® is a new surgical technique used in dentistry to section hard tissues without damaging adjacent soft tissues. We hypothesized that such a device could also be useful in craniofacial and orthognathic surgery.

**Material and method:** An ultrasonic device (Mectron) was employed in different craniofacial surgical procedures: a) to perform 144 Le Fort I osteotomies, 140 palatal expansions after Le Fort I osteotomies, and 134 bilateral sagittal osteotomies; b) to perform a Le Fort III osteotomy for treatment of Crouzon syndrome in 2 patients; c) to perform 5 segmental osteotomies and 3 osteotomies of the inferior edge of the mandible for facial asymmetry; d) to perform 12 cases of unicortical calvarial bone grafting; to remove the superior orbital roof in 20 cases of craniofaciostenosis and the frontal bone in 5 cases; e) to remove the external wall of the orbit or the anterior and posterior wall of the frontal sinuses in 10 cases of orbital cavity tumors; f) to approach the skull base through the frontal sinuses in 4 cases.

**Results:** Analysis of the results showed that PIEZOSURGERY®: a) allows very precise cutting; avoids bone cutting using an osteotome; b) spares soft tissue such as brain, dura-mater, palatal mucosa; c) preserves the nerves, especially those of the inferior alveolar; d) increases the time of bone cutting but not the overall operative time because of the absence of soft tissue protection.

**Discussion:** PIEZOSURGERY® is a new technical procedure, which can be particularly advantageous for bone cutting in maxillo facial surgery, sparing adjacent soft tissues such as brain, palatal mucosa, and the inferior alveolar nerve from any damage. The device's lack of power appears to be a minor problem compared with the advantages.

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**Guo ZZ, Liu X, Li Y, Deng YF, Wang Y.**

## The use of PIEZOSURGERY® osteotomy in treatment of longstanding maxillary fractures: report of 12 consecutive patients.

**EN** Shanghai Kou Qiang Yi Xue. 2007 Feb;16(1):97-9.

⇒ The use of PIEZOSURGERY® osteotomy in treatment of longstanding maxillary fractures: report of 12 consecutive patients.

**Purpose:** To evaluate the clinical applicability of PIEZOSURGERY® osteotomy: a new safe technique in managing long standing maxillary fractures.

**Methods:** 12 patients with long-standing maxillary fractures were surgically treated using Le Fort I osteotomy. During operation, PIEZOSURGERY® osteotomy was used for bone cutting and splitting. After repositioning, the bone segments were rigidly fixed with micro Ti-plate, Ti-mesh. All the patients were followed up for 6 to 12 months, and the functional and esthetic results were evaluated.

**Results:** Ultrasonic microvibrations allow accurate bone cutting without oscillating injuries to the soft tissue. All the wounds healed primarily without complications. The postoperative occlusion and appearance were satisfactory.

**Conclusions:** Maximal recovery of mastication and appearance can be achieved by using PIEZOSURGERY® osteotomy with fixation materials such as Ti-plates and Ti-meshes in selected patients with long-standing maxillary fractures.

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Robiony M, Polini F, Costa F, Vercellotti T, Politi M.

## Piezoelectric Bone Cutting in multipiece maxillary osteotomies. Technical Note.

EN J Oral Maxillofac Surg. 2004; 62:759-761.

### ⇒ Piezoelectric Bone Cutting in multipiece maxillary osteotomies. Technical Note.

PIEZOSURGERY® (patented by Mectron Medical Technology, Carasco [Genova], Italy.) is a new innovating technique used to perform safe and effective osteotomies using piezoelectric ultrasonic vibrations. It was first reported for preprosthetic surgery, alveolar crest expansion, and sinus grafting. We introduce and report the use of PIEZOSURGERY® for multipiece maxillary osteotomies, to overcome many of the complications of this delicate surgery on hard and soft tissues. Because of its micrometric and selective cut, the piezoelectric device produces safe and precise osteotomies without any osteonecrosis damage. This device works only on mineralized tissues, sparing soft tissues and their blood supply. Among various surgical phases, osteotomy is one of the most technique-sensitive procedures in maxillo-facial surgery. Osteotomies are usually conducted close to delicate anatomic structures, such as vestibular and lingual/palatal soft tissues that provide bone vascularization through the periosteum. Furthermore, bone is a hard tissue and many cutting or drilling osteotomes are very crude tools. In particular, rotating instruments are potentially injurious, due to the production of excessively high temperatures during osseous drilling, which can produce marginal osteonecrosis and impair bony regeneration. It is widely and strongly recommended to use a careful surgical technique and to reduce the amount of frictional heating with saline solution irrigation. Among cutting techniques, PIEZOSURGERY® is a new and innovating method that uses piezoelectric ultrasonic vibrations to perform precise and safe osteotomies, due to its characteristics of a micrometric and selective cut. It was first invented by Tomaso Vercellotti, MD, DDS, to overcome the limits of traditional instruments in oral bone surgery. The purpose of this report was to present the use of the piezoelectric cut in segmental maxillary Le Fort I osteotomy, a field in which effectiveness, precision, and safety of osteotomies are of paramount importance.

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Eggers G, Klein J, Blank J, Hassfeld S

## PIEZOSURGERY®: an ultrasound device for cutting bone and its use and limitations in maxillofacial surgery.

**EN** Br J Oral Maxillofac Surg. 2004 Apr;42(5):451-3.

⇒ **PIEZOSURGERY®: an ultrasound device for cutting bone and its use and limitations in maxillofacial surgery.**

PIEZOSURGERY® uses modulated ultrasonic vibration to allow controlled cutting of bony structures. Delicate bony structures can be cut easily and with great precision, without destruction of soft tissue. We have found this device useful when exact cutting of thin bones is essential. However, it is of only limited use in cutting thick bones and in regions with limited access.

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