EN

IT


DE

FR

ES

- Cleaning and sterilization of reusable inserts
- Pulizia e sterilizzazione di inserti pluriuso
- Reinigungs- und Sterilisationsanleitung für wiederverwendbare Instrumente
- Nettoyage et stérilisation des inserts réutilisables
- Limpieza y esterilización de los insertos reusables

PIEZOSURGERY®

→ Cleaning and sterilization of reusable inserts

PIEZOSURGERY®

→ TABLE OF CONTENTS

01 → INTRODUCTION	2
-------------------------	---

EN

02 → CLEANING AND STERILIZATION	2
---------------------------------------	---

02.1 Pre-cleaning with enzymatic detergent	3
--	---

02.2 Ultrasonic bath or Automated cleaning and thermodisinfection ...	5
---	---

02.3 Inspection and drying	6
----------------------------------	---

02.4 Sterilization	7
--------------------------	---

01 → INTRODUCTION

This manual describes the cleaning and sterilization procedures to be applied on reusable bone surgery inserts produced by Mectron.

Carefully read this manual before proceeding with the cleaning and sterilization operations, and always keep it within reach.

IMPORTANT: to prevent harm to persons or damages to objects, read all the “Safety precautions” present in the manual with special care. Depending on their degree of seriousness, the safety precautions are classified with the following indications:

WARNING
(always refers to personal injury)

CAUTION
(refers to possible damage to property)

02 → CLEANING AND STERILIZATION

The operator must perform a cleaning cycle with enzymatic detergent (PRE-CLEANING) as first step. After that it can be decided to perform the:

- Ultrasonic bath
- Automated cleaning and thermodisinfection.

The procedures applied on PIEZOSURGERY® accessories have been validated by an independent institute.

WARNING: First use: The inserts are delivered in NON STERILE condition and must be prepared prior the use by applying the procedures described in this manual.

Every use: Once used, each reusable insert must be thoroughly reprocessed prior to reuse, according to the procedures described in this manual.

WARNING: Cleaning processes must immediately be carried out after each use. Do not allow contaminated inserts to dry prior to reprocessing.

Limitations on reprocessing: Repeated processing has minimal effect on the instruments examined in this manual. End of life is normally determined by wear and damage due to use.

PREPARATION

- Unscrew the insert from the handpiece with the torque wrench.


⚠ **CAUTION:** Always unscrew the insert from the handpiece before proceeding to its cleaning and sterilization.

MATERIALS NECESSARY

- Enzymatic detergent with pH 6-9
- Water
- Container for immersion in enzymatic solution
- Brush with soft nylon bristles
- Compressed air
- Ultrasound tank
- Distilled water
- Single-use sterilization pouches
- Steam autoclave
- Syringe
- Thermodisinfectant

⚠ **CAUTION:** Do not use hydrogen peroxide. If you intend to disinfect the inserts, do not use hydrogen peroxide, but only disinfectants with neutral pH (pH7); always rinse with sterile water to preserve the disinfection.


⚠ **CAUTION:** Once used, correctly dispose of the enzymatic detergent solution, do not recycle.

02.1 → PRE-CLEANING WITH ENZYMATIC DETERGENT

1

Prepare a mild pH (6-9) enzymatic detergent solution* according to the manufacturer's specifications;

* Enzymatic detergent Enzymec (Mectron), process validated by an independent institute.


2

Lay the inserts horizontally in a clean container;
Add a sufficient amount of the prepared enzymatic detergent solution to the container to cover the inserts;

02 → PRE-CLEANING


Soak the inserts for 10 minutes at 40°C ±2°C. This procedure reduces the quantity of blood, protein, and mucous present in the insert;

3


During immersion into the enzymatic solution, delicately brush the surface of the insert with a brush with soft nylon bristles to eliminate all traces of visible dirt. Accurately clean the difficult areas such as sharp edges and especially the interstices between the cutting cusps;

4


Thoroughly rinse and scrub the insert with a with soft bristled nylon brush under running tap water;

5


By using a disposable syringe inject and aspirate the enzymatic detergent solution into hard-to-reach areas (lumen/cannulae); Repeat this operation three times for effective removal of soil from the inner lumen surface;

6


Flush the internal channel of the insert with distilled water injected at a pressure of 3.8 bars for at least 10 seconds so as to remove any residues;

7


02.2 → ULTRASONIC BATH OR AUTOMATED CLEANING AND THERMODISINFECTION

ULTRASONIC BATH


1

Place the insert tip/s in a mesh bottom stainless steel instrument tray;
Place the tray into the ultrasonic cleaner and ensure that the instruments in the bath, at 40°C, are submerged by the enzymatic solution;
Sonicate for at least 10 minutes;


2

Remove the insert from the ultrasonic tank and rinse with distilled water;

3

Scrub the inner and outer surfaces of the insert with a clean soft bristled nylon brush under running tap water;
Rinse with distilled water;

4

Flush the internal channel of the insert with distilled water injected under pressure (3,8 bar) for at least 10 seconds so as to remove any residues;

EN

AUTOMATED CLEANING AND THERMODISINFECTION

① **CAUTION:** The automated cleaning cycle and the thermodisinfection must always be preceded by the pre-cleaning steps with enzymatic detergent.

Lay the reusable inserts in a metallic tray and place it in thermodisinfector;

NOTE: Place the instruments in the washing machine so that dead zones do not arise and the water can properly drain. Also, make sure that the reusable inserts are properly held in place in the washing basket and cannot move during the washing process, as shocks could damage them.

⚠ **WARNING:** Avoid overloading. Failure to comply overloading compromises cleaning effectiveness.

Sequence and parameters applicable to the cycle:

- 1 min, Pre-washing with cold water;
- 5 min, Washing with detergent at 55°C ±2°C;
- 2 min, Neutralization with the appropriate agent;
- 1 min, Rinsing with cold water;
- 5 min, Thermodisinflect at 93°C with demineralized water following the national requirements about A0 values (thermal disinfection at a temperature of 90°C for 5min results in a A0 value of 3000).


Validated with Miele G78820 or PG8535, Program Vario TD; Detergent: Neodisher FA (cleaning agent 0.2% v/v) and Neodisher Z (neutralization agent 0.1% v/v).

02.3 → INSPECTION AND DRYING


Once you've finished the cleaning operations, check all the inserts under an adequate source of light, paying attention to the details that could conceal dirt residues (threading, cavities, grooves, cutting edges) and, if need be, repeat the cleaning cycle;

① **CAUTION:** Before starting the sterilization cycle, make sure that the insert is thoroughly dry, both externally and internally. For this purpose, blow compressed air both externally and through the internal thru hole; this will prevent the onset of stains and haloes on the surface or rust inside the insert.


1

2

02.4 → STERILIZATION

1

Individually seal the inserts in single-use sterilization pouches.
Sterilize the inserts in a steam autoclave.

The sterilization process in steam autoclave must be performed setting the parameters indicated below, which are validated by Mectron and reach the sterilization safety levels (SAL 10⁻⁶):

- **Cycle type:** 3 times Pre-vacuum (pressure min 60 mBar).
- **Sterilization temperature:** 132°C (range 0°C ÷ +3°C).
- **Minimum sterilization time:** 4 min.
- **Minimum drying time:** 20 min.

All the sterilization phases must be performed by the operator in compliance to the norms UNI EN ISO 17665-1:2007, UNI EN ISO 556-1:2002 e ANSI/AAMI ST 46:2002.

⚠ **CAUTION:** Do not sterilize the insert screwed onto the handpiece.

⚠ **WARNING: Checking infections - Parts that can be sterilized** - Diligently remove all residues of organic dirt before the sterilization.

⚠ **CAUTION:** Perform the sterilization using a water steam autoclave only. Do not use any other sterilization procedure (dry heat, irradiation, ethylene oxide, gas, low temperature plasma, etc.).

Do not use the following substances/ systems to sterilize the inserts: hydrogen peroxide, peracetic acid systems, formaldehyde sterilization, glutaraldehyde sterilization or comparable sterilization solutions/system

⚠ **CAUTION:** Do not exceed the allowed load of the steam sterilizer.

⚠ **CAUTION: Risk of contamination.** Do not use a gravity displacement autoclave to sterilize the inserts. The operating cycle of a gravity displacement autoclave do not ensure adequate sterilization of the lumens, cavities, recessed features of the instruments.

The indications that appear in this publication are not binding and can be modified without fore-notice.
The Italian version of this manual is the original document from which its translations have been obtained. In case of any discrepancy, the Italian version will have pertinence.

Texts, images, and graphics of this manual are property of Mectron S.p.A., Carasco, Italy. All rights reserved.
The contents cannot be copied, distributed, changed, or made available to third parties without the written approval of Mectron S.p.A.

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