EN

IT

DE

FR

ES

- Cleaning and sterilization manual
- Manuale di pulizia e sterilizzazione
- Reinigungs- und Sterilisationsanleitung
- Manuel de nettoyage et de stérilisation
- Manual de limpieza y esterilización

PIEZOSURGERY®

→ Cleaning and sterilization manual

PIEZOSURGERY®

→ TABLE OF CONTENTS

EN

01 →	INTRODUCTION	3
02 →	CLEAN/FLUSH FUNCTION	4
03 →	DISASSEMBLY OF PARTS FOR CLEANING AND STERILIZATION	7
04 →	CHART OF REFERENCE - CLEANING TYPOLOGY.....	10
05 →	CLEANING OF THE CONSOLE AND OF THE FOOTPEDAL	11
06 →	CLEANING OF THE ACCESSORIES.....	12
06.1	MANUAL/AUTOMATED CYCLE.....	13
06.1.1	Manual prewashing and cleaning with enzymaticm detergent	13
06.1.2	Thermodisinfection (optional).....	20
06.2	ONLY AUTOMATED CYCLE Complete cleaning cycle in thermodisinfectior	21
07 →	CLEANING CHECK	22
08 →	DRYING AND LUBRICATION	22
09 →	STERILIZATION	24

01 → INTRODUCTION

Carefully read this manual before proceeding with the cleaning and sterilization operations, and always keep it within reach.

IMPORTANT: to prevent harm to persons or damages to objects, read all the "Safety precautions" present in the manual with special care. Depending on their degree of seriousness, the safety precautions are classified with the following indications:

⚠ WARNING
(always refers to personal injury)

❗ CAUTION
(refers to possible damage to property)

The purpose of this manual is to make the operator knowledgeable of the safety precautions, the installation procedures, and the instructions for a correct use and maintenance of the device and its accessories. Use of this manual for purposes other than those strictly tied to cleaning and sterilization of the device is forbidden.

The information and illustrations in this manual are updated as of the date of issue reported on the last page.

Mectron is committed to continuously update its products with possible modifications to device components.

In case you uncover discrepancies between what described in this manual and the device in your possession, contact your Retailer or the After-Sales Service of Mectron for clarification and support.

Limitations on the number of cleaning and sterilization cycles: Cleaning and sterilization cycles repeated over time have minimum effects on the instruments examined in this Manual. The life expectancy of the products is usually determined by their consumption and by damages caused from their use.

⚠ WARNING: The operators who perform the cleaning and sterilization operations must be adequately protected and trained

⚠ WARNING: Infections control.

First use: All new and repaired accessories are supplied in NON STERILE conditions. Before use, and after each treatment, they must be cleaned and sterilised in strict compliance with the instructions given in the Cleaning and Sterilization Manual.

Subsequent uses: After every treatment, clean and sterilize all the reusable parts and accessories, following the instructions provided in the Cleaning and Sterilization Manual.

⚠ WARNING: The cleaning processes must be commenced immediately after each use. Do not allow that contaminated instruments dry up before starting the cleaning and sterilization process. To eliminate organic remains such as blood, bone, and other, use an enzymatic detergent with neutral pH (pH7), immediately after use.

❗ CAUTION: Do not use metallic brushes or abrasive sponges during the cleaning processes, because they damage the surfaces of the treated parts and the finishes of the inserts. Only use brushes with soft nylon bristles.

❗ CAUTION: The detergent solution must be completely removed from the treated parts to prevent accumulation of chemical residues.

❗ CAUTION: During cleaning procedures only use detergents with pH 6-9.

❗ CAUTION: If you intend to disinfect, we recommend that you use water-based disinfectant solutions with a neutral pH (pH 7). Alcohol-based disinfectant solutions and hydrogen peroxide are contra-indicated, because they can fade the color and/or damage the plastic materials. This also holds true for chemical products such as acetone and alcohol. Always rinse with sterile water to preserve the disinfection.

⚠ WARNING: Do not use tap water, unless where explicitly indicated.

⚠ WARNING: Reusable accessories that must be forwarded to an Authorized Mectron Service Center.

All the parts must be cleaned and sterilized in accordance to the procedures described in this Manual, before they are sent to an Authorized Mectron Service Center. Incorrectly treated and bio-contaminated parts will not be accepted by the Authorized Mectron Service Centers.

02 → CLEAN/FLUSH FUNCTION

The CLEAN/FLUSH function (according to the device used) enables to run a cleaning cycle on the irrigation circuit of the handpiece used during the treatment, by following the steps described in this chapter.

① CAUTION: Clean/flush function.

The CLEAN/FLUSH function must be used after every treatment, before starting the cleaning and sterilization procedures.

① CAUTION: If the tubes are not cleaned, the crystallization of the salts may seriously damage the device.


① CAUTION: The handpiece and the cord cannot be separated.

Close the flow of the adjuster.

Disconnect the flow adjuster from the irrigation bag and from the peristaltic pump tube or from the tube of the bone grafting kit, if used.


Dispose of the flow adjuster correctly, following the norms in force locally on the disposal of hospital wastes;

1


Prepare an enzymatic detergent solution with pH 6-9, as per the instructions of the manufacturer;


2

Immerse the end part of the peristaltic pump tube or the tube of the bone grafting kit, if used, in a container containing the enzymatic detergent solution prepared;

3


4

Position the handpiece, with or without insert, above a container to contain the liquid that will outflow during the cleaning cycle;

EN


PIEZOSURGERY® touch

5a

To enter the cleaning mode, select PUMP/ CLEAN on the keyboard: all the other selection options present on the display are disabled, and the PUMP/CLEAN writing flashes;


PLEASE NOTE You can exit the cleaning mode by pressing on PUMP/CLEAN again or, alternatively, by waiting for at least 15 seconds. The CLEAN function is disabled and the keyboard is enabled again, and displays the last setting used;

5b

To make the CLEAN function start while the PUMP/CLEAN writing is flashing, press the foot pedal once and release it: PUMP/ CLEAN stops flashing and the cleaning cycle starts;

PLEASE NOTE The Clean function can be interrupted at any time by pressing the PUMP/CLEAN writing again or, alternatively, by pressing on foot pedal. The keyboard is enabled again, and displays the last setting used;


5c

As soon as the peristaltic pump starts, the entire scale of values of the "irrigation" section lights up and during the liquid passage, the value of the irrigation shifts from 6 a 0 . The cycle lasts 20 seconds. Once it has ended, the keyboard is enabled again, and displays the last setting used; Go to point 6

PIEZOSURGERY® *white*


To enter the cleaning mode keep the PUMP/FLUSH button on the touch keyboard pressed; the LED of the level 6 of the function Irrigation flashes; all other buttons are disabled. To stop the FLUSH function release the button; Go to point 6;

5a


To perform a cleaning cycle of every handpiece, tube, or connector used during the same treatment, repeat the operations described above, from point 5 onward;

6

IMPORTANT After the cleaning cycle with the enzymatic detergent solution, perform a new cleaning cycle with distilled water to remove detergent residues.

7


Once you have finished cleaning the irrigation circuit, move on to the cleaning and sterilization of each single part.

03 → DISASSEMBLY OF PARTS FOR CLEANING AND STERILIZATION

Before moving on to the cleaning procedures described in the following chapters, disconnect all the accessories and components of the PIEZOSURGERY®.


EN

1


⚠ WARNING: Switch the device off. Always turn the device off via the switch, and disconnect the power supply cable from the power outlet and from the device console, before performing the cleaning and sterilization interventions.

⚠ WARNING: Supplementary grounding cable. If present, disconnect the supplementary grounding cable before you perform the cleaning and sterilization interventions.

2

Remove the handpiece fixed holder and the irrigation bag support rod from the device console;

3

Disconnect the foot pedal from the device: Grab the connector of the foot pedal, press the release flap, and pull the connector back;

⚠ CAUTION: Do not attempt to unscrew or turn the connector during the disconnection: the connector could get damaged.

⚠ CAUTION: When disconnecting the foot pedal, always and only hold the connector of the cord. Never actually pull on the cord itself.

Disconnect the handpiece from the device;

⚠ **CAUTION:** To avoid damaging the handpiece, disconnect it by always holding onto the connector only. Never pull the cord.




⚠ **CAUTION:** Do not attempt to unscrew or turn the connector when disconnecting the handpiece. The connector could get damaged.

⚠ **CAUTION:** The handpiece and cord cannot be separated.


Completely open the peristaltic pump and remove the peristaltic pump tube or the tube of the bone grafting kit (if present). Close the pump again;

Disconnect the peristaltic pump tube from the connector of the handpiece;


4


5


6

7


If the bone grafting kit has been used, disconnect the tube-to-tube connector from the peristaltic pump tube and from the tube of the bone grafting kit;

8

If present, unscrew the insert from the handpiece by using the torque wrench;


⚠ WARNING: During the tightening and removing operations, the user must pay particular attention to avoiding injury from inserts with sharp points and cutting edges.

9

Unscrew the front cone of the handpiece;

NOTE - USE OF THE HANDPIECE WITH LIGHT: on the cone with LED light the metal ring cannot be separated from the cone in plastic;

10

Remove the protective foil from the keyboard, if present.

04 → CHART OF REFERENCE CLEANING TYPOLOGY

EN

This chart serves merely as a quick reference. For the complete cleaning and sterilization procedures of the individual parts, refer to the following chapters.

Part	Enzymatic cleaning	Wash with ultrasounds	Thermo-disinfectant	Steam sterilization
Device Console	YES	NO	NO	NO
Foot pedal	YES	NO	NO	NO
Handpiece with cord	YES	NO	YES	YES
Inserts	YES	YES	YES	YES
Torque wrench	YES	YES	YES	YES
Peristaltic pump tube	YES	NO	YES	YES
Tube-to-tube connector	YES	YES	YES	YES
Irrigation bag support rod and handpiece fixed support	YES	NO	YES	YES
Handpiece mobile support	YES	NO	YES	YES

05 → CLEANING OF THE CONSOLE AND OF THE FOOT PEDAL

EN

PREPARATION

- Make sure that all the following accessories have been removed/disconnected from the device console (see Chapter 03):
 - electrical power supply cable
 - foot pedal
 - peristaltic pump tube
 - handpiece
 - handpiece fixed support
 - irrigation bag support rod
 - protective foil (if present)
 - supplementary grounding connector (if present)

⚠ WARNING: Do not sterilize the device console and the foot pedal: it could stop working and cause harm to person and/or damages to things.

⚠ WARNING: Always turn the device off via the switch and disconnect it from the electrical network before performing the cleaning and disinfection interventions.

⚠ WARNING: The device is not protected against penetration of liquids. Do not spray liquids directly onto the surface of the device.

⚠ CAUTION: Do not clean the device console and the foot pedal under running tap water.


⚠ CAUTION: Do not immerse the device console and the foot pedal in liquids and/or solutions of various nature.

MATERIALS NECESSARY

- Clean, soft cloths with low fiber-release
- Detergent solution (pH 6-9) and, if need be, non-aggressive disinfectant solution with neutral pH (pH7)

⚠ CAUTION: If you wish to disinfect, we recommend that you use water-based disinfectant solutions with neutral pH (pH7). Alcohol-based disinfectant solutions and hydrogen peroxide are contraindicated, because they can fade the color and/or damage the plastic materials. This also holds true for chemical products such as acetone and alcohol.

CLEANING METHOD


1

Clean the surface of the device with a clean, soft cloth with low fiber-release, dampened with a detergent solution (pH 6-9). If you wish to disinfect, use a non-aggressive disinfectant solution with neutral pH (pH7), following the instructions provided by the manufacturer of the solution;

Dry the device console and the foot pedal with a clean, non-abrasive cloth with low fiber-release.

2

06 → CLEANING OF THE ACCESSORIES

PREPARATION

- Perform the CLEAN/FLUSH function (see Chapter 02)
- Disconnect the handpiece from the device console
- Disconnect the peristaltic pump tube from the connector of the handpiece
- Remove the tube from the peristaltic pump. If the bone grafting kit has been used, disconnect the tube-to-tube connector from both tubes
- If present, disconnect the insert from the handpiece with the torque wrench
- Unscrew the front terminal of the handpiece
- Disconnect the irrigation bag support rod from the device console

Refer to Chapter 3 for the disconnection of all accessories.

- ⚠ **CAUTION:** The handpiece and the cord cannot be separated.
- ⚠ **CAUTION:** Disconnect the handpiece from the device exclusively acting on the connector.
- ⚠ **CAUTION:** Always disconnect the insert from the handpiece before proceeding to its cleaning and sterilization.
- ⚠ **CAUTION:** Do not immerse the handpiece in disinfectant solutions or other liquids, because it could get damaged.
- ⚠ **CAUTION:** Do not immerse the handpiece in an ultrasonic tank.

MATERIALS NECESSARY

- Enzymatic detergent with pH 6-9
- Detergent solution (pH 6-9) and, if need be, non-aggressive disinfectant solution with neutral pH (pH7)
- Water
- Container for the immersion in the enzymatic solution
- Ultrasound tank
- Clean, soft cloths with low fiber-release
- Brush with soft nylon bristles
- Syringe
- Distilled water
- Thermodisinfectant
- Compressed air
- Medical lubricant


- ⚠ **CAUTION:** If you wish to disinfect, we recommend that you use water-based disinfectant solutions with a neutral pH (pH7). Alcohol-based disinfectant solutions and hydrogen peroxide are contra-indicated, because they can fade the color and/or damage the plastic materials. This also holds true for chemical products such as acetone and alcohol. Always rinse with sterile water to preserve the disinfection.
- ⚠ **CAUTION:** Do not use hydrogen peroxide. If you wish to disinfect the accessories, do not use hydrogen peroxide, but only disinfectants with neutral pH (pH7); always rinse with sterile water to preserve the disinfection.

NOTE: Choose between the following cleaning methods.

In case the practice is not equipped with a thermoisinfector, the cleaning process can begin with the manual pre-washing (paragraph 06.1.1) and proceed with the sterilization (Chapter 09).

In case the practice is equipped with a thermoisinfector the pre-washing can be performed manually (paragraph 06.1.1) or can be automated (Chapter 06.2). The final step is always the sterilization (Chapter 09).


⚠ **CAUTION:** The following instructions have been validated by the medical device manufacturer as ABLE to reprocess a medical device for being reused. It is up to the process responsible to assure that the repeated processes are actually performed by means of equipment, materials and personnel structured as repeated processes in order to obtain the desired result. Generally this requires the validation and the systematic monitoring of the process. Similarly, in case the person in charge of the process deviates from the supplied instructions, he shall properly evaluate any action and their efficiency and the possible undesired consequences.

06.1 → MANUAL/AUTOMATED CYCLE**06.1.1 → MANUAL PRE-WASHING AND CLEANING WITH ENZYMATIC DETERGENT**

1

Prepare an enzymatic* detergent solution with pH 6-9, as per the instructions of the manufacturer;

⚠ **CAUTION:** Once used, correctly dispose of the enzymatic detergent solution, do not recycle.

HANDPIECE WITH CORD



2

Clean the surface of the handpiece and of its cord with a clean, soft cloth with low fiber-release, dampened with a detergent solution (pH 6-9) and, if you wish to disinfect use a non-aggressive disinfectant solution with neutral pH (pH7), following the instructions provided by the manufacturer of the solution;

Delicately clean the surface of the handpiece using the enzymatic detergent solution with the aid of a brush with soft nylon bristles, placing special care to the following areas:


- threading of the handpiece
- titanium stem
- front terminal in its external and internal parts

3

Accurately rinse with tap water to eliminate any detergent residues and perform the last rinse with distilled water;

4


Process validated by an independent institute with enzymatic detergent Enzymec (Mectron).


INSERTS

Place the insert horizontally in a clean container, and add a quantity of enzymatic detergent solution sufficient to completely cover the insert;

5


06 → CLEANING OF THE ACCESSORIES


6

Leave to soak for 10 minutes at $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$. This procedure reduces the quantity of blood, protein, and mucous present in the insert;


7

During immersion into the enzymatic solution, delicately brush the surface of the insert with a brush with soft nylon bristles to eliminate all traces of visible dirt. Accurately clean the difficult areas such as sharp edges and especially the interstices between the cutting cusps;

8

Delicately brush the surface of the insert again with the brush with soft nylon bristles;


9

Rinse the internal canal of the insert with distilled water injected under pressure (3,8 bar) for at least 10 seconds to eliminate all residues;

EN


Rinse the internal canal of the insert with distilled water injected under pressure (3,8 bar) for at least 10 seconds to eliminate all residues;

10


Place the insert in the ultrasound tank submerged by the enzymatic detergent solution at $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$, for at least 10 minutes;

11


Remove the insert from the ultrasonic tank and rinse with distilled water;

12


Delicately brush the surface of the insert again with the brush with soft nylon bristles;

13

Rinse the internal canal of the insert with distilled water injected under pressure (3,8 bar) for at least 10 seconds to eliminate all residues;


14

Process validated by an independent institute with enzymatic detergent Enzymec (Mectron).


INSERTS TORQUE WRENCH

EN


15

Put the wrench to soak in the enzymatic detergent solution for 10 minutes at $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$;


16

During immersion in the enzymatic solution, delicately brush the surface of the wrench with a brush with soft nylon bristles to eliminate all traces of visible dirt, both in the internal part and the external one;


17

Delicately brush the surface of the wrench with a brush with soft nylon bristles under tap water;

18

Rinse the wrench with distilled water until it is clearly visible that the contaminants have all been removed;


19

Place the insert in the ultrasound tank submerged by the enzymatic detergent solution at $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$, for at least 10 minutes;


Rinse with tap water to eliminate all residues;

20

Delicately brush the surface of the wrench with a brush with soft nylon bristles using distilled water;

21


Process validated by an independent institute with enzymatic detergent Enzymec (Mectron).

PERISTALTIC PUMP TUBES AND HANDPIECE MOBILE SUPPORT


Immerse the peristaltic pump tube in the enzymatic detergent solution prepared;

22


Rinse with abundant distilled water;

23


TUBE-TO-TUBE CONNECTOR

24

Place the connector in the ultrasound tank submerged by the enzymatic detergent solution at $40^{\circ}\text{C} \pm 2^{\circ}\text{C}$, for at least 10 minutes;


25

Remove the connector from the ultrasonic tank, brush it with the brush with soft nylon bristles, and rinse it with distilled water to eliminate all residues of detergent;


EN

IRRIGATION BAG SUPPORT ROD AND HANDPIECE FIXED SUPPORT

26

Clean the surfaces of the rod and of the handpiece holder with a soft cloth with low fiber-release, dampened with an enzymatic detergent solution;

27

Remove the residues of detergent with a soft cloth with low fiber-release dampened with distilled water;

06.1.2 → THERMODISINFECTION

❗ **CAUTION:** The following automated cleaning cycle and the thermodisinfection must always be preceded by the pre-cleaning steps with enzymatic detergent.

NOTE: Place the accessories in the washing machine so that dead zones do not arise and the water can properly drain. Also, make sure that the reusable inserts are properly held in place in the washing basket and cannot move during the washing process, as shocks could damage them.

⚠ **WARNING:** Avoid overloading. Failure to comply overloading compromises cleaning effectiveness.

❗ **CAUTION:** Do not kink the handpiece cord; this may damage it.

⚠ **WARNING:** Once the cleaning cycle in thermodisinfectors has been completed, the handpiece keeps the washing temperature for a long time. When removing the handpiece from the thermodisinfectors take proper measures to avoid damages for the user.

❗ **CAUTION:** Due to its shape, the handpiece can roll. For this reason, when it is not in use, it shall always be placed on its support.

Lay the accessories in a metallic tray and place it in thermodisinfectors;
Sequence and parameters applicable to the cycle:

CLEANING


- 1 min, Pre-washing with cold water;
- 5 min, Washing with alkaline detergent at $55^{\circ}\text{C} \pm 2^{\circ}\text{C}$;
- 1 min, Neutralization with proper solution (1/3 cold water, 2/3 warm water);
- 1 min, Rinsing with water (1/3 cold water, 2/3 warm water);

DISINFECTION

- 5 min, Thermodisinfection at 93°C with demineralized water.

The automated thermodisinfection was not experimentally tested. According to standard ISO 15883-1, Table B.1 [4] the thermodisinfection at a temperature of 90°C for 5min determines an A0 value of 3000.

1

Procedure validated with Miele PG8535 washer/disinfectors.
Programm Miele DES-VAR-TD.
Alkaline detergent: neodisher® FA (0.2 % v/v).
Neutralizer liquid: neodisher® Z (0.1 % v/v)

06.2 → ONLY AUTOMATED CYCLE



NOTE: Place the accessories in the washing machine so that dead zones do not arise and the water can properly drain. Also, make sure that the reusable inserts are properly held in place in the washing basket and cannot move during the washing process, as shocks could damage them.

⚠ **WARNING:** Avoid overloading. Failure to comply overloading compromises cleaning effectiveness.

⚠ **CAUTION:** Do not kink the handpiece cord; this may damage it.

⚠ **WARNING:** Once the cleaning cycle in thermidisinfector has been completed, the handpiece keeps the washing temperature for a long time. When removing the handpiece from the thermidisinfector take proper measures to avoid damages for the user.

⚠ **CAUTION:** Due to its shape, the handpiece can roll. For this reason, when it is not in use, it shall always be placed on its support.

1

Lay the accessories in a metallic tray. Connect the suitable adapter (supplied as optional) to the handpiece cord connector and then to the cleaning water jet connections of the washer/disinfector. Repeat the same step for the inserts and for the peristaltic pump tubes, connecting them to the suitable adapters. (supplied as optional).

2

Sequence and parameters applicable to the cycle:

CLEANING

- 1 min, Pre-washing with cold water;
- 5 min, Washing with alkaline detergent at $55^{\circ}\text{C} \pm 2^{\circ}\text{C}$;
- 1 min, Neutralization with proper solution (1/3 cold water, 2/3 warm water);
- 1 min, Rinsing with water (1/3 cold water, 2/3 warm water);

DISINFECTION

- 5 min, Thermidisinfection at 93°C with demineralized water.



The automated thermidisinfection was not experimentally tested. According to standard ISO 15883-1, Table B.1 [4] the thermidisinfection at a temperature of 90°C for 5min determines an A0 value of 3000.

Procedure validated with Miele PG8535 washer/disinfector.
 Programm Miele DES-VAR-TD.
 Alkaline detergent: neodisher® FA (0.2 % v/v).
 Neutralizer liquid: neodisher® Z (0.1 % v/v)

07 → CLEANING CHECK

Once you've finished the cleaning operations, check all the objects under an adequate source of light, if need be also with a magnifying glass, paying attention to the details that could conceal dirt residues (threading, cavities, grooves) and, if need be, repeat the cleaning cycle by following the manual/automated cycle as described in Paragraph 06.1. Finally, check the integrity of those parts and elements that may have worn out during their use;



1


08 → DRYING AND LUBRICATION

Thoroughly dry all parts of the handpiece with cord, especially the electrical contacts, blowing on them with compressed air;

1


⚠ CAUTION: The electrical contacts of the handpiece cord connector must be dry, before and after the sterilization cycle, before connecting the cord to the device. Always make sure that the connector electrical contacts are perfectly dry; if need be, dry them by blowing compressed air.


2

⚠ **CAUTION:** Before starting the sterilization cycle, make sure that the insert is thoroughly dry, both externally and internally. For this purpose, blow compressed air both externally and through the internal thru hole; this will prevent the onset of stains and haloes on the surface or rust inside the insert.


3

Thoroughly dry the accessories on all their surfaces, internal and external, using a clean and soft cloth with low fiber-release or by blowing compressed air, depending on the case.

4

Dry the torque wrench using a clean and soft cloth with low fiber-release;

5


Lubricate with medical degree lubricant in the indicated point;

⚠ **CAUTION:** do not use lubricants based on oil or silicon.

09 → STERILIZATION


Individually seal the handpiece (without inserts) and the front terminal in single-use packets for separate sterilization.

1


Individually seal the inserts in single-use sterilization packets.

2

Individually seal all the accessories, separately, in single-use sterilization packets.

3

STERILIZATION METHOD

The handpiece with cord is built with materials that resist to a maximum temperature of 135°C for a maximum time of 20 minutes.

Once the handpiece and all accessories have been individually packed, proceed with the sterilization cycle in steam autoclave.

The sterilization process in steam autoclave must be performed setting the parameters indicated below, which are validated by Mectron and reach the sterilization safety levels (SAL 10⁻⁶):

Cycle type: 3 times Pre-vacuum (minimum pressure 60 mBar).

Minimum sterilization temperature: 132°C (interval 0°C ÷ +3°C).

Minimum sterilization time: 4 minutes.

Minimum drying time: 20 minutes.

All the sterilization phases must be performed by the operator in compliance to the norms UNI EN ISO 17665-1:2007, UNI EN ISO 556-1:2002 and ANSI/AAMI ST 46:2002.

⚠ **CAUTION:** Do not sterilize the handpiece with the insert screwed onto it.

⚠ **WARNING: Checking infections - Parts that can be sterilized** - Diligently remove all residues of organic dirt before the sterilization.

⚠ **CAUTION:** Perform the sterilization using a water steam autoclave only. Do not use any other sterilization procedure (dry heat, irradiation, ethylene oxide, gas, low temperature plasma, etc.).

⚠ **CAUTION:** Do not exceed the allowed load of the steam sterilizer.

⚠ **WARNING:** Once the cleaning cycle in thermodisinfectant has been completed, the handpiece keeps the washing temperature for a long time. When removing the handpiece from the thermodisinfectant take proper measures to avoid damages for the user.

⚠ **CAUTION:** At the end of the sterilization cycle, let the handpiece cool down completely prior to using it.

The Italian version of this manual is the original document from which its translations have been obtained. In case of any discrepancy, the Italian version will have pertinence.

Texts, images, and graphics of this manual are property of Mectron S.p.A., Carasco, Italy. All rights reserved. The contents cannot be copied, distributed, changed, or made available to third parties without the written approval of Mectron S.p.A.

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