



Chris O'Brien
Lifehouse

RestorabiteTM

ARTG ID: 407358



Integrated
Prosthetics &
Reconstruction

the 1990s, the number of people in the world who are undernourished has increased from 600 million to 800 million (FAO 2001).

There are a number of reasons for this increase. First, the world population has increased from 5 billion in 1987 to 6 billion in 2000. Second, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000. Third, the number of people who are undernourished has increased from 600 million in 1987 to 800 million in 2000. Fourth, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000. Fifth, the number of people who are undernourished has increased from 600 million in 1987 to 800 million in 2000.

There are a number of reasons for this increase. First, the world population has increased from 5 billion in 1987 to 6 billion in 2000.

Second, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000.

Third, the number of people who are undernourished has increased from 600 million in 1987 to 800 million in 2000.

Fourth, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000.

Fifth, the number of people who are undernourished has increased from 600 million in 1987 to 800 million in 2000.

Sixth, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000.

Seventh, the number of people who are undernourished has increased from 600 million in 1987 to 800 million in 2000.

Eighth, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000.

Ninth, the number of people who are undernourished has increased from 600 million in 1987 to 800 million in 2000.

Tenth, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000.

Eleventh, the number of people who are undernourished has increased from 600 million in 1987 to 800 million in 2000.

Twelfth, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000.

Thirteenth, the number of people who are undernourished has increased from 600 million in 1987 to 800 million in 2000.

Fourteenth, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000.

Fifteenth, the number of people who are undernourished has increased from 600 million in 1987 to 800 million in 2000.

Sixteenth, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000.

Seventeenth, the number of people who are undernourished has increased from 600 million in 1987 to 800 million in 2000.

Eighteenth, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000.

Nineteenth, the number of people who are undernourished has increased from 600 million in 1987 to 800 million in 2000.

Twentieth, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000.

Twenty-first, the number of people who are undernourished has increased from 600 million in 1987 to 800 million in 2000.

Twenty-second, the number of people who are undernourished has increased from 15% of the world population in 1987 to 13% in 2000.

Table of contents

- 01 **Caution**
- 01 **Information**
- 02 **Indications**
- 02 **Contraindications**
- 02 **Risk/Potential benefit**
- 03 **Warnings and Precautions**
- 04 **Operating information**
- 05 **When to contact your doctor**
- 05 **Cleaning**
- 05 **Hygiene**
- 06 **Expected device life span**
- 07 **References**

- Contact information**

CAUTION

The Restorabite™ is to be used according to the recommendations and instructions for use. Please read all information in this Product Brochure carefully and use the contact details below if you require further information.

This Product Brochure contains general safety, operating, maintenance and care instructions for Restorabite™. This information should be reviewed by both prescribers and patients using this device. This Product Brochure does not replace the personalised guidance from a treating clinician.

INFORMATION

Purpose of the device

Restorabite™ is a device intended to improve mouth opening in patients with trismus (restricted mouth opening defined by an interincisal distance less than 35mm) (Dijkstra et al 2006), or to prevent trismus in patients who are at risk of developing trismus.

Intended Use

Restorabite™ is intended to assist with jaw stretching and strengthening for treatment of trismus (restriction of jaw opening)

Device Features

Its' unique multiple exercise level design allows for a varied force to be selected based on a patient's specific needs, to ensure that:

- Adequate force and effort are being used for maximal therapeutic gains.
- Excessive force is not applied to the teeth and jaw, to reduce or prevent the risk of complications and injuries relating to jaw stretching and strengthening exercises.
- The force applied to the jaw can be modified by user needs.

Device Description

The Restorabite™ consists of:

- Frame (Figures 1)
- Insert (Figure 2)
- Instruction For Use

The Restorabite™ outer frame can house an insert at different locations providing varying resistance.

The Restorabite™ is portable and hand-held device. The outer wings of the frame are squeezed together by the user and placed between the upper and lower teeth or gums. The user then gradually releases the pressure applied by their fingers to allow the wings of the frame to deploy force to the jaw, stretching it open (passive range of motion exercise). The user can then bite down to activate the jaw muscles for a strengthening exercise (active range of motion exercise).




Figure 1



Figure 2

INDICATIONS

The Restorabite™ is indicated for use in individuals whose mouth opening is limited. Restricted mouth opening is known as ‘trismus’. Trismus can arise because of many medical conditions. Commencing rehabilitation for trismus should always be undertaken under the guidance of the user’s physician and experienced clinician.

The Restorabite™ can be used for both rehabilitation and prevention of trismus. Rehabilitation exercises are undertaken when a user has trismus to improve their condition. Prevention or prophylactic use is for those who are at risk of developing trismus.

The minimum mouth opening distance Restorabite™ is able to accommodate for is 8.0 mm.

POTENTIAL BENEFIT and RISK

Potential Benefits

The purpose of the Restorabite™ is to increase mouth opening. Reduced ability to open the mouth can lead to several health-related problems such as pain, difficulty swallowing, speaking, conducting oral/dental care and accessing the airway.

Addressing trismus using Restorabite™ has led to improvements in mouth opening, speech and swallowing related quality of life and improved pain scores. The amount of improvement depends on the individual and their personal circumstances.

The use of Restorabite™ for trismus therapy provides a personalised exercise program for patients. Patients who used Restorabite™ report increased adherence to exercise program and motivation for exercise program, improved quality of life, and a decreased pain (Charters et al 2022)

The report has shown the following benefits:

- Improved mouth opening
- Increased adherence to exercise program and motivation for exercise program
- Improved quality of life

Risks

Failure to follow the instructions for use may result in complications. Incorrect device handling (such as overstretching of the device), and device cleaning (such as using the wrong cleaning agent) will increase the risks.

The Restorabite™ should not be used without an individualised exercise program provided by your physician and/or trained clinician. Application of excessive force to the teeth or jaw can cause complications such as fractures. A patient’s physician can select the maximum force suitable for a patient based on their individual circumstances to mitigate this risk.

The following risks were associated with the Restorabite™:

- Tension headache
- Neck pain
- Temporomandibular joint (TMJ) pain
- Tooth loss

If any of these complications occur, the user should contact their prescribing clinician.

CONTRAINDICATIONS

Restorabite™ is not recommended for patients with any of the following conditions:

- Osteomyelitis or osteoradionecrosis of their jaw (maxilla or mandible)
- Fracture of their jaw (maxilla or mandible)
- Unstable dentition
- Have active disease (e.g. untreated cancer of the jaw)
- The Restorabite has not been trialed in patients under the age of 18 years.



WARNINGS AND PRECAUTIONS

Warnings

- If you feel unexpected pain at any time while exercising, **DO NOT** continue to exercise. You should immediately stop using the device and contact your clinician.
- If you have weakened teeth, gum disease, weakened bones or joints in jaw, dental implants or dentures, **DO** exercise with great caution. Injury may occur if excessive force/inappropriate force selection is used.
- **DO** use appropriate force level to avoid injury.
- **DO NOT** share the Restorabite™. It is intended for single patient use only.
- **DO NOT** use the Restorabite™ if the package seal is broken before use. Contact your supplier.

Precautions

- Use of Restorabite™ must be done according to clinician led individualized exercise program to avoid injury.
- Do not sterilise the Restorabite™ or put the Restorabite™ in dishwasher. The sterilisation and dishwashing process may weaken the Restorabite™ and it may break during use.
- Store the Restorabite™ according to the Instruction for Use.

Care regimen

For the Restorabite™ to perform as intended, it is critical that you understand the warnings and precautions (as outlined in the section above) and follow the directions given by your prescribing clinician. It is important that you have regular follow-up appointments with your treating clinician. When using the Restorabite™ it is important to observe good hygiene (cleanliness) to minimise the risk of infection. For additional information, please see Hygiene Recommendations on page 9 of this manual.

OPERATING INFORMATION

Pre-exercise

The Restorabite™ allows a specified force to be exerted on the jaw to facilitate stretching (passive) and strengthening (active) exercises.

Your clinical team will conduct a thorough assessment to determine what level of force is best suited to your needs. This will involve selection of a frame location for the insert to commence your program and ensuring the correct position for the Restorabite™ within your mouth.

Exercise Program

Prepare

- Measure the distance between your top and bottom teeth or gums before and after an exercise set.
- Warm-up by placing heat packs on the jaws and/or manual jaw stretches.
- Set up the Restorabite™ with the insert in the position as per Figure 3.




Figure 3

Exercise

- Exercises should be incorporated into a daily routine.
- Stretch duration, number of repetitions and frequency should be tailored to the individual user.
- General exercise programs are outlined below:
 - 5-10 minutes of warm up, passive and active stretching, 5 times a day.
 - 30 minutes of stretching at a lower intensity comprised of warm up, passive and active stretching.
 - A combination of the two
- **Passive** range of motion exercises
 - Squeeze the wings of the Restorabite™ together.
 - Insert the Restorabite™ between upper and lower teeth or gums, keep your fingers in position as per Figure 4 .




Figure 4

- Gradually release the pressure from your fingers, allowing the Restorabite™ to open and press against your upper and lower teeth or gums.
- Relax and let the Restorabite™ stretch open your mouth.
- Hold the stretch for up to 30 seconds. Do not hold the stretch over 30 seconds.
- Release the Restorabite™.
- Rest briefly before the next repetition.
- It is normal for the stretch to be somewhat uncomfortable, however it should not be painful. Cease the stretch if you experience pain.

- **Active** range of motion exercises
 - Squeeze the wings of the Restorabite™ together.
 - Insert the Restorabite™ between upper and lower teeth or gum, keeping your fingers in position.
 - Gradually release the pressure from your fingers, allowing the Restorabite™ to open and press against your upper and lower teeth or gums.
 - Bite down on the Restorabite™ for up to 30 seconds. Do not hold the exercise over 30 seconds.
 - Release the Restorabite™.
 - Rest briefly before the next repetition.
 - It is normal for the stretch to be somewhat uncomfortable, however it should not be painful. Cease the stretch if you experience pain.

CLEANING

It is important to inspect and clean the Restorabite™ after each use.

- The Restorabite™ should be cleaned thoroughly following each use by rinsing in warm water, wash with a PH neutral detergent removing any soiling or debris with a soft toothbrush.
- Alternatively, the Restorabite™ should be wiped with a clean wet cloth and PH neutral detergent.
- Let all the components of Restorabite™ air-dry completely after each clean.

DO NOT sterilise or use a dishwasher to clean the Restorabite™, this can damage the device.

WHEN TO CONTACT YOUR DOCTOR

For your safety, comfort and best therapeutic gain, it is important to follow these recommendations:

YOU MUST IMMEDIATELY CONTACT YOUR TREATING DOCTOR IF:

- You experience sudden pain.
- There is a swelling in your jaw.
- The location where your Restorabite™ contacts your teeth or gum becomes inflamed, swollen or painful.
- There is a break in the tooth, or loss of tooth.

Contact your physician or clinician if you have any questions about your exercise program.

HYGIENE

It is important to inspect the area where the Restorabite™ contacts your dentition or gums prior to commencing your exercises. Good oral care involving twice daily cleaning (or as guided by your dentist) is recommended.

You should not use the Restorabite™ on any areas of broken gum or exposed bone.

IN CASE OF IRRITATION OR INFECTION

- If you have signs of infection stop using the device and notify your clinician.
- If you have any signs of irritation or infection such as redness, or mild pain stop using the device and notify your clinician
- **IN CASE OF ANAPHYLAXIS, HIGH FEVER AND/OR SEVERE PAIN, YOU SHOULD IMMEDIATELY GO TO THE HOSPITAL EMERGENCY ROOM.**

EXPECTED DEVICE LIFE SPAN

When used every day as part of a patient specific program, the Restorabite™ is tested for life-expectancy of 10 weeks. This duration may be extended with less frequent or less intensive use. Should the user fail to comply with the recommended use, care or cleaning guidelines, or if the device withstands trauma such as a fall or compression, the Restorabite™ may require replacement at an earlier date.

It is the responsibility of the user to ensure familiarity and compliance with the recommendations outlined in this document. It is the responsibility of the prescribing health care professional to select the appropriate insert and treatment regimen.

The Integrated Prosthetics & Reconstruction at Chris O'Brien Lifehouse is only responsible for the Restorabite™ performance when it is used in accordance with the Product Brochure and Instructions for Use.

It is important that you understand the recommendations, warnings and precautions outlined in this document and follow the directions of your doctor and/or health care professional.


REFERENCES

Dijkstra, P. U., Huisman, P. M., & Roodenburg, J. L. (2006). Criteria for trismus in head and neck oncology. *International journal of oral and maxillofacial surgery*, 35(4), 337–342.
<https://doi.org/10.1016/j.ijom.2005.08.001>

Charters, E., Cheng, K., Dunn, M., Wu, R., Palme, C., Howes, D., Low, T. H., Heng, C., Ricketts, V., Kneebone, K., Loy, J., & Clark, J. R. (2022). A pilot study of intensive intervention using a novel trismus device. *International journal of speech-language pathology*, 1–8. Advance online publication.
<https://doi.org/10.1080/17549507.2022.2130429>



Integrated Prosthetics & Reconstruction



Integrated Prosthetics & Reconstruction

Chris O'Brien Lifehouse

119-143 Missenden Rd

Camperdown NSW 2050

Australia

Tel: +61 2 8514 0268

Email: IntegratedPR@lh.org.au

Website: <http://mylifehouse.org.au/IntegratedPR>

