

UBrush!

Smile, because UBrush!






*How do you clean
between the teeth?*



UBrush!

The electric interdental brush
with LED Light

-  makes it so easy to clean
between the teeth
-  leads to better oral hygiene
-  produces an even more
attractive smile



Optimal oral hygiene – Now it's in your hands!

Even the best toothbrushes can't reach the interdental areas. That's why there are alternatives: piks, floss or interdental brushes. These usually end up being difficult to use, and patients give up.

UBrush! is the answer. This user-friendly electric interdental toothbrush helps you clean between the teeth, both easily and efficiently.

Ask your hygienist or dental therapist to show you how to use the UBrush! It really is easy. It will convince you how simple it is to clean between the teeth.

Our Tips:

- Teeth should be cleaned once or twice a day, in the morning and at night. After brushing with your regular toothbrush (manual or electric), continue with the UBrush! to clean interdentally.
- Choose the correct brush size. (You may need 2 or even 3 sizes). Feed carefully into the interdental area. It becomes easier when you get used to it.
- Increase the effect by inserting the brush 3-4 times in each interdental space.
- Rinse your mouth afterwards.

UBrush! available directly through your dentist or in the local pharmacy.

UBrush! Set

Contents: UBrush! handpiece, stand, battery, mirror, 11 brushes (assortment), instruction for use)

Ref.: 35000

UBrush! Refills

Content: 10 brushes / package

Diameter	Colour	Form	Ref.
0,5 mm	●	cylindric	35105
0,6 mm	●	cylindric	35106
0,7 mm	●	tapered	35207
0,8 mm	●	tapered	35208
0,9 mm	●	tapered	35209
1,0 mm	●	tapered	35210
1,2 mm	●	tapered	35212



UBrush! Enterprises B.V.

LOSER & CO
öfter mal was Gutes...

