



INNOVATIVE SOLUTIONS
FOR DENTAL AND MEDICAL
PROFESSIONALS

ORAL HEALTH CARE
P R E V E N T I V E

SUMMARY

PREVENTIVE DENTISTRY	AGENCIES	PAGE
About us	Biomedent	2
Cavity Prevention	CAO	3
Dry mouth solutions	Denta-Med	4
Dental Cleaning Device & Sports Mouth Guard	Dr Marks Hygenie	5-6
Tongue Cleaner	Embi	7
Tongue Retractors	Armor Dental	8-11
Floss	Manfloss & Queen of Clean	12
Tongue Posture	MyoSpots	13
Mouthwash & Gum Serum	OhCare!	14
Toothbrushes & Interdental Brushes	TePe	15
Post-Surgery Care	TePe / Vitis	16
Rehabilitation	Restorabite	17-18

ABOUT US



Biomedent specialises in supplying high-quality medical and dental products. To ensure our products remain accessible only to genuine healthcare professionals, purchases through our website are available exclusively to account holders.

A banner for the 'Online Account' featuring a background image of dental instruments. On the left, a teal button says 'Available Now'. To its right, the text 'Online Account' is displayed. Further right, six benefits are listed in two columns, each preceded by a teal checkmark icon. At the bottom, a teal bar contains a registration process flow: a person icon followed by 'REGISTER', a document icon followed by 'FILL FORM', a checkmark icon followed by 'APPROVED', and an arrow pointing to 'STAR ORDERING'.

Available Now

Online Account

- ✓ Order now, pay later
- ✓ Member pricing
- ✓ Pre-launch access
- ✓ View and track recent orders
- ✓ Wish list for fast reordering
- ✓ Manage billing and shipping addresses

REGISTER FILL FORM APPROVED STAR ORDERING

Creating an account is quick, easy, and comes with exclusive benefits:

- ✓ Order now, pay later – enjoy 7-day payment terms
- ✓ Member-only pricing – access special rates on selected educational courses
- ✓ Pre-launch access – be among the first to explore new technology
- ✓ Simplified ordering – view and track orders with ease
- ✓ Flexible account management – update billing and shipping details anytime
- ✓ Quick reordering – save favourites to your wish list for fast, hassle-free orders

Setting up your account takes only a few seconds. Once registered, you'll have everything you need at your fingertips.



CAVITY PHIGHTER

Test Strips:

Created with dental professionals in mind, Cavity pHighter test strips measure the oral pH balance to determine cavity risk. Mouth bacteria can lead to increased acidity, resulting in dental erosion and painful cavities.

Cavity pHighter strips measure the pH and buffer capacity in your saliva. Simply wet the pH strip with the patient's saliva, wait 10 seconds, and measure their color reading against the key on the bottle.

The easy-to-use color key will indicate their mouth pH level and associated risk.

Mints:

Created for at-home and professional dental use, Cavity pHighter Mints protect your patients' teeth after they've left your office. CAO's cavity-fighting mints protect your mouth against acid-producing bacteria by reducing acidity and balancing oral pH levels. The result? Decreased risk of cavities and minty-fresh-smelling breath!





DRY MOUTH SOLUTIONS

Dry Mouth Oral Hygiene Gel 65mL

- Provides up to 8 Hours relief from Dry Mouth
- Moisturises, hydrates, soothes & comforts
- Fluoridated formulation, replacing toothpaste

Dry Mouth Oral Hygiene Spray 50mL

- Moisturising, fluoridated oral hygiene spray
- Convenient oral hygiene spray with 360 degree nozzle
- Freshens breath & reduces halitosis (mouth odour)

Denture Adhesive Gel 50mL

- Secures ALL day (naturally) Seals, Cushions & Comforts
- PVM/MA Copolymer "glue" (free) formulation
- Neutral taste, Zinc free, Easy to clean

Moisturising Mouth Spray 60mL

- Convenient and immediate relief from Dry Mouth
- Enhanced moisture delivery with (prop.) chitin
- A 'nod' to the Australian Bush with a hint of lemon myrtle

Moisturising Mouthwash 500mL

- Enriched with all natural moisturising Malic Acid
- Leaves the mouth with a soft, velvety mouth-feel
- pH balanced, moisturising, calming mouth rinse

pH Neutralising Mouthwash 250mL

- Helps prevent enamel erosion
- Naturally assists to neutralise pH levels
- Refreshing flavour, safe to use ongoing





DR MARK'S HYGENIE



Dr Mark's HyGenie® is the simple solution for cleaning your dentures, night guard or Invisalign®. No more soaking in Polident, Steradent, lemon or vinegar solutions, and no more risky hand-held brushing. Simply lock your dental appliance safely inside HyGenie, add DentalFresh, twist and rinse. The high-quality rotating brushes inside do the work for you.

DR MARK'S DENTALFRESH

Mark's DentalFresh® is an easy-to-use, fast acting and fresh scented cleansing gel. Think of it as SOAP for washing your Sports Mouthguard (or dentures, night guard and clear aligners like Invisalign®). DentalFresh cleans dental appliances by washing away biofilms and attacking all the other nasties that collect on your mouthguard in practice or competition, giving you protection against harmful germs and bacteria. And remember, you just can't be 100% with a dirty mouthguard.



DR MARK'S HYGENIE REPLACEMENT BRUSH SET



With daily use, you should replace the brush set in your Dr Mark's HyGenie every 6 months to ensure optimal performance and good oral hygiene. There's no need for an entirely new device, simply replace your old brushes with the new brushes in just 90 seconds.



HYGENIE SPORT



Gone are the days of leaving your mouthguard in the bottom of a sweaty sports bag where it loves to collect dirt, bacteria, germs and who knows what else! HyGenie Sport is the simple solution for cleaning, storing and protecting your mouthguard. After every practice or competition, simply place your mouthguard inside HyGenie Sport, twist, and let the high quality rotating brushes do the work for you. You couldn't ask for an easier way to keep your mouthguard clean.

HYGENIE SPORT CAGE MOUTHGUARD CARRIER

Stop rummaging around the bottom of your bag for your mouthguard. With HyGenie Sport 'Cage' Mouthguard Carrier, you easily protect, store and find your mouthguard. The CAGE is bomb proof mouthguard protection, helping to keep it clean and safe. Clip The CAGE to the outside of your sports bag with the Carabiner Clips included, easy to see, easy to find and combat ready.



HYGENIE SPORT QUICKCLEAN SPRAY 100M



HyGenie Sport QuickClean Spray is a mouthguard sanitiser that quickly and effectively protects you and your mouthguard from the nasties that collect on them during practice, comps or floating around the bottom of your sports bag. Spray QuickClean onto your mouthguard after every use, then store in the CAGE mouthguard carrier and be confident that it's hygienic, protected and combat ready (even if it has been floating around the bottom of your sports bag). It's as easy as that.



EMBI FIRST STEP TONGUE SCRAPER

Tongue cleaning with a quality tongue scraper is the First Step towards Optimal Oral Care.

It's clinically proven to be the best method of removing the harmful bacteria that causes bad breath, tooth decay, gum disease and other health issues including poor gut health.

EMBI EDGE TONGUE CLEANER

Designed with dentists it's made from dental grade stainless steel and has a unique scoop that gives a deeper clean for your tastebuds. It has uniquely weighted forged stainless steel handles for ergonomic agility and is made by embi – known for their fine craftsman ship and quality.



EMBI EDGE CU29 COPPER TONGUE CLEANER

It's made from the purest quality copper and has a unique scoop that gives a deeper clean for your tastebuds. Copper is known for it's anti-microbial and anti-bacterial properties, which makes it a perfect and effective material for tongue cleaning as it promotes good bacteria and better oral hygiene.






Why a soft retractor?

Harsh Metal Is Harsh to Tissue

We understand there are common metal instruments used every day in your practice, but what if you had a **gentler choice**?

You do now! We want to give you an additional product that is more comfortable and gentler on tissue and better for your patient when they are healing.



1

Our soft retractors can be used chairside and can then immediately be given to the patient or caregiver to use and improve hygiene at home.

2

Developed to give clinicians a softer more comfortable option when retracting tissue compared to hard metal or plastic deal for any procedures.

3

You now have a comfortable choice to add to your dental tray to improve tissue management and home care.



ARMORDENTAL

Soft Retractors



This product can be autoclaved up to 30 cycles. If not designated as a single-patient-use (SPU) item, it should be replaced after repeated use to ensure optimal performance and hygiene.



HOME
Tongue Guide

Assessment, surgical and after care retractor for tongue evaluation and surgeries



SPU

SCAN·MATE


Your intraoral scanning partner. Faster and easier scans.



SPU

Mouth·Mate

Post procedure retractor to improve home care



SPU

ComfortSoft

Dual purpose retractor used chairside and at home.

ARMORDENTAL

MiROSoft™

**Soft.
Bendable.
Comfortable.**

HD (HIGH DEFINITION)
WATER REPELLENT
ANTI-FOG
SCRATCH-RESISTANT

SOFT BUMPER

NON-REFLECTIVE


SOFT
MATERIAL

BENDABLE &
CUSTOMIZABLE
STAINLESS STEEL
NECK TO ENSURE
THE PERFECT
FIT IN ACCESS

SLIP RESISTANT
GRIP

ERGONOMIC
HANDLE

Available Sizes



Bendable Mirror

- 4 soft mirror sizes available: #4, #5, Teardrop and #5 double-sided mirror. The unique teardrop mirror is a combination of the #4 and #5 mirrors, offering an elongated and improved view.
- Soft mirror lens is high definition, water repellent, anti-fog, and scratch-resistant.
- All soft mirrors have a soft bumper around the perimeter of the lens, which provides comfort to patients. Designed for silent retraction-no clanking against teeth!
- The first dental mirror featuring a soft, bendable, non-reflective back side and neck, allowing safe retraction for any procedure



ARMORDENTAL



TissuePro™ Minnesota Soft Sleeves

This product can be autoclaved up to 30 cycles. If not designated as a single-patient-use (SPU) item, it should be replaced after repeated use to ensure optimal performance and hygiene.

Soft Intraoral



A



B

Soft External



C

Soft All-In-One



D

TissuePro - to be placed over the Minnesota retractor to enhance both patient comfort and clinical control during retraction.

The first soft, non-reflective sleeve portfolio for the Minnesota retractor.

- A- Intraoral Soft Sleeve - Autoclavable
- B- Intraoral Soft Sleeve - SPU
- C- External Soft Sleeve - Autoclavable
- D- All-in-One Soft Sleeve - Autoclavable

Improve control and comfort, and reduce hand fatigue

Minimise trauma when retracting tissue, enhance blood supply for faster healing, and reduce slippage on bone.





MANFLOSS & QUEEN OF CLEAN has been designed to improve your oral hygiene by explaining it simply, in a fun way. Developed by a dental health professional who listened (for 20 years) to the everyday persons' frustrations with flossing, and decided to do something about it.

Finally, there is a luxurious wafer thin but also soft and wide tapes designed to stop you cursing and get you smiling! For the first time ever, precise instructions are enclosed to dispense the frustrations about flossing...learn how length is everything....





A new patented technology for tongue exercise.

Myospots are small biodegradable adhesive pads that are used to exercise tongue muscles and induce nasal breathing.

Myospots are small, circular pads, made of biodegradable natural polymers that are safe for internal use.

They have nice colours and flavours that give them a pleasant look and taste and are Vegan friendly. Their strong adhesive power enables them to adhere easily to the surface of the palate.

Upon adhering to palate, they involuntarily induce the tongue to elevate and touch them. The repetitive induced movement of the tongue to touch and dissolve the adhered pad provides good exercise for tongue muscles to increase tongue muscles tone.

Using MYOSPOTS daily for several months can help change the tongue posture to rest against the palate.



How to use MYOSPOTS?

Place a MYOSPOT with your thumb against the anterior part of the palate with gentle force for 10 seconds.

Remove your thumb and keep pressing the SPOT using your tongue until it dissolves.

oh!CARE



OH! Care Mouthwash, naturally improves your oral health whilst freshening your breath with the power of hyaluronic acid and xylitol as well as cavity preventing fluoride

Vegan friendly
100% Australian made
Recyclable packaging



The ultimate package to care for teeth and gums. Add HA-X Gel and Mouthwash to your daily routine to improve your oral Health.

Used in conjunction with twice daily toothbrushing and flossing, our mouthwash & gel may assist in:

- Reducing harmful bacteria
- Hydrating soft tissue
- Promoting wound healing
- Reducing inflammation
- Increasing saliva flow
- Reducing risk of tooth decay
- Reducing risk of gum disease
- Improving oral health
- Freshening breath

****If symptoms persist, it is important to seek advice from your dental health professional**



TEPE COMPACT TUFT TOOTHBRUSH

TePe Compact Tuft™ is a single tuft brush with a dome-shaped brush head for precision cleaning of hard-to-reach areas, such as implants, attachments for overdentures, fixed braces or along the gum line.

TEPE UNIVERSAL CARE TOOTHBRUSH

TePe Universal Care™ has a unique angle to facilitate access in hard-to-reach areas. Optimal for cleaning on inside surfaces and along the gum line of natural teeth, implants, lingual braces, or retainers.



TEPE SPECIAL CARE

TePe Special Care is an ultra-soft toothbrush recommended after oral surgery and for sore or delicate tissues.

TEPE INTERDENTAL BRUSHES – EXTRA SOFT

An interdental brush is the most efficient, user-friendly device for plaque removal between teeth. TePe Interdental brushes are designed to clean the gaps between the teeth, where a regular toothbrush is unable to reach. This Interdental brush has extra soft filaments for sensitive gums and teeth.





TEPE IMPLANT CARE KIT

The essential kit for effective plaque removal in and around implants for healthy smile. TePe Implant kit is a collection of carefully selected specialist products, developed together with dental experts.

THE ULTIMATE POST-SURGERY CARE BUNDLE!

This thoughtfully curated kit includes everything your patients need for a comfortable recovery.

Featuring ultra-soft TePe brushes kit for gentle cleaning and, soothing Oh Care products.



VITIS SURGICAL TOOTHBRUSH

The VITIS Implant Toothbrush is a daily brush specifically designed for patients with dental implants, to provide a delicate and deep care within the whole mouth.

TEPE ORTHODONTIC KIT

The essential brace kit for a healthy smile. TePe Orthodontic Kit features carefully selected products to help you keep your teeth and gums clean and healthy while wearing braces.



RestorabiteTM

ARTG ID: 407358



Purpose of the device

RestorabiteTM is a device intended to improve mouth opening in patients with trismus (restricted mouth opening defined by an interincisal distance less than 35mm) (Dijkstra et al 2006), or to prevent trismus in patients who are at risk of developing trismus




RestorabiteTM

ARTG ID: 407358



Restorabite is compact, easy to use and backed by clinical research. In medical studies, users' mouth opening improved on average 1 mm per week. However, results can vary from person to person.



Please use the RestorabiteTM according to the provided recommendations and instructions. For full details, refer to the Product Brochure available on the product page of our website.

Device Features

Its' unique multiple exercise level design allows for a varied force to be selected based on a patient's specific needs, to ensure that:

Adequate force and effort are being used for maximal therapeutic gains. Excessive force is not applied to the teeth and jaw, to reduce or prevent the risk of complications and injuries relating to jaw stretching and strengthening exercises. The force applied to the jaw can be modified by user needs.

This Product Brochure contains general safety, operating, maintenance and care instructions for RestorabiteTM. This information should be reviewed by both prescribers and patients using this device. This Product Brochure does not replace the personalised guidance from a treating clinician.

