

→ MECTRON PROPHYLAXIS POWDER AND MECTRON GLYCINE POWDER FOR BRIGHT SHINY TEETH – SUPRA- AND SUBGINGIVAL

→ reorder: find your local mectron dealer under www.mectron.com
mectron s.p.a., via Loreto 15/A, 16042 Carasco (Ge), Italia, tel +39 0185 35361, fax +39 0185 351374

→ A QUESTION OF CHOICE the right powder for the intended use

→ INDICATIONS	→ MECTRON PROPHYLAXIS POWDER	→ MECTRON GLYCINE POWDER
removal of plaque	😊	😊
removal of stain	😊	😐
use on enamel	😊	😊
use on restoration materials	😞	😊
use on dentine	😞	😊
use on root surfaces	😞	😊
subgingival use	😞	😊
cleaning of fissures	😊	😐
cleaning of cavity margins before etching	😊	😊
cleaning of implant surfaces	😞	😊