

BioMeDent Promotions

May - June 2022

The ultimate weapon against harmful
airborne viruses, bacteria and mould.



Yours Free With Purchases Over \$6,000



Offers ends 30th June 2022

Our Agencies

Bionnovation
biomedical

mectron
medical technology

deve
med. PARTNER FOR DENTAL INSTRUMENTS

PRF PROCESS
CHOUKROUN

TOOTH TRANSFORMER

GINGISAGE

oh!CARE

hismile

BREATHE

NOUVAG

VITIS

RESORBA
REPAIR AND REGENERATE

CARE DENT

DentaLife
AUSTRALIA

MYQSPOTS
AUS

MANFLOSS
SERIOUSLY GOOD FLOSS

TePe

MJK
Instruments

kikkoman
Kikkoman Biochemifa Company

Dr Mark's
HyGenie

1300 792 624

info@biomedent.com.au

www.biomedent.com.au


PO BOX 25, NORTH RYDE BC NSW 1670



Intellegent and Portable

The ultimate weapon against harmful airborne viruses, bacteria and mould.

Designed, developed and assembled in Australia.



Viruses Pollen Bacteria


Mould Dust Mites

THE ONLY AIR PURIFIER THAT **DETOXIFIES** THE AIR AROUND YOU.



VBreathe combines HEPA Filtration and proprietary all-natural VActive Gel technology to deliver the world's most intellegent, portable and effective indoor air purifier and detoxifier in its class.

Fight airborne mould, viruses and bacteria
with our all-natural indoor air purifier and detoxifier



\$690
Free With Purchases
over \$6,000

NEW

Independently tested - 99.999% efficiency
By leading 3rd-party Australian laboratories and universities



NEW

MicroPiezo

Related Treatments

prophylaxis - supragingival scaling
periodontics - subgingival scaling
perimplantitis - implant cleaning
endodontics - cleaning of root canals
endodontics - revision of root canals
endodontics - retrograde preparation

\$3,700



Set Scaling

Bonus Product



New compact and essential design, maximum ergonomics: reduced fatigue for the operator

Mectron Promotions

MICROPIEZO touch - \$3,700

Purchase: One (1) Micropiezo touch unit

Get: One (1) Set Scaling (ME-02160001) FREE!

(equipped with: S1, S2, S3, S4, S5,5 x K10 wrench)

You save \$502

Includes:

- Slim ultrasonic LED hand-piece with circular light
- 1 scaling insert S1-S
- 1 perio insert P3
- 1 liquid container, 500 ml
- 1 footpedal

MicroPiezo The most convenient unit in the market with optimal patient pain management.

Exclusive multifunction touch keyboard with colour code according to the clinical indication

3 functions with 12 different power settings

Scaling and perio

Endodontic

Prosthetics



Thanks to its innovative soft mode the MicroPiezo also sets new standards in insert control and patient pain management.

More than just a scaler!

- calculus to caries,
- endodontics to periodontics,
- extractions to restorations.

mectron

medical technology

- 1 combi touch
- 1 LED ultrasonic handpiece slim
- 1 scaling insert S1-S
- 1 perio insert P3
- 1 bottle 500 ml
- 1 bottle safety cap
- 1 air-polishing handpiece 120°
- 1 air-polishing handpiece PERIO
- 1 box subgingival perio tips (40 pcs)
- 1 wrench K9
- 1 bottle of mectron prophylaxis powder soft+, contents 250 g


start all in one

The combi touch combines ultrasound and air polishing in one unit, it allows complete prophylaxis treatment from removal of supragingival and subgingival calculus to gentle removal of stain and biofilm

COMBI TOUCH



Set Maintenance



Bonus Product

Mectron Promotions

Combi Touch Perio - \$8,039.90

Purchase: One (1) Combi touch Perio Pack (ME-5020015)

Get: One (1) Set Maintenance (ME-02160035) FREE!

(equipped with: S1, S2, S5, P2, P10, 5 X K10)

You save \$562.10

Air Polishing Slim HP



3 different nozzles

Always best access - 90° and 120° angled spray nozzles for supra and subgingival use in periodontal pockets up to 5 mm depth. The Perio nozzle equipped with the disposable sterile Subgingival Perio Tip provides optimal access in pockets deeper more than 5 mm. All the nozzles need just a click to get connected and they are safely fixed to the air-polishing handpiece thanks to a dedicated security connection system.

\$473 each.
Buy any 3 and get the 4th Free



Subgingival Perio Tip

Flexible and soft - disposable Subgingival Perio Tip will come sterile and gets simply attached to the perio nozzle by hand.

PiezoSurgery touch



\$12,764

Inserts not included



PIEZOSURGERY HAS DEFINED A NEW DIMENSION IN BONE SURGERY

- It's Simple, stress-free & selective.
- It increases the quality and success for you and your patients

PiezoSurgery white



\$6,820

Handpiece and inserts not included

Mectron Promotions



PIEZOSURGERY touch Promo - \$12,764

Purchase: One (1) PiezoSurgery touch unit (ME-05120061) - Delivered with a LED handpiece (complete with cord), a K8 torque wrench and suitcase

Get: One (1) PiezoSurgery Basic Kit (ME-01520001) FREE!

(Delivered with an autoclavable Insert box, OT2, OT7, OP1, OP3 and EX1 inserts)

Plus: One (1) Maintenance Pack FREE!

(includes - 1x Enzymec, 10 x Irrigation sets, 10 x Sterile Protection foils)

You save
\$1,353



PIEZOSURGERY white + LED Handpiece Promo - \$8,970

Purchase: One (1) PiezoSurgery white unit (ME-05120063) - Delivered with a K8 torque wrench and suitcase **PLUS PiezoSurgery LED Handpiece complete with Cord**

Get: One (1) PiezoSurgery Basic Kit (ME-01520001) FREE!

(Delivered with an autoclavable Insert box, OT2, OT7, OP1, OP3 and EX1 inserts)

Plus: One (1) Bottle of Enzymec FREE!

You save
\$1,235



PIEZOSURGERY white + Handpiece (non LED) Promo - \$8,635

Purchase: One (1) PiezoSurgery white unit (ME-05120061) - Delivered with a K8 torque wrench and suitcase **PLUS PiezoSurgery White Handpiece complete with Cord**

Get: One (1) PiezoSurgery Basic Kit (ME-01520001) FREE!

(Delivered with an autoclavable Insert box, OT2, OT7, OP1, OP3 and EX1 inserts)

Plus: One (1) Bottle of Enzymec FREE!

You save
\$1,235

INSERT QUALITY



Mectron ultrasound inserts are built with medical grade, high quality, stainless steel.



Depending on the indication, the inserts are coated with specially selected diamonds, The granulometry of the diamond coating is adapted to the respective treatment.



A coating of titanium nitride, applied to inserts which treat bone, increases working life. Each insert is checked in detail before getting the OK for sales.

ME-01520022
RESECTIVE PERIO KIT


RRP: \$1,755.60
NOW: \$1,492.26



Inserts:
OT13, OT14, OP8, OP9, OP5A

ME-01520023
SINUS LIFT LATERAL KIT

RRP: \$1,173.70
NOW: 997.64



Inserts:
SLC, SLO-H, SLS, SLE1, SLE2

BEST SELLING INSERTS

Take
15% off



OT12
\$353.10



OT12S
\$353.10



OT7
\$305.80



OT7A
\$305.80



OT5A
\$305.80



OT5B
\$305.80



OT7-20
\$404.80



EXP4-L
\$355.30



EXP4-R
\$3355.30



EX3
\$218.90



Innovation in Bone Regeneration

TT is an advanced system in the area of tissue engineering. Simply with a push of a button, an extracted tooth can be processed and transformed into useful bone graft material in a few minutes.

The instrument is able to transform the tooth into autologous bone graft material. Stimulate the adhesion, the proliferation and the cellular differentiation to promote the bone regeneration. Great bio-compatibility and the absence of adverse immune reaction ensure safety.

The high wettability allows the ease of handling and promotes bone regeneration. The tooth like a bone, is made of collagen type 1 and minerals of hydroxyapatite (HA).

The tooth HA is of high crystallinity and it's not resorbable. The Tooth Transformer reduces the crystallinity of the HA eliminate the bacteria and transforms the dentin into an autologous graft material.



RRP - \$10,400
NOW - \$9,950
Includes: 3 x single use

The easiest and safest way
to turn teeth into
autologous graft material

Tooth Transformer –
Single use cartridge
\$150.00 each
Buy 5 Get 1 FREE.





PROCESS FOR PRF



Centrifuge MD Class Ila - € 2797

INDICATIONS :

- Sinus lift
- Bone grafting
- Periodontology
- Implants
- Healing
- Extractions

**i-PRF Purple tubes not included in PRF System



*Photos non contractuelles

PRF SYSTEM - STARTER

PRF System "STARTER" contains :

- x1 Centrifuge PRF DUO Quattro
- x1 Polysteribox L
- x1 Silicone Mat for Polysteribox L
- x1 Box of 100 tubes A-PRF X
- x2 Box of 24 tubes S-PRF
- x2 Box of 25 blood collectors
- x1 PRF BoX Autoclavable
- x1 PRF tweezer
- x1 PRF tweezer "GIRAFFE"
- x1 PRF scissors
- x1 PRF Pad
- x1 Bowl
- x1 mini tray
- x1 tube holder
- x1 Tourniquet
- x1 Compactor Big
- x1 Compactor Small
- x1 Double Spoon

PLUS
1 x PomPac
1 x PomCol
Bonus Free



Everything you need to get you started
\$6,600
(RRP \$7016.21)

ADVANCED PRF


Why Choukroun's A-PRF™ ?

Many recent studies have shown the interest and potential of white cells in the inflammatory cascade, as a corollary, a prominent action in the early days of stimulation of osseo-progenitor cells. (Omar & Thomsen, Biomaterials, 2012, Kawazoe T, Cell. Transplant 2012)

It was therefore natural to try to capture the whole amount of monocytes in the PRF, to make it more active in stimulating bone grafts, but also to turn to a more rapid transformation of monocytes into macrophages to increase the effect bone stimulation.

This is done with the A-PRF™ or Advanced-PRF™.

V.I.S.T.A. stands for Vestibular Incision Subperiosteal Tun-neling Access. It contains a collection of 6 instruments for innovative tunneling techniques and practical surgery.




**RRP - \$660
NOW: \$561**

Fast System contains:

- 1 Pack box
- 1 screwdriver on c/angle
- 1 manual screwdriver
- 1 drill bit
- 2 Kelly clamps
- 1 fast clamp
- 23 screws (from 3.5 mm to 11 mm)
- 2 meshes 5 holes
- 1 mesh 12 holes
- 5 Spike with 5 screws
- 1 PolySteribox® SH + 1 mat + 50 sterilization labels + 100 security studs.


FAST SYSTEM

Allows the simplification of R.O.G. and bone augmentation (vertical and horizontal).






**RRP - \$3,586
NOW: \$3,048**

| Description | Length |
|--|--------|
|  Tent Screw dm | 6 mm |
| | 9 mm |
| | 12 mm |
| | 15 mm |
| \$63.76 each BUY 6 GET 1 FREE MIX N MAX OK | |



| | surgitime ptf ^e | surgitime 3d | surgitime titanium | surgitime seal |
|---|----------------------------|--------------|--------------------|----------------|
| Sinus Elevation | ■ | | | ■ |
| Infra-bone Defect | ■ | | | |
| Horizontal edge Increase | ■ | ■ | ■ | |
| Vertical edge Increase | ■ | ■ | ■ | |
| Concomitant edge increase With implant | ■ | ■ | ■ | |
| Alveolar bone preservation After extraction | ■ | | | ■ |
| Exposed thread treatment of implants | ■ | ■ | ■ | |
| Interproxima Bone formation | ■ | | | |
| Alveolar bone preservation Along with the implant | | | | ■ |

**The data in this document is provided as a general example only. The selection of suitable products and the amount of use should be determined based on the professional's clinical assessment considering the patient's systemic conditions, bone type, severity of the defect, and additional parameters, when necessary. Reported data may vary depending on patient and clinical condition.

| Product | Description | Price |
|---|---|-----------------|
|  | Surgitime Titanium Mesh 34 x 25 mm – Thickness 0.04 mm – Round hole 0.15 mm | \$134.79 |
| | Surgitime Titanium Mesh 34 x 25 mm – Thickness 0.04 mm – Round hole 0.85 mm | \$123.87 |
| | Surgitime Titanium Mesh 34 x 25 mm – Thickness 0.08 mm – Round hole 0.85 mm | \$123.87 |
| | Surgitime Titanium Mesh 60 x 60 mm – Thickness 0.08 mm – Round hole 0.85 mm | \$346.10 |
|  | Surgitime PTFE Membrane Rectangular 30 x 20 mm – Thickness 0.10 mm | \$65.58 |
| | Surgitime PTFE Membrane Rectangular 30 x 20 mm – Thickness 0.25mm | \$65.58 |
|  | Surgitime Titanium SEAL 34 x 25 mm • Thickness 0,04 mm | \$109.30 |

Bionnovation Bone Grafting Screw Fixation Kit

Surgically-inserted medical and dental device used on bone graft surgery for the fixation of grafts and membranes on the maxilla or mandible.

Set Pieces

- 13066 - Driver Handle
- 13129 - Philips Connection for Manual Drive (Ratchet) Short
- 13130 - Philips Connection for Manual Drive (Ratchet) Long
- 05051 - Twist Drill – Ø 1,0 x 15 mm
- 05053 - Twist Drill – Ø 1,2 x 15 mm
- 05055 - Twist Drill – Ø 1,4 x 15 mm
- 05057 - Twist Drill – Ø 1,6 x 15 mm
- 13132 - Philips Connection for Contra Angle Short
- 13133 - Philips Connection for Contra Angle Long
- 13127 - Installation rod – 70 mm
- 13085 - Screwdriver handle

RRP - \$947.23

**Take
15% off**





| | | | |
|---------|---------|---------|---------|
| BW001 | BW002 | BW003 | BW004 |
| \$59.91 | \$59.91 | \$59.91 | \$59.91 |

Buy 4 or more boxes of 10 and receive a stainless steel handle Free



| | | | | | |
|----------|----------|----------|----------|----------|----------|
| BW004M | BW0064M | SB001 | SB002 | SB003 | SB004 |
| \$104.79 | \$104.79 | \$104.79 | \$104.79 | \$104.79 | \$123.27 |

Buy 4 or more boxes of 10 and receive a Titanium Alloy Handle Free

IH001 - Stainless steel 11cm - \$61.45



IH001T - Titanium Alloy 11cm - \$152.46



IH002 - Stainless Steel 14cm - \$70.69



IH002T - Titanium Alloy 14cm - \$161.71





NOUVAG plays a leading role especially in the field of motor driven instruments for Dental Surgery, Micro Surgery and Plastic Surgery.

Every year, around 20,000 new appliances and instruments are used by customers around the globe. Thanks to intensive research and a high innovative potential, our small family owned Swiss company has succeeded in occupying a leading market position in the specialized market niche of medical technology.

Nouvag Mobile Dental Surgery Unit 30 – Complete Unit

Consists of:

MD 30 Implant Motor plus 20:1 contra angle

- control unit with LCD Display (3330)
- 2 micromotor socket connection
- built-in irrigation pump
- electric motor 21, 50,000 rpm – 2m autoclaveable cable with metal plug
- Vario-foot control IP 68 (1510), electronic, suitable for operating theatre
- Contra-angle 32:1 with internal / external cooling system (5201)
- bottle holder
- 2 x spray nozzle attachment for NOU-Clean (1974 /1958)

Vacuson 40 Suction Pump


- control unit with manometer 0-1 bar, vacuum regulator, 2 holders for bottles
- 2 x suction jar 2lt, unbreakable polysulfon
- 2 x lids for 2/5 litre suction jar
- 10 x disposable bacterial filters (4246)
- On/ Off pedaal
- 5 m power cord

RRP - \$16,135.39
NOW: \$13,000

Mobile Stand to house electrical cables and tubing



... makes the difference
devedmed unites best comfort and well thought out
details with our exclusive F-Line Design.
Functional qualities are standard implementations in
our F-Line products.



9100-00 – GPM Wash Tray (275 x 180 x 35mm) with top cover with handhold

Implant & Surgery Set
RRP: \$2,855.51
NOW: \$2,499

Devedmed Implant and Surgery Set,
 by Dr Henry Ho –
 comprising of
 18 implantology and surgical instru-
 ments.

German manufactured



To download the full instru-
 ment catalogue go to
[https://biomedent.com.au/
 wp-content/uploads/2021/01/
 devedmed_BSC4_Catalog.pdf](https://biomedent.com.au/wp-content/uploads/2021/01/devedmed_BSC4_Catalog.pdf)



1086-41 F - Minimal invasive needle hold-
 ers with ratchet "Barraquer" TC, 0 0,8mm,
 17,3cm, curved



1151-30 F - Scalpel handle #3 straight



1170-65 F - Operating micro-scissors
 16,0cm, "Gomel"; curved, fine serrated



2302-66 F - Atraumatic tweezers, micro
 "Cooley" 17,3cm, 0,6mm, straight, serrated



2385-00 F - Mouth mirror handle 1M



2406-91 F - Mazor #7



2406-91 F - Depth Gauge Mazor #7



2500-15 F - Periosteal 2.0 mm straight



2500-18 F - Periosteal 2.0/2,0 mm an-
 gled/curved over surface



2500-21 F - Membrane Lift Mazor #6



2500-30 F - Applicator "Palti"



2500-42 F - Compactor Applicator Mazor



2500-51 F - Papilla elevator "Wachtel"



2500-81 F - Membrane positioning
 instrument "Buser"



888-B1 F - Micro periosteal elevator
 "Buser #1" 17.3cm



888-B2 F - Micro periosteal elevator
 "Buser # 2"



901-25 F - Bone curettes Miller-Colburn
 #10

GINGISAGE




Kereden Founded in 1986 by Dr Patrick Lesage is the original inventor founder of the patented gingival retraction pastes technique. This new technique had revolutionized tissue management in dentistry, as it was the first paste that made it possible to manage sulcular opening and hemostasis simultaneously without the use of retraction cord. 20 years from inception KEREDEN laboratory improves the technique, formulation and delivery system to fit clinical demands.

- New formula means no need for specific injector
- Master dissolution in saliva and gingival fluid
- Optical impressions are made easier
- Remove the controversial use of aluminum.

TODAY HERE IS GINGISAGE™regular
The only aluminum free retraction paste
Composition is kaolin with ferric chloride 6%

Single dose capsules fit composite dispensers.



The First
aluminum chloride free
gingival retraction paste



INDICATIONS


- All procedures requiring the sulcus to be opened and/or dried
- Traditional crown and bridge impressions
- Digital impressions
- Cementation of prosthesis
- Composite or Glass Ionomer restorations
- In conjunction with a single cord
- After gingivectomy / gingival eviction

CONTRAINDICATIONS

- Do not use in the case of poor periodontal condition : furcations, deep pockets.

RRP - \$132
Buy 3 get 1 free

Short absorption time
Highly flexible with outstanding knotting properties
Controlled knot run-down with reliable knot security
Monofilament suture structure: no capillarity, no sawing action
Very good passage through tissue without additional trauma
Optimum tissue compatibility




Absorption properties

Absorbable suture materials approximate the tissues during the healing phase and progressively lose their tensile strength by standardised absorption. GLYCOLON® is a rapidly absorbed suture which loses half of its tensile strength after about 9 days.



Biocompatibility

The special processing of GLYCOLON® leads to rapid absorption in the body. GLYCOLON® with its high degree of chemical purity is broken down in the tissues via the harmless intermediate products polyglycolic acid and caprolactone. In turn, these are enzymatically cleaved and metabolised. Potential tissue reactions are further minimised by the monofilament suture structure.



| REFERENCE | PRODUCT | DESCRIPTION |
|-----------|--------------------------------|--|
| OZ-52160 | Glycolon 4/0 18mm DS18-PB41420 | 4/0 18 mm 3/8 Circle R/Cut Violet 45 cm DS18 - Box of 24 |
| OZ-52202 | Glycolon 3/0 18mm DS18-PB41413 | 3/0 18mm 3/8 Circle R/Cut Violet 45cm DS18 - Box of 24 |
| OZ-52507 | Glycolon 5/0 18mm DS18-PB41421 | 5/0 18mm 3/8 Circle R/Cut Violet 45cm DS18 - Box of 24 |
| OZ-52510 | Glycolon 5/0 16mm DSM16PB41509 | 5/0 16mm 3/8 Circle R/Cut Violet 45cm DSM16 - Box of 24 |
| OZ-52511 | Glycolon 5/0 13mm DSM13PB41503 | 5/0 13mm 3/8 Circ R/Cut Violet 45cm DSM13 - Box of 24 |
| OZ-52512 | Glycolon 6/0 11mm DSM11PB41512 | 6/0 11mm 3/8 Circ R/Cut Violet 45cm DSM11- Box of 24 |
| OZ-52514 | Glycolon 6/0 13mm DSM13PB41527 | 6/0 13mm 3/8 Circle R/Cut Violet 45cm DSM13 - Box of 24 |
| OZ-52515 | Glycolon 4/0 18mm DSM18PB41515 | 4/0 18 mm 3/8 Circle R/Cut Violet 70cm DSM18 - Box of 24 |

Dr Mark's HyGenie®

HyGenie - Demonstration Kit

Complete demonstration and presentation package of all Dr Mark's HyGenie products for your dental appliance wearing patients.

This is the right kit to help show your patients how they can best care for their important dental appliances, at home, quickly, effectively and safely. Saves the Practice or Clinic team time and helps build on your already strong patient relationships.

INCLUDES:

1 Dr Mark's HyGenie + 1 Dr Mark's SureGrip + 1 DentalFresh (50ml) + 1 DentalFresh (125ml)



\$42.95



HyGenie Sport- Demonstration Kit

Complete demonstration and presentation package of all HyGenie Sport products for your mouthguard wearing patients

INCLUDES

1 HyGenie SPORT + 1 50ml DentalFresh + 1 100ml QuickClean + 1 CAGE sports mouthguard carrier + 1 'Clip It' Carabiner Pack

\$42.95

HYGENIE SPORT - 5 PACK



\$199.95

Dental Professionals Practice Bundles

HYGENIE - 10 PACK



\$235.95

MANFLOSS™
SERIOUSLY GOOD FLOSS

QUEEN OF CLEAN™
DENTAL TAPE



SERIOUSLY GOOD FLOSS
AUSTRALIAN OWNED AND DESIGNED

- Refillable dispenser with 2 x 50m rolls of tape
- Super wide design makes it seriously effective
- Wafer thin to slip between teeth
- Soft to land on gums with a freshmint aftertaste
- sticky for grip and to absorb what's been dislodged
- Shred resistant

FOR OPTIMUM ORAL HYGIENE

Choice of two luxurions wafer thin but also soft and wide tapes designed to stop patients cursing about flossing! *PFAS FREE*

MANFLOSS - made black to show white plaque being removed.

QUEEN OF CLEAN - white for the squeaky clean amongst us who prefer not to know...

Dispensers -
Single pack - \$28 each
10 pack - \$246.40
Mix N Match OK

Refills -
Single Pack - \$11
10 pack - \$95.81
Mix N Match OK



1 x 50m Rolls



X-Soft Interdenal Brushes
and Angle Interdenal Brushes
Packs 25 - \$24.47 6 sizes to choose from.
Buy 3 packs of 25 get 1 free. (Mix N Match OK)

Orthodontic Kit




RRP \$22.28
Professionals
\$15.60



Implant / Ortho
TP-442100
Box 25 - \$70.13



Interspace Soft or Interspace
Medium
Box 14 - \$97



Compact Tuft
Box 25 - \$67.38



Universal Care
Box 14 - \$60.06



EasyPicks 36pc
TP-242310 Orange XS/S
TP-242310 Blue M/L
Box 10 - \$82.50
Mix N Max OK

- 1 x Supreme Compact Toothbrush
- 1 x Compact Tuft Brush
- 1 x Angle Yellow - sample (1 pc)
- 1 x EasyPicks XS/S - sample (2pc)
- 1 x EasyPicks M/L - sample (2pc)
- 1 x IDB Green and Yellow - sample (2pc)
- 1 x Orthodontic wax



Caredent Orthobrush

Designed for orthodontics. This unique toothbrush comes with a compact head together with special manufactured soft "V" shaped Dupont Tynex bristles. Stylish design also incorporates a non-slip TPE handle.

**Box of 6
\$7.04**



CareDent OrthoWax

Caredent Ortho Wax helps relieve wire and bracket irritation. Wax comes in a convenient carry case A natural wax product Easy to use Cut off the required amount of wax, roll it into a ball and apply it on the bracket or brace which is irritating the soft tissue.

**6 pks
\$6.99**

CareDent Interbrush In A Box

The soft rubber-like bristles are wire-free, gentle and flexible, making them ideal for hard to reach spaces in between teeth, around bridges, implants, crowns and orthodontic appliances. Latex-free. 30pc packs

**60 pks
\$237.60**



Caredent Orthodontic Kit

Each kit is designed to establish good oral care habits and contains everything needed for a beautiful smile just add your favourite toothpaste. The kit has been assembled to cater for all essential orthodontic needs. The all-in-one kit is the perfect item to carry when travelling and the sturdy zippered carry bag can be kept for future use.

Contains:

- Orthodontic V-Brush
- Disposable Dental Mirror
- Regular OrthoWax
- InterBrush Small (5 pack)
- Disclosing Tablets (2 tablets)

**10 packs
\$69.80**

hismile



\$44.10

NHpro Enamel Care is a concentrated serum manufactured from nani-hydroxyapatite that is intended to restore lost minerals and reduce sensitivity.

Clinically proven to improve enamel and alleviate tooth sensitivity.



Free Gift



\$31.20

HA5 Hyaluronic Gum Serum 30ml is a hydrating serum made from hyaluronic acid, designed to improve oral hygiene.

- hydrates soft tissue
- promotes healthy oral hygiene
- relief for sore gums
- anti-inflammatory
- antibacterial



\$28.00

HA5 Hyaluronic Mouthwash 290ml
A hydrating mouthwash made from hyaluronic acid, designed to improve oral hygiene.

- freshens Breath
- promotes healthy oral hygiene
- anti-inflammatory
- antibacterial

**Buy 2 x NHPro, 2 x ha5 Gum Serum, 2 x ha5 Mouthwash
\$206.60**

**receive 1 x Take Home Whitening Kit
Valued at \$149**



Worlds Smallest Mini Water Flosser Portable

Buy 1 - \$70.61
Buy 5 - 297.33 -
(\$59.47 each)

Mini Water Flosser is ideal for removing plaque and debris from deep between teeth and below the gum line, where standard brushing and flossing can't reach.

Dental professionals recommend flossing every day to help avoid gum disease. Many of us, however, do not use string floss since it is difficult and time consuming to use.

Water Flossing has been clinically demonstrated to be up to 50% more effective than traditional floss.

- Budget-friendly, portable, lightweight, collapsible design
- Ideal for daily use, travel and business trips
- 20-35 seconds of flossing time
- Two modes: regular floss mode and hydro-pulse massage mode for gum stimulation.
- Controlled 360* rotational nozzle
- Safe to use on implants, crowns, veneers, braces, and other dental work.

Product Parameter

- Water Tank capacity: 150ml
- Voltage: 100-240VAC
- Battery: 500mAh polymer battery, can be staked to 1000mAh
- Power: 5W 0.5A
- Operation mode: standard/pulse mode
- Outlet pressure: 80PSI/10g
- Use time: 15-20 days (full charge)
- Water jet pulsation frequency: 1400 - 1650 pulse/min
- Noise standard: less than 65 decibels
- Waterproof rating: IPX7

ADDRESS

PO Box 25, North Ryde BC NSW 1670

ONLINE

www.biomedent.com.au

Info@biomedent.com.au

CALL

1300 792 624