

Bionnovation's GRAFT Sample Data & Procedures Table

Bionnovation
implants and biomaterials
Engineering for the perfect smile!

Product Group	Bone Substitutes						Membranes • Barriers • Mesh			
Source	Xenograft			Synthetic			No reabsorbable			
Product Image										
Product Name	Bonefill (Bovine Bone)	Bonefill Porous (Bovine Bone)	Bonefill Porous Block	Hydroxyapatite	Hydroxyapatite Block	Beta TCP	Surgitime PTFE	Surgitime E-PTFE	Surgitime Titanium	Surgitime Titanium SEAL
Sinus Floor Elevation		X								
Socket/Ridge Preservation	X		X			X				X
Horizontal Augmentation	X		X	X			X	X	X	
Vertical Augmentation	X		X	X			X	X	X	
Periodontal Defects	X	X		X		X	X	X		
Covering Palatal Gingival Donor Sites										
Estimated Re-entry time (by months)	9	6	6	9	9	4				
Estimated Main Integration (by months)	9	6	6	9	9	4				