



Products Catalogue

 1300 792 624

 www.biomedent.com.au

Summary

Agency	Page
• Bionnovation	1 - 3
• Cervico by VPI	4
• DentaMed	5
• DeveMed	6
• Embi	7
• Gingisage	8
• Glycolon - Resorba	9
• Kulzer	10 - 11
• Mectron - Rex Implants	12 - 14
• MJK	15
• Nouvag	16
• Oh! Care	17
• Numedico	18
• PRF	19 - 21
• Tad Kit and Mini Plates	22 - 24



Triple Blade Drill - Bone Collector





- Easy and quick insertion.
- Auto cutting power even at low speed;
 - Collects and particles the bone quickly.




Cortical Drill



Decortication Bur with drill stop.
Perforation of the cortical bone layer advocated in GBR, because it was postulated that this increases vascularity of wound and releases growth factors and cells with angiogenic and osteogenic potential.



Description	Length
Tent Screw 	6 mm
	9 mm
	12 mm
	15 mm

Product	Description
	Surgitime Titanium Mesh 34 x 25 mm – Thickness 0.04 mm – Round hole 0.15 mm
	Surgitime Titanium Mesh 34 x 25 mm – Thickness 0.04 mm – Round hole 0.85 mm
	Surgitime Titanium Mesh 34 x 25 mm – Thickness 0.08 mm – Round hole 0.85 mm
	Surgitime Titanium Mesh 60 x 60 mm – Thickness 0.08 mm – Round hole 0.85 mm
	Surgitime PTFE Membrane Rectangular 30 x 20 mm – Thickness 0.10 mm
	Surgitime PTFE Membrane Rectangular 30 x 20 mm – Thickness 0.25mm
	Surgitime Titanium SEAL 34 x 25 mm • Thickness 0,04 mm



Bone Expander Kit



Used for bone expansion & condensation, executing preparation and instrumentation for installation of dental implants

Kit Includes:

- 13086 Container with lid
- 13066 – Digital Driver Adaptor
- 13113 – Expansor self-tapping 1,7/2,4 mm
- 13114 – Expansor self-tapping 2,0/3,1 mm
- 13115 – Expansor self-tapping 2,5/3,4 mm
- 13116 – Expansor self-tapping 3,0/3,9 mm

Bone Screws

Bionnovation Bone Screws	Diameter			
Length	1.2 mm	1.4 mm	1.6 mm	1.8 mm
3 mm	BN-7097			
4 mm	BN-7098	BN-7145	BN-7191	BN-7236
6 mm	BN-7092	BN-7090	BN-7093	BN-7238
8 mm	BN-7101	BN-7148	BN-7194	BN-7240
10 mm	BN-7103	BN-7150	BN-7094	BN-7242
12 mm	BN-7105	BN-7152	BN-7095	BN-7244
14 mm	BN-7107	BN-7091	BN-7198	BN-7246

Bionnovation Graft Screws: Are used for the fixation and stabilization of bone grafts and nonabsorbable membranes (barriers) used in guided bone regeneration (GBR).

Drilling speed: 200 rpm;
Placement speed: 30 rpm;





Bone Grafting Screw Fixation Kit



Surgically-inserted medical and dental device used on bone graft surgery for the fixation of grafts and membranes on the maxilla or mandible.

*Items also sold individually

Kit Includes:

- Driver Handle 13066
- Philips Connection for Manual Drive (Ratchet) Short 13129
- Philips Connection for Manual Drive (Ratchet) Long 13130
- Twist Drill – Ø 1,0 x 15 mm 05051
- Twist Drill – Ø 1,2 x 15 mm 05053
- Twist Drill – Ø 1,4 x 15 mm 05055
- Twist Drill – Ø 1,6 x 15 mm 05057
- Philips Connection for Contra Angle Short 13132
- Philips Connection for Contra Angle Long 13133
- Installation rod – 70 mm 13127
- Screwdriver handle 13085



Cervico Premium Kit



The Cervico Mold is an in-office tool used for the fabrication of the cervico healing abutments and impression posts.

It can also be used for the duplication of a modified healing abutment or temporary prosthesis to a duplicate impression post. It features 17 wells that can be used for all of the above purposes.

What's included:

- Base screw driver
- VPI prosthetic connection inserts*
- Centered Silicone Insert
- Mold Base

Cervico Essencial Kit



Essential Cervico Mold is a simplified version of the Cervico Premium Mold that works with the stock implant analogs available from the different implant companies. As per the premium mold you can fabricate Cervico healing abutments and impression posts by utilising one or two pieces prosthetic components, made out of titanium, or peek, or other suitable material as they come available by the implant companies.



Dry Mouth Oral Hygiene Gel 65mL

- Provides up to 8 Hours relief from Dry Mouth
- Moisturises, hydrates, soothes & comforts
- Fluoridated formulation, replacing toothpaste

Dry Mouth Oral Hygiene Spray 50mL

- Moisturising, fluoridated oral hygiene spray
- Convenient oral hygiene spray with 360 degree nozzle
- Freshens breath & reduces halitosis (mouth odour)

Moisturising Mouthwash 500mL

- Enriched with all natural moisturising Malic Acid
- Leaves the mouth with a soft, velvety mouth-feel
- pH balanced, moisturising, calming mouth rinse

Moisturising Mouth Spray 60mL

- Convenient and immediate relief from Dry Mouth
- Enhanced moisture delivery with (prop.) chitin
- A 'nod' to the Australian Bush with a hint of lemon myrtle

Denture Adhesive Gel 50mL

- Secures ALL day (naturally) Seals, Cushions & Comforts
- PVM/MA Copolymer "glue" (free) formulation
- Neutral taste, Zinc free, Easy to clean

pH Neutralising Mouthwash 250mL

- Helps prevent enamel erosion
- Naturally assists to neutralise pH levels
- Refreshing flavour, safe to use ongoing






Implant & Surgery Set



Devemed Implant and Surgery Set, by Dr Henry Ho – comprising of 18 implantology and surgical instruments. German manufactured

Sinus Lift Set – Small



Sinus instrument set, comprising of 8 sinus lift and grafting instruments.

Instruments



DM-1151-30F Scalpel Blade holder



DM-2525-10 F Periodontal Probe



DM-2500-30 F ISS Applicator Palti



DM-2302-66F Atraumatic Micro-Forceps 0.6mm



DM-1101P Extracting Forceps Upper Jaw



DM-1099-15 Needle holder



embi Edge Tongue Cleaner

Our designers were inspired by input from dentists and experienced consumer scraper users.

The Edge was developed to deliver a deeper clean than any other tongue cleaner. It is smooth and safe on the tongue. The design has a deeper scoop and ergonomically weighted handles to offer the most premium scraping experience.

The Edge is the ultimate all stainless steel tongue cleaner and with minimum care, can last you a lifetime. Use your embi every day and take it with you anywhere you go in your **FREE REUSABLE TRAVEL CASE.**



WHAT MAKES embi THE BEST YOU CAN BUY!

1. We are 100% focused on Tongue Scrapers and Tongue Cleaning.
2. embi is the leading brand in Australia, recommended by dentists and oral hygienists.
3. The unique design of our tongue scrapers is the best in class, giving a safe and deep clean every time, for life.
4. Each embi is flexible in movement, which helps fit the mouth and reduces the gag reflex.
5. You can be assured of the quality medical grade stainless steel and pure copper.
6. We are proudly Australian, 100% Earth Friendly, Sustainable and Reusable.




GINGISAGE™ regular

The only aluminum free retraction paste

Composition is kaolin with ferric chloride 6%

Single dose capsules fit composite dispensers.



The First
aluminum chloride free
gingival retraction paste



Indications:

- All procedures requiring the sulcus to be opened and/or dried
- Traditional crown and bridge impressions
- Digital impressions
- Cementation of prosthesis
- Composite or Glass Ionomer restorations
- In conjunction with a single cord
- After gingivectomy / gingival eviction

Contraindications:


- Do not use in the case of poor periodontal condition : furcations, deep pockets.



Biocompatibility

The special processing of GLYCOLON® leads to rapid absorption in the body. GLYCOLON® with its high degree of chemical purity is broken down in the tissues via the harmless intermediate products polyglycolic acid and caprolactone. In turn, these are enzymatically cleaved and metabolized. Potential tissue reactions are further minimized by the monofilament suture structure.

- Short absorption time
- Highly flexible with outstanding knotting properties
- Controlled knot run-down with reliable knot security
- Monofilament suture structure: no capillarity, no sawing action
- Very good passage through tissue without additional trauma
- Optimum tissue compatibility



REFERENCE	PRODUCT	DESCRIPTION
OZ-52160	Glycolon 4/0 18mm DS18-PB41420	4/0 18 mm 3/8 Circle R/Cut Violet 45 cm DS18 - Box of 24
OZ-52202	Glycolon 3/0 18mm DS18-PB41413	3/0 18mm 3/8 Circle R/Cut Violet 45cm DS18 - Box of 24
OZ-52507	Glycolon 5/0 18mm DS18-PB41421	5/0 18mm 3/8 Circle R/Cut Violet 45cm DS18 - Box of 24
OZ-52510	Glycolon 5/0 16mm DSM16PB41509	5/0 16mm 3/8 Circle R/Cut Violet 45cm DSM16 - Box of 24
OZ-52511	Glycolon 5/0 13mm DSM13PB41503	5/0 13mm 3/8 Circ R/Cut Violet 45cm DSM13 - Box of 24
OZ-52512	Glycolon 6/0 11mm DSM11PB41512	6/0 11mm 3/8 Circ R/Cut Violet 45cm DSM11- Box of 24
OZ-52514	Glycolon 6/0 13mm DSM13PB41527	6/0 13mm 3/8 Circle R/Cut Violet 45cm DSM13 - Box of 24
OZ-52515	Glycolon 4/0 18mm DSM18PB41515	4/0 18 mm 3/8 Circle R/Cut Violet 70cm DSM18 - Box of 24










Flexitime

With Flexitime you work with high precision, time flexibility, and according to your personal treatment style.



- **Intelligent time concept:**
The special formula controls the polymerization kinetics by oral temperature and sets rapidly only after insertion into the mouth.
- **Automatic mixing:** High-speed mixing with our Dynamix speed system is easy and saves precious time..
- **Wide range of indications:** Flexitime can be used for a wide range of indications from single-tooth preparation over multiple-unit prosthetics to a full denture. The A-silicone is available in 3 different delivery forms and 6 different viscosities.

Two-step, two-phase impression		Single-step, two-phase impression		Single-step, one-phase impression
Impression technique	Correction technique	Double-mixing technique	Sandwich technique	Monophase impression
Viscosity of tray material	Highly viscous, kneadable	Highly viscous	Highly viscous, kneadable	Medium viscous
Product recommendation	 Step 1 Heavy tray / Putty  Step 2 Light flow	 Heavy tray +  Medium flow/ Correct flow	 Heavy tray / Putty +  Medium flow / Correct flow	 Monophase
Intended use	<ul style="list-style-type: none"> • Crowns & Bridges 	<ul style="list-style-type: none"> • Transfer & Implant impression 	<ul style="list-style-type: none"> • Crowns • Bridges • Inlays • Onlays 	<ul style="list-style-type: none"> • Crowns • Bridges • Inlays • Onlays • Transfer & Implant impression • Functional impression



Venus Bulk Flow One

Always grab the right shade for the posterior region. The new Venus Bulk Flow ONE offers a one-shade solution within the Venus family for basic posterior restorations. With its incredible shade-matching properties, Venus Bulk Flow ONE, blends seamlessly into the surrounding dentition, independent of the tooth shade. Even when used for bulk-filling cavities of up to 4 mm.

The ONE shade has been optimised to deliver the perfect balance of chroma and translucency that allows the restoration to become virtually invisible. Therefore, the ONE shade can be used with confidence in any Class I and Class II restoration. And this colour adaptation is what we call Adaptive Light Matching



Streamlined restorative workflow in everyday posterior restorations. Venus Bulk Flow ONE combines the self-levelling effect of a flowable material with a high curing depth and good stability after application. This allows to place increments of up to 4 mm inside a cavity. With Venus Bulk Flow ONE dentists can reduce treatment time without the need for shade selection or the placement of a universal composite as an additional capping layer. With its range of indications from bulk filling through cavity lining to fillings in primary teeth, the efficient restoration procedure is further enhanced by the equally reliable and versatile bonding solution iBOND Universal.



PIEZOSURGERY® touch



What's included:

- 1 PIEZOSURGERY® touch unit
- 1 LED-handpiece, complete with cord
- 1 torque wrench K8
- 1 suitcase

METHOD PIEZOSURGERY®: EVIDENCED-BASED

With over 12 years of successful clinical research, the original PIEZOSURGERY® method is the only one supported by more than 120 clinical and scientific studies

PIEZOSURGERY white



Mectron PiezoSurgery White is as easy as it gets.

If you are looking for a device which provides simplicity and clarity in every detail – here you are: 100% simple handling for utmost treatment security.

Especially designed for easy cleaning, disinfection and sterilization. Cost-effective standard parts for perfect economy.

You call that perfect for daily use.

***Handpiece and inserts not included**



Rex Piezoimplants

AN IMPLANT DESIGNED FOR THE ANATOMY

The REX PiezoImplant System was designed specifically to manage reduced crestal width. Its rectangular section mimics the anatomy of the residual crestal bone, allowing for adequate bone thickness to be preserved on the vestibular and lingual sides. This level of bone preservation is not possible with traditional, screw-form implants because their circular section poses a geometrical discrepancy with the anatomy.

The morphological advantages of REX PiezoImplants allow preserving a crestal thickness greater than 1.5 mm buccally and lingually even when the residual crest is only 3.5 mm

Site Preparation Rexpander



Rex Piezoimplant Surgical Kit
and Surgical Instruments



Rex Implant Placement Device
(IPD)



Rex
Piezoimplants
Unique



Piezosurgery Inserts for Rex Implants

Implant site preparation is performed using dedicated PIEZOSURGERY® inserts with a rectangular section and increasing thickness.

Dedicated depth markings allow for the optimal degree of under-preparation in relation to bone density.



Piezosurgery Versatility

In addition to the inserts dedicated to the REX PiezoImplant (patent pending), PIEZOSURGERY® features a wide range of inserts from extractions to implant site preparation for all types of surgical needs.

In particular for soft tissue management: once the flap has been outlined with the scalpel, it can be lifted with PR1 and PR2 inserts that elevate the periosteum in a highly efficient and non-traumatic manner.





Micro-blades – Periodontology



Scalpel Blades – Implantology



Bone Block Applicator



Micro-blades Handles





Nouvag Mobile Dental Surgery Unit 30 – Complete Unit

Consists of:

MD 30 Implant Motor plus 20:1 contra angle

- control unit with LCD Display (3330)
- 2 micromotor socket connection
- built-in irrigation pump
- electric motor 21, 50,000 rpm – 2m autoclaveable cable with metal plug
- Vario-foot control IP 68 (1510), electronic, suitable for operating theatre
- Contra-angle 32:1 with internal / external cooling system (5201)
- bottle holder
- 2 x spray nozzle attachment for NOU-Clean (1974 /1958)

Vacuson 40 Suction Pump

- control unit with manometer 0-1 bar, vacuum regulator, 2 holders for bottles
- 2 x suction jar 2lt, unbreakable polysulfon
- 2 x lids for 2/5 litre suction jar
- 10 x disposable bacterial filters (4246)
- On/ Off pedel
- 5 m power cord

Mobile Stand to house electrical cables and tubing.





oh!CARE

OH! Care® is an exclusive, premium oral health care range.

Without alcohol, surfactants, titanium dioxide or animal by-products.

The focus relies on innovative, natural ingredients such as hyaluronic acid and xylitol to care for gums and protect teeth.



- Mouthwash
- HA-X gel



**Blood collection butterfly with
retractable needle**



Needle 30G 0.5'



Syringe Luer Slip 3ml with needle



**3ml Luer Lock Syringe without
needle**





PRF SYSTEM - STARTER



PRF System "STARTER" contains :

- 1 Centrifuge PRF DUO Quattro
- 1 Polysteribox L
- 1 Silicone Mat for Polysteribox L
- 1 Box of 100 tubes A-PRF X
- 2 Box of 24 tubes S-PRF
- 2 Box of 25 blood collectors
- 1 PRF BoX Autoclavable
- 1 PRF tweezer
- 1 PRF tweezer "GIRAFFE"
- 1 PRF scissors
- 1 PRF Pad
- 1 Bowl
- 1 mini tray
- 1 tube holder
- 1 Tourniquet
- 1 Compactor Big
- 1 Compactor Small
- 1 Double Spoon
- 1 Box of 100 tubes A-PRF
- 1 Box of 24 tubes S-PRF
- 1 Box of 25 Blood collectors
- 1 Box 25 Syringe 18G1.1/2"
- 1 Box of 100 alcohol swabs
- 1 PomPac – cooling block
- 1 PomCol – cooling cuff
- 1 PomSwing – Express Production of Sticky Bone


Indications:

- Sinus lift
- Implants
- Implants
- Healing
- Bone grafting
- Extractions

**Includes PRF
Online Course
Free of Charge**



Indications




A-PRF™
CHOUKROUN Platelet Rich Fibrin

At the end of the spin :

- Remove the tubes from the centrifuge, remove the caps & put them in the tube-holder
- Wait 5 min if the clots are not yet solid
- Scratch the red clot to only keep the fibrin


Don't leave the tube in the machine, neither the clots in the tubes

CLOT




- Sinus lift with PRF alone
- + Simultaneous implant placement
- Periodontal pocket < 6mm
- Cyst cavity filling

MEMBRANE




- Site coverage
- Schneiderian membrane protection
- Equivalent connective graft
- Poncho healing abutment

PLUG



Socket extraction filling




S-PRF™
CHOUKROUN Platelet Rich Fibrin

Green tubes can be launched *together* with the red tubes

At the end of the spin :

- Remove the tubes from the centrifuge, remove the caps & put them in the tube-holder
- Use a *2mL* syringe to take the supernatant, directly into the tube

STICKY BONE




- Socket extraction filling
- Guided Bone Regeneration
- Periodontal pocket > 6mm

Express option :
To accelerate the clotting, add a clot in the container

LARGE MEMBRANE


Extended site coverage




Express option :
To accelerate the clotting, add a membrane in the container

INJECTIONS

Injectable PRF stays liquid fro 20 to 30 minutes
Be careful not to shake the tube, the 2 phases could blend




PRE OP. FLAP INJECTION




PAPILLA / TMJ / ENDO

i-PRF™ **i-PRF M**
Women - 3 min. Men - 4 min.



AESTHETICS / ORTHOPEDICS

i-PRF™ +



Facial aesthetics / Alopecia
Cartilage regeneration
Articular pain



Fast System Kit

Allows the simplification of R.O.G. and bone augmentation (vertical and horizontal).



Soft Brushing Kit

Has been designed to extend the flap without incision and it also reduces the stiffness of the flap.



Veinlite EMS PRO

The world's bestselling vein access device just got better – now with an integrated white exam light.




Soft Brushing Kit

V.I.S.T.A. stands for Vestibular Incision Subperiosteal Tunneling Access. It contains a collection of 6 instruments for innovative tunneling techniques and practical surgery.





BioMeDent TAD Kit (Temporary Anchor Devices)




- 1 - Screwdriver for 1.5mm Screw (Square type)
- 2 - Screwdriver for 2.0mm Screw (Square type)
- 3 - Screwdriver for 1.5mm Screw (AO type)
- 4 - Screwdriver for 2.0mm Screw (AO type)
- 5 - Limited Drill Bit 1.1*6mm (common type)
- 6 - Limited Drill Bit 1.1*8mm (common type)
- 7 - Limited Drill Bit 1.6*6mm (common type)
- 8 - Limited Drill Bit 1.6*8mm (common type)
- 9 - Quick Coupling Handle (AO type)
- 10 - Limited Drill Bit 1.1*6mm (RA/Latch type)
- 11 - Limited Drill Bit 1.1*8mm (RA/Latch type)
- 12 - Limited Drill Bit 1.6*6mm (RA/Latch type)
- 13 - Limited Drill Bit 1.6*8mm (RA/Latch type)
- 14 - Aluminium Box




BioMeDent Plate and Screw Kit



- 1 - Hollow Reamer 2.0mm
- 2 - Hollow Reamer 1.5mm
- 3 - Screwdriver for 2.0mm Screw (AO type)
- 4 - Screwdriver for 1.5mm Screw (AO type)
- 5 - Depth Gauge 0-30mm
- 6 - Limited DrillBit (common type)
- 7 - Limited DrillBit (common type)
- 8 - Quick Coupling Handle (AO type)
- 9 - Guider 1.1mm& 1.6mm
- 10 - Limited DrillBit (RA/Latch type)
- 11 - Limited DrillBit 1.1*4mm (RA/Latch type)
- 12 - Screwdriver for 1.5mm Screw(RA/Latch type)
- 13 - Screwdriver for 2.0mm Screw(RA/Latch type)
- 14 - Aluminium Box





 1300 792 624

 info@bioment.com.au

 www.biomedent.com.au

