



2 0 2 4  
E O F Y  
EDITION

**INNOVATIVE SOLUTIONS**  
FOR DENTAL AND MEDICAL  
PROFESSIONALS



Your trusted source for Dental and Medical solutions, serving health professionals in Australia and New Zealand. Empowering through education

### AGENCIES WE REPRESENT



### INSTANT ASSET WRITE-OFF

A REMINDER: Businesses with an aggregated turnover of less than \$10 million can instantly write-off assets costing less than \$20,000 purchased between 1 July 2023 and 30 June 2024.

# SUMMARY

CATEGORIES	AGENCIES	PAGE
<b>ORAL SURGERY</b>		
Implant Motors	Bien Air	4
Bone Fixation Screws & Titanium Mesh	Bionnovation	5-7
Instruments	Devemed	8
Piezosurgery	Mectron	9
Micro Surgery Blades	MJK	10
Syringes, Needles & Blood Collectors	Numedico	11
Centrifuges & Tubes	PRF	12
Absorbable Sutures	Monoglyc/PRF	13
Absorbable & NonAbsorbable Sutures	Resorba	14-15
Piezol Implant System	Rex Implant	16
Autologous Bone Grafting	Tooth Transformer	17
Surgical Brushes	Vitis	18
TAD & Plates Surgical kit	BioMeDent	19-20
<b>GENERAL DENTISTRY</b>		
Gingival retraction paste	Gingisage	21
Curing Light & Digital Vision	CAO	22-23
Composite & Impression Material	Kulzer	24
Systematic Tool for Restorative Dentistry	Prosmate	25
Prophylaxis Powder	QueenDental	26
Customise Implant Abutment System	VP Innovato - Cervico	27

# SUMMARY

<b>CATEGORIES</b>	<b>AGENCIES</b>	<b>PAGE</b>
<b>FACIAL AESTHETICS</b>		
Aesthetic System / Dermapen	PRF/BioMeDent	28-29
<b>PREVENTIVE DENTISTRY</b>		
Cavity Prevention	CAO	30-31
Dry mouth solutions	Denta-Med	32
Dental Cleaning Device & Sports Mouth Guard	Dr Marks Hygenie	33-34
Tongue Cleaner	Embi	35
Toothpaste	Lucalut	36
Floss	Manfloss	37
Tongue Posture	MyoSpots	38
Mouthwash & Gum Serum	OhCare!	39
Toothbrushes & Interdental Brushes	TePe	40
Electric Interdental Brush	UBrush!	41
Whitening	CAO	42
<b>MAXILLOFACIAL PROSTHETICS</b>		
Silicones, Adhesives, Ocular & Magnets	Technovent	43
<b>BOOKS</b>		
Dental & Medical	Quintessence	44
Terms and Conditions		45-46



## Chiropro PLUS+3rd Gen set with Micro-Series CA 20:1 L



**20% OFF**

**WAS \$10,404.35**  
**NOW \$8,323.48**

Simplicity. The cornerstone of the new range of implant and oral surgery motors developed by Bien-Air Dental.

General dental practitioner or experienced surgeon, are you seeking a fast, simple and efficient solution for your implant dentistry procedures? Look no further than the Chiropro.

Your perfect solution for precise and smooth placement of implants.



### Bien Air PM 1:1 Ext Spray (Slowspeed) Handpiece

PM 1:1 is designed to withstand the stresses resulting from repeated lateral pressure, while maintaining exceptional precision.

**15% OFF**

**WAS \$1,352.49**  
**NOW \$1,149.61**



### Bien Air CA 20:1 L KM Handpiece

This contra-angle has an irrigation system which is internal to the instrument – a world first. So your grip is not impeded by the irrigation line.

**15% OFF**

**WAS \$4,261.07**  
**NOW \$3,621.90**

All prices are GST inc.



## Bone Grafting Screw Fixation Kit




### Kit Includes:

- Driver Handle 13066
- Philips Connection for Manual Drive (Ratchet) Short 13129
- Philips Connection for Manual Drive (Ratchet) Long 13130
- Twist Drill – Ø 1,0 x 15 mm 05051
- Twist Drill – Ø 1,2 x 15 mm 05053
- Twist Drill – Ø 1,4 x 15 mm 05055
- Twist Drill – Ø 1,6 x 15 mm 05057
- Philips Connection for Contra Angle Short 13132
- Philips Connection for Contra Angle Long 13133
- Installation rod – 70 mm 13127
- Screwdriver handle 13085

Surgically-inserted medical and dental device used on bone graft surgery for the fixation of grafts and membranes on the maxilla or mandible.

\*Items also sold individually

**\$1,089.00**





**Bionnovation Surgitime  
Titanium 3DF**

Biocompatible  
Flexible & Strong

12 x 18 mm | 0,08 mm | 0,15 mm

Built to promote bone formation  
Customized for all degrees of bone defect  
is categorized as 3 different shapes:



**\$123.87**



Buccal  
**1 wall**

2 customized shapes  
Cutting Guidelines **Thicker line**



Buccal | Proximal  
medial - distal  
**2 wall**



Buccal | Proximal  
medial - distal  
Lingual  
**3 wall**



## Bionnovation Bone Screws

Used for the fixation and stabilization of bone grafts and nonabsorbable membranes (barriers) used in guided bone regeneration (GBR).

Drilling speed: 200 rpm;  
Placement speed: 30 rpm;



**From \$25.39**

## Triple Blade Drill - Bone Collector



**\$210.65**

Easy and quick insertion.

- Auto cutting power even at low speed;
- Collects and particles the bone quickly.

Tent Screw DM is designed with a self-drilling tip, polished neck.

## Tent Screw DM



**From \$61.60**



**\$114.91**

## Cortical Drill

Perforation of the cortical bone layer has been advocated in GBR, because it was postulated that this increases the vascularity of the wound and releases growth factors and cells with angiogenic and osteogenic potential.



**10% OFF**



**Lucas Bone Curettes**



**WAS \$122.55  
NOW \$110.29**

**Periosteal Elevator Buser #1**



**WAS \$130.74  
NOW \$117.66**

**Watchel Papilla Elevator 173mm**




**WAS \$134.63  
NOW \$121.16**

**Micro Needle Holder – Barraquer**




**WAS \$768.06  
NOW \$691.25**

**Micro-Suturing Kit**



**WAS \$1,688.81  
NOW \$1,519.92**

**Sinus Lift Set – Small**



**WAS \$1,263.20  
NOW \$1,136.88**



## PIEZOSURGERY® white + Handpiece



**\$8,859.40**

When it comes to cutting bone, you can of course use traditional burs and saws. They do cut bone, too – but they do not differentiate: any soft tissue getting in their way will also be cut.


## PIEZOSURGERY® touch



**\$12,764.40**

**SAVE  
\$2,595.51**

Receive a FREE second LED Handpiece and Piezoelectric Bone Surgery book with the purchase of a Piezosurgery White or Touch.



All prices are GST inc.

Page 9



## Micro-blades Periodontology



Specific angulations of the blades allow to operate easily and precisely in the least accessible areas, and in most cases, with a single line of incision.

**From \$104.79**



## Scalpel Blades Implantology



- Approximal incisions (mesial and distal) and interdental spaces
- Papillary incisions
- Palatal incisions (flaps and connective grafts samples)
- Delicate incisions for mucogingival surgery
- Crestal incisions (implantology)
- Sulcus incisions in vestibular or lingual areas
- Vertical incisions in the anterior area (upper or lower)
- Sulcus and papillary incisions in narrow areas, or inaccessible to the blades 15 or 15c

**From \$59.51**



## Blood collection butterfly with retractable needle



Blood collection butterfly with retractable needle and holder.

- 21G x 3/4"
- Box contains 25 collectors
- Single use
- Non toxic

**\$152.42**

## Needle 30G 0.5' Box of 100

Made of Polypropylene/Polyisoprene  
Latex-free

Blister paper packed | EO Sterilization  
Box of 100 pcs

Certification: TGA Approved, CE Marked, FDA  
Approved, US 510k approved, International ISO  
standards 13485, ISO 7864-4



**\$13.72**




## PRF SYSTEM STARTER

Our carefully selected PRF Starter kit contains all the necessary components for your PRF journey



**WAS \$8,665.37**  
**NOW \$7,365.56**

Receive a **FREE** Venipuncture book with the purchase of any PRF System



## PRF SYSTEM PROFESSIONAL

Similar concept to the Starter pack, but with twice the number of instruments and consumables.



**WAS \$11,479.46**  
**NOW \$9,757.54**

All prices are GST inc.



## Monoglyc Absorbable Monofilament Sutures

Monoglyc PGA-PCL Suture is a highly pliable monofilament suture, that also absorbs predictably. With the smooth monofilament surface structures, it handles and ties easily, and is virtually inert in tissue. Basically, PGA-PCL sutures behave like PTFE sutures, but do not need to be removed!

24 Sutures Per Box



\$346.78

- Prevents bacterial wicking
- Fully Absorbable in approximately 90 Days
- Compares to Monocryl
- Sharp & high quality premium needle allow for atraumatic passage through tissue.

**BUY 3 GET  
ONE FREE**



## GLYCOLON SUTURES

GLYCOLON<sup>®</sup> is a monofilament suture material, whose degradation rate in comparison with PGA is clearly reduced by the specially selected ratio of polyglycolic acid to caprolactone.

- Short absorption time
  - Highly flexible with outstanding knotting properties
  - Controlled knot run-down with reliable knot security
  - Monofilament suture structure: no capillarity, no sawing action
  - Very good passage through tissue without additional
  - traumatisation
  - Optimum tissue compatibility
  - GLYCOLON<sup>®</sup> – specially developed for indications
  - with short wound healing:
  - Intracutaneous skin closure, Plastic surgery
- Gynaecology and obstetrics, Urology

**\$279.99**





**BUY 3 GET  
ONE FREE**

## RESOLON SUTURES

Long-term durability with simultaneous high biocompatibility

- Optimal tissue compatibility in the body due to the careful selection of materials and refining process
- Easy removal
- Non-absorbable suture materials remain virtually unchanged in body tissues. Once the scar tissue of the wound edges has become sufficiently strong to hold the wound together, the doctor removes the suture material by simply pulling it out.

Resorba Monofilament Polyamide  
Blue 36/box

- Non irritating
- Softer and more subtly
- Easy passage through tissue
- Perfect handling and for consistency knotting



**Bleu**

**\$217.80**





## REX PIEZOIMPLANT

The REX PiezoImplant System was designed specifically to manage reduced crestal width.

Its rectangular section mimics the anatomy of the residual crestal bone, allowing for adequate bone thickness to be preserved on the vestibular and lingual sides. This level of bone preservation is not possible with traditional, screw-form implants because their circular section poses a geometrical discrepancy with the anatomy.

The morphological advantages of REX PiezoImplants allow preserving a crestal thickness greater than 1.5 mm buccally and lingually even when the residual crest is only 3.5 mm

Includes:

- Handpiece
- Handpiece support
- Footpedal
- Terminals
- Wrench WR1

**TALK TO US  
ABOUT OUR SITE  
SPECIFIC  
SOLUTION FOR  
PATIENTS WITH  
NARROW RIDGES**





## TOOTH TRANSFORMER

### What's the TT Transformer?

It is a Medical Device that facilitates the work of the Dentist and ensures maximum safety for the patient, using autologous material (the patient's own tooth).



**An innovative patented system able to transform an extracted tooth into an autologous material for bone regeneration.**

TT Transformer® reduces the crystallinity of the hydroxyapatite, the mineral of which teeth are composed, and through its unique procedure makes available the morphogenetic proteins and growth factors present in dentine.

The particulate matter obtained is, therefore, able to determine an osteoinduction mechanism capable of stimulating cellular adhesion, proliferation and differentiation during bone regeneration.



## SURGICAL TOOTHBRUSHES



The **VITIS Implant Sulcular** is specifically indicated for cleaning the gingival and/ or peri-implant sulcus.

It can also be used in implant carriers and orthodontic patients.

**\$4.95**

The **VITIS Implant Toothbrush** is a daily brush specifically designed for patients with dental implants, to provide a delicate and deep care within the whole mouth.



**\$5.45**




The **VITIS Implant Monotip** is designed specifically to access very small spaces that require maximum cleaning precision of implant-supported prosthesis.

It can also be used in implant carriers and orthodontic patients.

**\$4.95**



## TAD KIT (Temporary Anchor Devices)



- 1 - Screwdriver for 1.5mm Screw (Square type)
- 2 - Screwdriver for 2.0mm Screw (Square type)
- 3 - Screwdriver for 1.5mm Screw (AO type)
- 4 - Screwdriver for 2.0mm Screw (AO type)
- 5 - Limited Drill Bit 1.1\*6mm (common type)
- 6 - Limited Drill Bit 1.1\*8mm (common type)
- 7 - Limited Drill Bit 1.6\*6mm (common type)
- 8 - Limited Drill Bit 1.6\*8mm (common type)
- 9 - Quick Coupling Handle (AO type)
- 10 - Limited Drill Bit 1.1\*6mm (RA/Latch type)
- 11 - Limited Drill Bit 1.1\*8mm (RA/Latch type)
- 12 - Limited Drill Bit 1.6\*6mm (RA/Latch type)
- 13 - Limited Drill Bit 1.6\*8mm (RA/Latch type)
- 14 - Aluminium Box


**\$1,083.50**



## PLATE AND SCREW KIT



- 1 - Hollow Reamer 2.0mm
- 2 - Hollow Reamer 1.5mm
- 3 - Screwdriver for 2.0mm Screw (AO type)
- 4 - Screwdriver for 1.5mm Screw (AO type)
- 5 - Depth Gauge 0-30mm
- 6 - Limited DrillBit (common type)
- 7 - Limited DrillBit (common type)
- 8 - Quick Coupling Handle (AO type)
- 9 - Guider 1.1mm& 1.6mm
- 10 - Limited DrillBit (RA/Latch type)
- 11 - Limited DrillBit 1.1\*4mm (RA/Latch type)
- 12 - Screwdriver for 1.5mm Screw(RA/Latch type)
- 13 - Screwdriver for 2.0mm Screw(RA/Latch type)
- 14 - Aluminium Box



**\$1,870.00**




**GINGISAGE regular** – Gingival retraction paste made by KEREDEN laboratory is aluminium free in order to avoid possible drawbacks:

The only aluminum free retraction paste, with a composition of kaolin and ferric chloride 6%

GINGISAGE regular, presents rheological properties that are based on Dr. Lesages studies regarding the mechanical behavior of the marginal gingiva and the process of obtaining sulcular opening without causing any damage to the epithelial attachment.

The capsule has a specially designed nozzle, which has been clinically tested to deliver a thin gingival retraction cord for improved ergonomics



\$132.00

The First  
aluminum chloride free  
gingival retraction paste



15% OFF

# MONET<sup>®</sup>

laser curing light



Cures composites, luting cements, adhesives, and sealants in **1 Second**

Laser curing has major advantages over LED curing. With the Monet laser cure, you will never have to worry about incomplete bonds or soft cures. The Monet laser produces a deeper cure than LED or any other light sources during the same exposure time.

TALK TO US  
ABOUT OUR  
SPECIAL  
TRADE-IN DEAL

All prices are GST inc.

Say goodbye to well known LED related problems like:

- LED technology was a significant advancement, but LED curing lights present their own set of challenges.
- LED lights create a light pattern that can be dispersed away from the curing target when the energy from the LED light waves disperse.
- LED is a dispersed light source emitting at random times and random directions.
- LED needs a complicated optical system to collimate and homogenize the light beam

Page 22

**WAS \$3,850.00**  
**NOW \$3,272.50**



Introducing a new idea on how we see patient treatment. Hubble combines magnification, illumination, & documentation into one piece of equipment designed for more effective patient care.



# HUBBLE™

## DIGITAL VISION HEAD-SET

Hubble allows the user to control magnification from 1 up to 20x. Hubble takes you beyond the fixed magnification and enhancing limitations that come with a pair of loupes, providing a feature of continuous zoom for fully adjustable magnification all while providing true-color LED illumination. With up to 20x magnification, Hubble eliminates the need for costly and obtrusive microscopes and replaces both loupes and headlights.



Due to a large number of requests for the Hubble Digital Vision Head-Set, we have setup a Pre-Order Deposit to save your place in line. The deposit amount of \$550 is fully refundable, will be applied towards the retail price of \$17,050.00 and will save your place in line.

All prices are GST inc.

Page 23

**\$17,050.00**





**KULZER**  
MITSUI CHEMICALS GROUP



## FLEXITIME

Time flexibility paired with the highest precision.  
Every time.



Offer

Buy any 4 Flexitime Cartridge Refills  
Get 1 Flexitime Cartridge Refill FREE



## VENUS® PEARL PURE SHADES

Do you know that only 4 shades can cover up to 90% of anterior & posterior restorations? Introducing Venus Pearl Pure Shades.



Offer

CHOOSE IT. USE IT. LOVE IT. Risk-Free Trial

Buy a Venus Pearl Pure Kit with bonus Venus Pearl Pure Medium sample. Try the sample and if you don't love it (why wouldn't you), return the unopened Venus Pearl Pure Kit within 2 weeks for a full credit.



Before



After

PROSMATE



15% OFF



### ProsMate System

Invented by Prosthodontist Dr Alan Yap, the ProsMate is an effective tool to aid in the efficient and systematic insertion of fixed prostheses.

The Baton allows the simultaneous etching, washing, drying and silanating of porcelain veneers.

With the ProsMate you will be able arrange prostheses in order of insertion, and in addition, you can handle silanted porcelain veneers without touching the bonding surfaces.

Milled out of metal, the ProsMate is autoclavable and can last a lifetime. The ProsMate is a simple yet brilliant system to add to your armamentarium for daily practice.

**WAS \$770.00**  
**NOW \$654.50**

All prices are GST inc.

Page 25



# QUEEN DENTAL



## PROPHYLAXIS POWDER – RHAPSODY FLOW


Designed for air-abrasive professional teeth cleaning.

Made on the basis of finely grained sodium bicarbonate particles

Rhapsody Flow is a sodium bicarbonate based powder specifically designed for professional oral hygiene. Our product has a mild and selective action, making it an ideal choice for professional tooth cleaning.



PROPHYLAXIS POWDER



WAS \$50.84  
NOW \$37.86

All prices are GST inc.

Page 26



**VP INNOVATO HOLDINGS**

MASTERING THE EMERGENCE PROFILE



### **CERVICO PREMIUM KIT**



The Cervico Mold is an in-office tool used for the fabrication of the cervico healing abutments and impression posts.

It can also be used for the duplication of a modified healing abutment or temporary prosthesis to a duplicate impression post. It features 17 wells that can be used for all of the above purposes.

**WAS \$3,882.29**

**NOW \$3,299.94**

### **CERVICO ESSENTIAL KIT**

Essential Cervico Mold is a simplified version of the Cervico Premium Mold that works with the stock implant analogs available from the different implant companies.

**WAS \$\$\$,137.20**

**NOW \$2,666.62**



15% OFF


PROCESS FOR PRF™



## COMPLETE HEAT AND COOLER PRF SYSTEM



Receive a **FREE**  
**Venipuncture**  
**book with the**  
**purchase of any**  
**PRF System**



The Aesthetic Heat and Cooling kit includes everything you need to get started.

The heating and cooling system is meant to extend the resorption properties of platelet rich fibrin

- 100% autologous
- No chemical additives

Indications:

- Biological filler in facial esthetics (Bio-Filler)
- Microneedling with PRF

**WAS \$9,585.01**  
**NOW \$8,147.25**

All prices are GST inc.

Page 28


15% OFF

PROCESS FOR PRF™



PRF SYSTEM AESTHETICS

Receive a **FREE**  
Venipuncture book  
with the purchase  
of any PRF System



**WAS \$4,708.93**  
**NOW \$4,002.59**

The Aesthetic kit includes everything you need to get started.

The concept of PRF (Platelet Rich Fibrin) is based on the centrifugation of whole blood without anticoagulants. (J. Choukroun et al. 2001). At the end of the spin, a fibrin clot containing the majority of the platelets and white blood cells is obtained.

This fibrin clot called Platelet Riche Fibrin or PRF will release gradually and growth factors or cytokines in the site (VEGF, PDGF, TGF Beta, Thrombospondin)

The expected objective of these growth factors is to accelerate the soft tissue and bone healing.

All prices are GST inc.

Page 29



## CAVITY PHIGHTER TEST STRIPS

### Test Strips:

Created with dental professionals in mind, Cavity pHighter test strips measure the oral pH balance to determine cavity risk. Mouth bacteria can lead to increased acidity, resulting in dental erosion and painful cavities.

Cavity pHighter strips measure the pH and buffer capacity in your saliva. Simply wet the pH strip with the patient's saliva, wait 10 seconds, and measure their color reading against the key on the bottle.

The easy-to-use color key will indicate their mouth pH level and associated risk.

### Features:

- Easy-to-use strips
- Quick results
- Includes pH risk guide & acidity level comparison

50 Strips per bottle



**BUY 3 GET  
ONE FREE**



## CAVITY PHIGHTER MINTS

### Mints:

Created for at-home and professional dental use, Cavity pHighter Mints protect your patients' teeth after they've left your office. CAO's cavity-fighting mints protect your mouth against acid-producing bacteria by reducing acidity and balancing oral pH levels.

The result? Decreased risk of cavities and minty-fresh-smelling breath!

Don't know if you are at risk for cavities? Cavity pHighter Test Strips quickly and accurately assess your risk level by testing your mouth pH.

### Features:

- Peppermint flavor for fresh breath
- Cavity protection
- Take one after each meal for optimal efficacy

Cavity pHighter Mints help protect your patients' mouths against acids produced by hard-to-kill bacteria. Oral bacteria multiply and produce acids that dissolve minerals in teeth, making them more porous and susceptible to cavities.



60 pieces per bottle





## DRY MOUTH SOLUTIONS

### Dry Mouth Oral Hygiene Gel 65mL

- Provides up to 8 Hours relief from Dry Mouth
- Moisturises, hydrates, soothes & comforts
- Fluoridated formulation, replacing toothpaste

### Dry Mouth Oral Hygiene Spray 50mL

- Moisturising, fluoridated oral hygiene spray
- Convenient oral hygiene spray with 360 degree nozzle
- Freshens breath & reduces halitosis (mouth odour)

### Denture Adhesive Gel 50mL

- Secures ALL day (naturally) Seals, Cushions & Comforts
- PVM/MA Copolymer "glue" ( free) formulation
- Neutral taste, Zinc free, Easy to clean

### Moisturising Mouth Spray 60mL

- Convenient and immediate relief from Dry Mouth
- Enhanced moisture delivery with (prop.) chitin
- A 'nod' to the Australian Bush with a hint of lemon myrtle

### Moisturising Mouthwash 500mL

- Enriched with all natural moisturising Malic Acid
- Leaves the mouth with a soft, velvety mouth-feel
- pH balanced, moisturising, calming mouth rinse

### pH Neutralising Mouthwash 250mL

- Helps prevent enamel erosion
- Naturally assists to neutralise pH levels
- Refreshing flavour, safe to use ongoing



# Dr Mark's HyGenie®



## DR MARK'S HYGENIE



Dr Mark's HyGenie® is the simple solution for cleaning your dentures, night guard or Invisalign®. No more soaking in Polident, Steradent, lemon or vinegar solutions, and no more risky hand-held brushing. Simply lock your dental appliance safely inside HyGenie, add DentalFresh, twist and rinse. The high-quality rotating brushes inside do the work for you.

## DR MARK'S DENTALFRESH

Mark's DentalFresh® is an easy-to-use, fast acting and fresh scented cleansing gel. Think of it as SOAP for washing your Sports Mouthguard (or dentures, night guard and clear aligners like Invisalign®). DentalFresh cleans dental appliances by washing away biofilms and attacking all the other nasties that collect on your mouthguard in practice or competition, giving you protection against harmful germs and bacteria. And remember, you just can't be 100% with a dirty mouthguard.



## DR MARK'S HYGENIE REPLACEMENT BRUSH SET



With daily use, you should replace the brush set in your Dr Mark's HyGenie every 6 months to ensure optimal performance and good oral hygiene. There's no need for an entirely new device, simply replace your old brushes with the new brushes in just 90 seconds.



### **HYGENIE SPORT**



Gone are the days of leaving your mouthguard in the bottom of a sweaty sports bag where it loves to collect dirt, bacteria, germs and who knows what else! HyGenie Sport is the simple solution for cleaning, storing and protecting your mouthguard. After every practice or competition, simply place your mouthguard inside HyGenie Sport, twist, and let the high quality rotating brushes do the work for you. You couldn't ask for an easier way to keep your mouthguard clean.

### **HYGENIE SPORT CAGE MOUTHGUARD CARRIER**

Stop rummaging around the bottom of your bag for your mouthguard. With HyGenie Sport 'Cage' Mouthguard Carrier, you easily protect, store and find your mouthguard. The CAGE is bomb proof mouthguard protection, helping to keep it clean and safe. Clip The CAGE to the outside of your sports bag with the Carabiner Clips included, easy to see, easy to find and combat ready.



### **HYGENIE SPORT QUICKCLEAN SPRAY 100M**



HyGenie Sport QuickClean Spray is a mouthguard sanitiser that quickly and effectively protects you and your mouthguard from the nasties that collect on them during practice, comps or floating around the bottom of your sports bag. Spray QuickClean onto your mouthguard after every use, then store in the CAGE mouthguard carrier and be confident that it's hygienic, protected and combat ready (even if it has been floating around the bottom of your sports bag). It's as easy as that.



## EMBI DENTAL GIFT PACK



The perfect bundle to kick-start your patients new oral hygiene routine!

### Features

- 100% Wheat Straw
- Eco Friendly
- Good for Travel
- Gentle
- Disposable

## EMBI EDGE TONGUE CLEANER

Designed with dentists it's made from dental grade stainless steel and has a unique scoop that gives a deeper clean for your tastebuds. It has uniquely weighted forged stainless steel handles for ergonomic agility and is made by embi – known for their fine craftsmanship and quality.



## EMBI EDGE CU29 COPPER TONGUE CLEANER

It's made from the purest quality copper and has a unique scoop that gives a deeper clean for your tastebuds. Copper is known for its anti-microbial and anti-bacterial properties, which makes it a perfect and effective material for tongue cleaning as it promotes good bacteria and better oral hygiene.



# **LACALUT**<sup>®</sup>



## **LACALUT AKTIV TOOTHPASTE**

Medical toothpaste – Tightens and strengthens your gums noticeably.  
Protects against gum bleeding, periodontitis and caries.



Especially for the prevention of periodontitis and bleeding gums. Sensitivity and sensitivity to pain is reduced, the gums are strengthened and gum disease is prevented.

The LACALUT<sup>®</sup> result:

- Reduces the teeth's sensitivity to pain
- Prevents gum disease
- Hardens the tooth enamel



MANFLOSS & QUEEN OF CLEAN has been designed to improve your oral hygiene by explaining it simply, in a fun way. Developed by a dental health professional who listened (for 20 years) to the everyday persons' frustrations with flossing, and decided to do something about it.

Finally, there is a luxurious wafer thin but also soft and wide tapes designed to stop you cursing and get you smiling! For the first time ever, precise instructions are enclosed to dispense the frustrations about flossing...learn how length is everything....





A new patented technology for tongue exercise.

Myospots are small biodegradable adhesive pads that are used to exercise tongue muscles and induce nasal breathing.

Myospots are small, circular pads, made of biodegradable natural polymers that are safe for internal use. They have nice colours and flavours that give them a pleasant look and taste and are Vegan friendly. Their strong adhesive power enables them to adhere easily to the surface of the palate. Upon adhering to palate, they involuntarily induce the tongue to elevate and touch them. The repetitive induced movement of the tongue to touch and dissolve the adhered pad provides good exercise for tongue muscles to increase tongue muscles tone.

Using MYOSPOTS daily for several months can help change the tongue posture to rest against the palate.



### **How to use MYOSPOTS?**

Place a MYOSPOT with your thumb against the anterior part of the palate with gentle force for 10 seconds.

Remove your thumb and keep pressing the SPOT using your tongue until it dissolves.

oh!CARE



OH! Care Mouthwash and Gum Serum are your new everyday essential!

This hydration-boosting combo has been formulated with hyaluronic acid and Xylitol to improve oral hygiene without alcohol and harsh chemicals.

Vegan friendly  
100% Australian made  
Recyclable packaging





### TEPE COMPACT TUFT TOOTHBRUSH

TePe Compact Tuft™ is a single tuft brush with a dome-shaped brush head for precision cleaning of hard-to-reach areas, such as implants, attachments for overdentures, fixed braces or along the gum line.

### TEPE UNIVERSAL CARE TOOTHBRUSH

TePe Universal Care™ has a unique angle to facilitate access in hard-to-reach areas. Optimal for cleaning on inside surfaces and along the gum line of natural teeth, implants, lingual braces, or retainers.



### TEPE SPECIAL CARE

TePe Special Care is an ultra-soft toothbrush recommended after oral surgery and for sore or delicate tissues.

### TEPE INTERDENTAL BRUSHES – EXTRA SOFT

An interdental brush is the most efficient, user-friendly device for plaque removal between teeth. TePe Interdental brushes are designed to clean the gaps between the teeth, where a regular toothbrush is unable to reach. This Interdental brush has extra soft filaments for sensitive gums and teeth.



25% OFF

**UBrush!**  
Smile, because UBrush!



## UBRUSH! ELECTRIC INTERDENTAL BRUSH



**WAS \$130.90**  
**NOW \$98.17**

UBrush Electric Interdental Brush is the answer.

This user-friendly electric interdental toothbrush helps you clean between the teeth, both easily and efficiently.

Even the best toothbrushes can't reach the interdental areas. That's why there are alternatives: piks, floss or interdental brushes.

These usually end up being difficult to use, and patients give up. Ask your hygienist or dental therapist to show you how to use the UBrush! It really is easy.

It will convince you how simple it is to clean between the teeth.

**BUY 3 GET  
ONE FREE**



**SHEER WHITE! TAKE-HOME WHITENING**



**\$66.75**

An advanced teeth whitening formula for dental professionals, Sheer White! Take-Home patient treatment kit for up to 5 more shades of whitening.  
Better Retention, Better Revenue, Better Whitening!

Features:

- SheerFilm® Advanced technology
- 20% carbamide peroxide
- Non-slip patented formula
- Reduced sensitivity
- Can be worn while drinking beverages
- 5-day treatment



## E-SKIN SPECTROCOLORIMETER



E-Skin is a “single measure and mix” system. The simplicity of the spectrophotometer and ease of use with minimal training make e-Skin the most cost effective digital skin matching solution for today’s clinicians.




E-Skin uses a digital library of nearly 22,000 skin tones to match to patient skin for prosthetic applications.

The e-Skin instrument measures skin colour and instantly retrieves and displays on its screen a matching colorant recipe from its database.


## SPECTROMATCH COLOURS

Individual color cartridges for use with the E-Skin Spectromatch System.







**Venipuncture for the Application of Blood Concentrates – Book**  
The benefits of PRF and other blood concentrates are well known in the dental profession, but venipuncture is not a skill dentists naturally possess. In this, as in every form of treatment, providing high-quality care is essential



**The Face – Book**  
For the first time, the highly complex topographic-anatomical relationships of facial anatomy are depicted layer by layer using extremely detailed anatomical illustrations with a three-dimensional aspect.



**Piezoelectric Bone Surgery**  
This book from the inventor of Piezosurgery provides the practical information clinicians need to successfully integrate piezoelectric bone surgery into their practices to improve outcomes in oral surgery and dental implantology.



**The Art of Clinical Anaplastology Book**  
Author: Keith F. Thomas  
Print Book, English, 2006



## Terms and Conditions

### Shipping: Delivery within Australia & New Zealand

BioMeDent will provide freight-free delivery for orders over the specified amount of \$500 BioMeDent reserves the right to make changes to the freight-free policy without notice. Freight Free terms are only applicable to online orders.

### Pricing Policy:

The prices advertised the BioMeDent website are in Australian Dollars (AUD) and include 10% GST ( Goods and services tax ). Prices are subject to change without notice.

### Return and Exchanges:

Refunds not accepted. We may from time to time and at our sole discretion agree to accept return of goods for a credit however we will only do so if the goods are unopened and are in a saleable condition when they are returned to us.

Please e-mail [info@biomedent.com.au](mailto:info@biomedent.com.au)

Please note following conditions:

- The exchange is made within 7 days of purchase; you are able to provide sufficient proof of purchase and the products are in perfect condition.
- No customized/special orders to be returned
- A 25% restocking fee maybe applicable in some instances
- All returns must be sent back at your own cost via postage/courier

**Please Scan the QR Code for full Terms & Conditions**



## **Terms and Conditions**

### **This Catalogue and Website Does not provide Medical Advice:**

The content on the catalogue and website, such as text, graphics, images, information obtained from licensors, and other material contained on the website are for informational purposes only. An individual's physical state, medical status and health requirements are not considered and thus Content is not intended to be a substitute for professional medical advice, diagnosis, or treatment.

We do not recommend or endorse any specific tests, products, procedures, opinions, or other information that may be mentioned on this website. Reliance on any information provided by us, our employees, others appearing on this website at their invitation, or other visitors to this website is solely at your own risk.

### **Product Description and Colours:**

We have done our best to display as accurately as possible the colours of the products shown on this web site. However, because the colours you see will depend on your monitor, we cannot guarantee that your monitor's display of any colour will be accurate.

### **Purchases, Charges and Billing**

Certain features offered through our website may be conditioned on entering into a financial transaction with us. We reserve the right to refuse any order you place with us or cancel your order for any reason. In the event we make a change to or cancel an order, we will attempt to notify you by contacting the e-mail and/or billing address/phone number provided at the time the order was made.

### **Please Scan the QR Code for full Terms & Conditions**



2 0 2 4  
E O F Y  
EDITION



 1300 792 624

 [www.biomedent.com.au](http://www.biomedent.com.au)