




Specialised Dental & Medical Solutions

# TAD Kit (Temporary Anchor Devices)

[www.biomedent.com.au](http://www.biomedent.com.au)


## BioMeDent TAD Kit



TAD Kit

REF CODE: 04A1035

## TAD Kit (Temporary Anchor Devices) Composition



- 1 - Screwdriver for 1.5mm Screw (Square type)
- 2 - Screwdriver for 2.0mm Screw (Square type)
- 3 - Screwdriver for 1.5mm Screw (AO type)
- 4 - Screwdriver for 2.0mm Screw (AO type)
- 5 - Limited Drill Bit 1.1\*6mm (common type)
- 6 - Limited Drill Bit 1.1\*8mm (common type)
- 7 - Limited Drill Bit 1.6\*6mm (common type)
- 8 - Limited Drill Bit 1.6\*8mm (common type)
- 9 - Quick Coupling Handle (AO type)
- 10 - Limited Drill Bit 1.1\*6mm (RA/Latch type)
- 11 - Limited Drill Bit 1.1\*8mm (RA/Latch type)
- 12 - Limited Drill Bit 1.6\*6mm (RA/Latch type)
- 13 - Limited Drill Bit 1.6\*8mm (RA/Latch type)
- 14 - Aluminium Box

## Orthodontic Screws



### SPECIFICATIONS

1.5 \*8mm Ortho S T Blue

1.5 \*6mm Ortho S T Blue




### SPECIFICATIONS

2.0 \*11mm Ortho S T Yellow

2.0 \*8mm Ortho S T Yellow

2.0 \*6mm Ortho S T Yellow



 1300 792 624

 [www.biomedent.com.au](http://www.biomedent.com.au)

 [info@biomedent.com.au](mailto:info@biomedent.com.au)