



# CATALOGUE 2025-2026

[www.nichrominox.fr](http://www.nichrominox.fr)



*Specialists for  
Instrument Organisation  
in the Dental Practice*

 **NICHROMINOX**

18, rue des Frères Lumière - F-69720 SAINT-BONNET-DE-MURE - FRANCE  
Tel. +33 (0)4 78 74 04 15 - Fax +33 (0)4 78 01 61 80 - [office@nichrominox.fr](mailto:office@nichrominox.fr)

An environmentally-friendly alternative to disposable sterilization pouches.

# Aluminium Sterilization Containers




These Containers are compatible with all Nichrominox Cassettes, including Classic and Galaxy models.

Available in colours : 1 2 3 4 5 7 13




**Lots of ergonomic accessories :**

- The 4 lugs on the cover make the Containers **stackable** : for easy storage and better circulation of the steam during sterilization.
- The Containers are equipped as standard with a **security sealing system** (seals are not included). The plastic security seal is broken when opening the box, so non-sterile containers can easily be identified.
- The Containers are supplied with an **identification label**. Delivered as neutral, you can affix a self-adhesive label or have it engraved.
  - Ref. 182022 - Label Engraving (15 digits max)



**HANDLE KIT** (see p.4)  
• Ref. 182015

View the video showing how to assemble :



**New filtration system, for paper filter or teflon filter !**

The Containers can be equipped either with a **single use paper filter**, or with a **teflon filter effective for up to 2000 cycles**. Changing the filter is quick and easy ! By using Nichrominox Teflon Filters, sterility is maintained for a minimum of 3 months.



**ALU CONTAINER 20X16**

• Ref. 300350 - Sizes :  
- Ext. 22,4 x 16,8 x 5,5 cm - Int. 20,2 x 15,9 x 4 cm

For Cassette 18x14 (Classic or Galaxy model)

- Ref. 187447 - Teflon Filter - per unit
- Ref. 187468 - Paper Filters - 100 pc



**ALU CONTAINER 30X20**

• Ref. 300300 - Sizes :  
- Ext. 32,2 x 20,8 x 5,5 cm - Int. 29,9 x 19,9 x 4 cm

For Cassette 28x18 (Classic or Galaxy model)

- Ref. 187446 - Teflon Filter - per unit
- Ref. 187458 - Paper Filters - 100 pc


The filters are not included within the Containers, they are sold separately.

Help your practice increase efficiency, move to cassette management !

The Galaxy range has a high quality lattice style, making the water action more efficient during cleaning in a washer disinfectant.

# GALAXY Cassettes


Available in colours : 2 3 4 5 7 13 **New colours !**



**CASSSETTE 28 X 18 GALAXY**

Size 28,5 x 18 x 3,4 cm  
**removable cover**

- Ref. 182049 - For 20 instruments
- Ref. 182027 - For 12 instruments + clear area



**CASSSETTE 18 X 20 GALAXY**

Size 18,5 x 20,5 x 3,4 cm  
**removable cover**

- Ref. 182048 - For 10 instruments + clear area



**CASSSETTE 18 X 14 GALAXY**

Size 18,5 x 14 x 3,4 cm

- Ref. 182046 - For 10 instruments



**CASSSETTE 18 X 10 GALAXY**

Size 18,5 x 10,8 x 3,4 cm

- Ref. 182025 - For 7 instruments



**CASSSETTE 18 X 7 GALAXY**

Size 18,5 x 7,5 x 3,4 cm

- Ref. 182044 - For 5 instruments

## «AUTO GRIP» SILICONE

**New !**

**Double the capacity of your cassettes !**

Thanks to their ingenious shape, these silicone holders allow you to **double the number of instruments** stored in the Cassette : instruments can be stored in the bottom and in the lid as well.

**CASSSETTE 18 x 14 GALAXY - AUTO GRIP SILICONE**


- Ref. 183046 - For 20 instruments

**CASSSETTE 18 x 10 GALAXY - AUTO GRIP SILICONE**

- Ref. 183025 - For 14 instruments

**CASSSETTE 18 x 7 GALAXY - AUTO GRIP SILICONE**

- Ref. 183044 - For 10 instruments




# Hand Instrument Holders

# Hand Instrument Holders

## Accessories for GALAXY Cassettes...

To be attached to the cover, using the fixing tool supplied with each kit.




2 models of Clip : for Pliers and for Scissors.

View the video showing how to assemble :



### SET OF 2 CLIPS

- Ref. 182053-B..... Pliers + Scissors
- Ref. 182053-P..... 2 Pliers
- Ref. 182053-C..... 2 Scissors



To easily remove the Galaxy Cassette out of the Container !

View the video showing how to assemble :



### HANDLE KIT


- Ref. 182015

## GALAXY Baskets

Expand your Cassettes !

The GALAXY Baskets are an extension to the 18x14\* and 28x18\*\* Cassettes : in addition to the instruments stored in the Cassettes, the Baskets will hold all kind of large instruments or accessories such as tweezers, hammers, elevators, syringes, bowls, etc...


Available in colours : 2 3 4 5 7 13



### GALAXY BASKET 28x18

- Ref. 182020 - Size 28,5 x 18 x 3 cm (without end clips)

\*\*Fits with Cassettes Ref. 182049 - 182027 - 182068 - 182058 - 182087



### GALAXY BASKET 18x14

- Ref. 182010 - Size 18,5 x 14 x 3 cm (without end clips)


\*Fits with Cassettes Ref. 182046 - 182066 - 182056



## Flat GALAXY Cassettes




AVAILABLE IN COLOURS :




### CASSETTE 18 x 13

- Ref. 183062
- Size 18 x 13 x 2,5 cm
- For 10 instruments



### CASSETTE 18 x 9

- Ref. 183055
- Size 18 x 9,5 x 2,5 cm
- For 7 instruments



### CASSETTE 18 x 6

- Ref. 183063
- Size 18 x 6,5 x 2,5 cm
- For 5 instruments

Help your practice increase efficiency, move to cassette management!

## Classic Cassettes

The essentials for reprocessing !

Available in colours : 2 3 4 5 7 13

New colours !



### CASSETTE 28 X 18

Size 28 x 18 x 3,4 cm

- Ref. 182068 - For 16 instruments
- Ref. 182058 - For 20 instruments
- Ref. 182087 - For 10 instruments
- + clear area



### CASSETTE 18 x 14

Size 18,5 x 14 x 3,4 cm

- Ref. 182066 - For 8 instruments
- Ref. 182056 - For 10 instruments

Holes in the edges for a more efficient water draining



### CASSETTE 18 X 20

Size 18 x 20,5 x 3,4 cm

- Ref. 182057 - For 8 instruments
- + clear area



### CASSETTE 18 X 10

Size 18,5 x 10,8 x 3,4 cm

- Ref. 182054 - For 6 instruments



### CASSETTE 18 X 7

Size 18,5 x 7 x 3,4 cm

- Ref. 182069 - For 4 instruments
- Ref. 182079 - For 2 handpieces



## Flat Cassettes



AVAILABLE IN COLOURS :



### CASSETTE 18 x 13

- Ref. 182062 - Size 18 x 13 x 2,5 cm
- For 10 instruments



### CASSETTE 18 x 9

- Ref. 182055
- Size 18 x 9,5 x 2,5 cm
- For 7 instruments



### CASSETTE 18 x 6

- Ref. 182063
- Size 18 x 6,5 x 2,5 cm
- For 5 instruments

Best Seller!

# Hand Instrument Holders



## Easy Clips

Retractable clip with a fixed hinge, so it can never be lost.


Stainless steel and silicone. Electropolished to give a bright shiny surface finish, no sharp edges. Very open holders for a better cleaning process.

Available in colours : 1 2 3 4 5 7 8 9 13



**EASY CLIP 3**

• Ref. 182903 - Size 18 x 5,5 x 3 cm  
Silicone holder for 3 instruments



**EASY CLIP 15**

• Ref. 182912 - Size 20,7 x 18 x 3 cm  
Silicone holder for 15 instruments



**EASY CLIP 5**

• Ref. 182900 - Size 18 x 8,5 x 3 cm  
Silicone holder for 5 instruments



**EASY CLIP 20**

• Ref. 182920 - Size 25 x 18 x 3 cm  
Silicone holder for 20 instruments

### Clip-on accessories!

Attach them to the edges of your cassettes...

new



**CLIP FOR SCALING TIP**

• Ref. 205030  
Size 1,5 x 1,6 x 1,5 cm



**ENDO PLUG 2**

• Ref. 182162  
Size 3,2 x 1,6 x 5,5 cm



**CLIP PLUG'IN 2**

• Ref. 205032  
Size 3,2 x 1,6 x 5,5 cm

**EASY CLIP 10**

• Ref. 182910 - Size 18 x 12,8 x 3 cm  
Silicone holder for 10 instruments

# Hand Instrument Holders

More than a tray, simpler than a cassette!

Electropolished for a bright shiny, smooth surface finish. Colours : 1 2 3 4 5 7 8 9 13

## Easy Trays

Retractable hinged closing clip!



**EASY TRAY 18 X 5**

Size 18 x 5 x 3 cm  
• Ref. 183703 - for 3 instruments  
• Ref. 183704 - for 4 instruments

**EASY TRAY 18 X 7**

• Ref. 183700 - Size 18 x 7 x 3 cm  
For 5 instruments

**EASY TRAY 18 X 10**

• Ref. 183705 - Size 18 x 10 x 3 cm  
For 5 instruments + clear area for accessories

**EASY TRAY 18 X 14**

• Ref. 183710 - Size 18 x 14 x 3 cm  
For 10 instruments + clear area for accessories

**EASY TRAY 21 X 18**

• Ref. 183715 - Size 21 x 18 x 3 cm  
For 15 instruments + clear area for accessories

**EASY TRAY 28 X 18**

• Ref. 183720 - Size 28 x 18 x 3 cm  
For 20 instruments + clear area for accessories

## Ergo Clips

Ergonomic rounded shapes!

**ERGO CLIP 3**

• Ref. 182703  
Size 18,4 x 8 x 2,5 cm  
For 3 instruments

**ERGO CLIP 5**

• Ref. 182700  
Size 18,4 x 11 x 2,5 cm  
For 5 instruments

**ERGO CLIP 10**

• Ref. 182710  
Size 18,4 x 15 x 2,5 cm  
For 10 instruments

**ERGO CLIP 15**

• Ref. 182712  
Size 23 x 18,4 x 2,5 cm  
For 10 instruments + clear area for accessories

**ERGO CLIP 20**

• Ref. 182720  
Size 28,8 x 18,4 x 2,5 cm  
For 20 instruments  
Also available For 10 instruments + waves - Ref. 182725

MODELS 5, 10, 15 AND 20 ARE ALSO AVAILABLE WITH FULL REMOVABLE COVER :

<b>Ergo Cassette 5</b>	18,7 x 11 x 2,5 cm - Ref. 182753
<b>Ergo Cassette 10</b>	18,7 x 15 x 2,5 cm - Ref. 182750
<b>Ergo Cassette 15</b>	22 x 18,8 x 2,5 cm - Ref. 182752
<b>Ergo Cassette 20</b>	28,4 x 18,8 x 2,5 cm - Ref. 182751



# Hand Instrument Holders

# Special Orthodontic Range



## Micro Cassettes

Stainless Steel and silicone. The ultimate cleaning and sterilizing container for all kind of small instruments and accessories. Complete with 4 interchangeable colour indicator strips 2 3 4 5



**MICRO CASSETTE 1**  
• Ref. 182060  
Size 7,4 x 4,5 x 2 cm



**MICRO CASSETTE 2**  
• Ref. 182070  
Size 7,4 x 5,9 x 2,4 cm



**MICRO CASSETTE 3**  
• Ref. 182040  
Size 4 x 5,5 x 2 cm



### Models for CLAMPS



• N°1 - Ref. 182090  
Size 8,4 x 6 x 2,5 cm



• N°2 - Ref. 182092  
Size 5,7 x 4,8 x 2,7 cm



## Divided Cassettes

The first cassettes for storage and sterilization of small dental items (eg, bands, brackets or implants components...) and also for longer instruments (eg. Endodontic instruments, implantology drills.) *Makes organisation easy.* Stainless steel and silicone. Available in 3 sizes with 2 types of compartments each.



**CASSETTES 21 X 12**  
Size 21 x 13 x 2 cm  
• Ref. 182076 ..... 12 compartments  
• Ref. 182075 ..... 24 compartments



**CASSETTES 21 X 6**  
Size 21 x 7 x 2 cm  
• Ref. 182096 ..... 6 compartments  
• Ref. 182095 ..... 12 compartments



**CASSETTES 11 X 6**  
Size 11 x 7 x 2 cm  
• Ref. 182094 ..... 3 compartments  
• Ref. 182093 ..... 6 compartments



## NEW! HOLDERS FOR ORTHODONTIC SETS



## Easy Trays ORTHO

Our solution for storage and management of orthodontic sets!

Stainless Steel and silicone. Cassettes are sold empty, instruments are not included.

Colours available : 1 2 3 4 5 7 8 9 13



**EASY TRAY ORTHO 18X9**  
• Ref. 183728 - Size 18 x 8,5 x 3 cm



**EASY TRAY ORTHO 18X14**  
• Ref. 183730 - Size 18 x 14 x 3 cm



**EASY TRAY ORTHO 21X18**  
• Ref. 183732 - Size 21 x 18 x 3 cm



**EASY TRAY ORTHO 28X18**  
• Ref. 183734 - Size 28 x 18 x 3 cm

The closing clip is hinged and retractable. When open, it can be used as a rack, making it easier to grab the pliers.



### Clip-on accessories!




**CLIP PLUG'IN 2**  
• Ref. 205032  
Size 3,2 x 1,6 x 5,5 cm

### CLIP FOR SCALING TIP

• Ref. 205030  
Size 1,5 x 1,6 x 1,5 cm



Attach your Plug'In and Cassettes together to form kits!



**S. STEEL CONNECTOR**  
• Ref. 190054 - Set of 2  
Size 1 x 2,2 x 0,8 cm



## ORTHO Cassette

Maximum storage!

Equipped with a special holder allowing vertical placement of a maximum of 10 pliers + 9 compartments for long instruments + 1 area for 4 ligature pliers.



**ORTHO CASSETTE**  
• Ref. 182030  
Size 12 x 8 x 17,5 cm

# Bur Stands



## Plug'In

The silicone plugs hold any type of bur securely!

Stainless Steel and silicone. Can be rotated to offer the burs flat or at an angle. Will hold FG or RA burs whatever their diameter. Keep the instruments safe and save assistant time. Available in colours **2** **3** **4** **5**.



**PLUG'IN 2 HOLES**  
• Ref. 190203 - Size 2,5 x 2 x 3 cm



**PLUG'IN 4 HOLES**  
• Ref. 190205 - Size 4,5 x 2 x 3 cm



**PLUG'IN 6 HOLES**  
• Ref. 190208 - Size 3,5 x 2,5 x 3 cm



**PLUG'IN 8 HOLES**  
• Ref. 190210 - Size 4,5 x 2,5 x 3 cm



**PLUG'IN 10 HOLES**  
• Ref. 190213 - Size 5,2 x 2,5 x 3 cm



**PLUG'IN 12 HOLES**  
• Ref. 190211 - Size 6,3 x 2,5 x 3 cm




**PLUG'IN 16 HOLES**  
• Ref. 190212 - Size 8,5 x 2,5 x 3 cm

### Attach your Plug'In!

The Connector enables Plug'In and Cassettes to be attached together to form kits!

**S. STEEL CONNECTOR**  
• Ref. 190054 - Set of 2  
Size 1 x 2,2 x 0,8 cm



# Bur Stands

## Ergo Plug

The 2 rows are inclined to give a better view and access

**ERGO PLUG 8**  
• Ref. 190310  
Size 5 x 3,3 x 2,8 cm

**ERGO PLUG 12**  
• Ref. 190311  
Size 6,5 x 3,3 x 2,8 cm

**ERGO PLUG 16**  
• Ref. 190312  
Size 8,5 x 3,3 x 2,8 cm

## Plug'In PLUS

Removable basin for collecting used burs!

**PLUG'IN 8 PLUS**  
• Ref. 190282  
Size 8,4 x 2,5 x 3 cm

**PLUG'IN 12 PLUS**  
• Ref. 190284  
Size 10,3 x 2,5 x 3 cm

**PLUG'IN 20 PLUS**  
• Ref. 190286  
Size 11 x 3,5 x 3,4 cm

## Plug'In « For Tall Burs »

**H.4 CM - 2 HOLES**  
• Ref. 206023  
Size 2,5 x 2 x 4 cm

**H.4 CM - 4 HOLES**  
• Ref. 206025  
Size 4,5 x 2 x 4 cm

**H.4 CM - 6 HOLES**  
• Ref. 206002  
Size 7,5 x 2 x 4 cm

**H.4 CM - 8 HOLES**  
• Ref. 206030  
Size 4,5 x 2,5 x 4 cm

**H.4 CM - 12 HOLES**  
• Ref. 206012  
Size 6,3 x 2,5 x 4 cm

**H.5 CM - 6 HOLES**  
• Ref. 206000  
Size 7,5 x 2 x 5 cm

**H.5 CM - 12 HOLES**  
• Ref. 206010  
Size 6,3 x 2,5 x 5 cm

**H.6 CM - 6 HOLES**  
• Ref. 206003  
Size 7,5 x 2 x 6 cm

**H.7 CM - 6 HOLES**  
• Ref. 206004  
Size 7,5 x 2 x 7 cm

# Bur Stands

## « Maxi » Plug'In

Maximum storage of burs!



### PLUG'IN 20 HOLES

• Ref. 190214 - Size 9,7 x 2,5 x 3,4 cm



### PLUG'IN 30 HOLES


• Ref. 190230 - Size 11 x 3,5 x 3,4 cm



### PLUG'IN 40 HOLES H.5 CM

• Ref. 206040 - Size 10,9 x 8 x 5 cm

New height : 5 cm!  
Convenient to all burs including taller ones.



### PLAQUE DE PLUG 40 HOLES and cassette

• Plug'In Plate - Ref. 190215 - Size 11 x 8 x 1,5 cm  
• Plug'In Plate + Cassette - Ref. 190217

The model for standard burs is still available :  
• Ref. 190218 - Size 10,9 x 8 x 3,5 cm

NEW

# Bur Stands

## UltraLight Holders

A very flat cassette, specifically designed for used in washer disinfector.

Available in colours 2 3 4 5



### 5 HOLES

• Ref. 190370  
Size 5,7 x 5,3 x 1,9 cm



### 9 HOLES

• Ref. 190380  
Size 9,3 x 5,3 x 1,9 cm

## Multi Plug'In



Storage of 4, 6 or 10 burs + 1 or 2 scaling tips or accessories.

- Plug'In 4+1 - Size 4 x 2,5 x 4,2 cm ..... Ref. 190408
- Plug'In 6+2 - Size 6 x 2,5 x 4,2 cm ..... Ref. 190410
- Plug'In 10+2 - Size 10 x 2,5 x 4,2 cm ..... Ref. 190412

## Multi Holders

S. Steel and silicone.



Size 6 x 4,5 cm

Ideal for surgical burs, implantology drills, scaling tips, and other accessories. For storage of 6 instruments.

- 6 holes Ø 2 à 3,5 mm ..... Ref. 205010
- 6 holes Ø 3,5 à 5 mm ..... Ref. 205012

## Multi Plugs

S. Steel and silicone.



Size 6,3x2,5x5,1 cm

Notches cut into the sides act as a wrench to make it easy to remove the tips from the hand piece.

- Multi Plug 2+2 - 2 scaler tips + 2 burs/brushes ..... Ref. 190421
- Multi Plug 5 - 5 scaler tips ..... Ref. 190420

## Maxi Bur Blocs

Removable cover, anodised aluminium.



### MAXI BUR BLOC 200

• Ref. 190260  
Size 18,5 x 9,5 x 4,5 cm  
160 holes FG + 40 holes RA



### MAXI BUR BLOC 100

• Ref. 190270  
Size 12 x 8,5 x 4,5 cm  
80 holes FG + 20 holes RA



### MAXI BUR BLOC 50

• Ref. 190240  
Size 10,6 x 4,7 x 4 cm  
40 holes FG + 10 holes RA



### MAXI BUR BLOC 25

• Ref. 190265  
Size 6,4 x 5 x 3,8 cm  
20 holes FG + 5 holes RA

## Multi Clips



- 2 holes - Size 3,5 x 2,5 x 4,1 cm ..... Ref. 205025
- 6 holes - Size 6 x 3,5 x 4,8 cm ..... Ref. 205022



Ideal for storage and sterilization of scaling tips and other accessories. For instruments Ø 3,5 à 5 mm.

## Bur Blocs

Size 4 x 3 x 2,3 cm - 12 holes.  
(MX model is composed of 6 FG + 6 RA)



- FG : Ref. 190250
- RA : Ref. 190252
- MX : Ref. 190253

Bur Blocs can be stored :

• per unit, in Box 5x4x3  
Alu (Ref. 180004) / S. Steel (Ref 182004)

• 4 pieces, in Box 13x5x3  
Alu (Ref. 180007) / S. Steel (Ref 182007)

(More information about boxes on pages 32-33)



# Bur Stands



## Bur Clip

Stainless steel, self locking cover.

Complete with 4 interchangeable colour indicator strips 2 3 4 5.

Bur Clip FG is fitted with a blue indicator and Bur Clip MX with a green indicator.



### BUR CLIP 8 HOLES

Size 3,5 x 2,8 x 2,8 cm

- Bur Clip 8 holes FG ..... Ref. 191266
- Bur Clip 8 holes MX (4 FG /4 RA)..... Ref. 191265

### BUR CLIP 17 HOLES

Size 6 x 4,5 x 2,8 cm

- Bur Clip 17 holes FG ..... Ref. 191261
- Bur Clip 17 holes MX (12 FG /5 RA)..... Ref. 191260



## Bur Clip Plus

*Solid, Reliable and Strong !*

Stainless steel, self locking cover. Complete with 4 interchangeable colour indicator strips 2 3 4 5.



### BUR CLIP PLUS 11 HOLES MX (8 FG + 3 RA)

• Ref. 191368 - Size 4,4 x 2,4 x 3,2 cm

### BUR CLIP PLUS 17 HOLES MX (12 FG + 5 RA)

• Ref. 191360 - Size 6 x 4 x 3,2 cm



## Color Blocs

*A low cost full colour range !*

Anodised aluminium - 5 colours available : 1 2 3 4 5



### COLOR BLOC 11 HOLES MX (8 FG + 3 RA)

• Ref. 190368 - Size 4,4 x 2,4 x 3 cm

### COLOR BLOC 17 HOLES MX (12 FG + 5 RA)

• Ref. 190360 - Size 6 x 4 x 3 cm

# Bur Stands

## Eco Blocs

Anodised aluminium - 3 sizes and 5 colours : 1 2 3 4 5



### ECO BLOC 4 HOLES

Size 3,5 x 1,5 x 3 cm

- Eco Bloc 4 holes FG ..... Ref. 190256
- Eco Bloc 4 holes MX (2 FG /2 RA) ..... Ref. 190255



### ECO BLOC 8 HOLES

Size 3 x 2,5 x 3 cm

- Eco Bloc 8 holes FG ..... Ref. 190236
- Eco Bloc 8 holes MX (4 FG /4 RA) ..... Ref. 190235

### ECO BLOC 14 HOLES

Size 5 x 2,5 x 3 cm

- Eco Bloc 14 holes FG ..... Ref. 190246
- Eco Bloc 14 holes MX (7 FG /7 RA) ..... Ref. 190245

## Eco Bloc PLUS



*Aluminium holders equipped with a decontamination bowl !*

Using those aluminium holders, soiled instruments can be isolated and put in a decontamination solution if desired, immediately after use. Available in colours : 1 2 3 4 5



### ECO BLOC PLUS 12 FG

• Ref. 190258  
Size 6 x 2,5 x 3 cm

### ECO BLOC PLUS 12 MX

• 8 FG /4 CA - Ref. 190257  
Size 6 x 2,5 x 3 cm



# Bur Stands



## Bur Flash *Simple and Practical!*

Stainless Steel, teflon and silicone.

An open and light holder making it easy to see, clean and sterilize all of your burs.

Bur Flash can be stored vertically or horizontally. Is supplied with 4 interchangeable silicone indicators in colours 2 3 4 5. Will support any type of cleaning and sterilization.



**A new, more functional design!** The stainless steel frame makes the item more resistant to multiple reprocessings. The teflon holder is completely removable for optimal hygiene.

**New Model!**



Silicone Refills p.42-43

### BUR FLASH 5 HOLES

Size 4,3 x 3,3 x 1,3 cm

- Bur Flash 5 holes FG.....Ref. 190226
- Bur Flash 5 holes MX (3FG + 2RA).....Ref. 190225

### BUR FLASH 10 HOLES

Size 7,3 x 3,3 x 1,3 cm


- Bur Flash 10 holes FG.....Ref. 190221
- Bur Flash 10 holes MX (5FG + 5RA).....Ref. 190220



## Bur Dispensers *Excellent Value!*

Stainless Steel and teflon - Size 4,5 x 3 x 1,5 cm - 5 Holes - Can be turned in any direction.

Bur Dispensers are available in colour white 8 : FG - RA or MX (3 holes FG + 2 RA) or colours blue 2 green 4 and yellow 5 : FG only.



**New colour shades!**

### WHITE BUR DISPENSER 5 HOLES

- Bur Dispenser - 5 holes FG 8 white.....Ref. 190290-8
- Bur Dispenser - 5 holes CA 8 white.....Ref. 190292-8
- Bur Dispenser - 5 holes MX 8 white.....Ref. 190293-8

### COLOURED BUR DISPENSER 5 HOLES

- Bur Dispenser - 5 holes FG Blue 2.....Ref. 190290-2
- Bur Dispenser - 5 holes FG Green 4.....Ref. 190290-4
- Bur Dispenser - 5 holes FG Yellow 5.....Ref. 190290-5

# Endodontic Stands



## Endo Flash *So Simple, So Successful!*

Stainless steel and teflon. An open and light holder making it easy to see the instruments. The top locks into position to hold the instruments safely during processing. Silicone cursor for counting the number of times processed.



**Special model for long files!**

Fits all endo files, including files with length 31 mm.

**NEW**



### ENDO FLASH 4

• Ref. 182104  
Size 4,2 x 1,3 x 5 cm

### ENDO FLASH 8

• Ref. 182108  
Size 7,3 x 1,3 x 5 cm



### ENDO FLASH 12

• Ref. 182112  
Size 9 x 1,3 x 5 cm



### ENDO FLASH 8 XL

• Ref. 182118  
Size 7,3 x 1,3 x 5 cm



## Endo Pro

**New Design!** More perforations, ideal for cleaning in washer disinfectant.

1 endometer, self locking cover system, silicone cursor on the cover for counting the number of times processed.

Silicone Refills p.42-43



**Model special for long files!**

### ENDO PRO 4 XL

Fits all endo files, including with length 31 mm!  
• Ref. 182340 - Size 4,8 x 6 x 2 cm

### ENDO PRO 8

• Ref. 182350  
Size 7,7 x 5,8 x 2 cm

### ENDO PRO 12


• Ref. 182360  
Size 9,2 x 5,8 x 2 cm



## Endo Dispenser

**Can be rotated in any direction!** Stainless steel and teflon. Colours 2 4 5 8.

Silicone cursor for counting the number of times processed.



**New colour shades!**

### ENDO DISPENSER 8 Holes

• Ref. 182210 - Size 6,8 x 4,6 x 1,7 cm

### ENDO DISPENSER 12 Holes

• Ref. 182260 - Size 9 x 5 x 1,7 cm

# Endodontic Stands



## Endo Cassettes

**Ergonomic and easy to use!** Stainless steel perforated box with teflon insert, 1 endometer, self locking cover system. Silicone cursor for counting the number of times processed.



**ENDO CASSETTE 1**  
16 holes for endo files

• Ref. 182151 - Size 3,5 x 3,9 x 5 cm



**ENDO CASSETTE 2**  
32 holes for endo files

• Ref. 182150 - Size 6,5 x 4,3 x 5 cm



**ENDO CASSETTE 3**

60 holes + 3 plastic cups + 1 free area  
• Ref. 182152 - Size 11,5 x 5,8 x 5 cm



## Endo Zippo <sup>new</sup>

**A new collection of ergonomic holders!** Compact in design, these holders can be stored in our Classic and Galaxy cassettes. More perforations, ideal for cleaning in washer disinfectant. Stainless Steel and Teflon with a bright shiny surface finish.



**ENDO ZIPPO 6**  
• Ref. 182153  
Size 2,3 x 3,5 x 5 cm



**ENDO ZIPPO 8**  
• Ref. 182154  
Size 2,3 x 4,5 x 5 cm



**ENDO ZIPPO 12**  
• Ref. 182155  
Size 2,3 x 5,5 x 5 cm



**ENDO ZIPPO 16**  
Ref. 182156  
Size 2,3 x 6,6 x 5 cm



# Endodontic Stands



## Endo Star The new Concept "All in one"

**NEW DESIGN, more perforations for a better cleaning process!**

*New Design!*



**ENDO STAR 5**

• Ref. 182160 - Size 5 x 2 x 5 cm  
5 endo instruments + 2 burs

**ENDO STAR 8**

• Ref. 182170 - Size 7,3 x 2,5 x 5 cm  
8 endo instruments + 3 burs

**ENDO STAR 12**


• Ref. 182180 - Size 9,4 x 2,5 x 5 cm  
12 endo instruments + 3 burs

Stainless steel and silicone. The Endo Star holds your standard or new Nickel Titanium instrumentation, as well as FG or RA burs. 1 ruler with stop to measure your instruments length, silicone cursor for counting the number of times processed, 1 non-slip silicone pad on the base (that can be removed during cleaning). **Bright shiny surface finish.**

Silicone Refills 0-42-43



## Endo Ultralight The top of the range!



**ENDO ULTRALIGHT 4**

• Ref. 182320 - Size 5 x 6 x 2 cm

**ENDO ULTRALIGHT 8**

• Ref. 182380 - Size 8,4 x 6 x 2 cm

**ENDO ULTRALIGHT 12**

• Ref. 182390 - Size 11,9 x 6 x 2 cm

Stainless steel and silicone. **Bright shiny surface finish.** Open holders, ideal for cleaning in washer disinfectant. When closed, they give a high level of protection for endodontic instruments. When opening the cassette you lift the internal holder to present the instruments. Silicone cursor on the cover for counting the number of times processed.

Stainless steel rings with a **graduated rotating ruler and stop attached**, designed to measure the working length of endodontic instruments. Suitable for the right or left hand.

### ENDO RING

• Ref. 182110 - Size 3,9 x 1,2 cm  
Set of 2



### ENDO RING PLUS

• Ref. 182115 - Size 4 x 3,2 cm  
Per unit



Foam Refill  
Set of 40  
Ref. 182116



# Endodontic Stands



## ENDO MICRO PLUS

Size 6,5 x 5 x 5 cm - Anodised aluminium, colours 1 2 3 4  
 • Ref. 180183 - 24 holes for endodontic instruments, 2 autoclavable plastic cups + 1 mini clean gripp

## ENDO MICRO PULPA

Size 6,5 x 5 x 5 cm - Anodised aluminium, colours 1 2 3 4  
 • Ref. 180180 - 30 holes for endodontic instruments



## MINI ENDO

Size 14 x 5,5 x 5 cm - Anodised aluminium, colours 1 2 3 4  
 • Ref. 180193  
 60 holes for endodontic instruments, 4 autoclavable plastic cups

## MINI ENDO with Clean Gripp

Size 14 x 5,5 x 5 cm - Anodised aluminium, colours 1 2 3 4  
 • Ref. 180192  
 36 holes for endodontic instruments, 3 autoclavable plastic cups, 1 clean gripp



## ENDO TOTAL

• Ref. 180196 - Colours 1 2  
 30 holes for endodontic instrum., 6 autoclavable plastic cups, 5 holes for burs (3 FG + 2 RA), 1 endometer, 1 clean gripp, 1 area for cotton pellets

## BASIC ENDO TOTAL

• Ref. 180197  
 Colour 1 only  
 96 holes for endodontic instruments

Size 12 x 9 x 6 cm - Anodised aluminium.

## ENDO ECO

• Ref. 180195 - Size 22 x 10,5 x 5 cm  
 Anodised aluminium, colours 1 2 3 4 5  
 180 holes for endodontic instruments, 8 autoclavable plastic cups.

*The ultimate holder !*

# Endodontic Stands

## MINI GRIPP

Gather your soiled instruments for pre-decontamination. Made of a cup for disinfecting liquid and a removable gauze.



## MINI GRIPP ALU

• Ref. 180250  
 Ø 4 - H. 3,5 cm - Anodised aluminium colours 1 2 3 4 5 7 13

## MINI GRIPP S. STEEL+TEFLON

• Ref. 182400  
 Ø 4 - H. 3 cm



*The Classic Style!*

## PULPA / GUTTA BOX

Size Ø 4,8 cm - H. 5,5 cm - Anodised aluminium, colours 1 2 3 4  
 • Pulpa Box - Ref. 180150 - 18 holes for endodontic instr.  
 • Gutta Box - Ref. 180160 - 6 autoclavable plastic cups

*Optimize and facilitate your root canal filling !*

Gutta Bath helps you ensure a safe and effective disinfection of your gutta-percha cones after calibration and before insertion into the root canal.



## GUTTA BATH

• Ref. 180135 - Size 7 x 5 x 1,2 cm



## COMPARTMENT BOXES

Size 22 x 10,5 x 2 cm  
 Anodised aluminium. Colour 1 only.  
 • Compartment Box 12 removable comp. - Ref. 180100  
 • Compartment Box 16 removable comp. - Ref. 180101



## ENDOMODULE

Size 13,5 x 5,5 x 1,8 cm  
 • Alu Endomodule - Ref. 180130  
 Anodised aluminium colours 1 2 3 4.  
 • S. Steel Endomodule - Ref. 182130

# Implantology

## Modulis Cassettes



### MODULIS CASSETTE 4

• Ref. 500914 - Size 10 x 7,5 x 2,3 cm  
1 ratchet holder module + 1 compartment module + 4 plug module.

For a customised cassette, or bigger size, please send a request to : office@nichrominox.fr



### MODULIS CASSETTE 8

• Ref. 500918 - Size 20 x 7,5 x 2,3 cm  
1 ratchet holder module + 1 compartment module + 12 tubes module + 4 plug module.

# Implantology

## Implant Cassettes

Ergonomic and compact !



• Ref. 500301 - Size 10 x 8 x 2 cm  
For 8 instruments and 1 ratchet (maximum length for ratchet : 9,5 cm).

• Ref. 500302 - Size 16 x 8 x 2 cm  
For 15 instruments and 1 ratchet (maximum length for ratchet : 10 cm) + 1 basin

## Implant Plug'In



• N°1 - Ref. 500401  
Size 10 x 4,5 x 5 cm  
6 silicone tubes - 1 ratchet holder (maximum length for ratchet 9 cm)

Silicone Refills p.42-43

• N°2 - Ref. 500501 - Size 15 x 4,5 x 5 cm  
13 silicone tubes - 1 ratchet holder (maximum length for ratchet 10 cm) - 1 stainless steel removable compartment with holes for an efficient cleaning process.

## Implantology Kits



### IMPLANTOLOGY KIT N°1

• Ref. 500571 - Size 10 x 9 x 5 cm  
14 silicone tubes - 1 ratchet holder (maximum length for ratchet 9 cm) - 1 stainless steel removable compar

### IMPLANTOLOGY KIT N°2

• Ref. 500631 - Size 15 x 8,5 x 5 cm  
24 silicone tubes - 1 ratchet holder (maximum length for ratchet 9,5 cm) - 1 stainless steel removable compartment.

## Easy Gauge

Place your implants with high accuracy



• Ref. 500550 - Size 12,7 x 6,6 x 0,7 cm  
Stainless Steel rulers (1 distal + 1 mesial). Perfect tool for very accurate first hole.

## Universal Implant Guide



• Ref. 500540  
Set of 2 (upper + lower jaw)  
The perfect tool for a high precision initial hole to be drilled.  
The «forks» can be slightly bent by hand and adapted to the shape of the dental arch.



Model T  
• Ref. 500850-T

### IMPLANTOLOGY KIT N°5

Size 8 x 5,5 x 4 cm



Model P  
• Ref. 500850-P

Silicone Refills p.42-43



new



Best Seller!

# Intra-Oral Photography



## Glass Mirrors

**High Clarity!** Our mirrors give a sharp image and are reversible. They are front surface coated so you do not get a double image. Thickness 3 mm. Please refer to instructions for use.



## Mirrors with Handle

**High Clarity!** These mirrors are made of Float Glass, with a thickness of 3 mm, and come with a malleable stainless steel handle attached with a bracket bonded to the rear of the mirror.



• Ref. 190792 - Size of the mirror 10 x 6,9 cm  
Palatal Mirror AE - Occlusal views


• Ref. 190791 - Size of the mirror 11,7 x 4,5 cm  
Lateral Mirror T1 - Buccal views

# Intra-Oral Photography



## Photo Kits

Stainless steel and silicone. Thanks to the silicone holders, mirrors and accessories are kept safe during cleaning and storage.




**Caution:** Plastic Retractors are autoclavable up to 121°C

### PHOTO KIT N°2

Size 29 x 18,3 x 3,7 cm  
 Ref. 182735 - includes 1 Large Palatal + 1 Small Palatal + 1 Lateral C + 1 Lingual P + 2 Plastic Retractors  
 Empty Cassette, Ref. 182730

Self locking and removable clip!



### PHOTO KIT N°1

Size 19,8 x 16,5 x 3,7 cm  
 Ref. 182734 - includes 1 Large Palatal + 1 Lateral C + 1 Lingual P  
 Empty Cassette, Ref. 182731

## Photo Cassettes


### PHOTO CASSETTE N°1

Size 18 x 15 x 4,7 cm  
 Ref. 182739 - includes 1 Large Palatal + 1 Small Palatal + 1 Lateral T1 + 1 blue Silicone Handle  
 Empty Cassette, Ref. 182737



Available in the same size:  
**PHOTO CASSETTE N°2**  
 Ref. 182769 (includes 1 Large Palatal + 1 Small Palatal + 1 Lateral C + 1 Lingual P)  
 Empty Cassette, Ref. 182767

New design: New twist lock system!



### PHOTO CASSETTE N°3

Size 22,5 x 12 x 4,7 cm  
 Ref. 182779 - includes 1 Universal L + 1 Double Ended Palatal Contraptor + 1 Double Ended Upper Contraptor  
 Empty Cassette, Ref. 182777

Silicone Refills on request



## Solo Cassettes

Protect your Photographic Mirrors from scratches during any cleaning and sterilization process.

Mirrors sold separately (see p.24)


Available in Colours: 2 3 4 5



**SOLO CASSETTE N°1**  
 • Ref. 182746 - Size 12 x 8,5 x 1,3 cm  
 For 1 Palatal or 1 Lateral T1



**SOLO CASSETTE N°2**  
 • Ref. 182747 - Size 16,2 x 7 x 1,3 cm  
 For 1 Lateral C or S or 1 Lingual



**SOLO CASSETTE N°3**  
 • Ref. 182748 - Size 18,6 x 9 x 1,3 cm  
 For 1 Universal or 1 Palatal XL



**SOLO CASSETTE N°4**  
 • Ref. 182749 - Size 12 x 8,5 x 1,9 cm  
 Special for AE or T1 mirrors

### Duo Cassette

• Ref. 183748 - Size 18,4 x 9,4 x 2,2 cm  
 For the storage of 2 mirrors (including XL models\*)  
 in an extra flat cassette!

\*except for mirror with brackets



# Intra-Oral Photography



**PALATAL CONTRASTORS**


## Contrastors

Made of Aluminium or Silicone, our Contrastors give excellent contrast to your images and eliminate those elements not required in the picture. Upper Contrastors are perfect for profiles. They allow you to hide the upper lip and the nose.



**UPPER CONTRASTORS**

### ALUMINIUM CONTRASTORS



**ALU PALATAL CONTRASTORS**

• Ref. 190780 - **Set of 2**  
(1 large + 1 small)



**ALU UPPER CONTRASTORS**

• Ref. 190776 - **Set of 2**  
(1 large + 1 small)

### SILICONE CONTRASTORS



**LARGE PALATAL CONTRASTOR**

• Ref. 190781 - **Per unit**

**SMALL PALATAL CONTRASTOR**

• Ref. 190782 - **Per unit**



**DOUBLE ENDED PALATAL CONTRASTOR**

• Ref. 190783 - **Per unit**



**UPPER LARGE CONTRASTOR**

• Ref. 190777 - **Per unit**

**UPPER SMALL CONTRASTOR**

• Ref. 190778 - **Per unit**



**DOUBLE ENDED UPPER CONTRASTOR**

• Ref. 190774 - **Per unit**

# Intra-Oral Photography



### SILICONE MIRROR HANDLE


• Ref. 190796 - Stainless Steel blade encased in silicone - Unbreakable, can be bent with your hands into any position. Available colour Blue 2 - Green 4 - Pink 11. Holds any mirror firmly in place without any damage thus avoiding finger prints and smudges on the mirrors.



### BLACK CHEEK RETRACTOR

• Ref. 190779 - Set of 2

Double ended retractor (small + large).



### METALLIC BLADES

• Ref. 190711 - N°1 Small


• Ref. 190712 - N°2 Large

Stainless Steel, Mirror polished effect.

Their specific shape gives very good results. For occlusal photos use the large end. For lingual photos use the small end.

## Plastic Retractors

Autoclavable at 121°C, please refer to instructions for use.



### PLASTIC DOUBLE CHEEK RETRACTOR

• Ref. 072750 - Small • Ref. 072752 - Medium • Ref. 072751 - Large

Self retaining retractor designed with tabs to assure a full view.

### PLASTIC RETRACTORS

• Ref. 072700 - Set of 2

Double ended retractor (small + large).

# Sterilization Trays



## Stainless Steel Trays



### STAINLESS STEEL TRAYS 28X18

- Stainless Steel
- Tray 28 x 18 x 1,9 cm ..... Ref. 182450
  - Cover 28 x 18 x 2,5 cm ..... Ref. 182451
  - Instrument holder (11) ..... Ref. 182452



### PERFORATED S. STEEL TRAYS 28X18

- Stainless Steel
- Tray 28 x 18 x 1,9 cm ..... Ref. 185450
  - Cover 28 x 18 x 2,5 cm ..... Ref. 185451
  - Instrument holder (11) ..... Ref. 182452



### STAINLESS STEEL TRAYS 18X14

- Stainless Steel
- Tray 18 x 14 x 1,9 cm ..... Ref. 182460
  - Cover 18 x 14 x 2,5 cm ..... Ref. 182461
  - Instrument holder (7) ..... Ref. 182462



### PERFORATED S. STEEL TRAYS 18X14

- Stainless Steel
- Tray 18 x 14 x 1,9 cm ..... Ref. 185460
  - Cover 18 x 14 x 2,5 cm ..... Ref. 185461
  - Instrument holder (7) ..... Ref. 182462

## Stainless Steel Racks



### STAINLESS STEEL RACK 28 X 18

- Size 30 x 19 x 27 cm ..... Ref. 182455 (empty)
- For 6 trays 28 x 18 with covers  
Can be wall mounted



### STAINLESS STEEL RACK 18 X 14

- Size 19 x 16 x 27 cm ..... Ref. 182465 (empty)
- For 6 trays 18 x 14 with covers  
Can be wall mounted

# Sterilization Trays




## Aluminium Trays

6 Colours 1 2 3 4 5 6 - Cover Colour 1 only



### ALUMINIUM TRAYS 28X18

- Anodised Aluminium
- Tray 28 x 18 x 1,9 cm ..... Ref. 181450
  - Cover 28 x 18 x 2,5 cm ..... Ref. 181451
  - Instrument Holder (16) ..... Ref. 181452




### PERFORATED ALUMINIUM TRAYS 28X18

- Anodised Aluminium
- Tray 28 x 18 x 1,9 cm ..... Ref. 184450
  - Cover 28 x 18 x 2,5 cm ..... Ref. 184451
  - Instrument Holder (16) ..... Ref. 181452



### ALUMINIUM TRAYS 18X14

- Anodised Aluminium
- Tray 18 x 14 x 1,9 cm ..... Ref. 181460
  - Cover 18 x 14 x 2,5 cm ..... Ref. 181461
  - Instrument Holder (10) ..... Ref. 181462



### PERFORATED ALUMINIUM TRAYS 18X14

- Anodised Aluminium
- Tray 18 x 14 x 1,9 cm ..... Ref. 184460
  - Cover 18 x 14 x 2,5 cm ..... Ref. 184461
  - Instrument Holder (10) ..... Ref. 181462

## Aluminium Racks



### ALUMINIUM RACKS 28 X 18

- ① Long side loading :
    - Size 30 x 19 x 27 cm ..... Ref. 181455 (empty)
  - ② Short side loading :
    - Size 19 x 30 x 27 cm ..... Ref. 181456 (empty)
- For 6 trays 28 x 18 with covers  
Can be wall mounted




### ALUMINIUM RACK 18 X 14

- Size 19 x 16 x 27 cm ..... Ref. 181465 (empty)
- For 6 trays 18 x 14 with covers  
Can be wall mounted

# Sterilization Trays


## Trays 20x10 and 20x15



### ALUMINIUM TRAYS

Anodised aluminium - Colours 1 2 3 4 5 6

- Size : 20 x15 x 2 cm ..... Ref. 181480
- Size : 20 x10 x 2 cm ..... Ref. 181470



### ALUMINIUM RACKS

- For 6 trays 20 x 15 :
  - Size 16 x 20 x 18,5 cm ..... Ref. 181485 (empty)
- For 6 trays 20 x 10 :
  - Size 10,3 x 20 x 18,5 cm ..... Ref. 181475 (empty)



### STAINLESS STEEL TRAYS

Stainless steel - High quality  
Perfect for your exam sets.

- Size : 20 x15 x 2 cm ..... Ref. 182480
- Size : 20 x10 x 2 cm ..... Ref. 182470



### STAINLESS STEEL RACKS

- For 6 trays 20 x 15
    - Size : 21 x 15,8 x 18,2 cm ..... Ref. 182485 (empty)
  - For 6 trays 20 x 10
    - Size : 20,6 x 10,2 x 18,2 cm ..... Ref. 182475 (empty)
- Can be wall mounted

## Aluminium Tubs



Ultra light, for use as a portable operatory drawer. One Tub is recommended for each different type of procedure. Can be fitted with a divider and arranged in an aluminium rack.

- Aluminium Tub - Colours 1 2 3 4 - Size 33,5 x 23,5 x 5,2 cm ..... Ref. 181510
  - Tub Divider ..... Ref. 181520
  - Empty Rack (for 4 Tubs) - Size 34 x 24 x 41 cm ..... Ref. 181525
- Can be wall mounted

# Sterilization Trays

## Steri Trays



## Containers

Sterilization boxes fitted with a paper filter (single use) - Anodised aluminium colour 1. Multiperforated covers with lockable system. Silicone heat resistant lockable seal (230°).



### STERI TRAY 18X14

- Ref. 187460
- Ext. size 19 x 14,5 x 4,6 cm - Int. size 18 x 13,5 x 2,8 cm
- Ref. 187468 - Refill of 100 Paper Filters

Delivered with an alu. instruments holder



### STERI TRAY 28X18

- Ref. 187450
- Ext size 29,5 x 19 x 4,6 cm - Int. size 28 x 18 x 2,8 cm
- Ref. 187458 - Refill of 100 Paper Filters



### CONTAINER 25X10X5

- Ref. 300100
- Ext. size 26 x 11 x 7,5 cm - Int. size 25 x 10 x 5,1 cm
- Ref. 300150 - Refill of 100 Paper Filters



### CONTAINER 25X15X6

- Ref. 300200
- Ext. size 26 x 16 x 8,5 cm - Int. size 25 x 15 x 5,8 cm
- Ref. 300250 - Refill of 100 Paper Filters



«High Quality» Aluminium

## Steri Boxes



Highly efficient and ecological sterilization boxes! Fully airtight, they are equipped with a teflon filter effective up to 2000 cycles, a holder for an identification label and a seal system. Anodised aluminium - Colour 1 only. For suitable cassettes, please contact us.



### STERI BOX 18X14

- Ref. 187445
- Ext. size 19 x 15,6 x 5,5 cm - Int. size 18 x 13,5 x 3,5 cm
- Ref. 187447 - Refill of Teflon Filter - per unit

Aluminium label + marking in addition & on demand :  
• Ref. 187448

Filtration system with 2 silicone gaskets and 1 Teflon filter.



### STERI BOX 28X18

- Ref. 187440
- Ext. size 30 x 19 x 6 cm - Int. size 28 x 18 x 4 cm
- Ref. 187446 - Refill of Teflon Filter - per unit



### STERI BOX 32X21

- Ref. 187470
- Ext. size 35 x 21,7 x 8 cm - Int. size 32 x 21 x 7 cm
- Ref. 187446 - Refill of Teflon Filter - per unit

NEW



# Sterilization Boxes

## Aluminium Boxes



Sizes (cm)*	Grey Colour
17 x 7 x 3	180014-1
18 x 9 x 3	180016-1
18 x 9 x 4	180017-1
21 x 10 x 3	180019-1
21 x 10 x 4	180020-1
21 x 10 x 5	180021-1



Sizes (cm)*	Blue Colour	Red Colour	Green Colour	Yellow Colour
17 x 7 x 3	180014-2	180014-3	180014-4	180014-5
18 x 9 x 3	180016-2	180016-3	180016-4	180016-5
21 x 10 x 3	180019-2	180019-3	180019-4	180019-5
21 x 10 x 5	180021-2	180021-3	180021-4	180021-5

## Aluminium Boxes Small Sizes



Suitable for Bur Blocs (p.13)

### BOXES 5 X 4 X 3 CM\*

Colour	Reference
Grey Colour	180004-1
Blue Colour	180004-2
Red Colour	180004-3
Green Colour	180004-4
Yellow Colour	180004-5



### BOXES 10 X 3 X 2 CM\*

Colour	Reference
Grey Colour	180005-1
Blue Colour	180005-2
Red Colour	180005-3
Green Colour	180005-4



### BOXES 13 X 5 X 3 CM\*

Colour	Reference
Grey Colour	180007-1
Blue Colour	180007-2
Red Colour	180007-3
Green Colour	180007-4

Suitable for Bur Blocs (p.13)

## Aluminium Perforated Boxes

SPECIAL AUTOCLAVE

Colour 1 only



Sizes (cm)*	Reference
17 x 7 x 3	184014
18 x 9 x 3	184016
21 x 10 x 3	184019
21 x 10 x 5	184021



## Aluminium Perforated Boxes

High quality anodised aluminium, 15/10 thickness.  
Colour 1 only.



Sizes (cm)*	Reference
25 x 10 x 5	180022-1
25 x 15 x 6	180024-1

## Stainless Steel Boxes

Sizes (cm)*	Reference
5 x 4 x 3	182004
13 x 5 x 3	182007

Suitable for Bur Blocs (p.13)



Sizes (cm)*	Reference
17 x 7 x 3	182014
18 x 9 x 3	182016
21 x 10 x 3	182019
21 x 10 x 5	182021



## S. Steel Boxes Perforated


SPECIAL AUTOCLAVE



Sizes (cm)*	Reference
17 x 7 x 3	185014
18 x 9 x 3	185016
21 x 10 x 3	185019
21 x 10 x 5	185021

\* Above mentioned sizes correspond to our trade descriptions. Precise dimensions are available on request.

\* Above mentioned sizes correspond to our trade descriptions. Precise dimensions are available on request.



## LOWER THIRD MOLAR RETRACTOR

- Ref. 074000 - The pair (right + left)

Specifically designed for removal of lower wisdom teeth. Claws for support. Improve visibility and access to the third molar region, by repelling the labial commissure both horizontally and vertically. Very comfortable handling.



## FLAP RETRACTOR

- Ref. 071108 - Per unit, width 8 mm
- Ref. 071112 - Per unit, width 12 mm

For removal of lower or upper wisdom teeth. Claws for support - Can be used for removing the upper or lower wisdom teeth, also for all types of oral surgery.



## DAUTREY TYPE RETRACTOR

- Ref. 072200 - Per unit, short (H. 6 et 4,5 cm)
- Ref. 072201 - Per unit, long (H. 9 et 7,5 cm)

Those retractors are curved, with very smooth edges, so they don't hurt commissures during surgeries on wisdom teeth.



## UPPER THIRD MOLAR RETRACTOR

- Ref. 071700 - The pair (right + left)

Specifically designed for removal of upper wisdom teeth. Claws for support. Improve visibility and access to the third molar region. The important factor is that the flap is pushed instead of being pulled.



## GINESTET TYPE RETRACTOR

- Ref. 072400 - The pair (right + left)

Keep the flap retracted during surgeries on impacted wisdom teeth.



## FARABEU TYPE RETRACTOR

- Ref. 072100 - Set of 2

This retractor is multi-purpose, it can be used for all types of oral surgery.



## MOUTH PROPS

- Ref. 070300 - Set of 2 (1 Small + 1 Large in 1 aluminium box)
- 2 sizes : 15 and 25 mm

Non reflective stainless steel. The patient bites on a soft silicone surface to keep the mouth open. Available in 2 sizes : 15 and 25 mm.




## MALLEABLE RETRACTORS

- Ref. 071000 - Set of 3

3 widths : 10, 15 and 20 mm

Can be easily bent with fingers or pliers to the shape you require.



## RETRACTOR FOR LINGUAL APPROACH ①


- Ref. 070800 - Set of 2

This retractor is designed for the lingual approach of the lower wisdom teeth.

## TONGUE RESTRAINER ②

- Ref. 072000 - Set of 2

This retractor allows you to gently keep the tongue down.



## REVERSIBLE RETRACTOR

- Ref. 072600 - Set of 2

One large + one small end.



## «VISIO-N» RETRACTOR

- Ref. 075100 - Per unit

Convenient to all types of daily procedures, creates the perfect gap without discomfort for the patient, between the cheek and jaw, from premolars to the distal surface of the rear molars.

Thanks to the ergonomic shape of the handle, it is easy-tiring to keep in place even for long procedures.



## UNIVERSAL RETRACTOR

- Ref. 075000 - Per unit

For all types of oral surgery.



## "S" SHAPE RETRACTOR

- Ref. 072500 - Per unit

Curved edges, retract both a flap and cheek without trauma. 2 ends of different widths : 12 and 20 mm.

# Various Accessories



## TONGUE PROTECTOR

- Ref. 190514 - The pair (right + left)


**Without mirror.** Tongue and cheek protection. Easy access to the whole arch. Allows you to use a handpiece at any angle safely.



## MIRROR TONGUE PROTECTOR

- Ref. 190515 - Per unit

**With a high quality mirror.** Tongue and cheek protection. Easy access to the whole arch. Allows you to use a handpiece at any angle safely.



## EXAMINATION MIRRORS

- Ref. 193510 - Set of 6

Anodised aluminium handles - Mirror Ø 2 cm. Sold as set of 6 with assorted coloured handles : 1 2 3 4 5 6.



## PARALLELISM MIRROR

- Ref. 190750 - Sold with handle

High Quality front surface mirror with laser marked lines on the surface. Parallel lines are spaced 3 mm, making it easy to check abutments parallelism. Size of the mirror 4,5 x 2,9 cm



## MIRRORS Ø4 & Ø5 CM

- Ref. 190772 - Mirror Ø 4 cm
- Ref. 190773 - Mirror Ø 5 cm

For Intra-Oral Photography or for Motivation. Sold with handle.



## SILICONE CODE RING

- Ref. 190202 - Box of 50 pieces

6 colours 2 3 4 5 7 8 - For identifying and organising instruments. Autoclavable 134°C.

# Various Accessories



## BOWLS

- Small Size - Ref. 182041  
Ø 5,3 cm - H. 2,5 cm - capacity 1,8 cl
  - Large Size - Ref. 182042  
Ø 6,1 cm - H. 2,7 cm - capacity 3,4 cl
- Stainless Steel.  
Small bowls for rinse, dilute and mix.



## GERMICIDE BOXES

- Ref. 182050 - Size Ø 7,5 cm - H. 3,5 cm

Stainless Steel. Composed of 3 parts : 1 bowl, 1 perforated insert and 1 opening cover with slot. Prevent airborne contamination.



## COTTON HOLDERS

- Size Ø 5,5 cm - H. 7,5 cm
- Ref. 190300-1 - Alu colour 1
- Ref. 190300-2 - Alu colour 2
- Ref. 190300-4 - Alu, colour 4
- Ref. 192300 - Stainless Steel



## WASTE BINS

- Size Ø 7 cm - H. 8 cm
  - Ref. 190400-1 Alu colour 1
  - Ref. 190400-2 Alu colour 2
  - Ref. 190400-4 Alu colour 4
  - Ref. 192400 - Stainless Steel
- Inside bowl



## BITE REGISTRATION TRAYS RIGID

- Ref. 190150 - Set of 2
- Stainless steel. To compress the impression material.

## WIRE BITE TRAYS

- Ref. 190110 - Box of 10
- Stainless steel. Makes the handling of the impression material easier.



## DENTAL FOX PLATE

- Ref. 119940 - Size 20 x 10,5 cm
- Anodized aluminium.  
Establish the occlusal plane to see its imperfections.



## STAINLESS STEEL GRIDS

- Ref. PM190650
- For resin reinforcement.\*
- Stainless Steel, 1 meter in 10 tubes of 10 cm.



## GROOVED PINS

- Nickel Chrome 18.10.
- Must be used as a temporary dental post.\*
- 6 diameters available :
- Ref. PM119709 to PM119714 - set of 12 Pins (per size)



## C4 SOLDERING


- Ref. 140100
- High flow, low bonding temperature (550°C), high mechanical strength, 35 to 45 kg/mm<sup>2</sup>.
- Delivered as a 0,35 mm filament, packed as coils of 5 meters (about 5 grams), with flux.

## S. STEEL STAPLES

For orthodontic, prosthetic equipment and for splinting.\*



from 1,75 to 4 mm



The S. Steel Staples are available in 10 sizes (mm) : 1,75 - 2,00 - 2,25 - 2,50 - 2,75 - 3,00 - 3,25 - 3,50 - 3,75 - 4,00



- Ref. PM119608 to PM119680 - Bag of 24 staples (per size)
- Ref. PM119699 - the Set (box + 6 staples each size)

## SPIREC DRILLS

*Back in Stock!*

- Ref. PN119851-M - Spirec Drills - Set of 5 Class IIa - Manufacturer : Hager&Meisinger GmbH
- Ref. PN119851-S - «Spirec type» Drills - Set of 6 Class IIa - Manufacturer : Stoner / Edenta



Allows you to measure the required spacing and - referring to the Set cover - to select the correct staple.

## INTRA ORAL PARALLELISOR

- Ref. 190530
- Essential for the correct placement of the staples in a splinting procedure.



## EXTRA SOFT S. STEEL WIRES NICHRO C4

- Ø 3/10 mm ..... Ref. PM119203
- Ø 4/10 mm ..... Ref. PM119204
- Ø 5/10 mm ..... Ref. PM119205
- Ø 2,5/10 mm ..... Ref. PM119325

Coils of 15 meters. Nickel Chrome 18.10 wire, particularly suitable to orthodontic appliances.

\* Intended for continuous use for not more than 30 days

## STAINLESS STEEL ROUND WIRES AND HALF ROUND WIRES



NICHROMINOX round wires are available in different strengths and diameters. They are sold in coils of 10 meters. Soft wires are packed in 10 coils of 1 meter.

- HARD ROUND WIRE** - Strength 160-170 kg
- Ø 4/10 to 9/10 ..... Ref. PM112204 to PM112209
- Ø 10 to 13/10 + 15/10 ..... Ref. PM112210 to PM112215

- MEDIUM ROUND WIRE** - Strength 140-150 kg
- Ø 5/10 to 7/10 + 9/10 ..... Ref. PM111205 to PM111209
- Ø 10/10 to 15/10 ..... Ref. PM111210 to PM111215

- SOFT ROUND WIRE** - Strength 60-80 kg
- Ø 3/10 ..... Ref. PM110103
- Ø 4/10 to 6/10 ..... Ref. PM110104 to PM110106

Wires composition : alloy NiCr 18.12

*Intended for continuous use for not more than 30 days*





NICHROMINOX half round wires are available in different strengths and diameters. They are sold in coils of 10 meters.

- HARD HALF ROUND WIRE** - Strength 160-170 kg
- Ø 13/10 to 15/10 ..... Ref. PM117213 to PM117215

- MEDIUM HALF ROUND WIRE** - Strength 140-150 kg
- Ø 13/10 to 15/10 ..... Ref. PM116213 to PM116215

## STRAIGHT B3 WIRES


Nickel Chrome 18.10  
Ø 6/10 to 16/10 + 20/10  
2 lengths : 15 cm and 30 cm




- Ref. PM119406 to PM119420 - length 15 cm
- Ref. PM119506 to PM119520 - length 30 cm

## TUBES FOR B3 WIRES

Nickel Chrome 18.10  
Ø 9/10 - 13/10 - 15/10  
2 lengths : 15 cm and 30 cm



- Ref. PM120109 to PM120115 - length 15 cm
- Ref. PM120209 to PM120215 - length 30 cm



## LINGUAL BAR

- Ref. PM119103
- Oval wire - Width 2,3 mm x thickness 1,4 mm. Strength 160 - 170 kg. Sold in 5 packs of 1 m.

*Intended for continuous use for not more than 30 days*



## STAINLESS STEEL PLATES

- Plate 2,5/10 mm ..... Ref. PM160105
- Plate 5/10 mm ..... Ref. PM160110
- Plate 6/10 mm ..... Ref. PM160112

Nickel Chrome 18.10. Size 150 x 85 mm.

# Surgery Accessories

**new** **DAUTREY ARCH BARS**



**Dautrey Arch Bars - Hard**  
• Ref. 010800  
Length 140 mm - set of 10



**Dautrey Arch Bars - Extra Soft**  
• Ref. 010700  
Length 150 mm - set of 10

## ARCH BARS

**Arch Bars with Hooks**  
• Ref. 010100  
Length 121 mm - set of 2



**Arch Bars with Gaps**  
• Ref. 010200  
Length 121 mm - set of 2



**Arch Bars Dr Duclos**  
• Ref. 010150  
(1 bar with hooks + 1 bar with gaps)  
Length 121 mm - set of 2




**Arch Bars with Hooks and Gaps**  
• Ref. 010300 (Dautrey type)  
Length 136 mm - set of 2



**Arch Bars with Double Hooks**  
• Ref. 010600 (Dr Champy)  
Length 121 mm - set of 2



**Arch Bars Flexible**  
(Dr Peri)  
• Ref. 010500  
Length 156,5 mm - set of 2



*Short term use : intended for continuous use for not more than 30 days*



*Intended for continuous use for not more than 30 days*

## EXTRA SOFT S. STEEL WIRE METALLIC COIL

Metallic coils with 5 meters wire. 3 diameters :

- Ø 3/10 mm - Ref. PM014003
- Ø 4/10 mm - Ref. PM014004
- Ø 5/10 mm - Ref. PM014005



## ELASTIC LIGATURES

Boxes of 100 pieces. 5 sizes available :

- n°0 to n°5 ..... Ref. 090800 to 090805
- Assorted ..... Ref. 090810



## STRAIGHT SHARP HEAD LIGATURE WIRE

To fix Arch Bars\*

Tubes of 28 straight wires. Length 135 mm

3 diameters available :

- Ø 3/10 mm - Ref. 012003
- Ø 4/10 mm - Ref. 012004
- Ø 5/10 mm - Ref. 012005

*\* Intended for continuous use for not more than 30 days*



## MOUTH GAGS

- Ref. 070000 - Adult standard size
- Ref. 070001 - Small size

Includes a slightly concave blunt blade fixed with a molt's to push the tongue aside.

# Surgery Accessories



**ORO PACKS** • Ref. 091500 - Box of 80 pc

Compress intended for oro-pharyngeal plugging during procedures under general anaesthetic (ENT, stomatology, maxillofacial surgery).

Oro Pack complete balloon protection. Placed at the aerodigestive junction, behind the tongue, it forms an essential barrier to liquids, blood, water, mucus.

The compress can easily be cut with scissors.

Autoclavable at 121°C.



## SUTURE SCISSORS

- Ref. 054000
- Length 14cm, lets you make all kinds of intra-oral sutures.



## SCISSORS FOR WIRE AND ARCH BARS

- Ref. 055000
- The serrated blade allows you to cut the wire very easily.



## SUTURE REMOVAL SCISSORS

- Ref. 055400
- The specially shaped cutting blade cuts the wire and nothing else.

# Silicone Refills

## PLUGS FOR BUR STANDS

• Ref. 190200  
Set of 8



Sold as sets of 8.  
Colours 1 2 3 4 5 7 8 9 10 11 12 13.  
Used for Plug In, Ultralight Bur Holders (p.10-11-12-13), Endo Star (p.19) and Implantology Kits (p.22-23).

## PLUGS FOR ENDODONTIC STANDS


• Ref. 182175  
Set of 8




Vendus par lots de 8.  
Colour 2 only.  
Used for Endo Ultralight and Endo Star (p.19).

## SILICONE TUBES


Set of 6. Central slot is available in 3 different sizes :



• Ref. 500600 - Set of 6 tubes  
For instruments Ø 2 to 3,5 mm  
Colours 1 2 3 4 5 7 8 9 11 13



• Ref. 500610 - Set of 6 tubes  
For instruments Ø 3,5 to 5 mm  
Colours 1 2 3 4 5 7 8 9 11 13



• Ref. 500620 - Set of 6 tubes  
For instruments Ø 5 to 6,5 mm  
Colours 1 2 3 4 5 7 8 9 11 13

For Implantology Kits (p.22-23), Multi Holders, Multi Clips and Multi Plug'In (p.13).

## CURSORS FOR STERIMETER


• Ref. 182370  
Set of 12



Colours 1 2 3 4 5 7 8 9 13  
For all endodontic stands (except Aluminium range) (p.17 à 19)

## INDICATOR STRIPS

• Ref. 190275  
Set of 12



Colours 1 2 3 4 5 7 8 9 13  
For Bur Flash (p.16), Multi Plugs (p.13), Bur Clips (p.14), Micro Cassettes and Divided Cassettes (p.8).

## FOR ERGO CLIPS / ERGO CASSETTES (page 7) and FLEXI CLIPS (off catalogue)

Colours : 1 2 3 4 5 7 8 9 13



Size and Cassette Ref.	Composition	Refill Ref.
18x8 cm 182703/182803	3 NOTCHES Ø 10 MM	182713
18x11 cm 182700/182753/182800	5 NOTCHES Ø 8,5 MM	182705
18x11 cm 182801	5 NOTCHES Ø 11 MM	182724
18x15 cm 182710/182750/182810	10 NOTCHES Ø 8,5 MM	182715
18x15 cm 182811	8 NOTCHES Ø 11 MM	182729
23x18 cm 182712/182752	10 NOTCHES Ø 8,5 MM + FLAT STRIP	182722
23x18 cm 182812	15 NOTCHES Ø 8,5 MM	182743
28x18 cm 182720/182751/182820	20 NOTCHES Ø 8,5 CM	182736
28x18 cm 182725	10 NOTCHES Ø 8,5 MM + WAVES	182721

## FOR FLAT CASSETTES (page 4 and 5) - Colours : 1 2 3 4 5 7 8 9 13



Size and Cassette Ref.	Composition	Refill Ref.
18x6 cm 182063/183063	5 NOTCHES Ø 8,5 MM	182073
18x9 cm 182055/183055	7 NOTCHES Ø 8,5 MM	182089
18x13 cm 182062/183062	10 NOTCHES Ø 8,5 MM	182072

# Silicone Refills

## FOR GALAXY CASSETTES (page 3), CLASSIC STYLE (page 5) AND ULTRALIGHT CASSETTES (off catalogue)


Colours : 2 3 4 5 7 13



**AUTO GRIP SILICONE :**  
doubles the cassette capacity! (see p.3)

Size and Cassette Ref.	Composition	Refill Ref.
18x7 cm 182069	4 NOTCHES Ø 11 MM	182078
18x7 cm 182039	5 NOTCHES Ø 8 MM	182052
18x7 cm 182079	2 HAND PIECES NOTCHES	182088
18x7 cm 182044	5 NOTCHES GALAXY	182045
18x7 cm 183044	10 «AUTO GRIP» NOTCHES	183045
18x10 cm 182054	6 NOTCHES Ø 11 MM	182084
18x10 cm 182025	7 NOTCHES GALAXY	182026
18x10 cm 183025	14 «AUTO GRIP» NOTCHES	183026
18x14 cm 182066	8 NOTCHES Ø 11 MM	182077
18x14 cm 182056/182036	10 NOTCHES Ø 8 MM	182083
18x14 cm 182046	10 NOTCHES GALAXY	182047
18x14 cm 183046	20 «AUTO GRIP» NOTCHES	183047
18x20 cm 182057	8 NOTCHES Ø 11 MM + CLEAR AREA	182082
18x20 cm 182048	10 NOTCHES GALAXY + CLEAR AREA	182028
28x18 cm 182068	16 NOTCHES Ø 11 MM	182071
28x18 cm 182058	20 NOTCHES Ø 8 CM	182081
28x18 cm 182049	20 NOTCHES GALAXY	182029
28x18 cm 182087	10 NOTCHES Ø 11 MM + CLEAR AREA	182074
28x18 cm 182037	12 NOTCHES Ø 8 MM + CLEAR AREA	182059
28x18 cm 182027	12 NOTCHES GALAXY + CLEAR AREA	182024

## FOR EASY CLIP CASSETTES (page 6) - Colours : 1 2 3 4 5 7 8 9 13



Size and Cassette Ref.	Composition	Refill Ref.
18x5 cm 182903	3 NOTCHES Ø 10 MM	182953
18x5 cm 182904	4 NOTCHES Ø 9 MM	182954
18x8 cm 182900	5 NOTCHES Ø 8,5 MM	182955
18x8 cm 182901	5 NOTCHES Ø 11 MM	182956
18x13 cm 182910	10 NOTCHES Ø 8,5 MM	182960
18x13 cm 182911	8 NOTCHES Ø 11 MM	182961
20x18 cm 182912	15 NOTCHES Ø 8,5 MM	182965
25x18 cm 182920	20 NOTCHES Ø 8,5 CM	182970

## FOR EASY TRAY CASSETTES (page 7) - Colours 1 2 3 4 5 7 8 9 13



Size and Cassette Ref.	Composition	Refill Ref.
18x5 cm 183703	3 NOTCHES Ø 10 MM	182953
18x5 cm 183704	4 NOTCHES Ø 9 MM	182954
18x7 cm 183700	5 NOTCHES Ø 8,5 MM	183701
18x10 cm 183705	5 NOTCHES Ø 8,5 MM + FLAT STRIP	183706
18x14 cm 183710	10 NOTCHES Ø 8,5 MM + FLAT STRIP	183711
21x18 cm 183715	15 NOTCHES Ø 8,5 CM + FLAT STRIP	183716
28x18 cm 183720	20 NOTCHES Ø 8,5 CM + FLAT STRIP	183721

## FOR EASY TRAY ORTHO CASSETTES (page 9) - Colours 1 2 3 4 5 7 8 9 13



Size and Cassette Ref.	Composition	Refill Ref.
18x9 cm 183728	WAVES	183729
18x14 cm 183730	2 NOTCHES Ø 8,5 MM + WAVES	183731
21x18 cm 183732	3 NOTCHES Ø 8,5 MM + WAVES	183733
28x18 cm 183734	3 NOTCHES Ø 8,5 MM + WAVES	183735

**THE EXPERTISE OF CRAFTSMEN**

**INDUSTRY OF EXCELLENCE**



Contact / Dealer Details :



18, rue des Frères Lumière - F-69720 SAINT-BONNET-DE-MURE - FRANCE  
Tel. +33 (0)4 78 74 04 15 - Fax +33 (0)4 78 01 61 80 - office@nichrominox.fr



PICTOGRAMMES / IZONS :	
Bain à Ultrasons Ultrasonic Bath	
Thermolaveur Thermoisinfecteur	
Autoclave Autoclave	