


Universal Prosthetic Kit



1 Universal Prosthetic Kit

- The Universal Prosthetic Kit is composed of various prosthetic drivers that are compatible to different brands and sizes of implant abutments.

2 Intended Use





- This Universal Prosthetic Kit is used when implanting a fixture in a place where a tooth is missing, or when connecting or detaching a cover screw or other upper structures.

Component

1. Driver


This tool is used to tighten an abutment or abutment screw.

1.2 Hex short, long / 1.25 Hex short, long / SCS short, long / Unigrip short, long driver

Name	Picture	Compatible implant or company	Code	Length(L)
1.2 Hex Driver (0.48")		Neobiotech Biomet 3i Dentsply-Friadent, Xive, Osstem, DIO, Megagen	HD1210S	10mm
			HD1215S	15mm
1.25 Hex Driver (0.5")		Biohorizons Zimmer, Astra Implant Direct Legacy Camlog, MIS, Dentis, Dentium	HD12510S	10mm
			HD12515S	15mm
SCS Driver		Straumann Implant Direct Replant	WSCSD10S	10mm
			WSCSD15S	15mm
Unigrip Driver		Nobel Biocare Implant Direct Swishplait Southern Implants Neoss	WUGD10S	10mm
			WUGD15S	15mm

2. Driver Holder Code : DH01

The Driver Holder connects to various prosthetic drivers to allow easy handling of the healing cover screws and abutment screws.



3. Torque Ratchet Code : TW60

The Torque Ratchet is used to tighten the hex driver and measure torque. By changing the direction of the tip, tightening and loosening can be controlled.



Neobiotech Co.,Ltd

10F, 36, Digitalro-27Gil, Guro-Gu, Seoul, Korea, 08381

Tel 82-2-582-2885, Fax 82-2-582-2883 www.neobiotech.com / Made in Korea

**Neo
Biotech**
Satisfaction to Dentists