

DISCOVER THE
BEST OF SWISS
TECHNOLOGY

BienAir⁺
Dental

iCHIROPRO



An implant device must first and foremost embody the very principle of reliability: An unwavering confidence in its ability to deliver the required torque and speed in the long term, and not an inch less. And this is where the iChiropro stands out: **it offers unrivalled precision in RPM and real-time Torque regulation achieved** by Bien-Air pioneering OBCS technology (Onboard Calibration Systems) which ensures **accuracy within 1% over the device's lifespan**.

Beyond its user-friendliness, you'll discover a comprehensive platform that completes your digital workflow. From an extensive range of data treatment **surpassing any other competitor** on the market to unrivalled advanced patient file management, to a vast library of fully illustrated **surgical and implant** procedures library from over **50 leading brands** constantly updated, iChiropro considerably simplifies your decision-making.

Thanks to its predefined motor settings, this unique device **allows you to concentrate solely on your procedure**, without having to worry about meticulous adjustments.


S W I S S  M A D E

+50 LEADING BRANDS

PRODUCT FEATURES

NEW

Implant loadability prediction algorithm based on cutting edge machine learning technology will be released in 2024. Stay tuned!



DIGITAL INTEGRATION

Implant barcode scanner, CoDiagnostix QR code import.

CUSTOMISABLE WORKFLOW

Implants brands, step sequences, fully personalised setting values.

OPERATION TRACKING

Downloadable operation reports with export functionality.

STEP-BY-STEP GUIDANCE

Over **900 OFFICIAL PROTOCOLS**, fully illustrated and easily accessible.

FREE ONLINE UPDATES

Regular updates of implant library and iOS versions.

USPs *The added benefits.*

- ✓ Preloaded surgical features, for a seamless digital workflow.
- ✓ iPad-controlled dental implantology motor.
- ✓ Direct or planned execution of your operation.
- ✓ End-to-end digital experience in surgery processes.
- ✓ Faster procedures using official protocols of **OVER 50** of the world's leading brands.
- ✓ Detailed operations files available and downloadable.
- ✓ Automatic patient data update after completion of implantation procedures.

iCHIROPRO *In a nutshell.*

- ✓ Cutting-edge technology to differentiate your practice.
- ✓ Reduced decision fatigue and accuracy concerns on the procedure day.
- ✓ Seamless integration into your digital workflow.
- ✓ Ultimate learning tool for beginners.



THIS ONE IS FOR YOU

The iChiropro is the perfect device for implantologists and oral surgeons who want to collect more data such as real-time insertion torque to determine primary stability, and secondary prognosis to plan the prosthesis. Moreover, this device is also made for professionals who want to trust the stability of RPM and torque at the bone site for long-term predictability. And above all, it is designed for anyone who wants to focus on their patient and not on adjusting the motor.

FREQUENTLY ASKED QUESTIONS

It's too complicated, I don't need all these features.

No worries, the «direct function» mode enables the insertion of an implant at the touch of a button. If you only need a «plug-in & operate» device, we invite you to discover the Chiropro.



I don't want to do regular SW updates.

Don't worry, no update is required to access all the available functions.

Do I need to calibrate my device?

You don't. The iCHIROPRO eliminates the need for external calibration thanks to its minimal torque transmission loss. In addition, our advanced software incorporates our OBCS technology compensation factor, which takes into account the gradual decrease in the efficiency of the contra-angle over time and ensures its accuracy with a margin of error of just 1% over the device's lifespan.

How do you handle the iPad obsolescence and its life cycle?

Every new iPad model is tested by Bien-Air engineers to ensure the widest compatibility with the iChiropro. To increase significantly your iPad lifespan once Apple is not supporting software updates, you can turn the device to airplane mode to continue using it comfortably.

How do you manage data protection?

The device and account are password-protected and comply with HIPPA security standards. The iCHIROPRO benefits from Apple's recognised excellence in data protection.

If I change my iPad model, will I lose all my data?

You can effortlessly manually export all operational data to the cloud. This operation is highly recommended before removing the application or changing iPad. By doing so, you can safely store all past operations in a protected location, such as the cloud, for future reference.

When a new app version is released, download it, and it will automatically replace the old version, guaranteeing that none of your data will be lost.

MORE INFORMATION

