



INNOVATIVE SOLUTIONS
FOR DENTAL AND MEDICAL
PROFESSIONALS

**RESTORATIVE
DENTISTRY**

SUMMARY



GENERAL DENTISTRY	AGENCIES	PAGE
Curing Light	CAO	2
Digital Vision	CAO	3
Composite & Impression Material	Kulzer	4
Systematic Tool for Restorative Dentistry	Prosmate	5
Customise Implant Abutment System	VP Innovato - Cervico	6

MONET[®]

laser curing light



Cures composites, luting cements, adhesives, and sealants in **1 Second**

Laser curing has major advantages over LED curing. With the Monet laser cure, you will never have to worry about incomplete bonds or soft cures. The Monet laser produces a deeper cure than LED or any other light sources during the same exposure time.

TALK TO US
ABOUT OUR
SPECIAL
TRADE-IN DEAL

Say goodbye to well known LED related problems like:

- LED technology was a significant advancement, but LED curing lights present their own set of challenges.
- LED lights create a light pattern that can be dispersed away from the curing target when the energy from the LED light waves disperse.
- LED is a dispersed light source emitting at random times and random directions.
- LED needs a complicated optical system to collimate and homogenize the light beam

Introducing a new idea on how we see patient treatment. Hubble combines magnification, illumination, & documentation into one piece of equipment designed for more effective patient care.



HUBBLE™

DIGITAL VISION HEAD-SET

Hubble allows the user to control magnification from 1 up to 20x. Hubble takes you beyond the fixed magnification and enhancing limitations that come with a pair of loupes, providing a feature of continuous zoom for fully adjustable magnification all while providing true-color LED illumination. With up to 20x magnification, Hubble eliminates the need for costly and obtrusive microscopes and replaces both loupes and headlights.



Due to a large number of requests for the Hubble Digital Vision Head-Set, we have setup a Pre-Order Deposit to save your place in line. The deposit amount of \$550 is fully refundable, will be applied towards the retail price of \$17,050.00 and will save your place in line.



KULZER
MITSUI CHEMICALS GROUP



FLEXITIME

Time flexibility paired with the highest precision. Every time.

- Intelligent time concept: its special formula controls the polymerisation kinetics by oral temperature. The setting (polymerisation) process accelerates when the tray is placed into the patient's mouth.
- Automatic mixing with Dynamix speed 2 system: homogeneously and void-free mixed impression materials are the basis for precisely fitting prosthetic restorations.



VENUS® PEARL PURE SHADES

Do you know that only 4 shades can cover up to 90% of anterior & posterior restorations? Introducing Venus Pearl Pure Shades.



Excellent results in monochromatic layering

A perfect balance of translucency, hue, and chroma reflects light in a way that our 4 shades can be used for both monochromatically placed anterior and posterior restorations.



Before



After


PROSMATE



ProsMate System

Invented by a Prosthodontist, the ProsMate is an effective tool to aid in the efficient and systematic insertion of fixed prostheses.

The Baton allows the simultaneous etching, washing, drying and silanating of porcelain veneers.



With the ProsMate you will be able arrange prostheses in order of insertion, and in addition, you can handle silanted porcelain veneers without touching the bonding surfaces..

Milled out of metal, the ProsMate is autoclavable and can last a lifetime. The ProsMate is a simple yet brilliant system to add to your armamentarium for daily practice



VP INNOVATO HOLDINGS

MASTERING THE EMERGENCE PROFILE



CERVICO PREMIUM KIT

The ultimate in-house custom healing abutment fabrication and ideal soft tissue profile generator set of tools. The Cervico Mold is an in-office tool used for the fabrication of the cervico healing abutments and impression posts. It can also be used for the duplication of a modified healing abutment or temporary prosthesis to a duplicate impression post. It features 17 wells that can be used for all of the above purposes.

What's included: • Base screw driver • VPI prosthetic connection inserts* • Centered Silicone Insert • Mold Base

CERVICO ESSENTIAL KIT

Essential Cervico Mold is a simplified version of the Cervico Premium Mold that works with the stock implant analogs available from the different implant companies.





2026
EDITION



 1300 792 624

 www.biomedent.com.au